



# Flavors and aromas

## Analysis of Essence de Térébenthine

### Application Note

Materials Testing & Research

#### Authors

Agilent Technologies, Inc.

#### Introduction

Gas chromatography with an Agilent CP-Wax 52 CB column identifies 14 components in turpentine in 20 minutes.

The most polar, bonded Carbowax phase, is the only column of this kind that separates  $\beta$ -phellandrene from 1,8-cineol and  $\alpha$ -terpinene from 1,4-cineol.




**Agilent Technologies**

## Conditions

Technique : GC-capillary  
Column : Agilent CP-Wax 57 CB, 0.32 mm x 50 m fused silica  
WCOT (df = 0.2  $\mu$ m) (Part no. CP97753)  
Temperature : 70 °C (15 min)  $\rightarrow$  210 °C, 30 °C/min  
Carrier Gas : N<sub>2</sub>, 90 kPa (0.9 bar, 13 psi)  
Injector : Split, 150 mL/min  
T = 210 °C  
Detector : FID  
T = 240 °C  
Sample Size : 0.03  $\mu$ L  
Concentration Range : undiluted sample

## Peak identification

1. tricyclene
2.  $\alpha$ -pinene
3.  $\alpha$ -fenchene
4. camphene
5.  $\beta$ -pinene
6. unknown
7.  $\beta$ -myrcene
8.  $\alpha$ -phellandrene
9.  $\alpha$ -terpinene
10. limonene (dipentene)
11.  $\beta$ -phellandrene
12.  $\gamma$ -terpinene
13. p-isopropyltoluene (p-cymene)
14. terpinolene
15. sesquiterpenes



[www.agilent.com/chem](http://www.agilent.com/chem)

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01460



**Agilent Technologies**