



Allergens

Application Note

BioPharma

Authors

Agilent Technologies, Inc.

Introduction


Analyze allergens in 400 seconds by GC using a 0.1 mm id Agilent FactorFour VF-5ms column.



Agilent Technologies

Conditions

Technique : GC
Column : Agilent FactorFour VF-5ms, 0.10 mm x 10 m fused silica (df = 0.4 μ m) (Part No. CP8934)
Temperature : 100 °C \rightarrow 320 °C, 40 °C/min
Carrier Gas : Helium, 70 kPa
Injector : Splitless
Detector : MS
Sample Size : 1.0 μ L
Concentration Range : ca. 100 ppm



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A02220



Agilent Technologies