

GC Tips and Tricks to Speed Up Your Analysis and Increase Your Throughput

Shannon Coleman
Application Scientist – Gas Chromatography

DE.3927662037



GC Method Translation Software

- decrease analysis time
- optimize GC parameters
- maintain peak patterns and resolution

GC Method Translation - EPACL.P.MXD

Criterion: Translate Only Best Efficiency Fast Analysis None **Speed gain: 1.17407**

	Original Method	Translated Method																																										
Column																																												
Length, m	30	<input checked="" type="checkbox"/> 30																																										
Internal Diameter, μm	320	<input checked="" type="checkbox"/> 320																																										
Film																																												
Thickness, μm	0.25	<input checked="" type="radio"/> Unlock <input checked="" type="radio"/> 0.25																																										
Phase Ratio	320.0	<input type="radio"/> 320.0																																										
Carrier Gas	Helium	<input type="checkbox"/> Helium																																										
Enter one Setpoint																																												
Head Pressure, psi	12.786	15.126																																										
Flow Rate, mL/min	2.0502	2.5600																																										
Outlet Velocity, cm/sec	56.20	70.17																																										
Average Velocity, cm/sec	38	44.61																																										
Hold-up Time, min	1.31579	1.12070																																										
Outlet Pressure (absolute), psi	14.696	<input type="checkbox"/> 14.696																																										
Ambient Pressure (absolute), psi	14.696	<input type="checkbox"/> 14.696																																										
Oven Temperature 3-ramp Program																																												
	<table border="1"> <thead> <tr> <th>Ramp Rate</th> <th>Final Temp.</th> <th>Final Time</th> </tr> <tr> <th>$^{\circ}\text{C}/\text{min}$</th> <th>$^{\circ}\text{C}$</th> <th>min</th> </tr> </thead> <tbody> <tr> <td>Initial</td> <td>120</td> <td>1.17</td> </tr> <tr> <td>Ramp 1</td> <td>25</td> <td>160</td> <td>0</td> </tr> <tr> <td>Ramp 2</td> <td>10</td> <td>260</td> <td>0</td> </tr> <tr> <td>Ramp 3</td> <td>15</td> <td>300</td> <td>4</td> </tr> </tbody> </table>	Ramp Rate	Final Temp.	Final Time	$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min	Initial	120	1.17	Ramp 1	25	160	0	Ramp 2	10	260	0	Ramp 3	15	300	4	<table border="1"> <thead> <tr> <th>Ramp Rate</th> <th>Final Temp.</th> <th>Final Time</th> </tr> <tr> <th>$^{\circ}\text{C}/\text{min}$</th> <th>$^{\circ}\text{C}$</th> <th>min</th> </tr> </thead> <tbody> <tr> <td>Initial</td> <td>120</td> <td>0.997</td> </tr> <tr> <td>Ramp 1</td> <td>29.352</td> <td>160</td> <td>0.000</td> </tr> <tr> <td>Ramp 2</td> <td>11.741</td> <td>260</td> <td>0.000</td> </tr> <tr> <td>Ramp 3</td> <td>17.611</td> <td>300</td> <td>3.407</td> </tr> </tbody> </table>	Ramp Rate	Final Temp.	Final Time	$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min	Initial	120	0.997	Ramp 1	29.352	160	0.000	Ramp 2	11.741	260	0.000	Ramp 3	17.611	300	3.407
Ramp Rate	Final Temp.	Final Time																																										
$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min																																										
Initial	120	1.17																																										
Ramp 1	25	160	0																																									
Ramp 2	10	260	0																																									
Ramp 3	15	300	4																																									
Ramp Rate	Final Temp.	Final Time																																										
$^{\circ}\text{C}/\text{min}$	$^{\circ}\text{C}$	min																																										
Initial	120	0.997																																										
Ramp 1	29.352	160	0.000																																									
Ramp 2	11.741	260	0.000																																									
Ramp 3	17.611	300	3.407																																									
Sample Information None																																												

Speed gain: 1.0000


Translate
 Best Efficiency

Last file imported:

	Original Method	Calculated Method																								
	Gas He	Gas He																								
Length (m)	80 m	80 m																								
Inner Diameter (μm)	250 μm	250 μm																								
Film Thickness (μm)	0.25 μm	0.25 μm																								
Phase Ratio	249.25	249.25																								
Inlet Pressure (gauge)	1.1739 psi	1.1739 psi																								
Outlet Flow (mL/min)	0.45587 mL/min	0.45587 mL/min																								
Average Velocity (cm/s)	24.844 cm/sec	24.844 cm/sec																								
Outlet Pressure (abs)	0 psi	0 psi																								
Holdup Time	2.0125 min	2.0125 min																								
Outlet Velocity (cm/s)	cm/sec	cm/sec																								
Isothermal																										
Ramps	<table border="1"> <thead> <tr> <th>#</th> <th>Ramp Rate ($^{\circ}\text{C}/\text{min}$)</th> <th>Final Temp ($^{\circ}\text{C}$)</th> <th>Final Time (min)</th> </tr> </thead> <tbody> <tr> <td>Init</td> <td></td> <td>70</td> <td>1.5</td> </tr> <tr> <td>1</td> <td>16</td> <td>200</td> <td>0.5</td> </tr> </tbody> </table>	#	Ramp Rate ($^{\circ}\text{C}/\text{min}$)	Final Temp ($^{\circ}\text{C}$)	Final Time (min)	Init		70	1.5	1	16	200	0.5	<table border="1"> <thead> <tr> <th>#</th> <th>Ramp Rate ($^{\circ}\text{C}/\text{min}$)</th> <th>Final Temp ($^{\circ}\text{C}$)</th> <th>Final Time (min)</th> </tr> </thead> <tbody> <tr> <td>Init</td> <td></td> <td>70</td> <td>1.5</td> </tr> <tr> <td>1</td> <td>16</td> <td>200</td> <td>0.5</td> </tr> </tbody> </table>	#	Ramp Rate ($^{\circ}\text{C}/\text{min}$)	Final Temp ($^{\circ}\text{C}$)	Final Time (min)	Init		70	1.5	1	16	200	0.5
#	Ramp Rate ($^{\circ}\text{C}/\text{min}$)	Final Temp ($^{\circ}\text{C}$)	Final Time (min)																							
Init		70	1.5																							
1	16	200	0.5																							
#	Ramp Rate ($^{\circ}\text{C}/\text{min}$)	Final Temp ($^{\circ}\text{C}$)	Final Time (min)																							
Init		70	1.5																							
1	16	200	0.5																							
Total Run Time	10.13 min	10.13 min																								

Method Translation Software (Not New)

Precise Time-Scaling of Gas Chromatographic Methods Using Method Translation and Retention Time Locking



Application

Gas Chromatography

May 1998

Authors

B. D. Quimby, L. M. Blumberg,
M. S. Klee, and P. L. Wylie
Agilent Technologies, Inc.
2850 Centerville Road
Wilmington, DE 19808-1610
USA

Abstract

Key Words

Pesticides, GC, GC-AED, retention
time locking, RTL, method transla-
tion, scalable RT libraries

Introduction

Interest in the analysis of pesticide

the observed element content of the
peak. The combination of time and
element content narrows rapidly the
possible compounds that could have
produced the heteroatom response to
a few pesticides.

The element-selective detection
is done with either gas

Acknowledgements & References


1. M. Klee and V. Giarrocco, "*Predictable Translation of Capillary GC Methods for Fast GC*," Publication 5965-7673E, March **1997**.
2. V. Giarrocco, B. D. Quimby, and M. S. Klee, "*Retention Time Locking: Concepts and Applications*," Publication 5966-2469E, December **1997**.
3. P. L. Wylie and B. D. Quimby, "*A Method Used to Screen for 567 Pesticides and Suspected Endocrine Disrupters*," Publication 5967-5860E, April **1998**.
4. L. Wool and D. Decker, "*Practical Fast Gas Chromatography for Contract Laboratory Program Pesticides*, *J. Chromat. Sci. Vol. 40*, September **2002**.
5. F. Bothe, K. Dettmer, and W. Enewald, "*Determination of Perfume Oil in Household Products by Headspace SPME and Fast Capillary Gas Chromatography*", *57*, **2003**.
6. M. Sinnott, S. Jones, "*Rapid Analysis of Food and fragrances Using High Efficiency Capillary GC Columns*", Publication 5989-7509EN, November **2007**.
7. "*Method Translation of HJ679-2013 for Intuvo*", August **2017**.

Can I speed up my runtime?

- What data do I really need for my analysis?
- Are there large amounts of baseline between peaks?
- Do I need to resolve all the peaks in my chromatogram?
- What carrier gas should I use?
- How much sample do I need to place on the column to detect my components of interest?

Variables for Shortening Runtime

- Stationary Phase
- Temperature Programming
- Carrier Gas: type and linear velocity
- Shorten Column Length
- Decrease Film Thickness
- Decrease Internal Diameter



Resolution

$$R_s = \frac{\sqrt{N}}{4} \left(\frac{k}{k+1} \right) \left(\frac{\alpha-1}{\alpha} \right)$$

Efficiency	$N = f$ (gas, L, r_c)	L = Length
Retention	$k = f$ (T, d_f , r_c)	r_c = column radius
Selectivity	$\alpha = f$ (T, phase)	d_f = film thickness
		T = temperature

Selectivity

Represents the separation power of a particular adsorbent to separate a mixture of components.

Stationary Phase - Common Types


Siloxane polymers

Poly(ethylene) glycols


Porous polymers

Start with the Right Phase

DB-1
15m x 0.32mm, 0.25 μ m
Oven:
40°C for 2 min
40-120°C at 5°C/min




DB-Wax
15m, 0.32mm, 0.25 μ m
Oven:
80-190°C at 20°C/min



Stationary Phase Selection

Utilize existing Information, Determine critical separations
Selectivity/Polarity, Temperature Limits, SP-Application specific



The screenshot shows the top navigation bar of the Agilent website with links for Home, News, Technology, Applications, Links, and Browse. Below the navigation bar is a breadcrumb trail: All Places > Technical Areas > Consumables > Documents. The main content area features a document titled "Software-Supported Method Development - the Scanview Program" with a "Version 13" badge. The document is attributed to don_gage and noreuter, with creation and modification dates. The title "The New Scanview Application Database, Version 2019" is prominently displayed, followed by the author information: Norbert Reuter*, Ingrid van der Meer, Agilent Technologies, Middelburg, The Netherlands. An "Introduction" section begins with a paragraph describing the new Scanview version 2019 as an application database containing full method descriptions and literature references for various analytical techniques.

<https://community.agilent.com/docs/DOC-2118-software-supported-method-development-the-scanview-program>

Choose the column phase that gives the best separation but not at the cost of robustness or ruggedness.

Helium Carrier Gas Alternatives

Important theoretical considerations relating to peak efficiency


Sharp, narrow peaks in a chromatogram is an indication of a high efficiency GC column.

- Remember that efficiency is represented mathematically by the symbol " N " called *Theoretical Plates*, and that the larger N is, the better the resolving power of the column (i.e., higher resolution).
- Resolution is described mathematically by the symbol R_s and its numeric value tells how well two adjacent peaks are separated from each other.

$$R_s = \frac{\sqrt{N}}{4} \left(\frac{k}{k+1} \right) \left(\frac{\alpha-1}{\alpha} \right)$$

A resolution value of 1.5 tells us that two peaks are baseline separated. The greater (higher) the R_s value, the more separation that has been achieved.

Calculating efficiency



We would like to know the actual time the component spends in the stationary phase.

$$t'_R = t_R - t_m \quad n = \frac{t'_R}{5.545 \left(\frac{t'_R}{W_h} \right)^2}$$

t'_R = corrected retention time.

n = effective theoretical plates.


Let's relate " n " to the length of the column.

$$\text{Plates per meter (N)} = \frac{n}{L} \quad \text{or}$$

$$\text{Height equivalent to a theoretical plate (HETP)} = \frac{L}{n}$$

Thus, the more efficient the column, the bigger the " N " the smaller the "HETP".

Efficiency and carrier gas linear velocity



Efficiency is a function of the carrier gas linear velocity or flow rate.

The minimum of the curve represents the smallest HETP (or largest plates per meter) and thus the best efficiency. "A" term is not present for capillary columns.

- Plot of HETP versus linear velocity is known as the **Van Deemter** plot.
- The linear velocity value at the minimum of the curve is the optimum value for achieving the best efficiency.

Helium Carrier Gas Alternatives

Let's make this easy


Sharp, narrow peaks in a chromatogram is an indication of a high efficiency GC column.

- Remember that **efficiency** is represented mathematically by the symbol "*N*" called *Theoretical Plates*, and that the larger *N* is, the better the resolving power of the column (i.e., higher resolution).
- Resolution** is described mathematically by the symbol R_s and its numeric value tells how well two adjacent peaks are separated from each other.

$$R_s = \frac{\sqrt{N}}{4} \left(\frac{k}{k+1} \right) \left(\frac{\alpha-1}{\alpha} \right)$$

A resolution value of 1.5 tells us that two peaks are baseline separated. The greater (higher) the R_s value, the more separation that has been achieved.

Calculating efficiency



Let's relate "n" to the length of the column.

Plates per meter (N) = $\frac{n}{L}$ or

Height equivalent to a theoretical plate (HETP) = $\frac{L}{n}$

Thus, the more efficient the column, the bigger the "N" the smaller the "HETP".


We would like to know the actual time the component spends in the stationary phase.

$$t'_R = t_R - t_m \quad n = \left(\frac{t'_R}{W_{0.5}} \right)^2$$

t'_R = corrected retention time.

n = effective theoretical plates.

Efficiency and carrier gas linear velocity



Efficiency is a function of the carrier gas linear velocity or flow rate.

The minimum of the curve represents the smallest HETP (or largest plates per meter) and thus the best efficiency. "A" term is not present for capillary columns.


- Plot of HETP versus linear velocity is known as the Van Deemter plot.
- The linear velocity value at the minimum of the curve is the optimum value for achieving the best efficiency.

Where does Helium come from?

Where does Helium come from?

Komperda


The **helium** is **formed** during the natural radioactive decay of elements such as uranium and thorium. These heavy elements were **formed** before the earth but they are not stable and very slowly, they decay. ... This alpha-particle is actually just the heart of a **helium** atom - its nucleus



https://s3.amazonaws.com/ai2-vision-textbook-dataset/dataset_releases/rc2/train/question_images/radioactive_decay_8172.png

Carrier Gas Selection



VAN DEEMTER CURVES



H = height equivalent of theoretical plates (goal is shorter height, to achieve more plates/meter)
u = velocity

Van Deemter

Why nitrogen gets a bad reputation for capillary GC





- N_2 actually provides the best efficiency, but at a slower speed
- Most helium methods have too much resolution
 - Lower N_2 efficiency at higher flows can still provide “good enough” resolution
- Most GC methods now use constant flow
 - N_2 efficiency losses with temperature programming are not as severe

Many Helium GC Have Excess Resolution


EN14103 – GC analysis of FAME content in biodiesel

HP-INNOWax, 30 m x 0.25 mm id x 0.25 μ m




Carrier Gas Selection

VAN DEEMTER CURVES



H = height equivalent of theoretical plates (goal is shorter height, to achieve more plates/meter)
ū = velocity

Same Holdup Time (T_r) Gives Consistent Retention Times Compared to Original Helium Method



Column Diameter - Theoretical Efficiency

	I.D. (mm)	n/m
	0.05	23,160
	0.10	11,980
	0.18	6,660
	0.20	5830
	0.25	4630
	0.32	3760
Retention Factor $k = 5$	0.45	2840
	0.53	2060

PHASE RATIO (β)

Column Dimensions

30 m x .53 mm x 3.0 μm

30 m x .32 mm x 1.8 μm

Phase Ratio β

44

44

$$\beta = \frac{r}{2d_f}$$

$$K_c = k\beta \quad (\text{same phase ratio gives same retention})$$

Method Translator

Method Translator interface showing parameters for Original Method and Calculated Method.

Speed gain: 0.7820

Translate:


Best Efficiency:

Last file imported: [] [] [] []


Parameter	Original Method	Calculated Method
Length (m)	30 m	30 m
Inner Diameter (μm)	530 μm	320 μm
Film Thickness (μm)	3.00 μm	1.80 μm
Phase Ratio	43.418	43.696

Different Column I. D.

Equal Phase Ratios



Column Diameter and Capacity



	I.D. (mm)	Capacity (ng)
	0.05	1-2
	0.10	6-13
	0.18	25-55
	0.20	35-70
	0.25	80-160
	0.32	110-220
	0.45	600-800
	0.53	1000-2000

Like Polarity
Phase/Solute
0.25 μm film thickness

Column Diameter and Carrier Gas Flow

Lower flow rates: Smaller diameter columns

Low flow rates : GC/MS
High flow rates: Headspace, purge & trap

Higher flow rates: Larger diameter columns


Diameter Summary

To increase	Diameter
Efficiency	Smaller
Resolution	Smaller
Pressure	Smaller
Capacity	Larger
Flow rate	Larger

Film Thickness and Capacity

Thickness (μm)	Capacity (ng)
0.10	50-100
0.25	125-250
0.50	250-300
1	500-1000
3	1500-3000
5	2500-5000


More stationary phase = More degradation products



0.32 mm I.D.
Like Polarity Phase/Solute

What column types/dimensions produce higher bleed?

- Polarity: More polar = higher bleed
- Low polarity = More thermally stable
 - Look at temperature limits as a general indicator of thermal stability
- The more total mass of polymer in the column the higher the bleed (within a given phase)
 - Larger diameters
 - Longer columns
 - Thicker films



Column Length and Efficiency (Theoretical Plates)

	Length (m)	N
	15	69,450
0.25 mm ID	30	138,900
n/m = 4630 (for k = 5)	60	277,800

More Meters = More Plates = More Resolution

Column Length and Resolution


$$R \propto \sqrt{N} \propto \sqrt{L}$$

Length X 4 = Resolution X 2

$$t \propto L$$

Upside = Cut a bunch off during routine inlet maintenance and not lose a lot of Resolution


Column Length VS Resolution and Retention: Isothermal




Double the plates, double the time but not double the resolution

DECREASE THE LENGTH

DB-5
30 m
0.53 mm I.D., 0.5 μ m




DB-5
15 m
0.53 mm I.D., 0.5 μ m




BTEX
Carrier: Helium, 36 cm/sec at 40°C
Oven : 40°C for 3 min, 5°/min to 100°C

Column Length and Cost



Length Summary

<u>To Increase</u>	<u>Length</u>
Efficiency	Longer
Resolution	Longer
Analysis Time	Longer
Pressure	Longer
Cost	Longer



Column Temperature

Developing Temperature Programs

Most powerful variable

Changes Selectivity and Retention

Natural log (ln) relationship between retention and temperature

Most difficult to predict and develop

Often involves trial and error

Changes in Column Dimensions, Gas Type or Velocity Require Changes in Temp Program Rates

Method Translator
38 min method to
19 minute method

Agilent Technologies Method Translator

Last file imported:

Speed gain
2.0000

Translate

Best Efficiency

Original Method Parameters

Gas: He

Length (m): 30 m

Inner Diameter (µm): 250 µm

Film Thickness (µm): 0.25 µm

Phase Ratio: 249.25

Inlet Pressure (gauge): 9.2554 psi

Outlet Flow (mL/min): 1.2109 mL/min

Average Velocity (cm/s): 39.903 cm/sec

Outlet Pressure (abs): 0 psi

Holdup Time: 1.253 min

Outlet Velocity (cm/s): Infinity

Calculated Method Parameters

Gas: He

Length (m): 20 m

Inner Diameter (µm): 180 µm

Film Thickness (µm): 0.18 µm

Phase Ratio: 249.25

Inlet Pressure (gauge): 26.373 psi

Outlet Flow (mL/min): 1.4351 mL/min

Average Velocity (cm/s): 53.204 cm/sec

Outlet Pressure (abs): 0 psi

Holdup Time: 0.62652 min

Outlet Velocity (cm/s): Infinity

Isothermal

Ramps

2

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		40	2.00
1	5.0000	120	0.00
2	12.0000	300	5.00

Total Run Time: 38.00 min

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		40	1.00
1	10.0000	120	0.00
2	24.0000	300	2.50

Total Run Time: 19.00 min


Pressure Units: PSI

Original Column Capacity: 1.71


Translated Column Capacity: 0.61

The column capacity of the translated method is 36% of the original column capacity. You may need to adjust your injection volume.

Conventional ASTM Mix E1618 for Fire Debris – 38mins.



Can you tell a difference?



Method Translator QC Mix

Speed Gain

Agilent Technologies Method Translator

Last file imported:

Speed gain
2.0000
 Translate
 Best Efficiency

Original Method Parameters		Calculated Method Parameters	
Gas: He		Gas: He	
Length (m)	30 m	20 m	
Inner Diameter (μm)	250 μm	180 μm	
Film Thickness (μm)	0.25 μm	0.18 μm	
Phase Ratio	249.25	249.25	
Inlet Pressure (gauge)	10.523 psi	28.546 psi	
Outlet Flow (mL/min)	1 mL/min	1.1852 mL/min	
Average Velocity (cm/s)	37.293 cm/sec	49.724 cm/sec	
Outlet Pressure (abs)	0 psi	0 psi	
Holdup Time	1.3407 min	0.67036 min	
Outlet Velocity (cm/s)	Infinity	Infinity	

Method Translator QC Mix

Speed Gain

Isothermal
 Ramps

1

#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		100	2.00
1	20.0000	300	18.00

Total Run Time 30.00 min

Pressure Units
PSI

Original Column Capacity: 1.71


#	Ramp Rate (°C/min)	Final Temp (°C)	Final Time (min)
Init		100	1.00
1	40.0000	300	9.00

Total Run Time 15.00 min


Translated Column Capacity: 0.61

The column capacity of the translated method is 36% of the original column capacity. You may need to adjust your injection volume.

QC Mix 16 minute



QC Mix 8 minute



In Summary

Diameter – Smaller allows shorter length but has less capacity

Small Change in ID Easier to Translate – Again think capacity

Length – Shorter might be possible without losing a lot of Resolution

Temperature Program – Use Method Translation Software to scale temperatures properly

Method Translation Software – FREE and it works

Hydrogen Carrier – Higher velocities for even faster analyses

Detector Selection – Have realistic expectations for what can be achieved

<https://www.agilent.com/en/support/gas-chromatography/gcmethodtranslation>