

## Petroleum/Chemicals Application Guide

*This guide contains up-to-date information about Supelco products and technology for separating hydrocarbons by chromatographic methods.*

### Key Words

- petroleum
- PIANO
- ASTM
- chemicals
- PONA
- Petrocol
- hydrocarbons
- SIMDIS

### Applications

<b>Hydrocarbons Analyses</b>	<b>1-10</b>
<b>Oxygenates</b>	<b>10</b>
<b>Simulated Distillation</b>	<b>11-14</b>
<b>Aliphatics / Aromatics</b>	<b>14-15</b>
<b>Polynuclear Aromatic Hydrocarbons</b>	<b>15</b>
<b>Trace Petroleum Components</b>	<b>16</b>
<b>Miscellaneous Petrochemicals Analyses</b>	<b>17-19</b>
<b>Phenols / Polymers</b>	<b>19</b>
<b>Solvents</b>	<b>20-21</b>
<b>Impurities in Solvents</b>	<b>21-22</b>
<b>High Speed / High Resolution Separations</b>	<b>22-23</b>
<b>Packed Column Applications</b>	<b>24-25</b>
<b>Gel Permeation of Polymers</b>	<b>26</b>
<b>Low Molecular Weight Compounds</b>	<b>26-27</b>

### Ordering Information

<b>Columns and Packings</b>	<b>28-32</b>
<b>Chemical Standards and Accessories</b>	<b>33-48</b>


## Detailed Hydrocarbons Analyses

### Nonpolar 50-Meter Capillary Column for High Resolution of Hydrocarbons (ASTM D5134)

Information about the paraffin, isoparaffin, olefin, aromatic, and naphthene (PIANO) content of many products is critical for process control and product quality testing. We developed the 50m x 0.20mm ID Petrocol™ DH 50.2 column specifically for this type of application (Figures 1-5). This high efficiency, bonded, nonpolar (methyl silicone) phase column typically offers 5000 plates/meter.


**Figure 1. Simple Hydrocarbon Mixtures**

Column: **Petrocol DH 50.2, 50m x 0.20mm ID, 0.5µm film**  
 Cat. No.: **24133-U**  
 Oven: 35°C (30 min) to 200°C at 2°C/min  
 Carrier: helium, 20cm/sec (set at 35°C)  
 Det.: FID, 225°C  
 Inj.: 225°C




**Figure 2. Unleaded Gasoline**


Column: **Petrocol DH 50.2, 50m x 0.20mm ID, 0.5µm film**  
 Cat. No.: **24133-U**  
 Oven: 35°C (30 min) to 200°C at 2°C/min  
 Carrier: helium, 20cm/sec (set at 35°C)  
 Det.: FID, 225°C  
 Inj.: 225°C




**Figure 3. Detailed Analysis of Qualitative Reference & Naptha**



**Figure 4. Quantitative PIANO Standard**



**Figure 5. Quantitative PIANO Standard, Using Rapid Temperature Program**




**Gasoline Range Petroleum Products**

Petrocol DH columns (Figures 6-8) can separate more than 300 gasoline components in less than 110 minutes — under ambient initial temperature conditions. They are ideal for PONA, PNA, and PIANO analyses.■

Nonpolar Petrocol DH columns provide greater separation of propane and C4 compounds, m- and p-xylene (Figure 6), and other light petroleum components. Because you can use Petrocol DH columns with modest split ratios (approximately 100:1), you can simultaneously monitor trace and concentrated sample components.

Our library of retention indices for more than 400 compounds is included with the column.

**Figure 6. 2,3-Dimethylbutane from Cyclopentane**




■ PONA: paraffins-olefins-naphthenes-aromatics  
 PNA: paraffins-naphthenes-aromatics  
 PIANO: paraffins-isoparaffins-aromatics-naphthenes-olefins

## Figure 7. Premium Unleaded Gasoline

Column: **Petrocol DH, 100m x 0.25mm ID, 0.50µm film**  
 Cat. No.: **24160-U**  
 Oven: 35°C (15 min) to 200°C at 2°C/min, hold 5 min  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.1µL gasoline, split (100:1)


- |                            |                             |
|----------------------------|-----------------------------|
| 1. Isobutane               | 17. 2,3,4-Trimethylpentane  |
| 2. n-Butane                | 18. Toluene                 |
| 3. Isopentane              | 19. 2,3-Dimethylhexane      |
| 4. n-Pentane               | 20. Ethylbenzene            |
| 5. 2,3-Dimethylbutane      | 21. m-Xylene                |
| 6. 2-Methylpentane         | 22. p-Xylene                |
| 7. 3-Methylpentane         | 23. o-Xylene                |
| 8. n-Hexane                | 24. 1-Methyl-3-ethylbenzene |
| 9. 2,4-Dimethylpentane     | 25. 1,3,5-Trimethylbenzene  |
| 10. Benzene                | 26. 1,2,4-Trimethylbenzene  |
| 11. 2-Methylhexane         | 27. 1,2,3-Trimethylbenzene  |
| 12. 3-Methylhexane         | 28. Naphthalene             |
| 13. 2,2,4-Trimethylpentane | 29. 2-Methylnaphthalene     |
| 14. n-Heptane              | 30. 1-Methylnaphthalene     |
| 15. 2,5-Dimethylhexane     | 31. Dimethylnaphthalenes    |
| 16. 2,4-Dimethylhexane     |                             |



713-1274, 0744


**Figure 8. A Difficult Separation: Isobutylene from Butene-1**

Column: **Petrocol DH, 100m x 0.25mm ID, 0.50µm film**  
 Cat. No.: **24160-U**  
 Oven: -20°C (10 min) to 200°C at 2°C/min  
 Carrier: helium, 19cm/sec  
 Det.: FID  
 Inj.: 25µL Phillips L6 blend, split (100:1)



**Figure 9. Light Hydrocarbons, Using Moderately Subambient Initial Temperature**

Column: **Petrocol DH 150, 150m x 0.25mm ID, 1.0µm film**  
 Cat. No.: **24155**  
 Oven: -20°C (30 min) to 75°C at 5°C/min  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 250µL nitrogen (0.1-1% each analyte), split (100:1)




**Figure 10. Gaseous Hydrocarbons from Liquid Hydrocarbons**

Column: **Petrocol DH 150, 150m x 0.25mm ID, 1.0µm film**  
 Cat. No.: **24155**  
 Oven: 35°C  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 25µL natural or liquified petroleum gas, split (100:1)


Natural Gas Standard

1. Methane
2. Ethane
3. Propane
4. Isobutane
5. n-Butane
6. Isopentane
7. n-Pentane



Liquified Petroleum Gas Standard

1. Ethane
2. Propylene
3. Propane
4. Isobutane
5. n-Butane
6. Isopentane



### Ambient Separations of Hydrocarbon Gases

Very high efficiency (typically 550,000 theoretical plates) and a low phase ratio (0.25mm ID, 1.0µm phase film) enable our 150-meter Petrocol DH 150 column to separate light hydrocarbons significantly better than other columns.

Use Petrocol DH columns to analyze hydrocarbon mixtures ranging from the light gases to the higher boiling components of gasoline and other samples of comparable boiling range (Figures 9-11). Perform most of these analyses using ambient or mildly subambient (0°C or -20°C) initial temperatures.

Perform separations that previously required special conditions, unusual columns, or prolonged analysis time.


Head pressure of approximately 75psig will provide the optimum linear velocity, 20cm/sec, for helium carrier gas. If your instrument cannot provide this pressure, consider using hydrogen as the carrier gas (40-50psig will provide the optimum 40cm/sec linear velocity), or consult your instrument manufacturer about alternative pressure regulators or the use of EPC.

Using an ambient initial temperature, gasoline analyzed on a Petrocol DH 150 column shows near baseline resolution of the light hydrocarbons cyclopentane and 2,3-dimethylbutane, and the heavier m-xylene and p-xylene isomers (Figure 12). Under these conditions, separations of the light hydrocarbons through toluene are optimized by using the isothermal hold period. The subsequent temperature program provides excellent detail for the heavier compounds. Conditions can be varied (e.g., analysis time shortened) to meet specific needs.

**Figure 11. Oxygenated Compounds in a C4 Gas Stream, Using Ambient, Isothermal Temperature**

Column: **Petrocol DH 150, 150m x 0.25mm ID, 1.0µm film**  
 Cat. No.: **24155**  
 Oven: 25°C  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 50µL nitrogen (0.1-2% each analyte), split (100:1)

1. Ethane
2. Propylene
3. Propane
4. Dimethylether
5. Isobutane
6. Methanol
7. Isobutylene
8. 1-Butene
9. n-Butane
10. trans-2-Butene
11. cis-2-Butene
12. Ethanol
13. n-Pentane
14. Methyl tert-butyl ether




713-0765

**Figure 12. Detailed Analysis of Gasoline, Using Ambient Initial Temperature**

Column: **Petrocol DH 150, 150m x 0.25mm ID, 1.0µm film**  
 Cat. No.: **24155**  
 Oven: 35°C (125 min) to 200°C at 2°C/min  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.1µL unleaded gasoline, split (100:1)

- |                         |                            |                                 |
|-------------------------|----------------------------|---------------------------------|
| 1. n-Butane             | 15. 2,2,4-Trimethylpentane | 37. 1-Methyl-3-ethylbenzene     |
| 2. Isopentane           | 16. n-Heptane              | 38. 1-Methyl-4-ethylbenzene     |
| 3. n-Pentane            | 17. 2,5-Dimethylhexane     | 39. 1,3,5-Trimethylbenzene      |
| 4. 2,2-Dimethylbutane   | 18. 2,4-Dimethylhexane     | 40. 3,3,4-Trimethylheptane      |
| 5. Cyclopentane         | 19. 2,3,4-Trimethylpentane | 41. 1-Methyl-2-ethylbenzene     |
| 6. 2,3-Dimethylbutane   | 20. Toluene                | 42. 1,2,4-Trimethylbenzene      |
| 7. 2-Methylpentane      | 21. 2,3,3-Trimethylpentane | 43. Isobutylbenzene             |
| 8. 3-Methylpentane      | 22. 2,3-Dimethylhexane     | 44. sec-Butylbenzene            |
| 9. n-Hexane             | 23. 2-Methylheptane        | 45. n-Decane                    |
| 10. 2,4-Dimethylpentane | 24. 3-Methylheptane        | 46. 1,2,3-Trimethylbenzene      |
| 11. Benzene             | 25. 2-Methyl-1-heptene     | 47. Indane                      |
| 12. 2-Methylhexane      | 26. n-Octane               | 48. 1,3-Diethylbenzene          |
| 13. 2,3-Dimethylpentane | 27. Ethylbenzene           | 49. n-Butylbenzene              |
| 14. 3-Methylhexane      | 28. m-Xylene               | 50. 1,4-Dimethyl-2-ethylbenzene |
|                         | 29. p-Xylene               | 51. 1,3-Dimethyl-4-ethylbenzene |
|                         | 30. o-Xylene               | 52. 1,2-Dimethyl-4-ethylbenzene |
|                         | 31. 1-Nonene               | 53. 1,2,4,5-Tetramethylbenzene  |
|                         | 32. n-Nonane               | 54. 1,2,3,5-Tetramethylbenzene  |
|                         | 33. Isopropylbenzene       | 55. Naphthalene                 |
|                         | 34. 3,3,5-Trimethylheptane | 56. 2-Methylnaphthalene         |
|                         | 35. 2,4,5-Trimethylheptane | 57. 1-Methylnaphthalene         |
|                         | 36. n-Propylbenzene        | 58. Dimethylnaphthalenes        |




713-1367,1368

## Reproducible Results

The highly reproducible Petrocol DH Octyl capillary column is capable of resolving compounds of different classes to the baseline (Figures 13 to 15). Table 1 summarizes some key shifts in elution order. As a general rule, polar compounds shift to shorter retention times relative to normal hydrocarbons.

### Figure 13. Key Hydrocarbon Pairs Completely Resolved


Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 110°C to 220°C at 8°C/min  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



94-0364

### Figure 14. PONA V Mix


Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C (15 min) to 200°C at 1°C/min (15 min)  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



94-0366

### Figure 15. Pyrolysis Gasoline (Py Gas)

Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C (15 min) to 200°C at 1°C/min (15 min)  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



94-0368,0369

### Table 1. Comparison of Selected Elutions Between Petrocol DH Octyl and SPB-1 Columns


Class	Compounds	Column	
		SPB-1	Petrocol DH Octyl
A	Benzene	Coelution	Baseline Resolution
O	1-Methylcyclopentene		
A	Toluene	Coelution	Baseline Resolution
I	2,3,3-Trimethylpentane		Resolution
A	p-Xylene	Partial Resolution	Baseline Resolution
A	m-Xylene		Resolution (reversed)
P	n-Paraffin	Baseline Resolution	Coelution
O	trans-2-Olefin		
OXY	Methanol	Partial Resolution	Baseline Resolution

## Difficult Petroleum Analyses

Use a Petrocol DH column for analyses ranging from resolving closely eluting isomers or light hydrocarbon gases to detailed separations of highly complex mixtures — such as petroleum naphthas and reformates and liquified coal fractions (Figures 16-20).

### Figure 16. Monitor Dicyclopentadiene for Breakdown Products


Column: **Petrocol DH, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24160-U**  
 Oven: 35°C (15 min) to 230°C at 2°C/min (15 min)  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.1µL DCPD, split 100:1



80-89

### Figure 17. Quantify ENB Isomers and VNB Isomers


Column: **Petrocol DH, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24160-U**  
 Oven: 35°C (15 min) to 200°C at 2°C/min  
 Det.: FID  
 Inj.: 0.05µL ENB/VNB mixture, split 100:1



796-0095

### Figure 18. Closely Eluting Hexene Isomerization Products


Column: **Petrocol DH, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24160-U**  
 Oven: 35°C (15 min) to 45°C at 2°C/min (30 min)  
 Det.: FID  
 Inj.: 3µL mixed hexene isomers, split 100:1



796-0096

### Figure 19. Light Hydrocarbon Gases, Using a Subambient Initial Temperature

Column: **Petrocol DH, 100m x 0.25mm, 0.5µm film**  
 Cat. No.: **24160-U**  
 Oven: 15 min, at -25°C, then to 50°C at 2°C/min, hold 5 min  
 Carrier: helium, 22cm/sec (set at 35°C)  
 Det.: FID  
 Inj.: 200µL, split (100:1)




80-88


**Figure 20. Petroleum Reformate**

Column: **Petrocol DH, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24160-U**  
 Oven: 35°C (15 min) to 200°C at 2°C/min (30 min)  
 Det.: FID  
 Inj.: 0.1µL reformate, split 100:1

1. Benzene
2. Toluene
3. Ethylbenzene
4. m-Xylene
5. p-Xylene
6. o-Xylene
7. 1-Methyl-3-ethylbenzene
8. 1,2,4-Trimethylbenzene
9. Naphthalene
10. 2-Methylnaphthalene
11. 1-Methylnaphthalene
12. Dimethylnaphthalenes



Higher Resolution of Polynuclear Aromatic Hydrocarbons



## Light Hydrocarbons

### Monitor Light Hydrocarbons Without Using Cryogenics



Because they do not require cooling to retain compounds with low (-40°C to 100°C) boiling temperatures, thick film SPB™-1 and SPB-5 capillary columns are more convenient to use than other types of columns.

Typical analyses are shown in Figures 21 and 22. When you need additional resolution, you still have the option of using subambient initial temperatures (Figure 21 inset).

**Figure 21. C1-C6 Hydrocarbons and C1-C4 Alcohols on a Thick Film Capillary Column**

Column: **SPB-1, 60m x 0.53mm ID, 5µm film**  
 Cat. No.: **25349**  
 Oven: 30°C (5 min) to 200°C at 20°C/min (inset: initial temp. 0°C)  
 Carrier: helium, 19-21cm/sec  
 Det.: FID, 220°C  
 Inj.: 250µL (balance N<sub>2</sub>), split 100:1, 200°C

Component	Concentration (Vol. %)
1. Methane	0.3
2. Ethylene	0.2
3. Ethane	0.2
4. Propylene	0.3
5. Propane	0.2
6. Propadiene	0.2
7. iso-Butane	0.2
8. Methanol	0.2
9. Butene-1	0.3
10. n-Butane	0.2
11. cis-2-Butene	0.3
12. Ethanol	0.2
13. iso-Propanol	0.1
14. n-Pentane	0.4
15. tert-Butanol	0.1
16. n-Propanol	0.1
17. Methyl tert-butyl ether	0.3
18. sec-Butanol	0.1
19. n-Hexane	0.4
20. iso-Butanol	0.1
21. n-Butanol	0.1
22. Benzene	0.2
Nitrogen	balance




796-0097,0098



## Figure 22. Natural Gas on a Thick Film Column

Column: **SPB-5, 30m x 0.53mm ID, 3 $\mu$ m film**  
 Cat. No.: **25343**  
 Oven: 30°C (5 min) to 200°C at 20°C/min  
 Carrier: helium, 19-21cm/sec  
 Det.: FID, 220°C  
 Inj.: 250 $\mu$ L (balance N<sub>2</sub>), split 100:1, 200°C



796-0099

## High Boiling Hydrocarbons

### Better Quantification of High Boiling Compounds


For analyses of paraffin mixtures, SPB-1 Thin Film bonded phase capillary columns (Figures 23-25) offer better sample resolution and more reliable quantification than packed columns.

When you use an SPB-1 Thin Film column with our cool on-column injection sleeve (designed for 1/4" packed column or capillary injection ports), you can virtually eliminate discrimination among sample components. Use a syringe with a 6" (15.24cm) needle to make injections.

## Figure 23. Polywax® 655


Column: **SPB-1 Thin Film, 15m x 0.53mm ID, 0.1 $\mu$ m film**  
 Cat. No.: **25360**  
 Oven: 50°C to 350°C at 15°C/min  
 Carrier: helium, 20mL/min (flow controlled)  
 Det.: FID, 370°C  
 Inj.: 1 $\mu$ L, cool on-column injection

713-0752



## Figure 24. Shoe Polish


Column: **SPB-1 Thin Film, 15m x 0.53mm ID, 0.1 $\mu$ m film**  
 Cat. No.: **25360**  
 Oven: 50°C to 350°C at 15°C/min  
 Carrier: helium, 5mL/min (flow controlled)  
 Det.: FID, 370°C  
 Inj.: 1 $\mu$ L, cool on-column injection



713-0751


**Figure 25. Paraffin Wax**

Column: **SPB-1 Thin Film, 30m x 0.53mm ID, 0.1µm film**  
 Cat. No.: **25361**  
 Oven: 40°C to 350°C at 5°C/min  
 Carrier: helium, 5mL/min (flow controlled)  
 Det.: FID, 370°C  
 Inj.: 1µL, cool on-column injection



**Figure 27. MTBE Contaminants (Petrocol DH Octyl Column)**

Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C (15 min)  
 Carrier: helium, 24cm/sec (set at 35°C)  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



**Oxygenates**


**Separate Oxygenated Compounds in Complex Hydrocarbon Matrices**

Because of the high toxicity of organic bound lead, the US EPA has mandated the reduction and eventual elimination of lead compounds in gasoline. The compounds substituted for lead are usually oxygen-containing substances such as alcohols and ethers.

Separation of these oxygenates from the matrix can be a problem. Our Petrocol DH Octyl and Petrocol DH 150 capillary columns (Figures 26 to 28) are excellent for this type of analysis.


**Figure 26. Oxygenates Mix (Petrocol DH Octyl Column)**

Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C (15 min) to 200°C at 1°C/min (15 min)  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



**Figure 28. Oxygenates (Petrocol DH 150 Column)**

Column: **Petrocol DH 150, 150m x 0.25mm ID, 1.0µm film**  
 Cat. No.: **24155**  
 Oven: 25°C  
 Carrier: helium, 20cm/sec (set at 35°C)  
 Det.: FID, 200°C  
 Inj.: approx. 0.1-2% each component in 50µL nitrogen, split 100:1, 200°C





## Simulated Distillation

### Packed Columns for More Reliable SIMDIS Analyses

Petrocol A and Petrocol B packed columns were developed specifically for Methods D3710 and D2887, respectively, of the American Society for Testing and Materials (ASTM). These 20" x 1/8" stainless steel columns minimize peak tailing, column bleed, and the deviation of apparent boiling point from true boiling point (Figures 29-33).

**Figure 29. Petrocol A Column Minimizes Hydrocarbon Peak Tailing**


Column: **Petrocol A, 20" x 1/8" SS**  
 Cat. No.: **12445**  
 Col. Temp.: -20°C to 200°C at 20°C/min  
 Carrier: helium, 25mL/min  
 Det.: FID  
 Inj.: 0.1µL of D3710 Qualitative Mix (Cat. No. 48884)



713-0772\_1395

**Figure 30. Gasoline by ASTM D3710**


Column: **Petrocol A, 20" x 1/8" SS**  
 Cat. No.: **12445**  
 Col. Temp.: -20°C to 250°C at 15°C/min  
 Carrier: helium, 25mL/min  
 Det.: FID  
 Inj.: 0.1µL of gasoline



713-0771

**Figure 31. Hydrocarbon Boiling Point and Retention on a Petrocol A Column Are Closely Correlated**

Column: **Petrocol A, 20" x 1/8" SS**  
 Cat. No.: **12445**  
 Col. Temp.: -20°C to 200°C at 20°C/min  
 Carrier: helium, 25mL/min  
 Det.: FID  
 Inj.: 0.1µL of D3710 Qualitative Mix (Cat. No. 48884)




713-1393


**Figure 32. Minimize Hydrocarbon Peak Tailing and Column Bleed for Method D2887**

Column: **Petrocol B, 20" x 1/8" SS**  
 Cat. No.: **12449**  
 Col. Temp.: -25°C to 350°C at 15°C/min  
 Carrier: helium, 30mL/min  
 Det.: FID  
 Inj.: 0.1µL mixed aliphatics and aromatics

Petrocol B Column: Negligible Bleed at 350°C




10% UC-W982 on 80/100 Chromosorb P AW:  
 High Bleed at 350°C



713-0770, 1394

**Figure 33. Negligible Column Bleed for Gas Oil Capillary Columns for SIMDIS Analysis**

Column: **Petrocol B, 20" x 1/8" SS**  
 Cat. No.: **12449**  
 Col. Temp.: -25°C to 350°C at 15°C/min  
 Carrier: helium, 30mL/min  
 Det.: FID  
 Inj.: 0.1µL reference gas oil



713-0769


**Capillary Columns for SIMDIS Analysis**

Petrocol 3710 and Petrocol 2887 capillary columns can improve your results for ASTM Methods D3710 and D2887, respectively (Figures 34-39 and Table 2). These columns offer important advantages over even the most reliable packed columns:

- Longer retention of early-eluting hydrocarbons — for a more linear boiling point/retention time curve (Figures 34-35)
- Low column bleed — for a sharp return to baseline and more reliable final boiling point data
- Symmetrical peaks
- Rapid conditioning (15 minutes to 1 hour)
- A potentially longer column lifetime, because the phase is bonded


You can connect these capillary columns to packed column injection ports and detectors by using our inexpensive adapter kits. Use the columns with packed column flow controllers and most autosamplers. Petrocol capillary columns satisfy all column performance criteria in the ASTM methods.

**Figure 34. Petrocol 3710 Capillary Column Provides Near-Linear Correlation Between Hydrocarbon Boiling Point and Retention Time**



713-1397, 1398, 1399


**Figure 35. Linear Boiling Point/Elution Time Relationship for n-Paraffins on a Petrocol 2887 Column**



713-1401

**Figure 36. D3710 Qualitative SIMDIS Standard**

Column: **Petrocol 3710, 10m x 0.75mm ID glass, 5.0µm film<sup>▲</sup>**  
 Col. Temp.: -20°C to 200°C at 20°C/min  
 Carrier: nitrogen, 15mL/min  
 Det.: FID  
 Inj.: 0.1µL D3710 Qualitative Mix (Cat. No. 48884)




713-0773

<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column.

**Figure 37. SIMDIS Analysis of Gasoline**


Column: **Petrocol 3710, 10m x 0.75mm ID glass, 5.0µm film<sup>▲</sup>**  
 Col. Temp.: -20°C to 200°C at 20°C/min  
 Carrier: nitrogen, 15mL/min  
 Det.: FID  
 Inj.: 0.1µL commercial gasoline



713-0775

**Figure 38. D2887 SIMDIS Calibration Blend**

Column: **Petrocol 2887, 5m x 0.53mm ID, 0.5µm film**  
 Cat. No.: **25323**  
 Col. Temp.: -20°C to 320°C at 20°C/min, hold 5 min  
 Carrier: nitrogen, 6mL/min  
 Det.: FID  
 Inj.: 0.1µL Quantitative Calibration Mix (Cat. No. 48882), direct



713-0776


**Table 2. SIMDIS Data from Petrocol Capillary Columns and Packed Columns Show Good Agreement<sup>■</sup>**

Sample Volume % Off <sup>○</sup>	Calc. Temp.(°F) <sup>▼</sup>		Relative Difference %	Boiling Pt. (°F) <sup>▲</sup>		Relative Difference %
	Petrocol 3710 Capillary Column	Packed Column		Petrocol 2887 Capillary Column	ASTM Packed Column Consensus	
IBP	20	34	37.0	231	238	2.9
5	76	87	9.5	280	289	3.1
10	95	94	0.7	325	336	3.3
15	138	115	13.0	370	384	3.6
20	163	150	5.9	412	429	4.0
25	186	164	8.9	451	—	—
30	197	185	4.4	485	496	2.2
35	227	203	7.9	515	—	—
40	232	215	5.4	542	548	1.1
45	260	243	4.8	568	—	—
50	278	250	7.5	590	594	0.7
55	283	278	1.3	611	—	—
60	311	295	3.7	630	629	0.2
65	322	308	3.1	650	—	—
70	332	336	0.8	670	669	0.1
75	345	347	0.4	690	—	—
80	364	365	0.2	711	709	0.3
85	387	392	0.9	734	732	0.3
90	405	416	1.9	761	759	0.3
95	440	458	2.8	794	797	0.4
FBP	624	662	4.2	868	887	2.1

<sup>■</sup> Means for 10 injections.  
<sup>○</sup> Each value represents 5% of the sample eluted from the column.  
<sup>▼</sup> Data obtained by using a commercial gasoline sample and Spectra-Physics D3710 SIMDIS software.  
<sup>▲</sup> Data obtained by using ASTM reference gas oil sample and Spectra-Physics D2887 SIMDIS software.

**Figure 39. Rapid Analysis of Reference Gas Oil, with Negligible Bleed**

Column: Petrocol 2887, 5m x 0.53mm ID, 0.5µm film  
 Cat. No.: 25323  
 Oven: -20°C to 320°C at 20°C/min, hold 5 min  
 Carrier: nitrogen, 6mL/min  
 Det.: FID, 350°C  
 Inj.: 0.1µL Reference Gas Oil (Cat. No. 48873), direct injection, 350°C



For more information, request Bulletin 864.


713-0774

## Aliphatics from Aromatics in Gasoline

Use a SUPELCOWAX™ 10 column to separate aliphatic and aromatic hydrocarbons. Nonane will elute before benzene, and the xylene isomers will be resolved (Figures 40 and 41).

**Figure 40. Aliphatics in Gasoline Elute Before Aromatics**

Column: SUPELCOWAX 10, 30m x 0.75mm ID, 1.0µm film<sup>▲</sup>  
 Oven: 60°C (8 min) to 230°C at 8°C/min (8 min hold)  
 Carrier: helium, 4mL/min  
 Det.: FID  
 Inj.: 0.1µL super unleaded gasoline




796-0101

<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column, or request Cat. No. 25301-1 (30m x 0.53mm ID fused silica, 1.0µm film).

**Figure 41. Nonane Elutes Before Benzene; Xylene Isomers Easily Resolved**

Column: **SUPELCO WAX 10, 30m x 0.75mm ID glass, 1.0µm film<sup>▲</sup>**  
 Oven: 50°C  
 Carrier: helium, 5.3mL/min, 50°C  
 Det.: FID  
 Inj.: 0.1µL CS<sub>2</sub> containing 1% each component, direct injection




<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column, or request Cat. No. 25301-U (30m x 0.53mm ID fused silica, 1.0µm film).

## Polynuclear Aromatic Hydrocarbons in Gasoline

SPB columns offer the efficiency and high thermal stability you need to isolate individual PAHs in gasoline (Figure 42).

**Figure 42. Gasoline Concentrated for Monitoring PAH Content**

Column: **SPB-5, 30m x 0.25mm ID, 0.25µm film**  
 Cat. No.: **24034**  
 Oven: 35°C to 320°C at 4°C/min  
 Carrier: helium, 22.2cm/sec  
 Det.: FID  
 Inj.: 0.5µL unweathered or 0.02µL weathered super unleaded gasoline



## Trace Compounds in Petroleum Products

Use a SUPELCOWAX 10 column to detect 10ppm or less of sulfolane in aromatic-rich or raffinate petroleum fractions, or in water washes (Figure 43). Sulfolane is eluted after the dimethylnaphthalenes, ensuring minimal interference.



A 0.75mm ID SUPELCOWAX 10 capillary column (Figure 44) can replace several packed columns for monitoring impurities in aromatic petroleum products (ASTM Method D4534). It ensures:

- Faster analyses
- Greater resolution
- More stable baseline for trace analyses
- Separation of xylene isomers
- Rapid recycle times
- Longer column life

The wide bore SUPELCOWAX 10 column can be used with packed column injection and detector systems.

### Figure 43. Detect Low Levels of Sulfolane in Petroleum Streams

Column: SUPELCOWAX 10, 30m x 0.75mm ID glass, 1.0µm film<sup>▲</sup>  
 Oven: 175°C to 225°C at 5°C/min  
 Carrier: helium, 22.2cm/sec (6mL/min)  
 Det.: FID  
 Inj.: 1µL gasoline (w/o or w/0.05% sulfolane), direct injection




796-0104

<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column, or request Cat. No. 25301-U (30m x 0.53mm ID fused silica, 1.0µm film).

### Figure 44. Impurities in Aromatic Petroleum Products by ASTM Method D4534

#### Major Component Overloaded to Detect Trace Components

Column: SUPELCOWAX 10, 30m x 0.75mm ID glass, 1.0µm film<sup>▲</sup>  
 Oven: 100°C (40°C for cyclohexane)  
 Carrier: helium, 19.3cm/sec (3.5mL/min)  
 Det.: FID  
 Inj.: 0.1µL each major aromatic compound containing 1000ppm benzene, direct injection



796-0105,0106,0107,0108,0109



## Miscellaneous Petrochemicals Analyses


### More Sample Information from Complex Hydrocarbon Mixtures

Our nonpolar Petrocol DH 50.2 and SPB-1 capillary columns provide excellent resolution of complex mixtures — such as charcoal lighter fluid, turpentine, and gasoline (Figures 45-47). These bonded phase columns offer a potentially long column life.

You can use 0.53mm ID SPB-1 columns in packed column instruments, with thermal conductivity or other less sensitive detectors. These wide bore capillary columns will give you more sample information than any packed column can provide.

### Figure 45. Charcoal Lighter Fluid


Column: **Petrocol DH 50.2, 50m x 0.20mm ID, 0.50µm film**  
 Cat. No.: **24133-U**  
 Oven: 35°C (5 min) to 200°C at 2°C/min  
 Carrier: helium, 19-21cm/sec  
 Det.: FID



### Figure 47. Gasoline on a Capillary Column in a Packed Column Chromatograph

Column: **SPB-1, 60m x 0.75mm ID glass, 1.0µm film<sup>▲</sup>**  
 Oven: 35°C (12 min) to 200°C at 4°C/min, hold  
 Carrier: helium, 5cc/min (flow controlled)  
 Det.: TCD  
 Inj.: 0.4µL gasoline, direct injection

- |                         |                        |
|-------------------------|------------------------|
| 1. Air                  | 17. 2,4-Dimethylhexane |
| 2. Propane              | 18. Toluene            |
| 3. Isobutane            | 19. 2-Methylheptane    |
| 4. n-Butane             | 20. 4-Methylheptane    |
| 5. Isopentane           | 21. 3-Methylheptane    |
| 6. n-Pentane            | 22. n-Octane           |
| 7. 2-Methylpentane      | 23. Ethylbenzene       |
| 8. 3-Methylpentane      | 24. m-Xylene           |
| 9. n-Hexane             | 25. p-Xylene           |
| 10. Methylcyclopentane  | 26. 2-Methyloctane     |
| 11. 2,4-Dimethylpentane | 27. 3-Methyloctane     |
| 12. Benzene             | 28. o-Xylene           |
| 13. Cyclohexane         | 29. n-Nonane           |
| 14. 2-Methylhexane      | 30. Propylbenzene      |
| 15. 3-Methylhexane      | 31. n-Decane           |
| 16. n-Heptane           | 32. n-Undecane         |




713-1351

<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column, or request Cat. No. 25418 (60m x 0.53mm ID fused silica, 1.0µm film).

### Figure 46. Turpentine


Column: **Petrocol DH 50.2, 50m x 0.20mm ID, 0.50µm film**  
 Cat. No.: **24133-U**  
 Oven: 35°C (5 min) to 200°C at 2°C/min  
 Carrier: helium, 19-21cm/sec  
 Det.: FID



Because our SPB-1 bonded phase columns are stable to 320°C, you can use them to resolve many high-boiling hydrocarbons (Figure 48). To minimize column bleed and avoid splitter discrimination, use a 0.5µm film, 0.53mm ID column and cool on-column injection. Stabilized, high polarity Nukol™ columns provide sharp separation of oxo-isooctyl isomers for accurate determination of 2-ethylhexanol (Figure 49). SPB-20 bonded phase columns resolve almost all alkylbenzene isomers to baseline (Figure 50). The Petrocol DH Octyl column completely separates methanol from isobutane (Figure 51).

### Figure 48. Petroleum Wax Components Resolved Without Splitter Discrimination


Column: **SPB-1, 30m x 0.53mm ID, 0.50µm film**  
 Cat. No.: **25315**  
 Oven: 200°C to 320°C at 16°C/min (hold 45 min)  
 Carrier: helium, 20mL/min  
 Det.: FID  
 Inj.: 10µg petroleum wax in 1µL n-dodecane



796-0104

### Figure 49. 2-Ethylhexyl Alcohol Isolated from Other Oxo-isooctyl Alcohol Isomers

Column: **Nukol, 30m x 0.53mm ID, 0.5µm film**  
 Cat. No.: **25327**  
 Oven: 90°C to 160°C at 4°C/min (hold 5 min)  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.2µL CH<sub>3</sub>OH, 2% isooctyl alcohols




713-0794

### Figure 50. Complex Mixtures of Alkylbenzenes

Column: **SPB-20, 30m x 0.25mm ID, 0.25µm film**  
 Cat. No.: **24086**  
 Oven: 170°  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.2µL mixed alkylbenzenes, split 100:1

- |                  |                   |
|------------------|-------------------|
| 1. 5øC10*        | 15. C11 tetralins |
| 2. 4øC10         | 16. 2øC12         |
| 3. 3øC10         | 17. 6øC13, 7øC13  |
| 4. 2øC10         | 18. 5øC13         |
| 5. 6øC11         | 19. 4øC13         |
| 6. 5øC11         | 20. 3øC13         |
| 7. 4øC11         | 21. C12 tetralins |
| 8. 3øC11         | 22. 2øC13         |
| 9. C10 tetralins | 23. 7øC14         |
| 10. 2øC11        | 24. 6øC14         |
| 11. 6øC12        | 25. 5øC14         |
| 12. 5øC12        | 26. 4øC14         |
| 13. 4øC12        | 27. 3øC14         |
| 14. 3øC12        |                   |


\* Position of phenyl group, chain length



796-0112

**Figure 51. Methanol and Isobutane**

Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



- 1. Methanol
- 2. Isobutane

94-0367


## Phenols / Polymers

### Cresol Isomers Separated, Using an Acidic, Bonded Phase Column

Phenols (and many other acidic polar compounds) are eluted from high efficiency Nukol bonded phase capillary columns as sharp, symmetrical peaks. A Nukol column will resolve p-cresol from the dimethylphenols (Figure 52).

**Figure 52. Resolve p-Cresol from Dimethylphenols**

Column: **Nukol, 30m x 0.25mm ID, 0.25µm film**  
 Cat. No.: **24107**  
 Oven: 200°C  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 1µL methylene chloride containing 1µg each phenol, split 100:1



2,4- & 2,5-Dimethylphenol

- 1. 2,6-Dimethylphenol
- 2. o-Cresol
- 3. Phenol
- 4. o-Ethylphenol
- 5. m-Cresol
- 6. 2,3-Dimethylphenol
- 7. p-Ethylphenol & 3,5-Dimethylphenol
- 8. m-Ethylphenol
- 9. 3,4-Dimethylphenol

796-0114


## Capillary Column Improves Analyses of Dodecylphenol Product

Our 0.53mm ID SPB-5 fused silica capillary column is ideal for this purpose. Compared to packed columns, the inert, bonded phase capillary column improves peak symmetry for the acidic sample components. Phenol elutes before the C12 olefins, rather than after them, for more reliable quantification of phenol and olefins (Figure 53).

This wide bore capillary column can be used in packed column instruments.

**Figure 53. Dodecylphenol Product Components Separate Well on a Bonded Phase Capillary Column**

Column: **SPB-5, 30m x 0.53mm ID, 1.5µm film**  
 Cat. No.: **25305-U**  
 Oven: 80°C (2 min) to 300°C at 8°C/min  
 Carrier: helium, 7mL/min  
 Det.: FID, 325°C  
 Inj.: 0.2µL carbon disulfide:dodecylphenol product (2:1), direct injection, 300°C



- 1. CS<sub>2</sub> (solvent)
- 2. Phenol
- 3. C12 olefins
- 4. Alkylphenol ethers
- 5. Dodecylphenol
- 6. Di-dodecylphenol

713-0839

## Solvents


### Analyze Almost Any Combination of Solvents

You can easily separate complex mixtures of common solvents on our nonpolar SPB-1 or polar SUPELCOWAX 10 capillary column. (SPB-1 columns separate the alcohols by boiling point.) Compare the separations shown in Figures 54-56 to determine which column is best for your samples.

Furthermore, you can use 0.53mm ID columns in systems designed for packed columns.

### Figure 54. Industrial Solvents

Column: **SPB-1, 30m x 0.32mm ID, 1.0µm film**  
 Cat. No.: **24045-U**  
 Oven: 30°C (8 min) to 125°C at 4°C/min, hold 5 min  
 Carrier: helium, 25cm/sec (set at 30°C)  
 Det.: FID, 250°C  
 Inj.: 0.2µL approximately equal proportions of each component, split 200:1, 250°C




### Solvents for Figures 54-56

- |                           |                            |
|---------------------------|----------------------------|
| 1. Methanol               | 20. 1,4-Dioxane            |
| 2. Methyl formate         | 21. Ethyl Cellosolve       |
| 3. Ethanol                | U Unknown                  |
| 4. Acetone                | 22. Methyl isobutyl ketone |
| 5. Isopropyl alcohol      | 23. Toluene                |
| 6. Ethyl formate          | 24. Isobutyl acetate       |
| 7. Methylene chloride     | 25. Mesityl oxide          |
| 8. n-Propyl alcohol       | 26. n-Butyl acetate        |
| 9. Methyl ethyl ketone    | 27. Diacetone alcohol      |
| 10. sec-Butyl alcohol     | 28. 5-Methyl-2-hexanone    |
| 11. Ethyl acetate         | 29. Ethylbenzene           |
| 12. Chloroform            | 30. m-Xylene               |
| 13. Tetrahydrofuran       | 31. p-Xylene               |
| 14. Methyl Cellosolve     | 32. Cyclohexanone          |
| 15. Isobutyl alcohol      | U Unknown                  |
| 16. 1,1,1-Trichloroethane | 33. o-Xylene               |
| 17. Isopropyl acetate     | 34. Cellosolve acetate     |
| 18. n-Butyl alcohol       | 35. Butyl Cellosolve       |
| 19. 2-Nitropropane        | 36. Isoamyl acetate        |


### Figure 55. Solvents on a Capillary Column in a Packed Column GC

Column: **SPB-1, 60m x 0.75mm ID, 1.0µm film<sup>▲</sup>**  
 Oven: 30°C (8 min) to 125°C at 4°C/min (hold)  
 Carrier: helium, 20cm/sec  
 Det.: FID  
 Inj.: 0.5µL, approx. equal proportions of each component, direct injection



### Figure 56. Solvents on a Capillary Column in a Packed Column GC

Column: **SUPELCOWAX 10, 30m x 0.53mm ID, 1.0µm film**  
 Cat. No.: **25301-U**  
 Oven: 50°C (4 min) to 175°C at 4°C/min (hold 4 min)  
 Carrier: helium, 5mL/min (flow controlled)  
 Det.: FID  
 Inj.: 0.1µL mixed neat solvents, direct injection



<sup>▲</sup>This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column, or request Cat. No. 25418 (60m x 0.53mm ID fused silica, 1.0µm film).

## Eliminate Glycol Peak Tailing and Sample Decomposition

Use 0.53mm ID Nukol columns to analyze these highly polar compounds with good resolution and peak symmetry, and at low temperatures that eliminate decomposition (Figure 57).

## Capillary Chromatography for Glycols in a Packed Column Chromatograph

You can have the inertness and high resolution of capillary GC without using expensive, specialized equipment. Use 0.53mm ID Nukol columns in packed column chromatographs, with injectors and detectors designed for packed columns.


## Analyze Trace Solvents in Pharmaceutical Products

Figure 58 shows the elution of six volatile organic solvents, commonly used in manufacturing pharmaceuticals, separated on a thick film SPB-1 capillary column. The low parts-per-million concentrations are representative of residual solvent levels in bulk drug preparations.

**Figure 57. Nukol Column Provides Sharp Glycol Peaks**

Column: **Nukol, 30m x 0.53mm ID, 0.5µm film**  
 Cat. No.: **25327**  
 Oven: 45°C (4 min) to 220°C at 10°C/min, hold 5 min  
 Carrier: helium, 20cm/sec (set at 155°C)  
 Det.: FID, 250°C  
 Inj.: 0.02µL ink solvents mixture, direct injection, 250°C


1. Heptane
2. Ethyl acetate
3. Isopropyl acetate
4. Methanol
5. Isopropanol
6. Ethanol
7. n-Propyl acetate
8. n-Propanol
9. Carbitol
- U Unknown
10. Ethylene glycol
11. Diethylene glycol



713-0793

**Figure 58. Solvents Commonly Used in Pharmaceuticals Manufacture**

Column: **SPB-1, 60m x 0.53mm ID, 5µm film**  
 Cat. No.: **25349**  
 Oven: 30°C (5 min) to 200°C at 20°C/min  
 Carrier: helium, 19-21cm/sec  
 Det.: FID, 220°C  
 Inj.: 2.0cc, split 100:1, 200°C



796-0159


<sup>▲</sup> Values are equivalent to headspace analysis for 10g of a pharmaceutical product.

## Impurities in Solvents

For monitoring impurities in ketones (ASTM Methods D2804, D3329, and D3893) and 2-ethylhexanol, SUPELCOWAX 10 and Nukol capillary columns improve peak resolution, shorten analysis times, and provide more repeatable results than packed columns (Figures 59 to 61).

**Figure 59. Trace Levels of Impurities in Methyl Ethyl Ketone, by Capillary GC**


Column: **SUPELCOWAX 10, 30m x 0.32mm ID, 0.5µm film**  
 Cat. No.: **24084**  
 Oven: 50°C (4 min) to 100°C at 8°C/min  
 Carrier: helium, 25cm/sec  
 Det.: FID, 220°C  
 Inj.: 0.2µL MEK, split 100:1




796-0118

\*Contaminant in sec-butanol, not spiked in sample.

**Figure 60. Choose Maximum Resolution or Save Time When Monitoring Purity of Methyl Isobutyl Ketone**



**Figure 61. Impurities from 2-Ethylhexanol**



**Table 3. Capillary Column Ensures Consistent Measurements of Impurities in Solvents**

Solvent/Impurity	Variation in Measurement (% absolute)				
	Injection* - Injection**		Column - Lab* - Lab**		Column
MIAK	0.008	0.15	0.005	1.07	0.011
Isopropyl alcohol	0.001	0.014	0.001	0.05	0.006
Mesityl oxide	0.001	0.04	0.001	0.05	0.002
Methyl isoamyl carbitol	0.008	0.04	0.003	0.11	0.005
Total ketones	0.004	0.15	0.002	0.62	0.009
MAK	0.046	0.24	0.050	0.84	0.116
Isopropyl alcohol	0.006	0.02	0.005	0.05	0.022
Methyl isobutyl carbitol	0.001	0.01	0.004	0.07	0.019
Total ketones	0.024	0.05	0.023	0.61	0.035


\* Value for SUPELCO WAX 10 capillary column.

\*\* ASTM specification for packed columns. Packed column specifications established at 0.1% each impurity, capillary column values determined at 0.05% each impurity.

### High Resolution or High Speed Analyses


When you need high resolution, a 0.20mm ID SUPELCO WAX 10 capillary column is ideal. Its high efficiency (approximately 25% greater than that of a 0.25mm ID column) lets you obtain the greatest resolution possible. When analysis time is most important, you can trade efficiency for more rapid separations (Figures 62 and 63).

**Figure 62. Trace Contaminants Completely Resolved from Major Sample Component**



**Figure 63. Fastest Separation of Aromatic Compounds in Balanced Proportions**

Column: **SUPELCO WAX 10, 30m x 0.20mm ID, 0.20µm film**  
 Cat. No.: **24169**  
 Oven: 75°C  
 Carrier: helium, 25cm/sec  
 Det.: FID  
 Inj.: 0.1µL hexane  
 (1.7% each analyte), split (150:1)




713-1355

Because aromatic and nonaromatic compounds have similar boiling points, they are hard to resolve on nonpolar columns. Using a narrow bore polar capillary column simplifies quantitative analysis of aromatic compounds. When gasoline is analyzed on a polar 30m x 0.20mm ID SUPELCO WAX 10 column, benzene and higher molecular weight aromatics are retained until after most nonpolar components (up to n-nonane) elute (Figure 64).

**Figure 64. Monitor Aromatic Components of Gasoline Rapidly and Efficiently**

Column: **SUPELCO WAX 10, 30m x 0.20mm ID, 0.20µm film**  
 Cat. No.: **24169**  
 Oven: 60°C (8 min) to 210°C at 8°C/min  
 Carrier: helium, 25cm/sec  
 Det.: FID  
 Inj.: 0.1µL gasoline, split (150:1)



713-1354

## Aromatics / Aliphatics by Packed Column GC

### Aliphatics from Aromatics

Most packed columns allow aliphatics and aromatics to elute together. Our high polarity 1,2,3-tris(2-cyanoethoxy propane (TCEP), BC-120, and BC-150 packings retard the aromatics relative to the aliphatics (Figure 65). The BC-150 phase has a higher temperature limit (240°C) than the BC-120 or TCEP phases (125°C).

### Xylene Isomers


5% SP-1200/5% Bentone® 34 packing separates o-, m-, and p-xylene and ethylbenzene in less than 7 minutes. If styrene, isopropylbenzene, or n-propylbenzene are present, use 5% SP-1200/1.75% Bentone 34 (Figure 66).

### Difficult-to-Resolve Aromatics

Both Carbpac™ F-TA and Carbpac C/0.1% SP-1000 columns will rapidly separate butylbenzenes and many other difficult-to-resolve isomers. Analyses on a Carbpac F-TA column require half as much time (Figure 67).


### Figure 65. Aromatics and Aliphatics

Packing: **GP 10% TCEP on 100/120 Chromosorb P AW**  
 Cat. No.: **12106-U** (packing, 20g/bottle)  
 Column: 8' x 1/8" stainless steel  
 Oven: 80°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.2µL, approx. equal volumes each analyte




### Figure 66. Xylene Isomers

Packing: **GP 5% SP-1200/1.75% Bentone 34 on 100/120 SUPELCOPORT**  
 Cat. No.: **12134** (20g/bottle)  
 Column: 6' x 1/8" stainless steel  
 Cat. No.: **12721** (general configuration stock column, other stock columns available)  
 Oven: 75°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 0.10µL, approx. equal volumes each component



### Figure 67. Analyses are Twice as Fast on a Carbpac F-TA Column

Cat. No.: 11820 (packing, 15g/bottle)  
 Column: 6' x 1/8" ID stainless steel  
 Cat. No.: 12495-U (general configuration stock column, other stock columns available)  
 Oven: 225°C  
 Carrier: nitrogen, 20mL/min  
 Det.: FID  
 Inj.: 1µL chloroform (~0.1% each analyte)






## Aromatics / Aliphatics by Capillary GC

The Petrocol DH Octyl capillary column (Figure 68) offers baseline resolution of the BTEX compounds, including the difficult-to-separate p- and m-xylene isomers.

**Figure 68. BTEX Compounds**

Column: **Petrocol DH Octyl, 100m x 0.25mm ID, 0.5µm film**  
 Cat. No.: **24282**  
 Oven: 35°C (15 min) to 200°C at 1°C/min (15 min)  
 Carrier: helium, 24cm/sec  
 Det.: FID, 250°C  
 Inj.: 0.1µL, split 215:1, 250°C



## C1-C6 Hydrocarbons by Packed Column GC

### Packings Specifically Prepared for Optimum Hydrocarbon Separations


23% SP-1700 columns offer excellent resolution of concentrated C1-C6 hydrocarbons (Figure 69). Use these extremely durable columns between 0°C and 110°C, with sample valves or direct injections.

Carbopack C/0.19% picric acid packing is prepared specifically for separating C3 compounds or C4 unsaturates (Figure 70). Use this packing at temperatures to 120°C.

**Figure 69. C1-C6 Saturated and Unsaturated Hydrocarbons**


Packing: **23% SP-1700 on 80/100 Chromosorb P AW**  
 Cat. No.: **11865** (25g) (packing only)  
 Column: 30' x 1/8" stainless steel  
 Cat. No.: **12809-U** (general configuration column)  
 Oven: 70°C  
 Carrier: helium, 25mL/min  
 Det.: FID  
 Inj.: 0.6µL ASTM Section L Blend No. 6, plus C5s

- |                         |                                       |
|-------------------------|---------------------------------------|
| 1. Ethane               | 11. n-Pentane                         |
| 2. Propane              | 12. 1-Pentene                         |
| 3. Propylene            | 13. 2-Methyl-1-butene/trans-2-Pentene |
| 4. Isobutane            | 14. cis-2-Pentene                     |
| 5. n-Butane             | 15. 2-Methyl-2-butene                 |
| 6. 1-Butene/Isobutylene | 16. 2-Methyl-1-pentene                |
| 7. trans-2-Butene       | 17. 3-Methyl-1-pentene                |
| 8. cis-2-Butene         | 18. n-Hexane                          |
| 9. Isopentane           | U Unknown                             |
| 10. 1,3-Butadiene       | 19. 3-Methylhexane                    |



**Figure 70. C4 Unsaturation**

Packing: **GP 80/100 Carbopack C/0.19% picric acid**  
 Cat. No.: **11824** (15g)  
 Column: 2m x 1/8" stainless steel  
 Cat. No.: **13867** (general configuration column)  
 Oven: 50°C  
 Carrier: nitrogen, 30mL/min  
 Det.: FID  
 Inj.: 5µL ASTM Section L Blend No. 6




## Polymers by Gel Permeation Chromatography

### Silica-Based Columns


By combining a 100Å pore and a 300Å pore SUPELCOSIL™ column in series, you can separate compounds having molecular weights from 500 to over 500,000 Dalton (Figure 71).

SUPELCOSIL LC-Diol columns are suitable for separating water-soluble synthetic polymers, or for normal phase analyses (Figure 72).

**Figure 71. Polystyrene Standards by Molecular Weight**



**Figure 72. Nonionic Surfactants Resolved Almost to Baseline (Normal Phase Analysis)**



## Low Molecular Weight Compounds by Gas-Solid Chromatography

### Porous Polymers for Highly Consistent Performance

HayeSep® porous polymers are thoroughly cleaned, then conditioned for 12 hours under oxygen-free nitrogen before they are packaged. This careful preparation ensures materials that are highly consistent, with no shrinkage and minimum bleed. And they will perform more consistently than any other commercial polymer you are currently using (Figures 73-77).

**HayeSep N, P, Q, R, S, and T** polymers are interchangeable with the corresponding Porapak® polymers for separating low molecular weight mixtures containing halogenated or sulfur-containing compounds, water, alcohols, glycols, free fatty acids, esters, ketones, or aldehydes.


**HayeSep A** polymer separates permanent gases (hydrogen, nitrogen, oxygen, argon, carbon monoxide, nitric oxide) at ambient temperatures. Use it at higher temperatures to analyze C2 hydrocarbons, hydrogen sulfide, or water.

**HayeSep B** polymer is designed to separate C1 and C2 amines, and trace levels of ammonia and water. It does not require caustic washes before it is packed into a column.

**HayeSep C** polymer is for analyses of polar molecules — such as hydrogen cyanide, ammonia, hydrogen sulfide, and water. Separation characteristics are similar to those of Chromosorb 104 polymer.


**HayeSep D** polymer exhibits superior separation characteristics for light gases. It can be used to isolate carbon monoxide and carbon dioxide from room air at ambient temperatures, or to elute acetylene before other C2 hydrocarbons. It is particularly useful for separations and analyses of water and hydrogen sulfide.

**Figure 73. Hydrocarbons and Sulfur Gases on HayeSep Q**




**Figure 74. Mapp® Gas on HayeSep R**

Packing: **HayeSep R**  
 Cat. No.: **10304** (packing)  
 Column: 10' x 1/8" stainless steel  
 Oven: 80°C  
 Carrier: helium, 30mL/min  
 Inj.: 15mL




**Figure 76. Trace Ammonia in Water on HayeSep C**

Packing: **80/100 HayeSep C**  
 Cat. No.: **10289**  
 Column: 5' x 1/8"  
 Oven: 115°C  
 Det.: TCD  
 Inj.: 1µL




**Figure 77. Scott Mix 216 on HayeSep D**

Packing: **HayeSep D**  
 Cat. No.: **10293** (packing)  
 Column: 20' x 1/8" Ni  
 Oven: 40°C (2 min) to 110°C at 24°C/min  
 Carrier: helium, 30mL/min  
 Det.: P.E. 900 T.C., 225ma, 140°C  
 Inj.: 100mL, Valco valve, ambient



**Figure 75. Permanent Gases on HayeSep A**

Packing: **80/100 HayeSep A**  
 Cat. No.: **10283** (packing)  
 Column: 36' x 1/8"  
 Oven: 25°C  
 Carrier: helium, 23cc/min  
 Det.: P.E. 900 T.C., 225ma  
 Inj.: 25mL air plus hydrogen and carbon monoxide



**Table 4. Characteristics of HayeSep Polymers**

HayeSep Polymer	Maximum Operating Temp. (°C)	Surface Area (m <sup>2</sup> /gram)	Tapped Bulk Density (gram/cc)	Polarity (1 = lowest; 10 = highest)
A	165	526	0.356	7
B	190	608	0.330	8
C	250	442	0.322	6
D	290	795	0.3311	1
N	165	405	0.355	9
P	250	165	0.420	3
Q	275	582	0.351	2
R	250	344	0.324	5
S	250	583	0.334	4
T	165	250	0.381	10

## Ordering Information:

### Capillary Columns for Petrochemical Analysis

#### Amines: Carbowax Amine

**For primary, secondary, and tertiary amines** - The Carbowax Amine column is a specially prepared, base-deactivated polyethylene glycol column designed for the analysis of primary, secondary, and tertiary amines and other volatile basic analytes.

**Phase:** nonbonded; base-modified poly(ethylene glycol)

**Temp. Limits:** 60°C to 200°C

#### Amines: PTA-5

This column is a specially prepared, base-deactivated poly(5% diphenyl/95% dimethylsiloxane) column designed for analyses of amines and other basic analytes.

**Phase:** bonded; base-modified poly(5% diphenyl/95% dimethylsiloxane)

**Temp. Limits:** -60°C to 320°C

Length (m)	D <sub>f</sub> (µm)	Beta	Cat. No.
<b>Carbowax Amine</b>			
<b>0.53mm ID Fused Silica</b>			
15	1.0	133	<b>25352</b>
30	1.0	133	<b>25353</b>
60	1.0	133	<b>25354</b>
<b>PTA-5</b>			
<b>0.25mm ID Fused Silica</b>			
30	0.50	125	<b>24277</b>
30	1.0	625	<b>24330</b>
<b>0.32mm ID Fused Silica</b>			
30	0.5	160	<b>24331</b>
30	1.0	80	<b>24332</b>
30	1.5	53	<b>24333</b>
<b>0.53mm ID Fused Silica</b>			
30	0.5	265	<b>25437</b>
30	1.5	88	<b>25438</b>
30	3.0	44	<b>25439</b>

#### Aromatics: TCEP

This highly polar phase offers unique polarity for certain separations, despite its relatively low temperature limit and the fact that it is not a bonded phase. Because many aromatic compounds have retention indices greater than 1100 on TCEP, it is used for analyses of aromatics in mineral spirits and impurities in individual aromatics and oxygenates. Maximum, continuous isothermal operation of this column should be limited to under 100°C.

**Phase:** nonbonded; 1,2,3-tris-2-cyano-ethoxypropane

**Temp. Limits:** subambient to 145°C

**McReynolds Nos.:** x'y' z' u' s' = 594 857 759 1031 917

Length (m)	D <sub>f</sub> (µm)	Beta	Cat. No.
<b>0.25mm ID Fused Silica</b>			
60	0.44	142	<b>24153</b>
<b>0.32mm ID Fused Silica</b>			
60	0.51	157	<b>24161</b>

#### Gases: Carboxen™-1006 PLOT▼

**For permanent gases and C1 – C3 light hydrocarbons** - The porous carbon molecular sieve (surface area ~ 750m<sup>2</sup>/gram) in Carboxen-1006 porous layer open tubular (PLOT) columns separates permanent gases and C1, C2, and C3 light hydrocarbons, using above-ambient initial temperatures. The columns also are ideal for resolving formaldehyde/water/methanol (formalin) mixtures and monitoring impurities in ethylene. Use Carboxen-1006 columns with high flow rates and rapid temperature programs, up to 250°C, to ensure excellent, fast separations.

#### Gases: Carboxen-1010 PLOT▼

**For permanent gases, C2 and C3 hydrocarbons** - A carbon molecular sieve column for separating hydrogen, nitrogen, carbon monoxide, methane, carbon dioxide, and C2 and C3 hydrocarbons. Oxygen is separated from nitrogen.

- 0.53mm ID column can be used in packed column chromatographs.

#### Gases: Mol Sieve 5A PLOT▼

**For permanent gases** - Oxygen, nitrogen, carbon monoxide and methane can be separated in less than 5 minutes. More difficult separations, such as argon from oxygen, can be achieved by using subambient temperatures (15°C or below).

- 0.53mm ID column can be used in packed column chromatographs.

#### Gases: Supel-Q™ PLOT▼

**For many hydrocarbon and other compounds** - Supel-Q PLOT columns contain a porous divinylbenzene polymer that effectively resolves carbon dioxide and C1-C4 hydrocarbons at above ambient temperatures. It also is suitable for analyses of other gases, such as sulfur gases, and alcohols, ketones, aldehydes, and many polar compounds. Gasoline and other petroleum fractions can be analyzed as well. These columns exhibit very little bleed, even at the maximum temperature. Relative to packed columns (e.g., Porapak-Q), Supel-Q PLOT columns offer better resolution in less time.

- 0.53mm ID column can be used in packed column chromatographs.

Dimensions (fused silica)	Max. Temp. (°C)	Cat. No.
<b>Carboxen-1006 PLOT</b>		
30m x 0.32mm ID	250	<b>24241-U</b>
30m x 0.53mm ID	250	<b>25461</b>
<b>Carboxen-1010 PLOT</b>		
30m x 0.32mm ID	250	<b>24246</b>
30m x 0.53mm ID	250	<b>25467</b>
<b>Mol Sieve 5A PLOT</b>		
30m x 0.32mm ID	300	<b>24243</b>
30m x 0.53mm ID	300	<b>25463</b>
<b>Supel-Q PLOT</b>		
30m x 0.32mm ID	250	<b>24242</b>
30m x 0.53mm ID	250	<b>25462</b>

▼A proprietary procedure fixes particles to the fused silica tubing and ensures they will not be dislodged in normal use. Manufactured under US patents 5,599,445, 5,607,580, 5,609,756, 5,620,603, and 5,630,937.

## Gases: SCOT Stainless Steel

Support-coated open tubular (SCOT) columns are prepared by depositing a layer of liquid phase-coated support particles on the inner wall of the tubing. This technology, developed by Perkin-Elmer, makes available many phases that cannot be obtained on conventional wall-coated open tubular capillary columns. SCOT columns combine the sensitivity and excellent sample resolution of capillary GC with the extensive stationary phase library of packed column GC.

50' x 1/32" OD x 0.02" ID with 1/16" connections

Liquid Phase	Max. Temp. (°C)	Cont.▲ Isothermal	Beta	Cat. No.
<b>Stainless Steel SCOT Columns</b>				
Bentone 34/DNDP <sup>Ⓟ</sup>	150	80	45	<b>25521</b>
BMEA	100	50	40	<b>25538</b>
Squalane	120	80	50	<b>25535</b>
TCEP	150	100	40	<b>25536</b>

<sup>Ⓟ</sup>Di-n-decylphthalate

▲This glass column is no longer available. Request the same dimensions in a fused silica column, as a custom column.

## High Temperature: HT-5

**For highest-temperature separations** - SGE aluminum-clad columns coated with a carborane phase, offering the highest maximum temperature of any commercially available column. They display low bleed for GC/MS and simulated distillation analyses.

**Phase:** bonded; siloxane-carborane  
(5% phenyl equivalent)

**Temp. Limits:** 10°C to 460°C

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>0.32mm ID Fused Silica</b>			
12	0.10	800	<b>25002</b>
25	0.10	800	<b>25003</b>
<b>0.53mm ID Fused Silica</b>			
6	0.10	1325	<b>25004</b>
12	0.15	883	<b>25005-U</b>

## Hydrocarbons: Petrocol DH 50.2

**For detailed hydrocarbon analyses** - A narrow bore column for detailed hydrocarbon analyses of naphthas, gasolines, and similar samples, according to ASTM Test Method D5134.

**Phase:** bonded; poly(dimethylsiloxane)

**Temp. Limits:** -60°C to 320°C

## Hydrocarbons: Petrocol DH

**For PNA, PONA, PIANO-type analyses** - A highly reproducible column displaying more than 400,000 theoretical plates, designed for detailed analyses of petroleum products. Includes an extensive retention index data sheet of 400+ analytes.

**Phase:** bonded; poly(dimethylsiloxane)

**Temp. Limits:** -60°C to 320°C

## Hydrocarbons: Petrocol DH 150

**For detailed hydrocarbon analyses** - The longest capillary column commercially available as a stock item. Columns typically display more than 600,000 theoretical plates. For detailed purity analyses of light hydrocarbon gases and petroleum products (oxygenates, solvents, naphthas, gasolines, etc.).

## Hydrocarbons: Petrocol DH Octyl

**For detailed analyses of petroleum products** - This highly reproducible column offers unique selectivity not obtainable with poly(dimethylsiloxane) columns, such as baseline separations of benzene/1-methylcyclopentene and toluene/2,3,3-trimethylpentane.

**Phase:** bonded; poly(50% n-octyl/50% methylsiloxane)

**Temp. Limits:** -60°C to 220°C

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>Petrocol DH 50.2</b>			
<b>0.20mm ID Fused Silica</b>			
50	0.50	100	<b>24133-U</b>
<b>Petrocol DH</b>			
<b>0.25mm ID Fused Silica</b>			
100	0.50	125	<b>24160-U</b>
<b>Petrocol DH 150</b>			
<b>0.25mm ID Fused Silica</b>			
150	1.0	63	<b>24155</b>
<b>Petrocol DH Octyl</b>			
<b>0.25mm ID Fused Silica</b>			
100	0.50	125	<b>24282</b>

## SIMDIS: Petrocol 2887

**For ASTM Test Method D2887** - Developed and tested to meet or exceed column performance requirements for simulated distillation of petroleum fractions having boiling points up to 1000°F, according to ASTM Test Method D2887.

**Phase:** bonded; poly(dimethylsiloxane)

**Temp. Limits:** subambient to 350°C

## SIMDIS: Petrocol EX2887

**For extended ASTM Test Method D2887** - A thin film version of the Petrocol 2887 column, developed for extended D2887 SIMDIS analysis of samples having final boiling points greater than 1000°F.

**Phase:** bonded; poly(dimethylsiloxane)

**Temp. Limits:** subambient to 380°C

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>Petrocol 2887</b>			
<b>0.53mm ID Fused Silica</b>			
5 <sup>■</sup>	0.50	265	<b>25323</b>
5 <sup>♣</sup>	0.50	265	<b>25380-U</b>
<b>Petrocol EX2887</b>			
<b>0.53mm ID Fused Silica</b>			
5	0.10	1325	<b>25337</b>

<sup>■</sup>5" cage  
<sup>♣</sup>8" cage

### SPB-1 Thin Film and SPB-1 Thick Film

Nonpolar methylsilicone columns that separate sample components according to boiling point. This bonded polymer matches the polarity of its nonbonded predecessors, SE-30 and SP-2100. The SPB-1 phase is used in many of our Petrocol specialty columns.

**Operating Conditions:** Chemically compatible with water and other injection solvents. Sensitive to strong inorganic acids and bases, stable to low levels of HCl in non-aqueous samples. Not damaged by organic acids or bases. Columns can be rinsed.

**Phase:** bonded; poly(dimethylsiloxane)

**Temp. Limits:** -60°C to 320°C

**McReynolds Nos.:** x' y' z' u' s' = 4 58 43 56 38

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>0.10mm ID Fused Silica</b>			
15	0.10	250	<b>24338</b>
<b>0.25mm ID Fused Silica</b>			
30	0.10	625	<b>24261</b>
<b>0.32mm ID Fused Silica</b>			
30	0.10	800	<b>24290</b>
30	5.0	16	<b>24296</b>
60	5.0	16	<b>24297</b>
<b>0.53mm ID Fused Silica</b>			
15	0.10	1325	<b>25360</b>
30	0.10	1325	<b>25361</b>
15	3.0	44	<b>25340</b>
30	3.0	44	<b>25341-U</b>
60	3.0	44	<b>25348</b>
15	5.0	27	<b>25344</b>
30	5.0	27	<b>25345-U</b>
60	5.0	27	<b>25349</b>

### SPB-5 Thin Film and SPB-5 Thick Film

The low phenyl content, 5%, improves thermal stability of the phase, while still providing essentially a boiling point elution order, and a slight increase in selectivity, especially for aromatic compounds.

**Operating Conditions:** Chemically compatible with water and other injection solvents. Sensitive to strong inorganic acids and bases, stable to low levels of HCl in non-aqueous samples. Not damaged by organic acids or bases. Columns can be rinsed.

**Phase:** bonded; poly(5% diphenyl 95% dimethylsiloxane)

**Temp. Limits:** -60°C to 320°C

**McReynolds Nos.:** x' y' z' u' s' = 19 74 64 93 62

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>0.10mm ID Fused Silica</b>			
15	0.10	250	<b>24341</b>
<b>0.53mm ID Fused Silica</b>			
15	3.0	44	<b>25342</b>
30	3.0	44	<b>25343</b>
60	3.0	44	<b>25350</b>
15	5.0	27	<b>25346</b>
30	5.0	27	<b>25347</b>
60	5.0	27	<b>25351</b>

### Solvents: OVI-G43

**For USP Analysis of organic volatile impurities (OVIs)** - This column is specially prepared and tested to meet the requirements of United States Pharmacopoeia (USP) Method 467 and the European Pharmacopoeia general method for determining residual organic solvents in pharmaceutical preparations. Use this column to separate OVIs for research purposes or qualitative analysis. The USP and European Pharmacopoeia methods also specify using a deactivated 5-meter guard column.

**Phase:** bonded; poly(6% cyanopropylphenyl/94% dimethylsiloxane)

**Temp. Limits:** -20°C to 260°C

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>0.53mm ID Fused Silica</b>			
30	3.0	44	<b>25396</b>
Description			Cat. No.
<b>Deactivated Guard Column for OVI-G43</b>			
5m x 0.53mm ID			<b>25339</b>
<b>Other Columns for Residual Solvents Analysis</b>			
G27 (SPB-5) 30m x 0.53mm ID, 5.0μm			<b>25347</b>
G16 (SUPELCO WAX 10) 30m x 0.53mm ID, 1.0μm			<b>25301-U</b>

### Sulfur Compounds (Volatile): SPB-1 SULFUR

A very thick film version of our SPB-1 columns, specially developed for analyses of sulfur gases and other volatile sulfur compounds. The column displays relatively low column bleed, even for the exceptionally thick film (4μm) of stationary phase, which makes it compatible for use with the Sievers Sulfur Chemiluminescence Detector (SCD) and other sulfur-specific detectors.

**Phase:** bonded; poly(dimethylpolysiloxane)

**Temp. Limits:** -60°C to 300°C

Length (m)	D <sub>f</sub> (μm)	Beta	Cat. No.
<b>0.32mm ID Fused Silica</b>			
30	4.0	20	<b>24158</b>

## Supelco Columns and Packings Used in Figures 1-75, Listed by Application

For columns of other dimensions, refer to pages 28-30, or see the current Supelco catalog.

Description	Cat. No.	Description	Cat. No.
<b>Detailed Hydrocarbons Analyses</b>		<b>Solvents</b>	
<b>Petrocol DH Fused Silica Capillary Column</b> 100m x 0.25mm ID, 0.50µm film	24160-U	<b>SPB-1 Capillary Columns</b> 30m x 0.32mm ID fused silica, 1.0µm film	24045-U
<b>Petrocol DH 50.2 Fused Silica Capillary Column</b> 50m x 0.20mm ID, 0.50µm film	24133-U	<b>SUPELCOWAX 10 Capillary Columns</b> 30m x 0.53mm ID fused silica, 1.0µm film	25301-U
<b>Petrocol DH 150 Fused Silica Capillary Column</b> 150m x 0.25mm ID, 1.0µm film	24155	<b>Nukol Wide Bore Capillary Column</b> 30m x 0.53mm ID, 0.5µm film	25327
<b>Petrocol DH Octyl Fused Silica Capillary Columns</b> 100m x 0.25mm ID, 0.5µm film	24282	<b>SPB-1 Thick Film Fused Silica Capillary Columns</b> 60m x 0.53mm ID, 5µm film 30m x 0.53mm ID, 3µm film	25349 25341-U
<b>Light Hydrocarbons</b>		<b>Impurities in Solvents</b>	
<b>SPB-1 Thick Film Fused Silica Capillary Column</b> 60m x 0.53mm ID, 5µm film	25349	<b>SUPELCOWAX 10 Fused Silica Capillary Columns</b> 30m x 0.32mm ID, 0.25µm film 30m x 0.32mm ID, 0.5µm film 15m x 0.53mm ID, 1µm film	24080-U 24084 25300-U
<b>SPB-5 Thick Film Fused Silica Capillary Column</b> 60m x 0.53mm ID, 3µm film	25343	<b>Nukol Fused Silica Capillary Column</b> 30m x 0.32mm ID, 0.25µm film	24131
<b>High Boiling Hydrocarbons</b>		<b>High Resolution or High Speed Analyses</b>	
<b>SPB-1 Thin Film Fused Silica Capillary Columns</b> 15m x 0.53mm ID, 0.10µm film 30m x 0.53mm ID, 0.10µm film	25360 25361	<b>SUPELCOWAX 10 Fused Silica Capillary Columns</b> 30m x 0.20mm ID, 0.20µm film 60m x 0.20mm ID, 0.20µm film	24169 24170
<b>Oxygenates</b>		<b>Aromatics / Aliphatics by Packed Column GC</b>	
<b>Petrocol DH Octyl Fused Silica Capillary Columns</b> 100m x 0.25mm ID, 0.5µm film	24282	<b>Packings (20g)</b> GP 10% TCEP on 100/120 Chromosorb® P AW GP 10% TCEP on 80/100 Chromosorb P AW GP 35% BC-120 on 100/120 Chromosorb P AW DMCS	12106-U 12122 11834
<b>Petrocol DH 150 Fused Silica Capillary Column</b> 150m x 0.25mm ID, 1.0µm film	24155	GP 35% BC-150 on 100/120 Chromosorb P AW DMCS 10% SP™-2100 on 100/120 SUPELCOPORT™ 10% SP-2100 on 80/100 SUPELCOPORT GP 5% SP-1200/1.75% Bentone® 34 on 100/120 SUPELCOPORT GP 5% SP-1200/5% Bentone 34 on 100/120 SUPELCOPORT	11840-U 11989 12140 12134 12133
<b>Simulated Distillation By Packed Column GC</b>		<b>Packings (Other Quantities)</b> 80/100 Carbo-pack™ C/0.1% SP-1000 15g 7g	
<b>Petrocol A Column</b> 20" x 1/8" stainless steel	12445	60/80 Carbo-pack F-TA packing is available in packed columns only.	11820 11908
<b>Petrocol B Column</b> 20" x 1/8" stainless steel	12449	<b>C1-C6 Hydrocarbons by Packed Column GC</b> 23% SP-1700 on 80/100 Chromosorb P AW packing in a 30' x 1/8" stainless steel column, general configuration (6" coil, 8" arms) 12809-U	
<b>Simulated Distillation By Capillary GC</b>		Inquire about other configurations or lengths.	
<b>Petrocol 2887 Capillary Column</b> 5m x 0.53mm ID fused silica, 0.5µm film	25323	<b>Packings</b> 23% SP-1700 on 80/100 Chromosorb P AW (25g) (21-22g fills a 30" x 1/8" column) GP 80/100 Carbo-pack C/0.19% picric acid (15g)	
<b>Polynuclear Aromatic Hydrocarbons in Gasoline</b>		11865 11824	
<b>SPB-5 Fused Silica Capillary Column</b> 30m x 0.25mm ID, 0.25µm film	24034		
<b>Miscellaneous Petrochemicals Analyses</b>			
<b>Petrocol DH 50.2 Fused Silica Capillary Column</b> 50m x 0.20mm ID, 0.50µm film	24133-U		
<b>SPB-1 Capillary Columns</b> 60m x 0.75mm ID glass, 1.0µm film▲ 30m x 0.53mm ID fused silica, 0.50µm film			
<b>SPB-20 Fused Silica Capillary Column</b> 30m x 0.25mm ID, 0.25µm film	25315		
<b>Nukol Fused Silica Capillary Column</b> 30m x 0.53mm ID, 0.50µm film	24086		
<b>Petrocol DH Octyl Fused Silica Capillary Columns</b> 100m x 0.25mm ID, 0.5µm film	25327		
<b>Phenols / Polymers</b>			
<b>Nukol Fused Silica Capillary Column</b> 30m x 0.25mm ID, 0.25µm film	24282		
<b>SPB-5 Fused Silica Capillary Column</b> 30m x 0.53mm ID, 1.5µm film	24107		
	25305-U		

▲This glass column is no longer available. Request the same dimensions in a fused silica column (custom column).

## Low Molecular Weight Compounds by Gas-Solid Chromatography

### HayeSep Polymers

HayeSep porous polymers are thoroughly cleaned and conditioned before they are packaged, which ensures high consistency, and minimal shrinkage and bleed. Uses include the analysis of low molecular weight halogenated and sulfur-containing compounds, water, alcohols, glycols, free fatty acids, esters, ketones and aldehydes. These polymers also are used for separating hydrocarbons, ammonia, and other solvents. (Request HayeSep Application Bulletin available from Supelco).

#### HayeSep A

Copolymer of high-purity divinylbenzene and ethylene glycol dimethacrylate. Separates permanent gases at ambient temperatures. Use at higher temperatures to analyze C2 hydrocarbons, hydrogen sulfide, and water.

#### HayeSep B

Copolymer of divinylbenzene and polyethyleneimine. Separates C1 and C2 amines and trace levels of ammonia and water.

#### HayeSep C

USP code [S10]. Copolymer of divinylbenzene and acrylonitrile. For the analysis of polar molecules. Separation characteristics are similar to those of Chromosorb 104.

#### HayeSep D

High-purity divinylbenzene polymer. Exhibits superior separation characteristics for light gases, carbon monoxide, carbon dioxide, and acetylene (ahead of other C2 hydrocarbons). For the analysis of water and hydrogen sulfide.

#### HayeSep N

Copolymer of divinylbenzene and ethylene glycol dimethylacrylate.

#### HayeSep P

Copolymer of divinylbenzene and styrene.

#### HayeSep Q

USP code [S3]. Divinylbenzene polymer.

#### HayeSep R

USP code [S4]. Copolymer of divinylbenzene and N-vinyl-2-pyrrolidone.

#### HayeSep S

USP code [S8]. Copolymer of divinylbenzene and 4-vinylpyridine.

#### HayeSep T

Polymer of ethylene glycol dimethacrylate.

### Accessories for Analyzing High Boiling Hydrocarbons (page 9)

Description	Cat. No.
Cool On-Column Injection Sleeve Kit	23630
Hamilton® 701 Syringe With 6" fixed needle	21574
Polywax® 655 6 x 1mL	48482

### HayeSep® Porous Polymers (75cc bottles)

HayeSep Series	Mesh Size	Use	Max. Temp. (°C)	Cat. No.
HayeSep A	60/80	permanent gases	165	10282
	80/100			10283
	100/120			10284
HayeSep B	60/80	C1-C2 amines, ammonia in H <sub>2</sub> O	190	10285
	80/100			10286
	100/120			10287
HayeSep C	60/80	polar compounds	250	10288
	80/100			10289
	100/120			10290
HayeSep D	60/80	light gases, CO, CO <sub>2</sub>	290	10291
	80/100			10292
	100/120			10293
HayeSep DB*	60/80	light gases, CO, CO <sub>2</sub>	290	10297-U
	80/100			10280-U
	100/120			10281-U
HayeSep N	60/80	acetylene, ethylene	165	10294
	80/100			10295
	100/120			10296
HayeSep P	60/80	ammonia, alcohols in water	250	10297
	80/100			10298
	100/120			10299
HayeSep Q	60/80	hydrocarbons, sulfur gases, general	275	10300-U
	80/100			10301-U
	100/120			10302
HayeSep R	60/80	light hydrocarbons, chlorinated compounds	250	10303
	80/100			10304
	100/120			10305-U
HayeSep S	60/80	C2-C3 hydrocarbons polar compounds	250	10306
	80/100			10307
	100/120			10308
HayeSep T	60/80	light hydrocarbons formaldehyde	165	10309
	80/100			10310
	100/120			10311

### Polymers by Gel Permeation Chromatography

#### SUPELCO SIL Silica-Based HPLC Columns

Phase Particle, Pore Size	Column Length (cm)	Cartridge with Fittings		
		4.6mm ID Cartridge Cat. No.	4.0mm ID Cartridge Cat. No.	3.0mm ID Cartridge Cat. No.
<b>LC-1</b> 5µm, 120Å	5	58237		58237C30
	15	58210U	58210C40	58210C30
	25	58296		
<b>LC-Diol</b> 5µm, 120Å	25	58201	58201C40	58201C30

#### SUPELCO SIL Preparative Scale HPLC Columns

Phase Particle, Pore Size	Column Length (cm)	Cartridge with Fittings
		6.2mm ID Cartridge Cat. No.
<b>LC-Si</b> 5µm, 120Å	25	58963
<b>LC-1</b> 5µm, 100Å	25	58961
<b>LC-Diol</b> 5µm, 300Å	25	58967
<b>LC-3Si</b> 5µm, 300Å	25	58965
<b>LC-301</b> 5µm, 300Å	25	58966
<b>LC-3Diol</b> 5µm, 300Å	25	58968

\*Deactivated for basic compounds.



## P-N-A, P-O-N-A, P-I-A-N-O Analyses

Use these mixes to determine retention times and indices and monitor response factors for components of complex mixtures of hydrocarbons. These standards are complex mixes of known quantities of hydrocarbons, accurately prepared by weight to **three decimal places**. Typical values are listed here; exact values will differ slightly from lot to lot.

- Formulations are weight percent.
- Each mix includes a detailed data sheet listing components by weight percent, mole percent, and liquid volume percent, retention times and retention indices for each component, and other information.
- A chromatogram from a 100-meter Petrocol DH capillary column (including conditions) is provided.
- Products are supplied in crimp-top vials with hole caps and septa. The shelf life of the unopened, refrigerated mixes is 1 year.

## Quantitative Reference Standards

Mixes and Solutions	Composition	Qty.	Cat. No.
n-Paraffins Mix <sup>n</sup>	<i>11 n-paraffins (typical values shown)</i> n-Pentane, 11.393% n-Hexane, 10.963% n-Heptane, 12.239% n-Octane, 10.088% n-Nonane, 13.204% n-Decane, 6.703%	0.1mL	<b>44585-U</b>
Isoparaffins Mix <sup>n</sup>	<i>37 isoparaffins (typical values shown)</i> Isopentane, 2.553% 2,3-Dimethylbutane, 1.438% 2-Methylpentane, 4.356% 3-Methylpentane, 3.147% 2,2-Dimethylpentane, 4.189% 2,4-Dimethylpentane, 3.116% 2,2,3-Trimethylbutane, 3.149% 3,3-Dimethylpentane, 1.308% 2-Methylhexane, 3.551% 2,3-Dimethylpentane, 1.612% 3-Methylhexane, 4.456% 3-Ethylpentane, 1.600% 2,2-Dimethylhexane, 1.439% 2,5-Dimethylhexane, 3.195% 2,2,3-Trimethylpentane, 3.185% 2,4-Dimethylhexane, 1.574% 2,3-Dimethylhexane, 1.365% 2-Methylheptane, 4.641% 4-Methylheptane, 3.377%	0.1mL	<b>44586-U</b>
Naphthenes Mix <sup>n</sup>	<i>30 naphthenes (typical values shown)</i> Cyclopentane, 5.805% Methylcyclopentane, 8.728% Cyclohexane, 5.758% 1,1-Dimethylcyclopentane, 3.370% cis-1,3-Dimethylcyclopentane, 1.962% trans-1,3-Dimethylcyclopentane, 2.689% trans-1,2-Dimethylcyclopentane, 3.087% Methylcyclohexane, 8.992% Ethylcyclopentane, 1.773% ctc-1,2,4-Trimethylcyclopentane, 3.721% ctc-1,2,3-Trimethylcyclopentane, 3.419% cct-1,2,4-Trimethylcyclopentane, 1.438% trans-1,4-Dimethylcyclohexane, 3.713% 1-Ethyl-1-methylcyclopentane, 1.634% trans-1,2-Dimethylcyclohexane, 3.644%	0.1mL	<b>44588</b>
Olefins Mix <sup>n</sup>	<i>25 olefins (typical values shown)</i> 3-Methyl-1-butene, 1.478% 1-Pentene, 5.640% 2-Methyl-1-butene, 1.880% 2-Methyl-1,3-butadiene, 2.389% trans-2-Pentene, 4.024% cis-2-Pentene, 2.305% 4-Methylpentene-1 1-Hexene, 7.693% trans-2-Hexene, 1.744% 2-Methylpentene-2, 4.054% cis-2-Hexene, 1.825% 1-Heptene, 8.716% 1-Decene, 7.665%	0.1mL	<b>44589</b>


■ For information on a kit containing these products, see the following page (Cat. No. 44594-U).

Mixes and Solutions	Composition	Qty.	Cat. No.
Aromatics Mix	<p><i>37 aromatics (typical values shown)</i></p> <p>Benzene, 3.139% Toluene, 4.563% Ethylbenzene, 1.570% m-Xylene, 2.939% p-Xylene, 2.982% o-Xylene, 3.208% Isopropylbenzene, 3.127% n-Propylbenzene, 1.705% 1-Methyl-3-ethylbenzene, 2.903% 1-Methyl-4-ethylbenzene, 3.166% 1,3,5-Trimethylbenzene, 3.427% 1-Methyl-2-ethylbenzene, 4.460% 1,2,4-Trimethylbenzene, 4.315% tert-Butylbenzene, 2.667% Isobutylbenzene, 0.845% sec-Butylbenzene, 1.469% 1-Methyl-3-isopropylbenzene, 1.403% 1-Methyl-4-isopropylbenzene, 2.789% 1-Methyl-2-isopropylbenzene, 2.631%</p>	0.1mL	<b>44587</b>
P-I-A-N-O Mix	<p><i>139 n-paraffins, isoparaffins, aromatics, naphthenes, and olefins</i></p> <p>A quantitative mix of the components in the five mixes described above and on the previous page. Data sheets for each class of compounds list weight percent, mole percent, and other information for each component. Approximate weight percentages:</p> <p>n-Paraffins, 18.9% Isoparaffins, 18.8% Aromatics, 23.3%</p>	0.1mL	<b>44593-U</b>
P-I-A-N-O Kit	<p><i>One of each of the following mixes</i></p> <p>n-Paraffins Mix (44585-U) Isoparaffins Mix (44586-U) Naphthenes Mix (44588)</p>	<p>Olefins Mix (44589) Aromatics Mix (44587) P-I-A-N-O Mix (44593-U)</p>	<b>44594-U</b>

**Included with P-I-A-N-O Mix and P-I-A-N-O Kit,  
Data Disc Makes it Easier to Organize P-I-A-N-O Data**

The P-I-A-N-O Certadisc Certification Worksheet and Database can greatly simplify the task of organizing certification data. A worksheet (Lotus 1-2-3 format) contains complete certification data for the n-paraffins, isoparaffins, aromatics, naphthenes, and olefins mixes, and for the P-I-A-N-O mixture. Calculations are easily performed with a few keystrokes. dBase III PLUS programs, designed for use with the Lotus worksheet, include information on all 139 compounds. The programs enable you to perform searches based on retention times, making it easier to identify peaks.

3½", 1.44 MEG disc. The disc is identified by product lot number.



P000165

## ASTM D5134

### Detailed Analysis of Petroleum Naphthas Through n-Nonane

Mixes and Solutions	Composition	Qty.	Cat. No.
ASTM D5134 Qualitative Column Evaluation Mix	<i>Nominal 0.5-1.0% by weight in 2-methylpentane</i> Use this mix of hydrocarbons for assessing column performance. Toluene, 0.5%    2,3,3-Trimethylpentane, 1%    4-Methylheptane, 1% n-Heptane, 1%    2-Methylheptane, 1%    n-Octane, 1%	1 mL	<b>502103</b>
ASTM D5134 Splitter Linearity Check Mix	<i>Neat at 10% by weight, each component</i> Use this quantitative mix of hydrocarbons to determine proper system performance and detector response factors. n-Hexane    Benzene    2,4-Dimethylhexane n-Heptane    Toluene    2-Methylheptane n-Octane    2-Methylhexane    2,4-Dimethylheptane n-Nonane	500mg	<b>506753</b>
2,3,3-Trimethylpentane	<i>Neat</i>	500mg	<b>502081</b>

### Highly Characterized Reference Materials

The following standards, taken from refinery process streams and exhaustively analyzed by GC/FID and GC/MS, are recommended for evaluating process performance, identifying sources of contamination, PIANO analysis, method development, and training. Each comes with a comprehensive data packet containing quantitative and qualitative data and chromatograms using a 100-meter Petrocol DH column.

Description	Qty.	Cat. No.
Petroleum Refinery Reformate	1 mL	<b>47489</b>
Petroleum Refinery Pyrolysis Gasoline (Py Gas)	1 mL	<b>47490-U</b>
Petroleum Refinery Heavy Straight Run Naphtha	1 mL	<b>47488</b>

### Qualitative Reference Materials

These refinery reference materials are the actual materials used in the ASTM D5134 round-robin method validation stage. They are referred to in the method, and are used for establishing component retention times for identification purposes. Each sample is accompanied by a comprehensive data booklet containing an expanded detailed chromatogram from a Petrocol DH 50.2 column, with identified peaks.

Mixes and Solutions	Composition	Qty.	Cat. No.
ASTM D5134 Qualitative Reference Alkylate Standard	<i>Approximately 30 identified components. Neat fraction.</i> Actual refinery alkylation product	6 x 1 mL	<b>48267-U</b>
ASTM D5134 Qualitative Reference Reformate Standard	<i>Approximately 70 identified components. Neat fraction.</i> Actual refinery reformate product	6 x 1 mL	<b>48266</b>
ASTM D5134 Qualitative Reference Naphtha Standard	<i>Approximately 100 identified components. Neat fraction.</i> Actual refinery naphtha product	6 x 1 mL	<b>48265-U</b>
ASTM D5134 Qualitative Reference Refinery Standards Kit	<i>2 x 1 mL each of the three qualitative reference standards (alkylate, naphtha, and reformate)</i>		<b>48268</b>

We can custom-manufacture petroleum standards. Please inquire.

## ASTM D5442

**Analysis of Petroleum Waxes by GC** - Qualitative and quantitative mixes of n-paraffins used for determining column resolution, retention times, response factors, and internal standard preparation.

Mixes and Solutions	Composition	Qty.	Cat. No.
ASTM D5442 Column Resolution Test Mix	<i>0.05 weight percent each of nC20 and nC24 in cyclohexane</i> Use to determine column resolution.	6 x 1mL	<b>506729</b>
<b>Retention Time Standards</b>			
ASTM D5442 C16-C44 Retention Time Mix	<i>Neat qualitative n-paraffins, each at 8.3 weight percent</i> Use to establish retention times from C16 to C44 C16                      C24                      C32 C18                      C26                      C36 C20                      C28                      C40 C22                      C30                      C44	500mg	<b>502251</b>
ASTM D5442 C12-C60 Retention Time Mix	<i>Neat qualitative n-paraffins, each at 6.25 weight percent</i> Use to establish retention times from C12 to C60 C12                      C24                      C36 C14                      C26                      C40 C16                      C28                      C44 C18                      C30                      C50 C20                      C32                      C60 C22	500mg	<b>500623</b>
<b>Quantitative Linearity Standards</b>			
ASTM D5442 C16-C44 Linearity Standard	<i>n-Paraffins, each at 0.01 weight percent</i>  Accurately prepared mix. Use to determine response factors from C16 to C44. C16                      C24                      C32 C18                      C26                      C36 C20                      C28                      C40 C22                      C30                      C44	1mL 6 x 1mL	<b>502278</b> <b>502286</b>
ASTM D5442 C12-C60 Linearity Standard	<i>n-Paraffins, each at 0.01 weight percent</i>  Accurately prepared mix. Use to determine response factors from C12 to C60. C12                      C24                      C36 C14                      C26                      C40 C16                      C28                      C44 C18                      C30                      C50 C20                      C32                      C60 C22	1mL 6 x 1mL	<b>502243</b> <b>502235</b>
<b>Internal Standards</b>			
ASTM D5442 n-Hexadecane Internal Standard	<i>n-Hexadecane (nC16), neat or in solution, added to the internal standard.</i>	5mL	<b>506699</b>

## Qualitative n-Paraffin Mixes

For determining retention indices and retention times.

Chemical Name	Concentration	Qty.	Cat. No.
n-Paraffin Mix C5, C6, C7, C8	<i>Neat, varied concentrations</i>		<b>47100</b>
n-Paraffin Mix C7, C8, C9, C10	<i>Neat, varied concentrations</i>		<b>47101</b>
n-Paraffin Mix C10, C12, C14, C16	<i>Neat, varied concentrations</i>		<b>47102</b>
n-Paraffin Mix C18, C20, C22, C24	<i>2% (wt./wt.) each n-paraffin in n-octane</i>	5mL	<b>47108</b>
n-Paraffin Mix C22, C24, C28, C32	<i>2% (wt./wt.) each n-paraffin in n-octane</i>	5mL	<b>47106</b>
n-Paraffin Mix C24, C28, C32, C36	<i>2% (wt./wt.) each n-paraffin in n-octane</i>	5mL	<b>47107</b>

## High Molecular Weight Hydrocarbon Standards

For high temperature SIMDIS or GC analyses. Polywax materials are polyethylene waxes having average molecular weights of 500 and 655 Dalton, respectively. Ethylene oligomers range in carbon number from approximately C20 to C100+ and are useful for establishing retention times.

Description	Concentration	Qty.	Cat. No.
Pentacontane (nC50)	Neat	50mg	<b>48595</b>
Hexacontane (nC60)	Neat	50mg	<b>48893</b>
Polywax 500	Neat	5000mg	<b>48475</b>
Polywax 500	10,000µg/mL in <i>p</i> -xylene	6 x 1mL	<b>48480-U</b>
Polywax 655	Neat	5000mg	<b>48477</b>
Polywax 655	10mg/mL in <i>p</i> -xylene	6 x 1mL	<b>48482</b>

## ASTM D3710

**Boiling Range Distribution of Gasoline (500°F maximum)** - These qualitative and quantitative hydrocarbon blends are prepared according to ASTM recommendations. Calibration mixes are either volume/volume or weight/weight formulations as indicated. Nominal concentration of actual values will differ from lot to lot. Qualitative calibration references are prepared to approximate weight/weight ( $\pm 10\%$ ) specifications and, because of the presence of gases, are not intended for quantitative use. All calibration references are accompanied by a data sheet.

Mixes and Solutions	Composition	Qty.	Cat. No.
ASTM D3710 Qualitative Calibration Mix (w/w) indicated	<i>Each of the following components in the approximate proportions</i>  n-Propane, 1.5% 2-Methylpropane, 1.5% n-Butane, 4.5% 2-Methylbutane, 9.7% n-Pentane, 7.6% 2-Methylpentane, 5.4% n-Hexane, 5.4%  2,4-Dimethylpentane, 5.4% n-Heptane, 9.7% Toluene, 10.8% n-Octane, 5.4% p-Xylene, 13.0% n-Propylbenzene, 4.3%  n-Decane, 3.2% n-Butylbenzene, 3.2% n-Dodecane, 3.2% n-Tridecane, 2.2% n-Tetradecane, 2.2% n-Pentadecane, 2.2%	1 x 1mL 6 x 1mL	<b>50642</b> <b>48884</b>
ASTM D3710 Quantitative Calibration Mix	<i>Each of the following components in the proportions (v/v) indicated</i>  2-Methylbutane, 10.5% n-Pentane, 8.1% 2-Methylpentane, 5.8% n-Hexane, 5.8% 2,4-Dimethylpentane, 5.8% n-Heptane, 10.5%  Toluene, 11.6% n-Octane, 5.8% p-Xylene, 14.0% n-Propylbenzene, 4.7% n-Decane, 3.5%  n-Butylbenzene, 3.5% n-Dodecane, 3.5% n-Tridecane, 2.3% n-Tetradecane, 2.3% n-Pentadecane, 2.3%	1 x 1mL 6 x 1mL	<b>506435</b> <b>48879</b>
ASTM D2887/D5307 Column Resolution Test Mix	<i>Each of the following components at 1% (w/v) in n-octane</i> n-Hexadecane n-Octadecane	6 x 1mL	<b>48889</b>
ASTM D2887 Quantitative Calibration Solution	<i>Prepared in carbon disulfide at 0.5 weight percent each component (except for nC16 and nC18, which are at 1 weight percent)</i>  Use for assessing column resolution as well as for quantitative analyses. n-Pentane n-Hexane n-Heptane n-Octane n-Nonane n-Decane n-Undecane  n-Dodecane n-Tetradecane n-Pentadecane n-Hexadecane n-Heptadecane n-Octadecane n-Eicosane  n-Tetracosane n-Octacosane n-Dotriacontane n-Hexatriacontane n-Tetracontane n-Tetratetracontane	1 x 1mL 6 x 1mL	<b>500631</b> <b>500658</b>

Mixes and Solutions	Qty.	Cat. No.
ASTM D2887 Reference Gas Oil Sample, Lot 2 <sup>■</sup>	1mL 6 x 1mL	<b>506419</b> <b>48873</b>

<sup>■</sup>This sample is a petroleum fraction having a boiling range from 250°F to 850°F, evaluated in round-robin studies by the ASTM. Use this sample and the supplied ASTM boiling range consensus values to evaluate system performance.



ASTM D5769

**Aromatics in Gasoline by GC/MS**

Mixes and Solutions	Composition	Qty.	Cat. No.
ASTM D5769/EPA Aromatics Internal Standard Mix	<i>Each of the following components at the weight percent indicated</i> Benzene-d <sub>6</sub> , 40.0      Ethylbenzene-d <sub>10</sub> , 40.0      Naphthalene-d <sub>8</sub> , 20.0	1 x 5.0mL	<b>47327</b>

ASTM D5599

**Oxygenates in Gasoline by GC and OFID Detection**

Mixes and Solutions	Composition	Qty.	Cat. No.
1,2-Dimethoxyethane (DME)	<i>Neat, Internal Standard</i>	5mL	<b>47214</b>

ASTM D5441

**Purity of Methyl tert-butyl ether (MTBE) by GC**

Mixes and Solutions	Composition	Qty.	Cat. No.
---------------------	-------------	------	----------

**MTBE Quantitative Solutions and Neat Materials**

	For analyses of oxygenates in gasoline.		
tert-Amyl methyl ether	<i>2000µg/mL in methanol</i>	1mL	<b>506737</b>
Methyl tert-butyl ether (MTBE)	<i>Neat</i>	1000mg	<b>48027</b>
Methyl tert-butyl ether (MTBE)	<i>2000µg/mL in methanol</i>	1mL	<b>48483</b>

**MTBE Contaminant Standards**

ASTM D5441 MTBE Contaminants (High) Mix A	<i>Each component at 1% w/w in MTBE</i> Use to identify and determine response factors for the common impurities in MTBE. tert-Amyl methyl ether      2-Methylbutane      trans-2-Pentene tert-Butanol      2-Methyl-2-butene      Triisobutylene (mixed isomers) tert-Butyl ethyl ether      Pentane      2,4,4-Trimethyl-1-pentene Methanol      cis-2-Pentene	1mL	<b>47942</b>
ASTM D5441 MTBE Contaminants (Low) Mix B	<i>Each component at 0.1% w/w in MTBE</i> Use to identify and determine response factors for the common impurities in MTBE. tert-Amyl methyl ether      2-Methylbutane      trans-2-Pentene tert-Butanol      2-Methyl-2-butene      Triisobutylene (mixed isomers) tert-Butyl ethyl ether      Pentane      2,4,4-Trimethyl-1-pentene Methanol      cis-2-Pentene	1mL	<b>47943</b>

## ASTM D4815

### Determination of Oxygenates (Ethers and Alcohols) in Gasoline by GC

Along with valve timing and peak identification mixes, two sets of quantitative calibration mixes are available, both with and without the method-specified internal standard. The internal standard-free mixes are packaged in quantities of approximately 9.5mL to allow the addition of the internal standard. All components used in preparing these standards have been analyzed for purity, water content, and the presence of other method components to 0.01%. Blends have been prepared using precise gravimetric techniques exceeding the requirements of ASTM Method D4815. Corrections are made for common impurities. All calibration blends are provided with a chromatogram and data verifying the purity and identity of the raw materials.

Mixes and Solutions	Composition	Qty.	Cat. No.
D4815 Quantitative Calibration Mix 1 (Without Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 3.00 tert-Butanol, 0.10 Methyl tert-butyl ether, 20.00	9.5mL	<b>47200</b>
D4815 Quantitative Calibration Mix 2 (Without Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 0.10 tert-Butanol, 3.00 Methyl tert-butyl ether, 15.00	9.5mL	<b>47201</b>
D4815 Quantitative Calibration Mix 3 (Without Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 6.00 tert-Butanol, 6.00 Methyl tert-butyl ether, 10.00	9.5mL	<b>47202</b>
D4815 Quantitative Calibration Mix 4 (Without Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 9.00 tert-Butanol, 8.00 Methyl tert-butyl ether, 5.00	9.5mL	<b>47203</b>
D4815 Quantitative Calibration Mix 5 (Without Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 12.00 tert-Butanol, 12.00 Methyl tert-butyl ether, 0.10	9.5mL	<b>47204</b>
1,2-Dimethoxyethane (DME)	<i>Neat, internal standard</i>	5mL	<b>47214</b>
ASTM D4815 Calibration Curve Verification Kit	<i>Each of the following products</i> Quantitative Calibration Mix 1 (47200) Quantitative Calibration Mix 2 (47201) Quantitative Calibration Mix 3 (47202)	Quantitative Calibration Mix 4 (47203) Quantitative Calibration Mix 5 (47204) 1,2-Dimethoxyethane (47214)	<b>47210</b>
D4815 Quantitative Calibration Mix 1 <sup>■</sup> (With Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 12.85 tert-Butanol, 0.095 Methyl tert-butyl ether, 19.00	1mL	<b>47205</b>
D4815 Quantitative Calibration Mix 2 <sup>■</sup> (With Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 0.095 tert-Butanol, 2.85 Methyl tert-butyl ether, 14.25	1mL	<b>47206</b>
D4815 Quantitative Calibration Mix 3 <sup>■</sup> (With Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 5.70 tert-Butanol, 5.70 Methyl tert-butyl ether, 9.50	1mL	<b>47207</b>

■ For information on a kit containing this product, see the following page (Cat. No. 47211).



Mixes and Solutions	Composition	Qty.	Cat. No.
D4815 Quantitative Calibration Mix 4 (With Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 8.55 tert-Butanol, 7.60 Methyl tert-butyl ether, 4.75	1 mL	<b>47208</b>
	tert-Pentanol, 3.56 isooctane:xylene (65:35), 70.54 DME (internal standard), 5.00		
D4815 Quantitative Calibration Mix 5 (With Internal Standard)	<i>Each component at the nominal weight percent indicated</i> Ethanol, 11.40 tert-Butanol, 11.40 Methyl tert-butyl ether, 0.095	1 mL	<b>47209</b>
	tert-Pentanol, 0.095 isooctane:xylene (65:35), 72.01 DME (internal standard), 5.00		
ASTM D4815 Valve Timing Mix	<i>Each component at the nominal weight percent indicated</i> Use for determining the proper backflush valve timing. Methylcyclopentane, 10.0 Diisopropyl ether, 10.0 Ethyl tert-butyl ether, 10.0	1 mL	<b>47212</b>
	Methyl tert-butyl ether, 10.0 n-Hexane, 60.0		
ASTM D4815 Qualitative Peak ID Mix	<i>Each component at the nominal weight percent indicated</i> Use for determining relative retention times of selected oxygenates. Methylcyclopentane, 4.0 Methanol, 7.3 Ethanol, 7.3 Isopropanol, 7.3 tert-Butanol, 7.3 n-Propanol, 7.3 Methyl tert-butyl ether, 4.0 sec-Butanol, 7.3	1 mL	<b>47213</b>
	Diisopropyl ether, 4.0 Isobutanol, 7.3 Ethyl tert-butyl ether, 4.0 tert-Pentanol, 7.3 1,2-Dimethoxyethane (internal standard), 6.0 n-Butanol, 7.3 Benzene, 5.0 tert-Amyl methyl ether, 7.3		
ASTM D4815 Quantitative Calibration Kit	<i>Each of the following products</i> Quantitative Calibration Mix 1 (47205) Quantitative Calibration Mix 2 (47206) Quantitative Calibration Mix 3 (47207) Quantitative Calibration Mix 4 (47208)		<b>47211</b>
	Quantitative Calibration Mix 5 (47209) Valve Timing Mix (47212) Qualitative Peak ID Mix (47213)		

## USP Solvent Standards

Described in US Pharmacopoeia methods for determining residual organic solvents in pharmaceutical preparations.

Description		Qty.	Cat. No.
International USP 467 Mix	<i>8 analytes at concentrations indicated in methanol</i> Acetonitrile (500µg/mL) Benzene (1000µg/mL) Chloroform (500µg/mL)	1 mL	<b>47632-U</b>
	1,2-Dichloroethane (1000µg/mL) Methylene chloride (5000µg/mL) 1,4-Dioxane (1000µg/mL) Pyridine (1000µg/mL) Trichloroethene (1000µg/mL)		
USP Organic Volatile Impurities Mix	<i>5 analytes at concentrations indicated in dimethyl sulfoxide</i> Benzene (1000µg/mL) Chloroform (500µg/mL)	4 x 1 mL	<b>47401</b>
	1,4-Dioxane (1000µg/mL) Methylene chloride (1000µg/mL) Trichloroethene (1000µg/mL)		
USP Modified Organic Volatile Impurities Mix	<i>5 analytes at concentrations indicated in methanol**</i> Benzene (1000µg/mL) Chloroform (500µg/mL) (5000µg/mL)	1 mL	<b>47398</b>
	1,4-Dioxane (1000µg/mL) Methylene chloride (1000µg/mL) Trichloroethene (1000µg/mL)		
USP 467 OVI Mix (High Concentration), 23rd Ed.	<i>5 analytes at indicated concentrations in DMSO</i> Benzene, 1000µg/mL Chloroform, 500µg/mL	1 mL 4 x 1 mL	<b>47538-U</b> <b>47539-U</b>
	1,4-Dioxane, 1000µg/mL Methylene chloride, 5000µg/mL Trichloroethylene, 1000µg/mL		
USP 467 OVI Mix (Low Concentration), 23rd Ed.	<i>5 analytes at indicated concentrations in DMSO</i> Benzene, 100µg/mL Chloroform, 50µg/mL	1 mL 4 x 1 mL	<b>47540-U</b> <b>47541-U</b>
	1,4-Dioxane, 100µg/mL Methylene chloride, 500µg/mL Trichloroethylene, 100µg/mL		

\*\*Alternate solutions are available. Please inquire.

## Transportable Gases

We offer an expanded line of pure gases and gas mixtures in convenient transportable cylinders, manufactured for Supelco by Scott Specialty Gases. Scott's popular SCOTTY® products are available in a variety of sizes and cylinder capacities. The portability of SCOTTY gases makes these products suitable for many types of applications, both in the laboratory and in the field.

Supelco warrants that the Scott calibration gas products listed below meet the analytical specifications for the period of time stated on the cylinder and/or the Certificate of Analysis.

- No cylinder demurrage or rental
- Long-term guaranteed stability for extended shelf life
- Sufficient product in every cylinder for many applications
- Save on shipping costs
- "Go anywhere" portability makes calibration jobs easier
- Quality accessory equipment available
- High accuracy increases dependability
- A large number of calibration gases are available

## Scott Specialty Gases

Shelf life 1 year, unless otherwise noted.

## Specifications

### SCOTTY I

Contents: 4 liters  
 Pressure: 120psig  
 Outlet Fitting: Aerosol-type push button with applicator tube  
 Weight: ~100g  
 Dimensions: 2.5 x 8 in.  
 D.O.T. Specs: 2Q

### SCOTTY II

Contents: 14 liters  
 Pressure: 240psig  
 Outlet Fitting: CGA-160-1/8" NPT F  
 Weight: 1.5 lb.  
 Dimensions: 3 x 11 in.  
 D.O.T. Specs: 4B240

### SCOTTY IV

Contents: 48 liters  
 Pressure: 300psig  
 Outlet Fitting: CGA-165  
 Weight: 1.75 lb.  
 Dimensions: 4 x 16.25 in.  
 D.O.T. Specs: 39 NRC

### SCOTTY HP

Contents: 104 liters  
 Pressure: 1800psig  
 Outlet Fitting: CGA-180  
 Weight: 2.2 lb.  
 Dimensions: 3.25 x 12.25 in.  
 D.O.T. Specs: 3AL1800

Gas Composition and Concentration	SCOTTY I Cat. No.	SCOTTY II Cat. No.	SCOTTY IV Cat. No.
<b>Pure Gases</b>			
Air, zero (THC < 1ppm)	501212	501220	501239
Argon 99.995%	501247	501255	_____
Carbon dioxide 99.8%*	501271	23402	501298
Hydrogen 99.99%	_____	300100	_____
Methane 99.0%	_____	22562	_____
Oxygen 99.6%	_____	300500	_____
<b>Two-Component Mixtures</b>			
Benzene in air (1ppm)	_____	_____	303402-U
Benzene in air (10ppm)	_____	_____	_____
Benzene in air (100ppm)	_____	_____	303404
1,3-Butadiene in nitrogen (10ppm)	_____	303405	303406
Carbon dioxide in helium (100ppm)*■	_____	308200	_____
Carbon dioxide in nitrogen (100ppm)	_____	308300	501301
Carbon dioxide in nitrogen (1000ppm)	_____	501336	501344
Chlorine in nitrogen (10ppm)** See page 439.	_____	_____	_____
Ethylene in air (10ppm)*	_____	501379	501387
Ethylene in helium (100ppm)	_____	22572	_____
Hydrogen in helium (100ppm)■	_____	301200	_____
Hydrogen in nitrogen (1%)	501409	501417	501425
Hydrogen in nitrogen (100ppm)	_____	301300	501433
Methane in helium (100ppm)■	501441	307200	501468
Methane in nitrogen (100ppm)■	_____	307300-U	_____
Methane in nitrogen (1%)■	501476	23443	501484
Nitrogen in helium (100ppm)	_____	303200	_____
Nitrous oxide in nitrogen (1ppm)■	_____	501514	501522

Gas Composition and Concentration	SCOTTY I Cat. No.	SCOTTY II Cat. No.	SCOTTY IV Cat. No.
<b>Two-Component Mixtures (contd.)</b>			
Oxygen in helium (100ppm)■	————	305200	————
Oxygen in nitrogen (2%)■	501530	501549	501557
Oxygen in nitrogen (6%)■	501565	501573	501581
1,1,1-Trichloroethane in nitrogen (10ppm)	————	————	303408
Trichloroethylene in nitrogen (10ppm)	————	303400	303401
Vinyl chloride in nitrogen (1ppm)	————	22554	501603
Vinyl chloride in nitrogen (10ppm)	————	22553	501611
Vinyl chloride in nitrogen (50ppm)	————	22555-U	————
Vinyl chloride in nitrogen (100ppm)	————	22552	————
Vinyl chloride in nitrogen (1000ppm)	————	22556	————
<b>Three-Component Mixtures</b>			
Carbon Dioxide and oxygen in nitrogen (1% / 20%)	————	23441	501638
<b>Multi-Component Mixtures</b>			
Carbon monoxide, carbon dioxide, hydrogen and oxygen each at 0.5% in nitrogen	501646	23438	501654
Carbon monoxide, carbon dioxide, methane, ethane, ethylene and acetylene each at 1% in nitrogen	501662	23437	23462
Methane, carbon monoxide, carbon dioxide, hydrogen and oxygen each at 1% in nitrogen	501670	22561	23463
Carbon monoxide, carbon dioxide, nitrogen, and oxygen each at 5%, methane and hydrogen each at 4% in helium	————	501697	————
Carbon monoxide at 7%, carbon dioxide at 15% and oxygen at 5% in nitrogen	501719	23442	————
Carbon monoxide and oxygen each at 7%, carbon dioxide at 15% and methane at 4.5% in nitrogen	————	501743	501751
Branched paraffins, each at 15ppm: iso-butane, 2-methylbutane, 2,2-dimethylpropane, 2-methylpentane, 3-methylpentane, 2,2-dimethylbutane in nitrogen	22501	23445	————
C1-C6 n-Paraffins, each at 15ppm: methane, ethane, propane, butane, pentane and hexane in nitrogen	501778	23444	501786
C1-C6 n-Paraffins, each at 100ppm: methane, ethane, propane, butane, pentane, hexane in nitrogen	501840	330300	501859
C1-C6 n-Paraffins, each at 100ppm: methane, ethane, propane, butane, pentane, hexane in helium	501794	330200	501808
C1-C6 n-Paraffins, each at 1000ppm: methane, ethane, propane, butane, pentane, hexane in helium	501816	501824	501832
C2-C4 Alkynes, each at 15ppm: acetylene, propyne, 1-butyne, 2-butyne in nitrogen	22508	————	————
C2-C6 Olefins, each at 100ppm: ethylene, propylene, 1-butene, 1-pentene, 1-hexene in nitrogen	————	332300-U	501875
C2-C6 Olefins, each at 100ppm: ethylene, propylene, 1-butene, 1-pentene, 1-hexene in helium	————	332200	501867
Methane, ethane, ethylene, acetylene, propane, propylene, propyne, n-butane each at 15ppm in nitrogen	22566	————	23470-U
n-Butane, iso-butane, cis-2-butene, trans-2-butene, 1-butene, iso-butylene, 1,3-butadiene, ethyl acetylene each at 15ppm in nitrogen	22567	————	23471
BTEX Mix**: benzene, ethylbenzene, toluene, m-xylene, o-xylene, p-xylene, each at 10ppm in nitrogen	————	————	501883
<hr/>			
<b>Gas Composition and Concentration</b>	<b>SCOTTY HP Cat. No.</b>		
Chlorine in nitrogen (10ppm)**	501352		

■ Custom concentrations available. Please inquire.

\*Pressure restricted due to either vapor pressure restriction or flammable oxidizer restriction.

\*\*Shelf life 6 months.



P000152

## Natural Gas Reference Standards

Prepared gravimetrically with weights traceable to the National Institute of Standards and Technology, then verified by analysis. In 14-liter SCOTTY II cylinders. Shelf life: 1 year.

Component (mole percent)	GPA Standard	Calorimetric Standard	High Ethane Standard	Helium-Enriched Standard
Helium	0.50	–	–	2.00
Nitrogen	5.00	2.50	9.00	1.60
Carbon dioxide	1.00	3.00	0.50	0.20
Methane	70.50	88.73	64.00	88.90
Ethane	9.00	3.50	12.50	3.00
Propane	6.00	1.00	7.00	1.70
Isobutane	3.00	0.40	3.00	1.00
n-Butane	3.00	0.40	3.00	1.00
Isopentane	1.00	0.15	0.50	0.30
n-Pentane	1.00	0.15	0.50	0.30
Neopentane	–	0.10	–	–
n-Hexane	–	0.05	–	–
n-Heptane	–	0.02	–	–

Standard	BTU	Qty.	Cat. No.
GPA	1298	785g	<b>303100-U</b>
Calorimetric	1028	790g	<b>303101</b>
High Ethane	1500	763g	<b>303102</b>
Helium-Enriched	1083	774g	<b>303103</b>

### Model 24 SCOTTY Single-Stage, General Purpose Regulator for SCOTTY II and IV Cylinders

Designed for noncorrosive service.

- Brass body with acetal resin bonnet/Viton diaphragm
- Tamper resistant, locking control knob
- Inlet connection 1/4 inch AN flare (CGA-165 or CGA-160)
- Maximum inlet pressure 300psig
- Delivery pressure range 1-60psig, can be preset
- Miniature

Description	Cat. No.
<b>Model 24 Single-Stage Regulator</b> with CGA-160 (for SCOTTY II)	<b>507911</b>
with CGA-165 (for SCOTTY IV)	<b>501395</b>



507911

P000154



501395

P000154

### Syringe Adapter for SCOTTY II and IV Cylinders•

Withdraw calibration gas into a syringe, through a silicone rubber septum. A vent at the septum permits purging prior to filling the syringe. Constructed of chromium-plated brass. Maximum pressure 240psi (16.9kg/cm<sup>2</sup>). 1/8" NPT male fitting.

Description	Cat. No.
Syringe Adapter	<b>609010</b>
Additional Septa (pk. of 10)	<b>608010</b>



P000160

### Pressure Regulator for SCOTTY II and IV Cylinders•

Reliable pressure regulation to 1psig (0.07kg/cm<sup>2</sup>), indicated on a 0-60psi (4.2kg/cm<sup>2</sup>) gauge. Easily connects to SCOTTY cylinders with the 1/8" NPT connector provided. You can also attach the syringe adapter (Cat. No. 609010) to the regulator for low pressure sample removal. Aluminum body with acetyl resin bonnet. Maximum inlet pressure: 400psi (28.1kg/cm<sup>2</sup>).

•**Note:** For safety reasons, this regulator does not fit higher pressure SCOTTY IV-EL cylinders.

Description	Cat. No.
Regulator	<b>513010</b>



P000155


### Stand for SCOTTY IV Cylinders

Stabilizes your cylinder on a benchtop or other surface.

Description	Cat. No.
Stand	<b>500410</b>



P000158




### Using SCOTTY IV Cylinders

SCOTTY IV gas cylinders should be used with a pressure regulator and SCOTTY syringe adapter. SCOTTY IV cylinders are fitted with two ports, one of which is a 1/8-inch NPT port, as shown at left. Fit the pressure regulator to the 1/8-inch NPT port, then connect the syringe adapter to the regulator.

**Custom Gases** - We can provide custom gas blends to your specifications. To inquire, contact our Technical Service department.

## Environmental Standards

Mixes and Solutions	Composition	Qty.	Cat. No.	
<b>BTEX Standards*</b>				
UST BTEX Mix 	6 analytes, 200µg/mL each in methanol		1 mL	<b>48026</b>
	Benzene	Toluene	1 mL	<b>458026</b>
HC BTEX Mix	Ethylbenzene	m-Xylene		
	6 analytes, 2000µg/mL each in methanol		1 mL	<b>47993</b>
HC BTEX/MTBE Mix	Benzene	Toluene		
	Ethylbenzene	m-Xylene		
	7 analytes, 2000µg/mL each in methanol		1 mL	<b>47505-U</b>
	Benzene	Toluene		
	Ethylbenzene	m-Xylene		
	Methyl tert-butyl ether			

## Gasoline Range Organics (GRO) Standards

Gasoline Additives Mix	4 analytes, 200µg/mL each in methanol		1 mL	<b>47905</b>
	1,2-Dibromoethane	1,2-Dichloroethane		
	Dibromomethane	Methyl tert-butyl ether		
Gasoline Range Organics (GRO) Mix	9 analytes, 2000µg/mL each in methanol		1 mL	<b>47576-U</b>
	Benzene	Naphthalene		
	Ethylbenzene	Toluene		
	3-Methylpentane	1,2,4-Trimethylbenzene		
PVOC Mix	8 analytes, 2000µg/mL each in methanol		1 mL	<b>47916</b>
	Benzene	Toluene		
	Ethylbenzene	1,2,4-Trimethylbenzene		
	Methyl tert-butyl ether	1,3,5-Trimethylbenzene		
EPA Gasoline Range Organics (GRO) Mix	9 analytes at indicated concentrations in methanol		1 mL	<b>47577-U</b>
	Benzene, 500µg/mL	2-Methylpentane, 1500µg/mL		
	Ethylbenzene, 500µg/mL	Toluene, 1500µg/mL		
	Heptane, 500µg/mL	1,2,4-Trimethylbenzene, 1000µg/mL		
		2,2,4-Trimethylpentane, 1500µg/mL		
		m-Xylene, 1000µg/mL		
		o-Xylene, 1000µg/mL		
UST Modified Gasoline Range Organics (GRO)	10 analytes, 1000µg/mL each in methanol		1 mL	<b>48167</b>
	Benzene	Toluene		
	Ethylbenzene	1,2,4-Trimethylbenzene		
	Methyl tert-butyl ether	1,3,5-Trimethylbenzene		
	Naphthalene			
Revised PVOC/GRO Mix	10 analytes, 2000µg/mL each in methanol		1 mL	<b>47578-U</b>
	Benzene	Toluene		
	Ethylbenzene	1,2,4-Trimethylbenzene		
	Methyl tert-butyl ether	1,3,5-Trimethylbenzene		
	Naphthalene			

Mixes and Solutions	Composition	Qty.	Cat. No.
<b>Surrogate Standards</b>			
α,α,α-Trifluorotoluene	10,000µg/mL in methanol	1 mL	<b>47582-U</b>
1-Chlorooctane	10,000µg/mL in methanol	1 mL	<b>47583-U</b>
<b>Internal Standard</b>			
1-Chloro-4-fluorobenzene	5000µg/mL in methanol	1 mL	<b>48194</b>
<b>Diesel Range Organics (DRO) Standards</b>			
UST Modified Diesel Range Organics (DRO)	10 analytes, 1000µg/mL each in hexane	1 mL	<b>48166</b>
	n-Decane	n-Octadecane	
	n-Dodecane	n-Eicosane	
	n-Tetradecane	n-Docosane	
	n-Hexadecane		
		n-Tetracosane	
		n-Hexacosane	
		n-Octacosane	

\*BTEX standards also are available in gas.

 Separate Source standards. See our current catalog.

Mixes and Solutions	Composition	Qty.	Cat. No.
<b>Internal Standard or Surrogate Standard for Diesel</b>			
1-Chlorooctadecane	10,000µg/mL in methylene chloride	1mL	47584-U
<b>Surrogate Standard for Diesel</b>			
2-Fluorobiphenyl	10,000µg/mL in methylene chloride	1mL	47581-U
<b>Surrogate Standard</b>			
o-Terphenyl	2000µg/mL in acetone	1mL	48169
	10,000µg/mL in methylene chloride	1mL	47580-U
<b>Internal Standard for DRO</b>			
5-α-Androstane	2000µg/mL in methylene chloride	1mL	48168
<b>Fuel Standards</b>			
Aviation Gasoline <sup>■</sup>	20,000µg/mL in methanol	1mL	47531-U
Gasoline <sup>■</sup>	20,000µg/mL in methanol	1mL	47516-U
Jet (Turbine) Fuel <sup>■</sup>	20,000µg/mL in methanol	1mL	47533-U
JP-4 Military Fuel Standard	10,000µg/mL in methylene chloride	1mL	47585-U
JP-5 Military Fuel Standard	10,000µg/mL in methylene chloride	1mL	47586-U
JP-8 Military Fuel Standard	10,000µg/mL in methylene chloride	1mL	47587-U
Kerosene Reference Standard	50,000µg/mL in hexane	1mL	47517-U
No. 1 Fuel Oil <sup>■</sup>	20,000µg/mL in methanol	1mL	47518-U
No. 2 Fuel Oil <sup>■</sup>	20,000µg/mL in methanol	1mL	47515-U
No. 3 Fuel Oil <sup>■</sup>	50,000µg/mL in hexane	1mL	47534-U
No. 4 Fuel Oil <sup>■</sup>	50,000µg/mL in hexane	1mL	47535-U
No. 6 Fuel Oil	20,000µg/mL in hexane:chloroform (50:50)	1mL	47536-U
<b>State-Specific Petroleum Method Standards</b>			
<b>Massachusetts</b>			
EPA Matrix	31 analytes, 25µg/mL each in hexane	1mL	46859-U
Spike Standard Mix	n-Nonane n-Decane n-Dodecane n-Tetradecane n-Hexadecane n-Octadecane n-Nonadecane n-Eicosane n-Docosane n-Tetracosane n-Hexacosane n-Octacosane n-Triacontane n-Hexatriacontane Acenaphthene Acenaphthylene Anthracene Benzo(a)anthracene Benzo(b)fluoranthene Benzo(k)fluoranthene Benzo(ghi)perylene Benzo(a)pyrene Chrysene Dibenz(a,h)anthracene Fluorene Fluoranthene Indeno(1,2,3-cd)pyrene 2-Methylnaphthalene Naphthalene Pyrene Phenanthrene		
<b>Florida</b>			
n-Hydrocarbons Mix	17 analytes, 1000µg/mL each in methylene chloride:carbon disulfide (50:50)	1mL	46855-U
	n-Octane n-Decane n-Dodecane n-Tetradecane n-Hexadecane n-Octadecane n-Eicosane n-Docosane n-Tetracosane n-Hexacosane n-Octacosane n-Triacontane n-Dotriacontane n-Tetracontane n-Hexatriacontane n-Octatriacontane n-Tetracontane		

<sup>■</sup>Refinery grade.

#### Trademarks

Apiezon – Biddle Instruments  
Bentone – National Lead Co., Baroid Sales Div.  
Carbowax – Union Carbide Corp.  
Cellosolve – Union Carbide Corp.  
Chromosorb – Celite Corp.  
dBase III Plus – Ashton-Tate Co.  
Dexsil – Dexsil Chemical Corp.  
Hamilton – Hamilton Co.  
HayeSep – Hayes Separations Inc.  
Lotus 1-2-3 – Lotus Development Co.

OV – Ohio Valley Specialty Chemical Co.  
Polywax – Petrolite Specialty Polymers Group  
Oxyblend, SCOTTY – Scott Specialty Gases, Inc.  
UCON – Union Carbide Corp.  
Viton – Dupont Dow Elastomers  
Carbopack, Carboxen, ChromFax, Nukol, Petrocol, Progel,  
PTE, SP, SPB, Supelco, SUPELCOGEL, SUPELCOPORT,  
SUPELCOSIL, SUPELCOWAX,  
Supel-Q, VOCOL – Sigma-Aldrich Co.  
Fused silica columns manufactured under HP US Pat. No.  
4,293,415.

## Custom Chemical Standards

We can prepare chemical standards for your environmental, food and beverage, petroleum, chemical, or pharmaceutical applications, using carefully monitored raw materials and solvents. We ship most custom products within 7-10 working days of receiving the order.\* Every custom standard includes a material safety data sheet (MSDS) and a certificate of composition; if testing is requested, we can provide a certificate of analysis.

### Quality Control

We offer three forms of quality control for custom-prepared standards:

#### **Gravimetric assurance (routine, no charge)**

Components are guaranteed to be within  $\pm 0.5\%$  of the amounts you specify. A Certificate of Composition is included.

#### **Qualitative testing (optional)\***

Typically the method will be a verification of the GC or LC elution and of the number of components. We will provide a copy of the chromatogram.

#### **Quantitative testing (optional)\***

Typically the standard will be compared to a separately prepared reference batch of the same formulation, using the internal standard method (GC), or the external standard method (LC). We will provide a copy of the chromatogram and a Quantitative Certificate of Analysis.

### Packaging

Typical packaging is 1mL, 2mL, 5mL, or 10mL flame-sealed amber ampuls. Other packaging is available (see below). Minimum order volume is 3mL, but this can be packaged in three 1mL ampuls or in a single ampul.

### To Obtain a Quote

**Telephone 814-359-5419 or 814-359-5752, or Fax 814-359-5750**

Within 24 hours we will confirm receipt of your request and provide a quote number.

## Materials Packaging Service

We can accept any materials repackaging job, large or small, simple or difficult. Our facilities and equipment enable us to package hygroscopic, light-sensitive, or other delicate materials under suitable conditions and into inert atmospheres. Our innovative ampul filling procedure eliminates product burning during the sealing process. We have experience repackaging biocides under ISO and cGMP requirements.

### Bulk Liquids or Solids

We can convert any volume of bulk liquid into smaller volumes, accurately and precisely dispensing your material into clear or amber glass ampuls, from 0.2mL to 20mL. We can weigh solids into glass vials or ampuls, with accuracy and precision to meet your most demanding needs.

### Labeling and Certificates

You choose the information to be included on your labels: ingredients, volume, company identification, etc. Labels can be sequentially numbered for inventory control or other needs. We can tailor certificates to your requirements—just tell us what you need.

### Testing

Testing can be performed by you or by Supelco (see above). We will provide evaluation samples before we ship your product.

### Storage and Distribution

We can inventory your product, then deliver it in a single shipment, or according to your schedule. We have more than 30 years of experience in shipping chemicals, domestically and internationally. We can ship anywhere, using the best and most cost-effective means.

### To Obtain a Quote

**Telephone 814-359-5450, Fax 814-359-5750, or email [custompackaging@sial.com](mailto:custompackaging@sial.com)**

Within 24 hours we will confirm receipt of your request and provide a quote number.

\*Qualitative and/or quantitative testing will add about one week to the delivery time for the product.

[sigma-aldrich.com/supelco](http://sigma-aldrich.com/supelco)

**Order/Customer Service** 800-247-6628, 800-325-3010 (US only) • Fax 800-325-5052 (US only) • E-mail [supelco@sial.com](mailto:supelco@sial.com)  
**Technical Service** 800-359-3041 (US only), 814-359-3041 • Fax 800-359-3044 (US only), 814-359-5468 • E-mail [techservice@sial.com](mailto:techservice@sial.com)  
**SUPELCO** • 595 North Harrison Road, Bellefonte, PA 16823-0048 • 814-359-3441

We are committed to the success of our Customers, Employees and Shareholders through leadership in **Life Science, High Technology and Service.**  
The **SIGMA-ALDRICH** Family  **SIGMA**  **ALDRICH**  **SUPELCO**  **ISOTEC**

©2004 Sigma-Aldrich Co. Printed in USA. Supelco brand products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.



T109858E  
AYD



**SIGMA-ALDRICH**