

# Advances in hardware design and function of the new timsOmni MS platform

**D. Papanastasiou<sup>1</sup>; A. Smyrnakis<sup>1</sup>; M. Kosmopoulou<sup>1</sup>; A. Grigoriadis<sup>1</sup>; I. Orfanopoulos<sup>1</sup>; N. Manolis<sup>1</sup>; I. Panagiotopoulos<sup>1</sup>; R. Gioves<sup>1</sup>; A. Lekkas<sup>1</sup>; F. Busch<sup>2</sup>; J.F. Greisch<sup>2</sup>; S. Pengelley<sup>3</sup>; C. Gebhardt<sup>3</sup>; A. Apalategui<sup>3</sup>; J. Koehling<sup>3</sup>; M. Krause<sup>3</sup>; J. Decker<sup>3</sup>; N. Goedecke<sup>3</sup>; E. Carrascosa<sup>3</sup>; O. Raether<sup>3</sup>**

<sup>1</sup>Fasnatech Science & Technology, Athens, Greece;  
<sup>2</sup>Bruker Switzerland AG, Fällanden, Switzerland  
<sup>3</sup>Bruker Daltonics, Bremen, Germany


## Introduction

Innovations in MS instrumentation are driven by the growing demand for comprehensive characterization of proteins and a wide array of molecular compounds. Central to meeting this demand is the advancement of complementary fragmentation techniques, which are critical for generating high-resolution structural and sequence information required for precise analyte identification. Here, we detail recent innovations in the design and method implementation on the timsOmni<sup>TM</sup> MS platform, highlighting its outstanding versatility and broad applicability.


## Instrumentation

Essential components of the timsOmni platform shown in Fig. 1 include:


- An array of ESI sources covering the entire range of flow rates, from analytical and low micro-flow ESI (VIP-HESI), to online nano-flow ESI (CaptiveSpray), and offline nano-flow ESI (NEOS) for native MS.
- 1.3 mm i.d. inlet capillary delivering nanoamperes of ion current to an RF ion funnel operated at 10 mbar.
- A stacked-ring RF ion guide configured for collisional activation enabling desolvation, Collision Induced Unfolding (CIU) and in-source Collision Induced Dissociation (isCID) upstream of the Trapped Ion Mobility Spectrometer (TIMS).
- A high-capacity dual TIMS analyzer with ion accumulation and ion mobility separation regions, further configured with a gate to allow only ions of a selected mobility to be transmitted.
- An additional collisional activation region at the exit of the TIMS operated at 0.5 mbar pressure.
- A quadrupole mass filter with isolation capabilities extended to 4500 m/z.
- The Omnitrap<sup>TM</sup> platform facilitating MS<sup>n</sup> functionality with ion enrichment and trapped electron-based fragmentation (eXd) with fine electron energy control.
- A new collision cell design supplied with a ramped AC signal for mass selective ejection of ions to improve the duty cycle of the OA TOF analyzer.



**Fig. 1** Schematic of the timsOmni MS platform. Essential components along the ion path are labelled. A magnified view shows the external electron source coupled to Q5. The eXd reaction landscape enabled by fine electron energy control is portrayed. Desolvation of native NIST mAb T ~3 mbar pressure is demonstrated upstream of the TIMS analyzer.




**Fig. 2** CIU-TIMS-ECD experiments performed with native carbonic anhydrase charge state 11+ and ECD sequence maps for different activation levels providing a unique unfolding fingerprint for this protein.




**Fig. 3** Dynamic precursor accumulation actively regulates ion flow in MS/MS scans. (a) LC separation of Adalimumab subunits with the MS2 signal controlled dynamically, (b) Fc/2 sequence map for a 2e7 target, and (c) sequence coverage vs number of scans for Fc/2 in offline experiments comparing different dissociation methods (~300ms/scan).

## Glycopeptide characterization with eXd

LC-MS analysis of a transferrin digest was performed, and O-/N-glycopeptides were targeted for fragmentation. Fig. 4 shows the LC trace (100ng/μL, 0.2μL) and traces corresponding to selected N-glycopeptides accumulated in the Omnitrap and subjected to eXd. Fig. 5 shows the annotated glycopeptide ECD and EID mass spectra.



**Fig. 4** LCMS analysis of transferrin digest and details highlighting the selection of the two low abundance N-glycopeptides and the corresponding MS1 and MS2 traces.



**Fig. 5** (a) ECD of C(Carbamidomethyl)GLVPVLAENY(NHexNAc4Hex5NeuAc)K highlighting the c<sub>10</sub><sup>+</sup> and c<sub>122</sub><sup>+</sup> fragments verifying the position of the glycan and (b) the corresponding EID mass spectrum with abundant fragments related to the glycan.

## Conclusions

- The timsOmni combines TIMS and omnitrap MS<sup>n</sup> eXd technologies for structural characterization of biomolecules.
- Dynamic accumulation of precursor ions in MS/MS scans provides high s/n mass spectra for proteins and peptides.

## New Instrumentation