

Unlocking the discriminative potential of MALDI-TOF profiling data with AI-enabled software for biomarker discovery, classification and typing



Luis Mancera¹, Jesús Jimenez¹, Manuel J. Arroyo¹, Arndt Asperger²

¹ Clover Bioanalytical Software, S.L. Av. del Conocimiento 41, 18016 Granada, Spain

² Bruker Daltonics GmbH & Co. KG, Bremen, Germany

Introduction

Because of its fast analysis speed, high sample throughput, and access to a wide analyte space, MALDI-TOF MS is an ideal tool for rapid and simple profiling analysis allowing for biomarker discovery, classification and typing in a broad range of applications.

We describe here the use of an AI-enabled, workflow-oriented software platform for efficient analysis of MALDI-TOF profiling data in various application examples:

- i) Identification of biomarkers in human blood serum protein profiles
- ii) Classification (i.e., Hierarchical Clustering / PCA) of research-grade biosimilars based on their N-glycan profiles
- iii) Verification of closely related fish species based on protein profiles from tissue extracts applying a supervised machine-learning based classifier model (SVM)

Methods

MALDI-TOF MS data were acquired on neoflex™, autoflex® and microflex® series instruments (Bruker, Bremen, Germany), respectively. Protein profiles were acquired in positive linear mode, N-glycan profiles in positive reflector mode.

Spectra pre-processing (i.e., smoothing, baseline correction, m/z alignment), m/z feature table generation and statistical data analysis were performed using Clover MSDA software (Clover Bioanalytical Software S.L., Granada, Spain), a cloud-based application accessible via web browser interface. Raw spectra were uploaded in Bruker proprietary format and were organized in a customizable study and experiment structure. Clover MSDA software offers three tailored analysis workflows:

Biomarker Analysis: Workflow aiming for detection of discriminative m/z features separating samples from known categories. Uses univariate statistical methods, such as T-test, Mann-Whitney U test, and ROC.

Classification: Workflow enabling classification of MALDI-TOF spectra employing various unsupervised and supervised methods (multivariate statistics; clustering; ML based classifiers), such as Principal Component Analysis (PCA), Hierarchical Clustering, K-Means Clustering, Partial Least Squares Discriminant Analysis (PLS-DA), Support Vector Machine (SVM), Random Forest (RF), K Nearest Neighbors (KNN), and Light Gradient-Boosting Machine (LightGBM). Results from PLS-DA, SVM, RF, KNN and LightGBM can be saved as prediction models for category assignment of samples of unknown origin (*Identification* workflow).

Identification: Workflow intended for category assignment of spectra from unknown samples using one or multiple prediction models generated in the *Classification* workflow

Reference

G.Huffman, L.Mancera, A.Asperger, Bruker Technical Note TN-63, www.bruker.com



Results

Workflow example I - Biomarker Analysis

Experiment: Identification of most discriminative m/z features in MALDI-TOF protein profiles to differentiate between human blood serum samples spiked with artificial peptides and non-spiked controls

MALDI-TOF dataset: 20 spectra per sample category (i.e., spiked, non-spiked), m/z range 1000-10,000, average-mass m/z features




Figure 1: Result of biomarker analysis spiked vs. non-spiked human blood serum

Top left: Butterfly plot of representative MALDI-TOF spectra from non-spiked and spiked human blood serum. Marked peaks represent spiked peptides. **Top right:** Outlier detection result (PCA reconstruction and spectra correlation plots). **Center left:** Peak heatmap. **Center right:** Pearson correlation heatmap. **Bottom left:** T-test result table listing m/z features with q < 0.05. **Bottom right:** ROC result table reporting m/z features with AUC ≥ 0.8.

All m/z features revealed by T-test and ROC analysis as discriminative between sample categories represent the spiked peptides and their oxidized versions.

Workflow example III – Identification of unknown samples using a trained & validated SVM prediction model

Experiment: Training and validation of an SVM prediction model for verification of closely related *Salmonidae* family fish species based on MALDI-TOF protein profiles generated from fish tissue extracts.

MALDI-TOF dataset: Training: 15-16 spectra per fish species, Validation: 6-8 spectra per fish species, m/z range 2000-30,000, average-mass m/z features




Figure 3: Verification of *Salmonidae* fish species applying a SVM prediction model

Left: MALDI-TOF protein profiles of *Salmonidae* fish species. **Center top:** Result from 10-fold cross validation of the SVM model. **Center bottom:** SVM loadings used for model generation. **Right:** Identification results obtained from the SVM model for 18 “unknown” samples. 17 spectra were correctly verified with high probability (>75%), one with medium probability (55%-75%).

Workflow example II - Classification (here: Hierarchical Clustering/PCA)

Experiment: Classification of 4 nivolumab research-grade biosimilar samples based on their MALDI-TOF N-glycan profile spectra

MALDI-TOF dataset: 30 spectra per sample, m/z range 1000-3500, monoisotopic m/z features




Figure 2: Classification result obtained for nivolumab research-grade biosimilars

Top: MALDI-TOF N-glycan spectra from nivolumab research-grade biosimilars. **Bottom left:** Hierarchical clustering result (dendrogram). **Bottom right:** PCA result (score plot).

Hierarchical clustering and PCA successfully classified the spectra by sample origin. Furthermore, both methods confirmed the fucosyltransferase (FT) knockout IgG1 biosimilar being particularly different in its N-glycan profile from the other three biosimilars, which is explained by the lack of core-fucosylated N-glycans in the FT knockout IgG1. Additionally, PCA separated the biosimilar samples by IgG subtype (IgG1 vs. IgG4 Fc regions) along PC2.

Conclusions

- Clover MSDA software provides a powerful statistical toolbox to efficiently leverage the discriminative power of MALDI-TOF MS profiling data.
- Clover MSDA provides access to a wide range of AI-enabled data analysis methods via an intuitive, workflow-oriented web interface. It features user management, project- and experiment-driven data management, as well as detailed reporting of results.
- Clover MSDA supports tailored data analysis workflows allowing for time-efficient analysis of large cohorts of MALDI-TOF spectra with minimal effort enabling simple and fast biomarker discovery, classification and typing across a broad range of life science and further applications.

MALDI-TOF; Classification & Typing