


Rapid-Turnaround Top-Down Sequence Verification and Characterization of Modified Proteins and Oligonucleotides Using a Benchtop MALDI-TOF/TOF Instrument



Sinduri Vuppala¹, Savannah Snyder¹, Arndt Asperger²
¹ Bruker Scientific LLC, Billerica, MA, USA
² Bruker Daltonics GmbH & Co. KG, Bremen, Germany

Introduction

MALDI Top-Down Sequencing (MALDI-TDS) is a proven method for rapid-turnaround verification and characterization of intact proteins, oligonucleotides, PEGylated peptides and other modalities. MALDI-TDS benefits from the formation of singly charged MALDI In-Source Decay (MALDI-ISD) fragments eliminating the need for demanding and potentially error-prone charge deconvolution routines and, thus, making interpretation of MALDI-TDS data simple and straightforward.

Here we highlight the unique MALDI-TDS capabilities of the neofleX™ benchtop MALDI-TOF/TOF for verification and characterization of modified proteins and oligonucleotides using the following samples for demonstration:

- i) a recombinant protein expression product featuring a disulfide bridge and a minor, terminally truncated proteoform
- ii) a heavily modified RNA 39-mer (2’O-methyl, 2’fluoro, phosphorothioate)
- iii) a 5’-DBCO labeled DNA 21-mer.

Methods

A 29.7 kDa recombinant protein of *Thermobifida alba estI* gene coding for cutinase F7IX06-9ACTN, expressed in L-arabinose inducible *E.coli* host cells (concentration 40 pmol/μl), was prepared on the MALDI plate using either 2,5-DHAP (for intact mass measurement) or SDHB (for MALDI-ISD top-down sequencing) as a matrix.

Two oligonucleotide samples (concentration 10 pmol/μl), a heavily modified RNA 39-mer (MW 12.5 kDa) and a 5’-DBCO-labeled DNA 21-mer (MW 7 kDa), were prepared on the MALDI plate using 3-HPA or 2,4,6-THAP as MALDI matrix.

All data were acquired on a Bruker neofleX benchtop MALDI-TOF/TOF instrument using default application methods in positive linear mode (intact mass analysis) and positive reflector mode (MALDI-ISD top-down sequencing).

Data analysis of MALDI-TDS spectra was performed using Bruker OmniScape™ and Biotools software.

References

Bruker Application Note MT-142, www.bruker.com

Acknowledgement

Sample courtesy of the recombinant protein sample by Prof. Bart Devreese and Dr. Isabel Vandenberghe, University of Gent, Belgium, is gratefully acknowledged.



Results

i) 29.7 kDa recombinant protein of *Thermobifida alba estI* gene coding for cutinase F7IX06-9ACTN, expressed in L-arabinose inducible *E.coli* host cells




Figure 1: Intact-mass MALDI-TOF spectrum

Intact-mass analysis confirmed the molecular weight of the expected recombinant expression product, but indicated a protein heterogeneity in form of an additional minor proteoform appearing approx. 203 Da below the mass of the expected expression product. Based on its measured mass, this minor proteoform could suggest a truncated sequence version lacking the N-terminal dipeptide MA.




Figure 2: MALDI-ISD top-down sequencing result

The MALDI-ISD spectrum acquired from the protein sample (Figure 2, left) provided the following information:

- i) The protein sequence of the **major proteoform** was confirmed by means of MALDI-ISD sequence readout at a depth of **80 amino acid residues from N- and C-terminus each**, resulting in a **MS/MS sequence confirmation percentage of >50% provided in a single spectrum** (Figure 2, left, sequence annotation above dashed line, Figure 2, top right, sequence map).
- ii) A **minor proteoform** detected in the intact-mass spectrum at a distance of approx. [-203 Da] was elucidated as **truncated sequence version lacking the N-terminal dipeptide MA** based on the occurrence of a low abundant matching series of N-terminal c-type MALDI-ISD fragment ions (Figure 2, left, sequence annotation below dashed line, Figure 2, bottom right, sequence map).
- iii) A gap observed in the series of matching C-terminal fragment ions (y, z+2) revealed the presence of a disulfide bridge between cysteine residues 242 and 260.

ii) Modified RNA 39-mer (2’O-Me, 2’fluoro, phosphorothio.)




Figure 3: Intact-mass MALDI-TOF spectrum (top) and MALDI-ISD spectrum (bottom)

The sequence of the RNA 39-mer, including its complex modification profile (2’fluoro, 2’O-methyl, phosphorothioate) was confirmed by the MALDI-ISD data with ~67% MS/MS sequence coverage based on matching d- and w-type fragment ions mainly.

iii) 5’-DBCO labeled DNA 21-mer




Figure 4: Intact-mass MALDI-TOF spectrum (top) and MALDI-ISD spectrum (bottom)

The sequence of the 5’-DBCO-labeled DNA 21-mer was successfully verified, including localization of the DBCO-label at its 5’-end, based on matching d-, y- and w-type fragment ions yielding ~85% MS/MS sequence coverage.

Conclusions

- MALDI-TDS performed on neofleX benchtop MALDI-TOF/TOF provides a simple but highly efficient LC-free tool for fast-turnaround top-down sequence verification and characterization of modified proteins, oligonucleotides and other biomolecular modalities.
- MALDI-TDS uses singly charged MALDI-ISD fragments making data analysis particularly easy and straightforward, without a need for charge deconvolution.
- MALDI-TDS allows for seamless switching between sequence analysis of proteins and oligonucleotides without any down-time as is typically the case in LCMS based workflows due to different eluent/modifier system requirements.
- Speed and flexibility of MALDI-TDS facilitate swift decision making in biopharma development as well as deeper insights in structural biology research.

MALDI-TOF/TOF; MALDI-TDS