

Definitive Protein Sequencing by Top-Down MS with timsOmni and OmniScape

Dalila Bensaddek¹; Miriam Escarlet Diaz Galicia²; Magnus Rueping²; Mariangela Kosmopoulou³; Athanasios Smyrnakis³; George Alevizos³; Dimitris Papanastasiou³; Maan Amad¹

¹Analytical Chemistry Core Lab, King Abdullah University of Science and Technology, KAUST, Thuwal, Saudi Arabia; ²KAUST Catalysis Center (KCC), Division of Physical Sciences & Engineering, King Abdullah University of Science and Technology, KAUST, Thuwal, Saudi Arabia; ³Fasmatech Science & Technology, Chalandri, Greece

Introduction

- Bottom-up (BU) proteomics is the dominant approach for protein identification and characterization; however, its reliance on peptide-centric measurements and protein inference can lead to incorrect or incomplete sequence assignments, particularly in the presence of sequence scrambling or proteolytically inaccessible regions.
- In contrast, top-down (TD) mass spectrometry directly analyzes intact proteoforms, preserving sequence connectivity that may be lost in bottom-up workflows. Despite this advantage, *de novo* sequence analysis by top-down MS remains challenging due to limited availability of dedicated software.
- Here, a protein sequence originally inferred from bottom-up data was re-examined using top-down experiments on a timsOmni™ platform with multiple fragmentation modes. Comprehensive OmniScape™ analysis enabled reconstruction of a corrected protein sequence.

Methods

- The ~9 kDa protein was initially dialyzed against water and lyophilized. Prior to analysis, the sample was reconstituted at 4 μM in 50:49.5:0.5 H₂O:ACN:FA and analyzed offline on a timsOmni platform using an Apollo II source at 3 μL/min. MS² experiments were performed on the 8+ precursor charge state using multiple fragmentation methods, including cCID, rCID, ECD, ECcID and EID.
- Data processing was performed in OmniScape using multiple workflows. Data analysis followed a stepwise approach (Results 1-4), starting from zero sequence coverage of the initially proposed, BU-derived protein sequence, progressing through *de novo* sequence reconstruction and clipping analysis, and concluding with mutation screening to confirm the revised sequence, ultimately achieving 100% sequence coverage (SC).


Results

1

Observation

TD-MS measured mass is higher than the BU-predicted mass, indicating sequence inconsistency.


Additionally, the confirmation workflow on MS² ECD data resulted only in random fragment matches (Fig. 1).



2

Method


De novo on MS² ECD data: the top-scored sequence (Fig. 2) tag includes a His tag connected to the N-terminal part of the protein. Additionally, the first fragment identified has a mass shift corresponding to MHH amino acids.



3

Revised Sequence

Confirmation result of the first sequence revision, based on the *de novo* analysis, is shown in Fig. 3.




4

Method

De novo on iSCID data: another sequence tag was consistent with the BU-derived sequence near C-terminus; however, the systematically lower y-ion masses further underlined that the BU-proposed sequence is erroneously longer.

Clipping analysis was subsequently applied to iteratively trim the C-terminal region of the first sequence revision.


Confirmation result (Fig. 4) of the second sequence revision provided great coverage near both protein termini.



5

Method


Coverage gap is observed between O45-O59 for the second sequence revision, while some high-abundance peaks remain unassigned in the MS² ECD.



6

Revised Sequence

Targeted amino acid substitution analysis for O45 and O59, based on the MS² ECD data, revealed the substitution of both O residues with Q (Fig. 5, 6).




7

Method

Despite the almost full coverage of the third sequence revision, further validation of the heavily revised N-terminal region was considered necessary.

Single-site Substitution Screening across residues M1–V20 via PTM Screening workflow, further validated the revised N-terminal region, with the *de novo*-derived amino acid series yielding the highest sequence and intensity coverage for the MS² ECD data.

Finally, application of the confirmation workflow on the MS² CID, MS² EID, and MS² ECcID data further confirmed the sequence identified in OmniScape, with MS² ECcID achieving SC=100% (Fig. 7).



Summary

Definitive sequence identification

- ✓ His-tag confirmed at N-terminus
- ✓ A new 12 amino acid sequence was identified *de novo* as part of the N-terminal region of the protein
- ✓ C-terminal over-length was identified *de novo* and corrected by Clipping Analysis
- ✓ O45Q and O59Q substitutions identified via PTM Screening
- ✓ Revised sequence was validated by Single-site Substitution Screening and further confirmed by multiple fragmentation methods with SC=100%

Conclusion

De novo sequencing & PTM Screening workflows enable protein sequence reconstruction using TD timsOmni data

De novo & PTM Screening in OmniScape