

# Countermeasures and Solutions for Helium Gas Supply Shortages - GCMS

Analytical & Measuring Instruments Division  
Shimadzu Corporation

# Introduction

## Countermeasures and Solutions for Helium Gas Supply Shortages

 [Price Inquiry](#) [Product Inquiry](#) [Technical Service / Support Inquiry](#) [Other Inquiry](#)

Due to helium gas supply shortages, soaring prices, and other factors, reducing the quantity of helium gas consumed and considering alternative carrier gases have become urgent issues. The following describes functionality for reducing helium gas consumption rates and indicates precautions for switching to a different carrier gas.



GC

### Reducing Helium Consumption

- [Reducing Consumption During Analysis](#)
- [Reducing Consumption After Analysis](#)
- [Shutting off the Gas Supply After Analysis](#)

### Changing the Carrier Gas



GCMS

### Reducing Helium Consumption

- [Reducing Consumption During Analysis](#)
- [Reducing Consumption During Standby or After Analysis](#)
- [Estimating Reductions by Ecology simulation](#)

### Changing the Carrier Gas



Because of the challenges seen by our customers, we have created a dedicated page for [Countermeasures and Solutions for Helium Gas Supply Shortages](#)

Please scan the above QR code to be taken to the webpage seen here.

# Contents

---

- **Approach to Helium Gas Supply Shortages**
- **Reducing Helium Gas Consumption**
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas
- **Countermeasures for Each Sample Introduction System**

# Contents

---

- **Approach to Helium Gas Supply Shortages**
- Reducing Helium Gas Consumption
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas
- Countermeasures for Each Sample Introduction System


# Approach to Helium Gas Supply Shortages

Reducing Helium Consumption

Changing the Carrier Gas

# Introduction

## Countermeasure Flowchart



\* Nitrogen carrier gas can also be used in GCMS systems. However, sensitivity will probably decrease to one-tenth or less of normal, so hydrogen is recommended if hydrogen usage is possible.

# Contents

---


- Approach to Helium Gas Supply Shortages
- **Reducing Helium Gas Consumption**
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas
- Countermeasures for Each Sample Introduction System

# Cases where Using an Alternative Carrier Gas would be Difficult



**Sensitivity is low near the limit of quantitation (older models, ultra-trace analysis, etc.)**

- Alternative carrier gases (hydrogen or nitrogen) generally decrease sensitivity.
- Using alternative gases might be difficult unless S/N is at least 50 to 100 at the limit of quantitation required for using helium.



**System is not compatible with alternative carrier gases (certain sample introduction systems, databases, etc.)**

- Some sample introduction systems or optional products are not compatible with alternative carrier gases.
- Databases for semi-quantitative analysis cannot be used with alternative carrier gases.




# Methods for Reducing Helium Gas Consumption

- 1. Reducing Consumption during Analysis**
  - Carrier gas saver mode
- 2. Reducing Consumption during Standby or after Analysis**
  - Ecology mode
  - Gas selector (optional)
- 3. Utilization of Correct Purity Gas is Key**
  - Using a Gas Filter (optional)
- 4. Reducing Consumption by Shutting Off System on Days Not Used**
  - GCMS systems with quick startup

# Reducing Consumption during Analysis

## - Carrier Gas Saver

### Helium Gas Consumption during Analysis



$$\text{Total flow} = \text{Column flow} + \text{Split flow} + \text{Purge flow}$$

Analysis is normally performed with the split ratio set at 30 to 50 mL/min immediately after injection. However . . .




**Required immediately after injection, but not during entire analysis**

**Reduce split flow rate during analysis. Reduce total carrier gas flow rate.**

# Reducing Consumption during Analysis

## - Carrier Gas Saver Mode

Helium gas can be saved during analysis by using the carrier gas saver mode to reduce the split ratio after injection.



Helium consumption per analysis can be reduced by **74 %**.

Simulation Example

Analysis time: 30 min  
Total flow: 50 mL/min

Injection mode: Splitless  
Split ratio during carrier gas save: 5

Sampling time: 1 min  
(1455 mL normally, 375 mL during carrier gas save)

Carrier gas save: 3 min

# Reducing Consumption during Standby / after Analysis - Ecology Mode

## Manual


One click to start ecology mode




During Consecutive Analyses  
Total carrier gas flowrate: 50 mL/min



## Automatic (After Batch Analysis is Finished)



| Sample ID | Sample Type    | Analysis Type | Method File | Data File            |
|-----------|----------------|---------------|-------------|----------------------|
| STD-0001  | 1:Standard (I) | IT QT         | SIM.qgm     | pest_std_001.qgd     |
| STD-0002  | 1:Standard     | IT QT         | SIM.qgm     | pest_std_002.qgd     |
| STD-0003  | 1:Standard     | IT QT         | SIM.qgm     | pest_std_003.qgd     |
| UNK-0004  | 0:Unknown      | IT QT         | SIM.qgm     | pest_unknown_001.qgd |



After Consecutive Analyses  
Total carrier gas flowrate: 20 mL/min




Reduces helium consumption during instrument standby by **60 %**.

The ecology mode is available on GCMS-QP2010 Ultra or higher-end QP-series models and all TQ-series models. It cannot be used with some pretreatment systems, such as TurboMatrix HS and Optic-4 units.

# Example of Carrier Gas Savings Using the Ecology and Carrier Gas Saver Modes

Number of days a 7000 L helium gas cylinder can be used\*




\* Simulation Parameters  
One system connected per helium gas cylinder  
Per day: 30 min × 10 analyses (with 5 hrs of analyses and 19 hrs of standby)

# Reducing Consumption during Standby - Gas Selector



## Gas Selector

Installing a gas selector in a GCMS-NX series system can significantly reduce helium gas consumption by automatically switching to an alternative gas after analysis.



ガスセレクタ模式図

This is a dedicated option for Nexis GC-2030 and GCMS-NX series.



# Utilization of Correct Purity is Key


## Using a Gas Purifier Filter



- Purifier filter for helium or nitrogen gases


Removes moisture, oxygen, and hydrocarbons (replaced about once a year)

### Effectiveness of Filtering Out Hydrocarbons



3 L of He gas trapped by Tenax TA and measured by TD-GC/MS

### Effectiveness of Filtering Out Water and Oxygen




SIM measurement of water and oxygen in carrier gas

- By installing a gas purifier filter, 99.995 % helium can be used.
- That also helps reduce costs.
- Gas cylinders are easier to obtain as well.


# Reducing Consumption by Shutting Off System on Days Not Used

## Shutting Off the System on Weekends or Other Days the System is Not Used

Shimadzu's high-end GC-MS model uses a high-capacity differential vacuum pump to achieve extremely fast instrument startup times (about 2 hrs for qualitative analysis or 4 hrs for quantitative analysis).



**Comparison of He Consumption Using the Ecology Mode and Shutting Down the System for 48 Hours on Weekends**





# Information on ECO Simulation


## Savings during Analysis


### ECO Simulation

ECO simulation can simulate the benefits of using the carrier gas saver function and increasing analysis speed.




To access the ECO simulation function, use this [link](#).


Go to GCMS top >>



## Ecology simulation

Comparing running costs and CO<sub>2</sub> emissions output between new and previous models



This tool simulates the reduction in running CO<sub>2</sub> emissions using the new GC-MS system as

Simul

\* Please utilize working  
\* Simulation results are  
\* This simulation is based on a standard GC-MS system consisting of operation for standby even when the GC-MS

| Your current GC-MS<br>(Eco-mode un-corresponding.) | GCMS-QP2010 Plus<br>/GCMS-QP2010 |
|--|----------------------------------|
| CO <sub>2</sub> emissions                          | <b>4,485</b> kg / year           |
| Electricity consumption                            | <b>7,994</b> kWh / year          |
| *Including a PC                                    | <b>1,839</b> US\$ / year         |
| Gas consumption                                    | <b>3.7</b> Cylinder / year       |
|  | <b>1,481</b> US\$ / year         |

➔

New GC-MS  
(Eco-mode corresponding.)

| GCMS-QP2020 NX<br>/GCMS-QP2020 |                             |
|--------------------------------|-----------------------------|
| CO <sub>2</sub> emissions      | <b>3,437</b> kg / year      |
| Electricity consumption        | <b>6,126</b> kWh / year     |
| *Including a PC                | <b>1,409</b> US\$ / year    |
| Gas consumption                | <b>1.98</b> Cylinder / year |
|                                | <b>790</b> US\$ / year      |

In comparison with your current GC-MS, there is a reduction of

**about 1,048 kg of CO<sub>2</sub> and**

**US\$ 1,121 in running costs per year.**

# Contents

---

- Approach to Helium Gas Supply Shortages
- Reducing Helium Gas Consumption
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas
- Countermeasures for Each Sample Introduction System

# Comparison of Various Carrier Gas Properties - GCMS

|                                | Helium  | Hydrogen  | Nitrogen  |
|--------------------------------|---|---|---|
| <b>Sensitivity</b>             | <b>A</b><br>Suitable for high-sensitivity analysis                  | <b>B</b><br>Sensitivity is about 1/3 to 1/5 of helium   | <b>C</b><br>Sensitivity is about 1/10 of helium due to high noise |
| <b>Price</b>                   | <b>B</b><br>Rapidly increasing prices and sometimes long lead times | <b>A</b><br>Available inexpensively. Installing a hydrogen generator eliminates need to manage gas cylinders. | <b>A</b><br>Available inexpensively                               |
| <b>Safety</b>                  | <b>A</b><br>Non-explosive inert gas                                 | <b>B</b><br>Can explode, but explosion risk is low at volumes used for GC-MS.                                 | <b>A</b><br>Non-explosive gas with high safety                    |
| <b>Column I.D.</b>             | <b>A</b><br>Supported up to 0.53 mm                                 | <b>B</b><br>0.18 mm or less recommended   | <b>B</b><br>0.18 mm or less recommended                           |
| <b>Compatible GC/MS Models</b> | All QP and TQ series  | All QP except QP2010 SE and TQ series   | QP2010 Ultra and QP2020 series (not supported by TQ series)       |

Note: Argon cannot be used as a GCMS carrier gas.

# Alternative Carrier Gas Compatibility Table

| GCMS Main Unit     | Discontinued Year   | He | H <sub>2</sub> | N <sub>2</sub> | N <sub>2</sub><br>(During Standby) |
|--------------------|---------------------|----|----------------|----------------|------------------------------------|
| GCMS-QP5050A       | 2004                | ✓  | —              | —              | —                                  |
| GCMS-QP2010        | 2006                | ✓  | ✓              | —              | ✓                                  |
| GCMS-QP2010 Plus   | 2010                | ✓  | ✓              | —              | ✓                                  |
| GCMS-QP2010 Ultra  | 2015                | ✓  | ✓              | ✓              | ✓                                  |
| GCMS-QP2010 SE     | Currently available | ✓  | —              | —              | ✓                                  |
| GCMS-QP2020 Series | Currently available | ✓  | ✓              | ✓              | ✓                                  |
| GCMS-TQ Series     | Currently available | ✓  | ✓              | —              | ✓                                  |


✓: Indicates supported gas

# Sensitivity and Separation with Alternative Carrier Gases

## Sensitivity

### Approximate Measurement Ranges for Various Carrier Gases (Quantity for SIM and On-Column)


These measurement ranges are only approximate and might not be suitable for some target compounds due to sensitivity or compound properties.



Hydrogen and nitrogen have lower ionization energies than helium, which tends to increase noise.


## GC Separation

### Separation Performance for Various Carrier Gases



- The smaller the height equivalent to Height equivalent to a theoretical plate (HETP), the better the separation.
- Hydrogen provides good separation even at higher linear velocities, which means analysis times can be shortened.
- Nitrogen has a steeper HETP curve, which results in a tendency for separation to become worse.

# Column Selection for Alternative Carrier Gases



Max. column flowrate  
(single GC-MS)  
He: 15 mL/min  
H<sub>2</sub>: 10 mL/min  
N<sub>2</sub>: 3 mL/min

Scale down length  
and internal diameter.



30 m × 0.25 mm



20 m × 0.18 mm

- Maximum column flowrate is lower for hydrogen or nitrogen, because the MS vacuum level decreases.
- For nitrogen, scaling down the column length and internal diameter is recommended due to the 3 mL/min maximum column flowrate.

# Contents

---

- Approach to Helium Gas Supply Shortages
- Reducing Helium Gas Consumption
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas**
  - Key Points for Using a Nitrogen Carrier Gas
- Countermeasures for Each Sample Introduction System

# Using a Hydrogen Carrier Gas

## Hydrogen Gas Cylinder



Note: Order gas cylinders from a gas supplier.

## Hydrogen Generator



- Contact a Shimadzu sales representative before changing the carrier gas.
- Using dichloromethane or carbon disulfide as a solvent could generate hazardous gases. Prepare a suitable exhaust ventilation system in advance.





# Hydrogen Safety Measures – Optional GC Accessory

## Using Hydrogen Safely—Optional GC Safety Countermeasures

### Hydrogen Sensor Helps Ensure Hydrogen Carrier Gas is Used Safely

Shimadzu GC's newly designed hydrogen sensor monitors the hydrogen concentration in the GC oven and can detect potential leaks early. When leaks has been detected, it lowers the temperature and automatically switches to a safe standby mode. If the hydrogen concentration rises continuously, the main power is turned off to prevent accidents. (This is an option for Nexis GC-2030)



External Appearance of Hydrogen Sensor


Hydrogen sensor monitors the atmosphere inside the GC oven

# Safety Countermeasures Using a Hydrogen Gas Generator

## Using an Alternative Carrier Gas—Precautions for Using Hydrogen

A hydrogen gas generator can be a good choice when installing a hydrogen cylinder is difficult. For more details, refer to the following link.

<https://www.peakscientific.com/products/hydrogen/>




### Features

- Suitable for flame gas and carrier gas at trace detection limits
- 99.99999%\* Purity
- Internal leak detection with automatic shutdown features
- Proven PEM technology to generate hydrogen safely and reliably
- Regenerative PSA dryers to ensure highest level of purity
- Automatic loading pump as standard
- Maintenance limited to replacing de-ionizer cartridge
- Compact, space-saving modular design
- Creates hydrogen on demand, minimal storage of hydrogen in the system
- Combine multiple units for higher flow requirements
- GC in-oven hydrogen leak detector available as an optional extra
- Peak offers a 3 year cell warranty with this generator as standard.

\*based on O2 content independently verified by National Physical Laboratory, UK

# Separation Using a Hydrogen Carrier Gas

The optimal linear velocity for hydrogen carrier gas is about the same as for helium carrier gas. With hydrogen, separation does not worsen if the linear velocity is increased, which enables high-speed analysis.



# Converting Analytical Conditions for a Hydrogen Carrier Gas

EZGC Method Translator from Restek Corporation

**EZGC™ Method Translator**

| Carrier Gas | Original | Translation |
|-------------|----------|-------------|
| Carrier Gas | Helium   | Hydrogen    |

| Column         | Original | Translation |
|----------------|----------|-------------|
| Length         | 30.00    | 20.00       |
| Inner Diameter | 0.25     | 0.18        |
| Film Thickness | 0.25     | 0.36        |
| Phase Ratio    | 2.00     | 4.00        |

| Control Parameters | Original | Translation |
|--------------------|----------|-------------|
| Flow Rate          | 1.42     | 1.28        |
| Carrier Gas        | 43.02    | 76.14       |
| Pressure           | 1.14     | 0.44        |
| Outlet Pressure    | 99.00    | 101.02      |
| Temperature        | 0.00     | 0.00        |

| Oven Program | Original | Translation |
|--------------|----------|-------------|
| Ramp         | 80       | 80          |
| Temp         | 180      | 180         |
| Hold         | 2        | 1.15        |
| Ramp         | 30.7     | 30.7        |
| Temp         | 180      | 180         |
| Hold         | 0        | 0           |
| Ramp         | 7.1      | 7.1         |
| Temp         | 280      | 280         |
| Hold         | 2.1      | 2.1         |

**Control Method**

Constant Pressure


**Results**

| Run Time | Original | Translation |
|----------|----------|-------------|
| Run Time | 30.00    | 20.59       |
| Speed    |          | 1.46 x      |


Use FC Values for Original | Use FC Values for Translation

The translator is available free via the website or by downloading the software.

Helium carrier  
(30 m × 0.25 mm)



Hydrogen carrier  
(20 m × 0.18 mm)



1. C17, 2. Pencycuron, 3. Dimethoate, 4. Simazine, 5. Atrazine, 6. C18

Note: For more information, refer to the Restek Corporation website.  
<http://www.restek.com/ezgc-mtfc>


Elution patterns similar to using helium can be obtained.

# Converting Analytical Conditions for a Hydrogen Carrier Gas


## Retention Indices can be Used Even if the Carrier Gas is Switched to Hydrogen

By using hydrogen in combination with the Automatic Adjustment of Retention Time (AART) function, retention times for target compounds in databases and existing methods can be adjusted to support identification.

Correlation between Retention Indices before and after Analytical Condition Conversion



Using AART Function to Automatically Adjust Retention Times




# Sensitivity Using a Hydrogen Carrier Gas

| Name of Compound | S/N Ratio (H <sub>2</sub> /He) | Name of Compound | S/N Ratio (H <sub>2</sub> /He) | Name of Compound | S/N Ratio (H <sub>2</sub> /He) |
|------------------|--------------------------------|------------------|--------------------------------|------------------|--------------------------------|
| Dichlorvos       | 0.16                           | Malaoxon         | 0.18                           | Flutolanil       | 0.29                           |
| Dichlobenil      | 0.68                           | Simetryn         | 0.18                           | Isoprothiolane   | 0.15                           |
| Etridiazol       | 0.22                           | Tolclofos-methyl | 0.40                           | Buprofezin       | 0.35                           |
| Chloroneb        | 1.80                           | Alachlor         | 0.23                           | Mepronil         | 0.26                           |
| Isoprocarb       | 0.84                           | Dithiopyr        | 0.44                           | Chlornitrofen    | 0.82                           |
| Molinate         | 1.41                           | Fenitrothion     | 0.39                           | Edifenphos       | 0.29                           |
| Fenobcarb        | 0.32                           | Esprocarb        | 0.35                           | Propiconazole-1  | 0.37                           |
| Trifluralin      | 0.39                           | Thiobencarb      | 0.48                           | Endosulfate      | 0.76                           |
| Pencycuron       | 0.47                           | Fenthion         | 0.27                           | Propiconazole-2  | 0.40                           |
| Dimethoate       | 0.53                           | Chlorpyrifos     | 0.45                           | Thenylchlor      | 0.51                           |
| Simazine         | 0.30                           | Fthalide         | 0.47                           | Pyributicarb     | 0.37                           |
| Atrazine         | 0.18                           | Dimethametryn    | 0.24                           | Iprodione        | 0.90                           |
| Propyzamide      | 0.65                           | Pendimethalin    | 0.31                           | Pyridafenthion   | 0.45                           |
| Pyroquilone      | 0.27                           | Methyl daimuron  | 0.46                           | EPN              | 0.60                           |
| Diazinon         | 0.19                           | Isofenphos       | 0.20                           | Piperophos       | 0.34                           |
| Ethylthiometon   | 0.17                           | Captan           | 0.37                           | Anilofos         | 0.32                           |
| Chlorothalonil   | 0.70                           | Phenthoate       | 0.39                           | Pyriproxyfen     | 0.43                           |
| Iprobenfos       | 0.16                           | Procymidone      | 0.65                           | Cafenstrole      | 0.32                           |
| Bromobutide      | 0.36                           | Methidathion     | 0.21                           | Ethofenprox      | 0.31                           |
| Terbucarb        | 0.47                           | Butamifos        | 0.56                           |                  |                                |

Mean: 0.43




Distribution of H<sub>2</sub>/He S/N Ratios






- Reproducibility is 10 % or lower for all compounds.
- For most compounds, sensitivity is about 1/3 to 1/5 as high as using helium.

# Mass Spectra Using a Hydrogen Carrier Gas

## 2-Methylisoborneol



## Isoxathion Oxon



- Though an existing mass spectral library can be used, similarity scores in library search results might be about 10 to 20 lower depending on the compound.
- For more accurate identification, measuring a standard sample to check quantitation values and reference ions in mass spectra or to check retention times is recommended.

# Contents


---

- Approach to Helium Gas Supply Shortages
- Reducing Helium Gas Consumption
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas**
- Countermeasures for Each Sample Introduction System



# Using a Nitrogen Carrier Gas

## Nitrogen Gas Cylinder



## Nitrogen Gas Generator




Nitrogen generators designed for LC-MS cannot be used as GC-MS carrier gas because the purity is too low (99.0 % or less). Be sure to prepare a nitrogen gas cylinder.

**Caution: Nitrogen leakage cannot be checked with a peak monitor. Use a leak detector or Snoop solution to check for leaks.**

# Separation Using a Nitrogen Carrier Gas



With a nitrogen carrier gas, separation is best at a low linear velocity so prioritizing separation will result in long analysis times.

Also, the range of linear velocities that result in good separation is narrow, so the nitrogen carrier gas is not suitable for conditions that cause linear velocity variation.



# Converting Analytical Conditions for a Nitrogen Carrier Gas

- If constant pressure or constant flowrate control is used for a temperature-programmed analysis, the column linear velocity will vary.
- Because nitrogen has a narrower range of linear velocities that result in optimal separation than helium or hydrogen, constant linear velocity control may provide the best results.



### EZGC™ Method Translator

| Carrier Gas | Original | Translation |
|-------------|----------|-------------|
|             | Helium   | Nitrogen    |

| Column         | Original | Translation |
|----------------|----------|-------------|
| Length         | 30.00    | 20.00 m     |
| Inner Diameter | 0.25     | 0.18 mm     |
| Film Thickness | 0.25     | 0.18 μm     |
| Phase Ratio    | 250      | 250         |

| Control Parameters     | Original | Translation  |
|------------------------|----------|--------------|
| Outlet Flow            | 1.40     | 0.32 mL/min  |
| Average Velocity       | 43.54    | 26.56 cm/sec |
| Holdup Time            | 1.15     | 1.25 min     |
| Inlet Pressure (gauge) | 98.04    | 40.49 kPa    |
| Outlet Pressure (abs)  | 0.00     | 0.00 kPa     |

| Program         | Original | Translation |
|-----------------|----------|-------------|
| He              |          |             |
| Ramp (°C/min)   | 10       | 9.2         |
| Temp (°C)       | 300      | 300         |
| Hold (min)      | 5        | 5.45        |
| Number of Ramps | 1 (1-4)  |             |

Control Method: **Constant Linear Velocity**



| Results  | Solve for  | Original | Translation |
|----------|------------|----------|-------------|
| Run Time | Efficiency | 29.00    | 31.56 min   |
| Speed    | Speed      |          | 0.92 x      |



Constant linear velocity

# Converting Analytical Conditions for a Nitrogen Carrier Gas

## Example of Using Py-GC/MS for Instantaneous Pyrolysis of a Circuit Board



## Library Search Results for Identifiable Peaks


| ID | Name of Compound    | He               | N <sub>2</sub>   |
|----|---------------------|------------------|------------------|
|    |                     | Similarity Score | Similarity Score |
| 1  | Phenol              | 99               | 98               |
| 2  | 2-Methylphenol      | 98               | 98               |
| 3  | Dimethylphenol      | 97               | 97               |
| 4  | MethylethylPhenol   | 95               | 97               |
| 5  | p-Isopropenylphenol | 95               | 95               |
| 6  | p-Cumylphenol       | 93               | 93               |
| 7  | Biphenol A          | 95               | 96               |
| 8  | Dibromobisphenol A  | 83               | 89               |
| 9  | Tribromobisphenol A | 83               | 87               |

Mass spectral patterns obtained with a nitrogen carrier gas tend to be almost equivalent to those with helium.



EZGC Method Translator can be used to align chromatogram patterns to some extent.

# Sensitivity Using a Nitrogen Carrier Gas

S/N Ratio Comparison Using Helium and Nitrogen Carrier Gases  
(27 Azo Compounds at 1 ppm)



1 ppm 2-methyl-5-nitroaniline



For most compounds, sensitivity was about 1/10 as high as using helium.

# Application Examples

## Using Alternative Carrier Gases

Examples of applications using alternative carrier gases are posted on the Shimadzu website.

Please refer to it.

<https://www.shimadzu.com/an/service-support/technical-support/technical-information/gas-chromatograph-mass-spectrometry/mesure/top/index.html#applihe>



### Application

- > [Measurement of VOCs in Vehicle Interiors Using Thermal Desorption GC-MS with Nitrogen as the Carrier Gas](#)
- > [Py-GC/MS Analysis of Electronic Circuit Board Parts Using Nitrogen Carrier Gas](#)
- > [High-Sensitivity Analysis of Phenols in Drinking Water Using Nitrogen Carrier Gas](#)
- > [Simultaneous Analysis of Pesticides by GC-MS Using Hydrogen Carrier Gas](#)
- > [Using a Method Translator Program for GC-MS Analysis with a Hydrogen Carrier Gas](#)
- > [Simultaneous Analysis of 66 Pesticides by GC-MS Using Hydrogen Carrier Gas](#)
- > [Analysis of Di\(2-ethylhexyl\)phthalate by GC-MS Using Hydrogen Carrier Gas](#)

NEW

LAN-XMS-E125

---

**Application Data Sheet**

No. 126

GC-MS  
Gas Chromatograph Mass Spectrometer

**Simultaneous Analysis of Pesticides by GC-MS Using Hydrogen Carrier Gas**

Helium gas is used as a carrier gas in GC/MS. However, in recent years, the use of hydrogen as an alternative gas has increased due to helium gas supply shortages and soaring prices. There are advantages when hydrogen is used as an alternative gas; for instance, sensitivity is similar to that with helium, and high-speed analysis can be performed. However, caution is necessary in handling hydrogen as it is flammable. Hydrogen is obtained from hydrogen gas generators through the electrolysis of water. As a result, smaller quantities of hydrogen need to be stored, which dramatically improves the degree of safety in comparison to when gas cylinders are used. In addition, there is no need to purchase gas cylinders consecutively, which reduces running costs. This Data Sheet presents an evaluation of the usefulness of the combination of a hydrogen gas generator with the GCMS-QP2020 for the simultaneous analysis of pesticides. The GCMS-QP2020 is equipped with a new type of turbomolecular pump, and is more than capable of accommodating hydrogen as the carrier gas.

---

**Experiment**

Mixed standard solutions were prepared at 0.005 mg/L, 0.01 mg/L, 0.05 mg/L, 0.1 mg/L, and 0.5 mg/L, by diluting a pesticide standard sample containing 59 pesticides. The method was created utilizing the EZGC Method Translator\*, which is provided online by Restek Corporation (<http://www.restek.com/ezgc-mtfc>). For details on EZGC Method Translator\*, refer to Application Data Sheet No. 120.

**Table 1: Analytical Conditions**

|   |  |   |
|---|--|---|
| <p>GC-MS: GCMS-QP2020<br/>                     Hydrogen Gas Generator: Precision H<sub>2</sub> Trace (PEAK Scientific Corp.)<br/>                     Column: SH-PeS-SM2 (20 m long, 0.18 mm I.D., φ = 0.36 μm) (P/N: 227-3617-01)<br/>                     Gas Inert: Sky Single Tapel Mini Liner w/ Wood (P/N: 23336-5)</p> | <p>MS<br/>                     Ionization Mode: EI<br/>                     Interface Temperature: 250 °C<br/>                     Ion Source Temperature: 230 °C<br/>                     Measurement mode: SIM mode<br/>                     Event Time: 0.3 sec</p> | <p>Fig. 1: Precision H<sub>2</sub> Trace (PEAK Scientific Corp.) Hydrogen Gas Generator and GCMS-QP2020</p> |
|---|--|---|

**Table 2: SIM Monitoring m/z**


| Compound Name  | Quantitation m/z | Reference m/z | Compound Name     | Quantitation m/z | Reference m/z | Compound Name   | Quantitation m/z | Reference m/z |
|----------------|------------------|---------------|-------------------|------------------|---------------|-----------------|------------------|---------------|
| Dichlorvos     | 165              | 100.0         | Malathion         | 127              | 99.0, 166.0   | Flutolanil      | 173              | 281.0         |
| Dichlorobenzil | 171              | 173.0         | Simetryn          | 213              | 170.0         | Isoprothiolane  | 189              | 118.0         |
| Ethionazole    | 213              | 211.0         | Tolalophos-methyl | 205              | 125.0         | Buprofezin      | 105              | 175.0         |
| Chlorobenz     | 163              | 161.0         | Alachlor          | 188              | 160.0         | Mipronil        | 119              | 269.0         |
| Isopropcarb    | 136              | 121.0         | Dithiopyr         | 354              | 308.0         | Chlorotolifen   | 319              | 317.0         |
| Molinate       | 126              | 86.0          | Fenitrothion      | 277              | 260.0         | Edifenphos      | 310              | 109.0         |
| Fenobucarb     | 150              | 121.0         | Espinozab         | 61               | 222.0         | Propiconazole-1 | 259              | 281.0         |
| Trifluralin    | 306              | 290.0         | Thiofenecarb      | 100              | 72.0          | Endosulfan      | 272              | 274.0         |
| Phencycluron   | 125              | 163.0         | Fenflon           | 278              | 125.0, 163.0  | Propiconazole-2 | 259              | 261.0         |
| Dimethoate     | 87               | 126.0         | Chlorpyrifos      | 314              | 197.0         | Thiophan-methyl | 127              | 238.0         |
| Simazine       | 201              | 186.0         | Fthalide          | 243              | 241.0         | Pyridofuticarb  | 165              | 108.0         |
| Azinphos       | 215              | 200.0         | Dimethamethyln    | 212              | 255.0         | prothione       | 314              | 316.0         |
| Propoxazine    | 175              | 173.0         | Permethrin        | 252              | 261.0         | Pyracloproprion | 340              | 199.0         |
| Pyoxolon       | 130              | 173.0         | Methylglymon      | 107              | 119.0         | EPN             | 157              | 169.0         |
| Diazinon       | 304              | 179.0         | Isofenphos        | 213              | 155.0         | Piperophos      | 122              | 140.0         |
| Ethionmethon   | 99               | 87.0          | Captaf            | 79               | 117.0, 149.0  | Anilox          | 226              | 125.0         |
| Chlorothalonil | 266              | 264.0         | Phenothoate       | 274              | 125.0         | Pyridopyfen     | 136              | 226.0         |
| Iprobenfos     | 204              | 91.0          | Procyimifone      | 149              | 263.0         | Calenstrole     | 100              | 188.0         |
| Bromobutide    | 120              | 119.0         | Methidathion      | 145              | 65.0          | Ethionphos      | 163              | 135.0         |
| Terbufcarb     | 220              | 205.0         | Butamifos         | 286              | 200.0         |                 |                  |               |

# Contents


---

- Approach to Helium Gas Supply Shortages
- Reducing Helium Gas Consumption
- **Using an Alternative Carrier Gas**
  - Key Points for Using a Hydrogen Carrier Gas
  - Key Points for Using a Nitrogen Carrier Gas
- **Countermeasures for Each Sample Introduction System**

# HS-20 NX Headspace Sampler




|               | Carrier Gas                         | Vial Pressurizing Gas | Purge Gas          |
|---------------|-------------------------------------|-----------------------|--------------------|
| Supported Gas | He, N <sub>2</sub> , H <sub>2</sub> | He, N <sub>2</sub>    | He, N <sub>2</sub> |



Note: Purge gas is required for HS -20 Trap models only.




# TD-30/20 Thermal Desorption Systems






|               | Carrier Gas        | Purge Gas          |
|---------------|--------------------|--------------------|
| Supported Gas | He, N <sub>2</sub> | He, N <sub>2</sub> |

# PT7000/Aqua PT6000 Purge & Trap System




|                             | Carrier Gas                         | Sample Purge Gas   | Pressurizing Pumping Gas |
|-----------------------------|-------------------------------------|--------------------|--------------------------|
| Supported Gas (PT7000)      | He, N <sub>2</sub> , H <sub>2</sub> | He, N <sub>2</sub> | He, N <sub>2</sub>       |
| Supported Gas (Aqua PT6000) | He, N <sub>2</sub>                  | He, N <sub>2</sub> | He, N <sub>2</sub>       |

-  Carrier gas
-  Sample purge gas
-  Pressurizing pumping gas


# AOC-6000/5000 Series Multi-Functional Autosamplers



Purging the HS Syringe




Purging the SPME Fiber



|               | Carrier Gas                         | Purge Gas for HS and SPME |
|---------------|-------------------------------------|---------------------------|
| Supported Gas | He, N <sub>2</sub> , H <sub>2</sub> | He, N <sub>2</sub>        |

# PY-3030D/PY-2020iD Multi-Shot Pyrolyzer



— Carrier gas


|               |                    |
|---------------|--------------------|
|               | <b>Carrier Gas</b> |
| Supported Gas | He, N <sub>2</sub> |

# Measures for Using Py-Screener

## Py-Screener Ver. 2 Phthalate Ester and Brominated Flame Retardant Screening System




Note: Identical as when using Py-Screener Ver. 1




# Measures for Using Py-Screener

- IEC standards specify using helium, so measures involve reducing helium consumption.
- Additional helium gas reductions can be achieved by the following measures.


Method (1): Change the Cup Collection Gas



Method (2): Change the Analysis Method



| Cup Collection Gas | He     | N <sub>2</sub> | N <sub>2</sub> |
|--------------------|--------|----------------|----------------|
| Method             | Normal | Normal         | Fast*          |
| Total He Reduction | —      | 36 %           | 73 %           |



\* If the cup collection gas is changed to N<sub>2</sub>, change the Auto-Shot sampler purge time to 3 minutes. Consequently, 0.2 L more helium will be used for Py-GCMS than when using helium.

\* Indicates the helium gas reduction when using the dedicated fast analysis method for phthalate esters. The helium gas reduction rate will be almost equivalent when the fast analysis method for simultaneous analysis of phthalate esters and brominated flame retardants.

# Summary

## Countermeasures-1: Reduce helium gas consumption

- Use the carrier gas saver mode and ecology mode.
- Use a gas selector.
- Use purity gas filter.
- Shut off the system if not being used for 2 to 3 days.

## Countermeasures-2: Use an alternative carrier gas

- Alternative carrier gas choices differ depending on the system and options used.
- As a first choice, consider using hydrogen as an alternative carrier gas.
- Implement safety measures whenever hydrogen is used.

Be sure to read “Gas Chromatograph Hydrogen Gas Safety” on the Shimadzu website before using hydrogen.

<https://www.shimadzu.com/an/service-support/technical-support/handling-precautions/gas-chromatography/index.html>



Thank you for your attention.

GCMS-QP, GCMS-TQ, and Nexis are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.

Third-party trademarks and trade names may be used in this presentation to refer to either the entities or their products/services, whether or not they are used with trademark symbol “TM” or “®”.