

Kit for Direct Probe Ionization Mass Spectrometer

# DPiMS QT



# DPiMS™ QT

## Quick, Simple Measurements with Simple Pretreatment

- Analysis can start immediately with just simple pretreatment.
- Capable of connections to quadrupole TOF mass spectrometers, enabling high-resolution mass spectrometry.
- The slight amount of the sample adhering to the probe is ionized, so the analysis is highly resistant to contamination of the MS unit.

### 3 Step Workflow

---

#### Analysis

Step 1  
Collect the sample.



Step 2  
Add it to the plate.




Step 3  
Start the analysis.




#### Ionization


Step 1  
Sample  
Probe  
Sampling



Step 2  
Probe  
Ions  
MS  
Ionization



Step 3  
Probe  
Ions  
MS  
Mass spectrometry





SHIMADZU

LCMS-9030  
LIQUID CHROMATOGRAPHY-MASS SPECTROMETER

POWER  
START  
STOP  
METHOD  
DATA  
SYSTEM

# DPIIMS™ Can be Used in All Sorts of Fields

## Measure Liquid and Solid Samples with Simple Pretreatment

- Reduces the time required to analyze blood, urine, and other biological samples with conventional methods by approximately 50 %.
- Foods, tissue sections, and other solid samples can be analyzed as it.

Pretreatment times are significantly reduced.

Note: Analytical precision can be improved even for solid samples by performing simple pretreatment such as dilution and centrifugation.



Biological sample



Solid sample

## The Functional Components of Foods Can be Detected Quickly with Simple Pretreatment

### Quick and Easy Detection of Isoflavones in Processed Soy Bean Foods

10 mg of each food product was permeated with 1 mL of 50 % ethanol. This mixture was agitated and centrifuged, after which the supernatant was measured, and each component was detected on-plate.

#### MS Measurement Conditions

Mass range: *m/z* 100–800 / IF Voltage: 3.5 kV  
 Desolvation line temp: 250°C  
 Heat blocktemp: 50°C  
 Polarity: Positive / Event Time: 0.02sec




Boiled soy beans




Chocolate made with processed soy beans

Daidzein *m/z* 255.0652




Genistein *m/z* 271.0601



Daidzein *m/z* 255.0652




Genistein *m/z* 271.0601




## Analysis System with No Carryover


Only a few dozen pL of the sample adheres to the probe after each sampling. There is no worry of contamination inside the mass spectrometer. Carryover can even be prevented by changing probes. Concentrated samples and samples of unknown concentration can be measured.



Nal (400 mg/L) Base Peak (*m/z* 472.6710) MS Chromatogram



Blank solvent measurement after 30 analyses of high concentration samples




Blank Solvent


## Speedy Qualitative Analysis

The qualitative analysis sequences for DPiMS QT and Q-TOF LC/MS were compared. Using DPiMS QT, it was possible to significantly reduce the measurement time in comparison to conventional methods using Q-TOF LC/MS.

### Comparison of Sample Pretreatment Times




### Comparison of Qualitative Analysis Times



## Comparison of Pretreatment Times

### Q-TOF LC/MS Measurement: 10 minutes



Add 300  $\mu\text{L}$  of acetonitrile, 200  $\mu\text{L}$  of water, and 100  $\mu\text{L}$  of whole blood to the Micro Volume QuEChERS kit.

Centrifugation

Instill the supernatant into an LC vial.

### DPiMS QT Measurement: 5 minutes




Mix 20  $\mu\text{L}$  of whole blood and 200  $\mu\text{L}$  of ethanol with 180  $\mu\text{L}$  of water.

Centrifugation


Instill the supernatant onto the sample plate.

## Speedy Qualitative Analysis

Only a few dozen seconds are spent on DPiMS QT qualitative analysis. When performing qualitative screening analyses, the time before results are checked can be reduced, improving procedural efficiency.



Single DPiMS QT Measurement: 0.3 minutes




Implementing Accurate Mass and Composition Estimation from the MS Spectrum




## Switching between DPiMS QT and Q-TOF LC/MS is Simple

Only about 15 seconds are required to remove the DPiMS QT unit, reconfiguring the system for Q-TOF LC/MS. The resources (solvents and columns) required for Q-TOF LC/MS analysis can be reduced by implementing qualitative analysis and primary screening via DPiMS QT, thereby reducing the number of samples requiring quantitative analysis.



## PESI MS Solution Software


PESI MS Solution software is used to control the probe, edit MS method files, and start measurements. Detailed MS analytical conditions are configured from the LabSolutions™ LCMS method editing window. Analysis can be started easily by selecting the probe control mode and MS method in PESI MS Solution. During analysis, the status can be confirmed in the LabSolutions LCMS window.



PESI MS Solution



Probe control method file




LabSolutions method file


## Sequence for Qualitative Analysis via MS and MS/MS

A qualitative analysis was performed on a whole blood sample mixed with trazodone (500 ng/mL). A cycle of MS and MS/MS measurements was performed simultaneously. LabSolutions Insight Explore was used for the data analysis. LabSolutions Insight Explore enables composition estimation, library searches, and structural analysis.

### MS Measurement



Trazodone C<sub>19</sub>H<sub>22</sub>ClN<sub>5</sub>O Chromatogram extracted at *m/z* 372.1586



MS spectrum in the vicinity of *m/z* 372.1586


### Checking for Peaks

The chromatogram extracted for the target is depicted by entering the composition formula or the applicable ion *m/z* into the compound table.


### Composition Estimation

From the measured spectrum, select a spectrum of any *m/z* and use the composition estimation function to list estimated composition formulas in order by score.

### MS/MS Measurement



Results for the Structural Estimation of Trazodone



Product Ion Spectrum Measured with DPiMS QT

### Assigning Fragments

Using the structural analysis assign function in LabSolutions Insight Explore, fragments are assigned with respect to the product ion spectrum.

### Scoring via the Spectral Library

Scoring is performed with respect to the spectrum measured with DPiMS QT by creating a spectral library using Q-TOF LC/MS.

## Specifications

Hardware	DPiMS QT controller	Control of the DPiMS QT unit (Installed in the mass spectrometer)
	DPiMS QT unit	Probe Voltage: $\pm 5$ kV max. (set voltage)
		Probe Stroke: 46.30 mm max.
		Number of Extraction Cycles: 300 max.
		Probe Speed: 300 mm/sec max.
		Probe Acceleration: 1 G max.
	Probe Stop Time: 60,000 msec max.	
Software	PESI MS Solution	Used to select analytical conditions and start analysis
	LabSolutions	Used to specify MS analysis method settings and analyze data
	(Probe Control Software)	Used to specify probe control method settings
Consumables	Probe	Set of 10 or 50 pcs.
	Plate	Plates available for both liquid and biological samples (set of 100 pcs. each)
Installation requirements*	Room temperature	18 to 28 °C
	Humidity	20 to 70 % (No-condensing and no-discharge due to static electricity)
	Others	Should be installed in an environment without dust, oscillation, electromagnetic noise, corrosive gas and electromagnetic interference.

\*Equivalent to the installation requirement of mass spectrometer.

DPiMS and LabSolutions are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.



Shimadzu Corporation  
[www.shimadzu.com/an/](http://www.shimadzu.com/an/)

**For Research Use Only. Not for use in diagnostic procedures.**

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.