

Identification of Unknown Polymer Additives by LCMS-IT-TOF

Technical Report vol.8



1. Introduction

It would be very difficult to imagine what our lifestyle would be like without plastics. We depend on plastics for just about everything. There are hundreds of polymers available for thousands of different applications. What is not so well-known is that synthetic polymers nearly always require polymer additives to achieve key performance properties. For example, when plastics are exposed to heat and light, or when they come into contact with oxygen, they require the inclusion of heat stabilizers or antioxidants to keep them from degrading in performance. There are thousands of common additives and mixtures used commercially, including ultraviolet light absorbing agents (UVA: Ultraviolet Absorbers) and stabilization agents (HALS: Hindered Amine Light Stabilizers). Even one polymer from a single manufacturer can contain different types of additives and blends of additives depending on the grade and intended use.

For this reason, identification of the additives used in the polymer material becomes important for adjusting product characteristics with respect to competitors' products, as well as for improving products within a company's product line.

This report introduces a unique method for identifying polymer additives.

2. Strategy and Pretreatment

This study was undertaken to identify the polymer additives used in a particular type of polymer. Starting with the assumption that the most likely additives were fairly common commercial materials, we recognized that if the composition formula candidates could be sufficiently narrowed, we could predict the composition based on formulae published in catalogs posted on the Internet by chemical additive manufacturers. Accurate mass measurement using MS^n with the LCMS-IT-TOF instrument can quickly narrow the composition formula candidates to just the most likely compounds. The predicted structures for the composition formulas were obtained by comparing them to known structures in different chemical databases. Sample preparation consisted of adding one of the beads from among those provided by the customer in 1 mL THF/methanol (50/50) solution, and then extracting the additives by placing the vial in an ultrasonic bath for 30 minutes. The analysis was conducted by injecting 10 μ L of the supernatant directly into the HPLC column.




Fig. 1: Polymer Beads Used for Measurement

The 3-D ion trap-time-of-flight mass spectrometer (LCMS-IT-TOF) is used to analyze the masses of the components eluting from the HPLC column, and for structural elucidation. The design of the instrument allows for high-speed MSⁿ analysis because the ion acceleration method used to remove ions from the trap ballistically ejects all of the ions from the ion trap towards the TOF

simultaneously. Taking advantage of the high-speed performance of this instrument, the precursor ion automatic selection measurement mode is utilized to conduct MS, MS² and MS³ analysis in order of intensity during peak elution. High-speed polarity switching was used to allow analysis of different types of ions in a single run.




Fig. 2: LCMS-IT-TOF Ion Path

3. Analytical Conditions

Column : Imtakt Cadenza CD-C18 3.0 mm I.D. x 75 mm L
 Mobile phase A : 5 mM ammonium acetate - water
 Mobile phase B : Acetonitrile
 Gradient program : 50%B (0 min) → 100%B (20-40min) → 50%B (40.01-50min)
 Flow rate : 0.4 mL/min
 Injection volume : 10 μL
 Column temp. : 40°C
 Ionization mode : ESI(+) ESI(-)
 Nebulizing gas : 1.5 L/min
 Drying gas pressure : 100 kPa
 Probe voltage : +4.5 kV -3.5 kV
 CDL temperature : 200°C
 BH temperature : 200°C

4. Results

4-1. Chromatogram and Mass Chromatograms




Fig. 3: UV Chromatogram and Mass Chromatogram Peaks A, B, C and E were detected by ESI+, peak D was detected by ESI-.

4-2. Analysis of Peak A




Fig. 4: MS, MS² and MS³ Spectra of Peak A ▼ indicates the selected precursor ion.




Fig. 5: Predicted Composition Results by Formula Predictor (Electron mass is not considered.)




Fig. 6: Database Search Results Using ChemFinder

(M+H)⁺ expected m/z 481.4005




Fig. 7: Predicted Structure and MS², MS³ Spectra Assignments

4-3. Identification Table

Other peaks were identified by this technique. For some of these compounds, the predicted composition formula and mass accuracy are summarized in the table below.

	Compound name	MW	Formula	Calc. (M+H) ⁺	Meas. (M+H) ⁺	ERR (ppm)	Calc. (M-H) ⁻	Meas. (M-H) ⁻	ERR (ppm)
A	Decanedioic acid bis (2,2,6,6-tetramethyl-4-piperidyl) ester	480.3927	C ₂₈ H ₅₂ N ₂ O ₄	481.4005	481.4014	1.87			
B	Triphenyl Phosphate	326.0708	C ₁₈ H ₁₅ O ₄ P	327.0786	327.0788	0.61			
C	Tinuvin P	225.0902	C ₁₃ H ₁₁ N ₃ O	226.0980	226.0981	0.44			
D	2,2-Bis (3-sec-butyl-4-hydroxyphenyl) propane	340.2402	C ₂₃ H ₃₂ O ₂				339.2324	339.2322	-0.59
	Compound name	MW	Formula	Calc. (M+NH ₄) ⁺	Meas. (M+NH ₄) ⁺	ERR (ppm)	Calc. (M-H) ⁻	Meas. (M-H) ⁻	ERR (ppm)
E	Irganox 1076	530.4699	C ₃₅ H ₆₂ O ₃	548.5043	548.5067	4.38			

5. Conclusion

The technique of determining the compositional formula by MSⁿ along with searching a compound database using that composition formula as a keyword is a realistic method for determining polymer additive candidates. Utilization of the predicted structure obtained from the compound database compared with the MSⁿ spectrum obtained from the accurate mass can be considered a unique

method enabled by the use of the LCMS-IT-TOF. It is the ability to conduct MSⁿ measurements along with accurate mass that makes identification fast and easy. This technique can be applied not only to identification of polymer additives, but also to identification of many compounds, including natural products, metabolites, impurities, and degradation products.

Founded in 1875, Shimadzu Corporation, a leader in the development of advanced technologies, has a distinguished history of innovation built on the foundation of contributing to society through science and technology. We maintain a global network of sales, service, technical support and applications centers on six continents, and have established long-term relationships with a host of highly trained distributors located in over 100 countries. For information about Shimadzu, and to contact your local office, please visit our Web site at www.shimadzu.com



SHIMADZU CORPORATION. International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan

Phone: 81(3)3219-5641 Fax: 81(3)3219-5710

URL <http://www.shimadzu.com>

The contents of this brochure are subject to change without notice.