

For LabSolutions LCMS

LC/MS/MS Method Package for D/L Amino Acids





LCMS-8060

Most important amino acids exist as stereoisomers. D- and L- forms of mirror image isomers, or enantiomers, are named according to their activity on polarized light. By using CROWNPAK CR-I(+) and CR-I(-) columns with chiral stationary phases, the D- and L-forms of amino acids can be analyzed separately. With CR-I(+) elution order is from D- to L-, and with CR-I(-) the elution order is reversed.

In Just Ten Minutes, Chiral Amino Acids Can Be Analyzed Simultaneously

With conventional chiral amino acid analysis, it is necessary to perform derivatization or use very long run times. With this method package, derivatization is not necessary, and high-sensitivity analysis can be performed in a short period of time, bringing efficiency to the chiral separations laboratory.



HPLC Conditions

Column	: CROWNPAK CR-I(+)/(-) (3 mmI.D. x 150 mmL., 5 μm)
Mobile Phase	: Acetonitrile/Ethanol/Water/TFA = 80/15/5/0.5
Flowrate	: 0.6 mL/min
Injection Volume	: 1 μL
Column Temp.	: 25 °C

MS Conditions


Nebulizer Gas Flowrate	: 3.0 L/min
Drying Gas Flowrate	: 15.0 L/min
Heating Gas Flowrate	: 5.0 L/min
Interface Temp.	: 250 °C
DL Temp.	: 250 °C
Heat Block Temp	: 300 °C

LC/MS/MS Method Package for D/L Amino Acids


All D/L Amino Acids Can Be Quantified by Column Switching

The physicochemical properties of Glutamine and Lysine, Isoleucine and allo-Isoleucine, Threonine and allo-Threonine are extremely similar. They have virtually the same MS/MS fragmentation patterns, and share many of their MRM transitions. Chromatographic separation is therefore required to analyze them individually by LC-MS/MS. Even if there is coelution on the CR-I(+) column, confirmation can be made by automated switching to a secondary CR-I(-) column.

Analysis of D/L-Isoleucine, D/L-*allo*-Isoleucine, and D/L-Leucine



Analysis of D/L-Threonine, and D/L-*allo*-Threonine



List of Registered Amino Acids

D/L-Alanine
D/L-Arginine
D/L-Asparagine
D/L-Aspartic acid

D/L-Cysteine
D/L-Glutamine
D/L-Glutamic acid
Glycine

D/L-Histidine
D/L-Isoleucine
D/L-*allo*-Isoleucine
D/L-Leucine

D/L-Lysine
D/L-Methionine
D/L-Phenylalanine
DL-Proline

D/L-Serine
D/L-Threonine
D/L-*allo*-Threonine
D/L-Tryptophane

D/L-Tyrosine
D/L-Valine

Precautions

- DL-Proline is a secondary amine, so it cannot be separated with these analysis conditions.
- LabSolutions LCMS Ver. 5.86 or later is required.

CROWNPAK CR-I(+) and CR-I(-) are products of Daicel Corporation.



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.