


Analytical to Preparative Scale-up System

Nexera ASAPrep



Unique Algorithm Increases Efficiency of Compound Purification

During compound synthesis, both the target product and impurities are generated. These impurities must be removed in a subsequent purification step to ensure sufficient purity for actual use. Preparative liquid chromatography (preparative LC) is commonly used for this process. However, because it requires specialized expertise and takes a significant amount of time, there is a clear need to improve efficiency.

The Nexera™ ASAPrep™ (automated scaling-up from

analytical to preparative chromatography) system is equipped with a unique algorithm that automatically evaluates the complexity of preparative chromatography and generates optimal conditions.

Furthermore, an intuitive UI design enables easy purification by anyone, regardless of experience level.

Because the system is able to help assign purification tasks depending on the difficulty level, it can enable a more efficient purification workflow.



Nexera ASAPrep Features

- The complexity of preparative separation can be evaluated automatically from screening results.
- Intuitive software enables screening and purification operations to be executed with ease.
- Unique algorithm ensures accurate scale-up.



Automatically Evaluates Difficulty Level of Preparative Separation

Specialized screening methods are executed to automatically evaluate the difficulty level of preparative separation. Based on that difficulty level, the software then automatically creates a preparative method.





Intuitive Analytical/Preparative Operation



Screening analysis can be executed by simply specifying the vial position, sample name, injection volume, and Target Mass.

The suggested optimal preparative method can be executed by simply importing the screening results obtained from the analytical system into the preparative system.



Unique Algorithm Ensures Accurate Scale-Up


Equipped with Shimadzu's unique chromatography scale-up algorithm, a focused gradient is used to automatically generate a preparative method from screening results.

In addition to using the standard algorithm, which supports a wide variety of compounds, users can create their own custom scale-up algorithm.


Application



Screening Analysis Results



Purification Analysis Results



System Contents

Analytical System: For Screening Analysis

This Nexera ASAPrep-compatible system is configured with Nexera XR and LCMS-2050 units.

The system uses built-in solvent and column-switching valves to screen samples based on either acidic or basic conditions.



Preparative System: For Purification Analysis

This Nexera ASAPrep-compatible system is configured with Nexera Prep and LCMS-2050 units.

The system uses built-in solvent and column-switching valves to switch between acidic and basic conditions. Optimized preparative condition are generated and executed based on screening results.



Shim-pack Scepter™

The columns used in the Nexera ASAPrep system must allow the use of acidic and basic mobile phases across a wide pH range, and enable seamless scale-up from analytical to preparative analysis.

Shim-pack Scepter C18 columns are packed with an organic silica hybrid base material that provides excellent pH resistance, allowing them to be used with the wide pH range required by the Nexera ASAPrep system.


The standard method included with Nexera ASAPrep systems has been optimized using the Shim-pack Scepter C18 column.



Open Solution

Open Solution software enables efficient analysis by multiple users sharing the same system in a research laboratory. It supports both LC and LC-MS analysis, allowing users to efficiently perform qualitative analysis and scale-up to purification simply by selecting from preset methods.

Remote data analysis using a personal computer helps minimize the exclusive use of analytical computers. Furthermore, built-in functionality helps maintain system performance.



Nexera, ASAPrep, and Shim-pack Scepter are registered trademarks of Shimadzu Corporation or its subsidiaries in Japan and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The contents of this publication are provided to you "as is" without warranty of any kind, and are subject to change without notice. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.