

In-Source Fragmentation Searches in NIST26 Chromatogram Window

Video/Handout

James Little

Mass Spec Interpretation Services

April 24, 2026

mzinterpretation.com

See **Full Course** on NIST26 with new **Integrated** Deconvolution/Library Searching for **EI GC-MS** and **LC-MS/MS**!

Mass Spec (m/z) Interpretation Services

Organic Mass Spectrometry



“In-source MS/MS”, spectra like these are typical when fragmentation occurs in the source

1. Mixed precursor + fragment ions

- You usually see the **protonated molecule** (e.g., $[M+H]^+$) still present
- Alongside **multiple fragment ions** produced before any collision cell


👉 That’s the key difference vs true MS/MS:

- **MS/MS**: isolated precursor → clean fragment spectrum
- **In-source**: everything fragments together → “mixed” spectrum

Other Names


- In-source fragment ions (*best, most explicit*)
- Source-induced fragments
- In-source fragmentation products
- Source fragmentation ions
- In-source CID products

If NIST saw these ions when infusing a reference standard, they would obtain MS2 spectra for the molecular ion adducts *and* fragment ions, of course *at high resolution*



Molecules that primarily lose water from their [M+H]⁺ in Positive Ion Electrospray


- For certain compounds such as cholesterol, there would be **no** MS/MS (MS2) Spectrum in Library if the [M+H]⁺ ion was necessary
- Essentially **no significant ion intensity** is noted for [M+H]⁺ ion adduct for cholesterol
- **Thus**, the major MS/MS spectra is for the [M+H-H₂O]⁺ peak at *m/z* 369.3506 in NIST MS/MS library



Chromatogram Results for Small Losses

- Thus, for such molecules as cholesterol, the chromatogram window search results are shown for $[M+H-H_2O]^+$ **with the correct structure**
- This is **only done** for small losses such as NH_3 , H_2O , etc.
- Or sodium adducts, $[M+Na]^+$
- Another example is shown below for Oxosorbicillinol

#	Scan	RT	Score	Score (Unfiltered)	Abund.Rel.	Prec. m/z	Prob	dPPM	XIC Num.	nSpec	Iso. Profile	Width	Prec. Type	Lib	Obs	Lib ID
548	959	6.2938	22	928	4.81	247.0994	0	-11.7	881	3	0.9999	16.4	$[M+H-H_2O]^+$	hr_msms_nist#2	2	Oxosorbicillinol



Name: TESTMIX2_180504_MAS011_06_shortened.959.959. File: TESTMIX2_180504_MAS0
 MW: N/A ID#: 1299 DB: Text File
 Comment: RTINSECONDS=377.629000 Parent=247.0994 Charge=1
 10 largest peaks:
 247.0993 999.00 | 191.0370 101.27 | 179.0370 69.93 | 205.0526 30.67 | 156.0682 10.1
 192.0448 6.20 | 178.0292 3.89 | 164.0261 3.65 | 155.0604 3.50 | 163.0182 2.1
 Synonyms:
 no synonyms.

Name: Oxosorbicillinol
 Precursor type: $[M+H-H_2O]^+$
 Instrument type: HCD
 Collision energy: NCE=25% 12eV
 Precursor m/z: 247.0965
 Formula: $C_{14}H_{16}O_5$
 MW: 264 Exact Mass: 264.099773 CAS#: 271795-00-7 NIST#: 4764506 ID#: 933068 DB: hr_
 Comment: NIST Mass Spectrometry Data Center
 Ion mode: P
 Instrument: Thermo Finnigan Elite Orbitrap
 Ionization: ESI
 Collision gas: N2
 Sample inlet: direct flow injection
 Spectrum type: MS2
 Notes: micromol/L in water/acetonitrile/formic acid (50/50/0.1); Spec=Consensus Nreps=
 InChIKey: REQRQKXUWRLQY-KGAWMYLDSA-N Non-stereo
 10 largest peaks:
 247.0959 999.0 | 179.0335 60.6 | 219.1010 44.7 | 197.0440 21.3 | 205.0855 19.1 |

Chromatogram Results for Small Losses

- For large fragment losses, the correct structure *is not noted* for the entry in the list
- The list only shows structures for M+H, M+Na, etc. ions and for small molecule loss as mentioned previously
- **Thus, to see** the proposed structure or structures, one **must send** the spectrum to NIST Lib Search Window
- Add Score(Unfiltered) as a property in the lower window
- Pick a reasonable value of >800 in top bar filter settings
- Also useful to sort by Score (unfiltered) with highest value at top of list
- Look for Entries whose Score is **MUCH** lower than Score (unfiltered)

#	Scan	RT	Score	▼ Score (Unfiltered)	Abund.
1	937	6.2258	978	978	39.4
2	3379	13.3839	975	975	7.89
3	1728	8.5602	970	970	10.2
4	1893	9.0431	968	968	15.0
5	951	6.2674	960	960	12.5
6	920	6.1776	960	960	22.0
7	1047	6.5469	958	958	3.78
8	924	6.1942	946	946	1.40
9	3194	12.8413	939	939	6.78
10	779	5.7726	937	937	34.8

Reload on Startup Best Hits OnlyIDs Highlight Score Bkgd Filter Scores < Average Ab. Threshold

#	Scan	RT	Score	▼ Score (Unfiltered)	Abund.Rel.
11	959	6.2938	22	928	4.81



Example in list, Big Difference!

In-Source Candidates for Sending to Lib Search for Identity

- Possible In-source spectra are noted with **blue dot** in Chromatogram TIC plot
- Reasonable structures **ONLY found** by sending spectrum to Lib search

PRE-RELEASE NIST MS Search 4.0 (Tandem mode) - [TESTMIX2_180504_MAS011_06_shortened, MS/MS, nHits=534]


File View Window Help

MS m/z

Reload on Startup Best Hits OnlyIDs Highlight Score $>=800$ Bkgd $<=90\%$ Filter Scores < 600 Average Ab. Threshold Auto


Abundance(Log)

TIC



Sending In-Source Spectrum to Lib Search




- When selected with mouse, turns from blue to green to show active spectrum linked to list
- Note peak detected has very large difference in Score of 22 and Score Unfiltered of 928!
- Right click on entry and select library search






Several *Candidates* Found in Lib Search Results from In-Source Spectrum with High Scores

#	Lib.	Score	DotProd	Prob. (%)	Prec. Type	Instr. Type	Energy	DBs	dPPM	Name
1	hr	928	951	51.6	[M+H] ⁺ =>[...]	IT-FT	35%	5	-1.2	Cletoquine
2	hr	924	953	48.3	[M+H] ⁺ =>[...]	IT-FT	35%	4	-1.2	Desethylchloroquine
3	hr	909	952	51.6	[M+H-C2H...]	HCD	44%	5	-1.2	Cletoquine
4	hr	901	946	51.6	[M+H-C2H...]	HCD	35%	5	-1.2	Cletoquine
5	hr	877	923	51.6	[M+H-C2H...]	IT-FT	35%	5	-1.2	Cletoquine
6	hr	832	916	51.6	[M+H-C2H...]	HCD	54%	5	-1.2	Cletoquine
7	hr	825	906	51.6	[M+H-C2H...]	HCD	27%	5	-1.2	Cletoquine
8	hr	748	803	48.3	[M+H-C2H...]	IT-FT	35%	4	-1.2	Desethylchloroquine
9	hr	737	854	51.6	[M+H-C2H...]	HCD	20%	5	-1.2	Cletoquine

Make Sure Search Options Correct



Unclicked. off!



MS/MS Hit Filter Settings