

Ultra Fast Analysis of Amino Acids in Cultured Cell Extracts Using UHPLC/MS/MS

IMSC 2012 PWe-007

¹Taku Tsukamoto, ²Yuki Sato, ¹Satoshi Yamaki
¹Shimadzu Corporation, Kyoto, Japan,
²Shimadzu GLC Ltd., Tokyo, Japan

Ultra Fast Analysis of Amino Acids in Cultured Cell Extracts Using UHPLC/MS/MS


Introduction

An essential aspect for the research of metabolic behavior involves the analysis of amino acids (AAs), which is typically carried out using ion-exchange HPLC with post-column derivatization. Recently, not only relatively high concentrations of AAs such as proteinogenic AAs but also relatively low concentrations of AAs and related compounds such as non-proteinogenic AAs and dipeptides are also monitored as significant targets for this purpose.

In case of biological samples, however, interferences by other compounds which involve high concentrations of AAs

significantly affect the sensitivity for low concentrations of AAs due to low selectivity of this traditional method. Moreover this method requires over two hours for the trace analysis of AAs including related compounds (Fig. 1).

Herein, we describe highly selective and ultra fast trace analysis method of AAs using reversed-phase UHPLC/MS/MS with pre-column derivatization specifically aimed at detection of low concentrations of AAs in biological samples, along with its application for the analysis of cultured cell extracts.



Interference to low concentrations of AAs by matrix
Fig. 1 Ion-exchange post-column derivatization method (Sample: Rat plasma)

Methods and Materials

Preparation

Samples were prepared using the EZ:faast amino acid kit (Phenomenex). After addition of 3 internal standards and solid phase extraction, samples were derivatized using alkylchloroformate which modifies N-terminal and C-terminal of AAs (Fig. 2). This preparation takes about 7 minutes.




Fig. 2 Procedure of EZ:faast AAs kit preparation


Ultra Fast Analysis of Amino Acids in Cultured Cell Extracts Using UHPLC/MS/MS

Instruments and analytical methods

Reversed-phase core-shell type column Kinetex C18 (Phenomenex) and a UHPLC system (Shimadzu Nexera) coupled with a triple quadrupole type mass spectrometer (Shimadzu LCMS-8030) was used for ultra fast analysis of AAs.

Using our analytical method, 36 AAs and related compounds were analyzed in 7 minutes (Fig. 3). This means total analysis time including preparation was 14 minutes, significantly shorter than that of traditional HPLC methods.

Column :Phenomenex Kinetex C18 2.1x150 mm, 1,7 μ m
 Mobile Phase :5 mmol/L Ammonium Formate - Water/Methanol
 Gradient Prog. :B conc. 50%(0 min)->65%(1-3 min)->90%(4-4.5 min)->50%(4.51-7 min)
 Flow rate: :0.4 mL/min
 Column Temp. :50°C
 Inj. Vol. :2 μ L



**Analysis time 7 min
(+preparation 7 min)
...total time 14 min**



**About 1/10
comparing Ion-
exchange method**

Fig. 3 MRM chromatograms of 36 AAs standard

Table 1 shows calibration curve ranges, correlation coefficient (R) and reproducibility (CV%) of 36 AAs.

Most of AAs were detected selectively with MRM and lower limits of calibration curves were ranged from 2 to 500 pmol/mL.

Table 1 Calibration curves and reproducibility of 36 AAs

MRM transition	calibration curve range (pmol/mL)	R	CV %
Arg	303.2>70.1	20 - 80000	0.9934 4.5
Gln	275.2>171.9	50 - 80000	0.9953 3.7
Ans	369.2>309.2	200 - 80000	0.9892 7.7
Cit	304.2>156.1	50 - 80000	0.9907 1.2
Ser	234.1>146.0	50 - 10000	0.9938 3.8
Asn	243.1>115.1	50 - 10000	0.9926 4.1
Pro-Hyp	357.2>156.0	5 - 40000	0.9991 1.9
Hyp	260.1>172.0	10 - 4000	0.9959 1.0
Gly	204.1>118.0	500 - 10000	0.9991 1.7
Thr	248.1>159.9	50 - 10000	0.9964 7.5
Gly-Pro	301.2>157.9	2 - 10000	0.9980 10.1
β Ala	218.1>98.0	5 - 10000	0.9948 7.1
Ala	218.1>130.0	100 - 80000	0.9941 2.1
GABA	232.2>130.0	100 - 10000	0.9780 1.2
Sar	218.1>115.9	100 - 80000	0.9932 1.8
β AIBA	232.2>129.9	2 - 10000	0.9976 6.9
ABA	232.2>143.9	2 - 10000	0.9963 8.3
Pro	244.2>156.0	50 - 10000	0.9993 1.8

MRM transition	calibration curve range (pmol/mL)	R	CV %
Met	278.1>190.0	2 - 4000	0.9934 11.4
Orn	347.2>287.0	50 - 80000	0.9925 4.2
Car	441.2>284.0	20 - 10000	0.9909 9.5
Tpr	262.1>173.9	5 - 10000	0.9997 3.7
Asp	304.2>216.0	50 - 80000	0.9954 4.1
Val	246.2>157.9	50 - 80000	0.9930 3.8
His	370.2>196.0	50 - 80000	0.9971 5.2
Lys	361.2>301.0	50 - 80000	0.9911 1.8
Glu	318.2>171.9	50 - 80000	0.9886 2.4
Trp	333.2>245.0	2 - 10000	0.9968 7.9
AAA	332.2>244.0	20 - 10000	0.9949 2.9
Leu	260.2>172.0	10 - 80000	0.9947 2.3
Phe	294.2>205.9	5 - 80000	0.9918 5.3
Ile	260.2>129.9	50 - 80000	0.9946 0.9
APA	346.2>198.0	10 - 80000	0.9941 1.9
Cth	479.2>230.0	2 - 80000	0.9963 4.6
C-C	497.2>248.0	2 - 10000	0.9984 4.0
Tyr	396.2>308.0	50 - 10000	0.9971 2.1

Ultra Fast Analysis of Amino Acids in Cultured Cell Extracts Using UHPLC/MS/MS

This result indicates the matrix in cultured cell extracts didn't cause severe suppression effect nor enhancement effect.




Fig. 6 Results of recovery test

Conclusions

- Ultra fast and highly selective analysis of amino acids using reversed-phase UHPLC/MS/MS with pre-column derivatization was investigated.
- 36 amino acids could be analyzed in 7 min for preparation and further 7min for instrumental analysis.
- In the cases of cultured cell extracts, low concentrations of amino acids were detected without severe interference from the matrix.

Acknowledgment

The authors would like to thank Dr. Koji Okamoto, National Cancer Center (Japan) for providing samples.