

## ASMS 2017 MP 448

Junichi Masuda<sup>1</sup>, Satoshi Yamaki<sup>1</sup>, Mami Okamoto<sup>2</sup>, Muneo Sato<sup>2</sup>, Yoshihiro Hayakawa<sup>1</sup>, Yuji Sawada<sup>2</sup>, Masaru Furuta<sup>1</sup>, Masami Y. Hirai<sup>2</sup>

<sup>1</sup> Shimadzu Corporation. 380-1 Horiyamashita, Hadano, Kanagawa 259–1304,

<sup>2</sup> RIKEN CSRS. . 1-7-22 Suehiro-cho, Tsurumi-ku, Yokohama, Kanagawa 230-0045, Japan



PO-CON1726E

# Introduction

Here we demonstrate the metabolomic analysis for carotenoids in Chrysanthemum flower as a model samples by a liquid chromatograph tandem mass spectrometry. Chrysanthemum flower (*Chrysanthemum morifolium*) is one of the most important horticulture plant and the color of flower gains control of the market value. To enhance the market vale, breed improvement has been performed by hybridizing or mutation since it is known that major determinant factors of the colors are several metabolites such as carotenoids or flavonoids.

Recently, systematical metabolomic approach has been applied to improve the breeding of plants with efficiency. Under this situation, a metabolomic approach of carotenoids was investigated to establish the metbolomic platform using liquid chromatograph tandem mass spectrometry in order to develop novel and efficient breeding process for horticulture plant.




Figure 1. Image of the Radiation breeding

# Experimental

### Standard sample and Analytical condition

Standard regents of major carotenoids were purchased from Wako Chemical (Japan). Specified carotenoids in Chrysanthemum such as epoxides compounds were obtained from the Institute of Vegetable and Floriculture Science, NARO (Tsukuba, Japan). The analysis of carotenoids was performed by a reversed phase chromatography followed by the detection using a triple quadrupole mass spectrometer, LCMS-8050 (Shimadzu, Japan) equipped with a APCI under SIM or MRM mode. Analytical condition is shown in Table 1.



HPLC : Nexera UHPL	_C system					
Column	: SunShell C18, 2.6um, 2.1 x 50mmL.					
Mobile phase	: A - Water					
	: B - Acetonitrile / 2-propanol (2/1)					
Gradient program	: 70%B(0 min) -100%B(50 - 9 min) - 100%B(9 - 14 min) -					
	70%B(14.01 - 20min)					
Flow rate	: 0.3 mL / min					
Column temperature	: 40 °C					
Injection volume	: 1.5 µL as loop injection					
MS : LCMS-8050 Triple quadrupole mass spectrometer						
Ionization	: APCI (Positive / Negative)					
lon spray voltage	: +4.5 kV, -3.5 kV					
Nebulizing Gas Flow	: 3 L/min					
Drying Gas Flow	: 5 L/min					
IF Temp.	: 350 °C					
DL Temp.	: 200 °C					
HB Temp.	: 300 °C					
SIM and MRM	:					
	Triple Quadrupole LC/MS/MS [LCMS-8050]					

Table 1. Analytical condition

SIM or MRM parameters were optimized by flow injection analysis of each carotenoid standard. Though both positive and negative ions were observed for some carotenoid, the SIM or MRM parameter was determined depending on the sensitivity and CID spectrum of each carotenoids. Optimized SIM or MRM parameter is in Table 2 with each standard retention time.

#### Metabolomic Approach for Carotenoids Analysis of Chrysanthemum flower by Liquid Chromatograph Tandem Mass Spectrometry

Compounds		+/-	Transition	CE	Retention	Start	End
ζ-carotene	met. 1	+	541.50>81.05	-39	11.069	10.158	12.158
Phytofluene	met. 2	+	543.45		11.687	10.700	12.700
Neurosporene	met. 3	+	539.50>81.00	-42	10.638	9.670	11.670
(3S,5S,6R,3'R,6'R)-5,6-Dihydro-5,6- Dihydroxylutein	met. 4	-	730.50>326.25	31	7.946	6.944	8.944
(all-E)-Lutein-5,6-epxide	met. 5	-	584.4		4.369	3.465	5.465
(9'Z)-Lutein-5,6-epxide	met. 6	-	584.35		5.338	4.338	6.338
(9Z)-Lutein	met. 7	-	568.40>550.35	23	6.354	5.357	7.357
(9'Z)-Lutein	met. 8	-	568.2		6.363	5.357	7.357
β-carotene	met. 9	+	537.40>177.35	-20	11.654	10.678	12.678
α-carotene	met. 10	+	537.40>123.15	-27	11.5	10.503	12.503
Zeinoxanthin	met. 11	-	552.35		8.989	7.990	9.990
Zeaxanthin	met. 12	-	568.35		5.327	4.322	6.322
Anthraxanthin	met. 13	-	584.35		4.371	3.305	5.305
Violaxanthin	met. 14	-	600.45		3.441	2.435	4.435
Neoxanthin	met. 15	-	600.4		3.104	2.104	4.104
Cantaxanthin	met. 16	-	564.20>472.20	22.0	6.531	5.532	7.532
Fukoxanthin	met. 17	-	658.45>59.15		2.964	1.970	3.970

Table 2. SIM or MRM parameter for each carotenoid

Although ESI or APCI can be used for carotenoids analysis, we found that APCI can give the stable result than ESI due to ion-suppression.

Thus, APCI was used for this study.

### Extraction of carotenoids

Carotenoids in each flower, which was provided by the Technology Center for Agriculture, Forestry and Fisheries in Hyogo, Nagasaki and Kagoshima prefectures, were extracted with n-Hexane following the procedure shown in Figure 2. The final supernatants obtained were injected onto the system.






Figure 2. Procedure of extraction for carotenoids

### Result and Discussion

### Difference of carotenoids in Wild type or Mutant

Obtained chromatograms of standard carotenoids and breeds (Wild type and mutants) are shown in Figure 3. Divergence in chromatograms among the breeds was found such as Zeinoxanthn, Lutain-related carotenoids as below.





Figure 3. MRM Chromatograms of carotenoids



The results suggest that radiation may have an influence for the pathway of carotenoids such as hydroxylation and epoxidation, especially for the pathway from Zeinoxanthin to Lutein-epxide in the map shown as Figure 4.




Figure 4. Metabolic pathway of major carotenoids

### Cluster analysis for the breed vs. carotenoids

Cluster analysis was performed for the series of result (breeds vs. carotenoids) by the self-developed program using the R language.

This heatmap shown in Figure 5 suggests that there are some mutants which may have several specific carotenoids. The result indicates the possibility for screening of mutant breeds using the developed analytical procedure.




Figure 5. Heatmap of correlation for breeds vs. carotenoids

# Conclusions

- Major carotenoids in Chrysanthemum flower were simultaneously identified using a LC-MS/MS system with a convenient procedure of extraction.
- Difference of the specific carotenoids in mutant suggested the influence for the metabolic pathway.
- The result of cluster analysis suggests the potential for evaluation of breeding such as mutants vs. flower color with metabolic approach of carotenoids identification.
- Further, we have been investigating the relation between color of flowers and difference of metabolic carotenoids to establish the screening frame work

## Reference


"Diversity of Carotenoid Composition in Flower Petals" A. Omiya; JARQ 45 (2), 163 – 171 (2011)



# Acknowledgment

Authors wish to thank Dr. Akemi Omiya, the Institute of Vegetable and Floriculture Science, NARO for providing the standard of specific carotenoids as well as the Technology Center for Agriculture, Forestry and Fisheries in Hyogo, Nagasaki and Kagoshima prefectures for supplying the Chrysanthemum flowers.

This research was conducted by Japan's Council for Science, Technology and Innovation, Cross-ministerial Strategic Innovation Promotion Program (SIP) "Technologies for creating next-generation agriculture, forestry and fisheries" (administrative body: NARO Bio-oriented Technology Research Advancement Institution (BRAIN).





Shimadzu Corporation

www.shimadzu.com/an/

#### For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "@". Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.