

EFFICIENT LC METHOD DEVELOPMENT

with Agilent InfinityLab Poroshell 120 columns

STEP 1 Select particle size

UHPLC

- Complex separations
- High-throughput samples

We recommend

InfinityLab Poroshell 120 1.9 μm

- Fast performance
- More runs per day
- High sensitivity
- Best resolution
- Solvent savings

UHPLC/HPLC

- Multi-purpose analysis

We recommend

InfinityLab Poroshell 120 2.7 μm

- Faster runs on UHPLC and existing HPLC systems
- For UHPLC: Higher performance, more samples analyzed per day, and solvent savings

HPLC

- Legacy methods

We recommend

InfinityLab Poroshell 120 4 μm

- Faster performance with legacy methods
- More samples analyzed per day

STEP 2 Select column chemistry

BEST ALL AROUND

InfinityLab Poroshell 120 EC-C18

1.9 μm , 2.7 μm , 4 μm

Pore Size 120 Å

InfinityLab Poroshell 120 EC-C8

1.9 μm , 2.7 μm , 4 μm

Temp Limit 60 °C

InfinityLab Poroshell 120 Phenyl-Hexyl

1.9 μm , 2.7 μm , 4 μm

pH Range 2.0-8.0

LOW-pH MOBILE PHASES

InfinityLab Poroshell 120 SB-C18, 2.7 μm

Pore Size 120 Å
Temp Limit 90 °C
pH Range 1.0-8.0

InfinityLab Poroshell 120 SB-C8, 2.7 μm

Pore Size 120 Å
Temp Limit 80 °C
pH Range 1.0-8.0

HIGH-pH MOBILE PHASES

InfinityLab Poroshell HPH-C18

1.9 μm , 2.7 μm , 4 μm

Pore Size 100 Å

Temp Limit 60 °C

InfinityLab Poroshell HPH-C8

2.7 μm , 4 μm

pH Range 3.0-11.0

ALTERNATIVE SELECTIVITY

InfinityLab Poroshell 120 Bonus-RP

2.7 μm

Pore Size 120 Å

Temp Limit 60 °C

InfinityLab Poroshell 120 PFP

1.9 μm , 2.7 μm , 4 μm

pH Range 2.0-8.0

POLAR ANALYTES

InfinityLab Poroshell 120 SB-Aq
2.7 μm

Pore Size 120 Å
Temp Limit 80 °C
pH Range 1.0-8.0

InfinityLab Poroshell 120 HILIC
1.9 μm , 2.7 μm , 4 μm

Pore Size 120 Å
Temp Limit 60 °C
pH Range 0.0-8.0

InfinityLab Poroshell 120 EC-CN
2.7 μm

Pore Size 120 Å
Temp Limit 60 °C
pH Range 2.0-8.0

Suggested columns to use as a starting point

Acidic compounds: InfinityLab Poroshell 120 SB-Aq or SB-C18

Acidic/neutral compounds: InfinityLab Poroshell 120 EC-C18

Basic compounds: InfinityLab Poroshell HPH-C18

Polar compounds: InfinityLab Poroshell 120 HILIC

Columns with Alternate Selectivity

PFP, Phenyl Hexyl, Bonus-RP, SB-CN, SB-Aq

For help with choosing a different selectivity, visit www.hplccolumns.org

Accelerate your analytical method development workflow

The Agilent InfinityLab LC Method Development Solution gives you automated access to hundreds of unique separation conditions—without changing a single column or solvent bottle.

For more information, visit www.agilent.com/chem/welivelcmetdev


Helpful tools and links

LC Column Navigator Tool: Get recommendations for a new column, based on method parameters. www.agilent.com/chem/navigator

Method Validation Kits: Convenient, cost-effective kits put columns from three separate lots right at your fingertips. www.agilent.com/chem/method_validation_kits


Find the right parts and supplies for your LC system. www.agilent.com/chem/agilentresources

Troubleshooting Chromatographic Issues




Poor peak shape, tailing, and poor resolution

- Avoid overloading the column; test different injection volumes/concentrations.
- Check for loose connections and leaks. **TIP:** Use InfinityLab QuickConnect fittings for leak-free connections.
- Basic compounds: Try a high pH mobile phase to reduce analyte interactions.
- Acidic compounds: Use ion pairing reagents to minimize analyte interactions.
- Increase the separation temperature.



Inconsistent results from run to run

- Ensure adequate column equilibration.
- Check lamp signal quality. **TIP:** Track lamp performance with an InfinityLab lamp with RFID.
- Run a standard test mixture to rule out instrument contributions to inconsistency.
- Be sure that your injection solvent is compatible with method starting conditions.



Ghost peaks

- Confirm mobile phase identity and purity. Perform a blank run to check for carryover; if carryover persists, change the needle seat and rotor seal.
- Make sure the flow cell is clean. **TIP:** Track usage with InfinityLab MaxLight flow cells with RFID.
- Know the properties of all injected components—even those not UV visible.
- Wash all components off the column with an appropriate % organic.
- Increase wash duration at the end of the run.
- Ensure the needle wash solvent strength is appropriate for sample components.
- Use an InfinityLab column ID tag to track column history; if necessary, try a new column.

STEP 3 Use this general protocol to develop fast, robust LC methods

- 1 Test initial conditions: 10-95% gradient with ACN plus 0.1% formic acid (low pH) at 30 °C.
- 2 Adjust slope (% organic over time) to incorporate the first and last peak. Optimize for analytes of interest.
- 3 If ACN gives poor peak shape or retention, change the organic modifier to MeOH.
- 4 Adjust column temperature (up to allowable limits).
- 5 Try alternate column chemistries to vary retention and selectivity.



Agilent Technologies