

My Chromatography Has Changed: Steps for Effective Troubleshooting

Jean Lane

Applications Engineer


LC Columns and Consumables Technical Support

June 18, 2024



Common Symptoms and Problems

- Pressure
- Peak shape
- Retention
- Baseline



- Increased pressure
- Low pressure
- Leaks
- Pressure fluctuations
- Tailing
- Peak splitting and doubling
- Fronting
- Broadening
- Changing retention time
- Loss of resolution
- Noisy baseline
- Drifting baseline
- Reduced intensity or sensitivity

Steps For Effective Troubleshooting

You've recognized that there is a problem

Ask questions:

- When did the system or chromatography last function properly?
- Has anything been changed?
- For the method, was the procedure followed correctly?
- Are the instrument settings correct?
- What exactly is the problem that is being seen?

Steps For Effective Troubleshooting

You've recognized that there is a problem

Considerations for where the problem might be?

- Pump
- Injector/autosampler
- Column
- Detector
- Data system
- Mobile phase
- Sample
- Tubing/fittings
- User

The problem could be one, some, many, or all of these independently, or together.

How's that for a challenge?

First Step in Troubleshooting


Understand your HPLC system

Detector

Column compartment

Autosampler

Pump



Know your flow path

Changes in System Pressure

Causes of Increases in Back Pressure


- Particles leading to blockage can come from sources located both *outside* and *inside* the LC system:
 - Solvent, buffer
 - Microbial growth in solvent reservoirs
 - The sample
 - Wear of LC components – piston seals, autosampler valve
- Debris will either be captured on the filter, frit, or inline filter (inexpensive replaceable frit), guard column, or a column frit (column = expensive)

Reduce LC problems by eliminating the most common sources of flow blockage that can cause increased pressure. Preventing this is the key.


Filter, filter, filter

Blockages and Clogging

| Characteristics | |
|--------------------------|--|
| Parts affected | Blockages: <ul style="list-style-type: none"> Capillaries, needle, and needle seat Detector flow cells Clogging: <ul style="list-style-type: none"> Filter frits (inline filter, column filter) |
| Characteristic | ● |
| Identification | <ul style="list-style-type: none"> Start by disconnecting the capillary at the column inlet Install a test setup with a restriction capillary Continue disconnecting capillaries, one-by-one, moving back toward the pump |
| Possible root cause | <ul style="list-style-type: none"> Debris from mechanically worn parts (needle seat material, rotor seal at injection valve) Coring of vial septa material |
| Instant action/first aid | <ul style="list-style-type: none"> Backflush affected part Replace part |
| Preventive measures | <ul style="list-style-type: none"> Replace wear parts in time; apply proper preventive maintenance schedules Use high-quality septa Install inline filters Use a guard column |



Blockages: instant pressure increase step



Clogging: constant pressure increase over time

Microbial Growth

- Potential problems
 - Increased system pressure or pressure fluctuations
 - Increased column pressure, premature column failure
 - Can mimic application problems
 - Gradient inaccuracies
 - Ghost peaks
- Prevent or reduce microbial growth
 - Use freshly prepared mobile phase
 - Filter
 - Do not leave mobile phase in the instrument for days without flow
 - Always discard “old” mobile phase
 - Do not add fresh mobile phase to old. No “topping off”
 - Use an amber solvent bottle for aqueous mobile phase
 - If possible, you can:
 - Add 5% organic to water – this can be used to reduce bacterial growth
 - Use a few mg/L of sodium azide





P/n 3150-0577
Solvent filter/degasser assembly



Glass solvent inlet filter (20 mm),
p/n 5041-2168
Stainless steel solvent inlet filter
p/n 01018-60028
Amber solvent bottle 1 L,
p/n 9301-6526
Clear solvent bottle 1 L
p/n 9301-6524


Why Use an Inline Filter?

Particles lead to blockage



Filter particles to prevent column clogging

Extend column lifetime and reduce cost per sample




Accelerated lifetime test shows how an inline filter removes particles

InfinityLab Quick Change Inline Filter and Filter Discs

Ultimate ease-of-use

- **Finger-tight, tool-free** replacement of filter disc
- **Click and seal:** A click alerts users when the filter is tight up to 1300 bar, assuring no risk of over- or under-tightening




[InfinityLab Quick Change Inline Filters](#)

High efficiency, easy-to-use filter discs


- **Various dimensions and porosities** – filter discs are available in 2.1 mm and 4.6 mm inner diameters with different pore sizes. The filter housing is compatible with all types of filter discs.
- **Touchless packaging to avoid potential contamination** – with specially designed packaging, you're able to insert the filter disc into the filter housing without touching it, to avoid potential contamination.
- **In situ replacement** of filter disc – no need to disconnect the inline filter from the system
- **Smart alert** to remind users when filter discs need replacing



Different dimensions and porosities of filter discs



Filter discs in touchless packaging





No-touch insertion of filter disc into filter housing

Why Filter the LC Sample

- Capillaries, frits, and the column inlet are less likely to end up with blockages
- Less wear and tear of injection and switching valves
- Less downtime

Agilent Syringe Filter Selector tool
[Captiva Syringe Filter Selector | Agilent](#)




| Applications | | |
|--|-------------------|--|
| Type of Filtration | Recommended | Alternatives |
| HPLC • UHPLC • LC/MS • GC | RC | PTFE or Nylon |
| ICP-MS | PTFE | Glass Fiber/PTFE (High Particle Samples) |
| CE | RC | Nylon |
| Undiluted Organic Solvents | PTFE | Nylon |
| Protein Analysis • Samples with Biomolecules – Buffers | PES | RC or CA |
| Tissue Culture Media | PES | RC or CA |
| High Particle-Load Samples – Organic Solvents | Glass Fiber/PTFE | - |
| High Particle-Load Samples – Aqueous Solutions | Glass Fiber/Nylon | - |

Filtration

Captiva premium syringe filters

Captiva syringe filters guide
Pub no: [5991-1230EN](#)

Column lifetime test



Sample filtration
improves column lifetime

The impact of filtering a 0.3 µm latex bead suspension on the lifetime of a sub-2 µm column

Filtering helps to reduce clogging, which can lead to increased pressure problems

Agilent technical note: [5994-1947EN](#)

Guard Columns



Agilent Guard Cartridge and holder



[Video: Fast Guards for HPLC and UHPLC | Agilent](#)

InfinityLab Fast Guard



Guards protect your column in many ways

Poroshell Column + Poroshell Guard Sample: Infant Formula* (1:300 in Water)



*Unfiltered infant formula including proteins and other precipitated ingredients.

Changes in System Pressure


Low pressure

| | Potential Cause | Recommended Action |
|---|---------------------------------|---|
| ● | Leak in high-pressure flow path | <ul style="list-style-type: none">• Visual inspection of the flow path• Instrument diagnostic tests LA |
| ● | Wrong mobile phase | <ul style="list-style-type: none">• Check for correct mobile phase• Check solvent reservoir and tube connections |

Helpful Troubleshooting Tool





LA: With its advanced diagnostic and maintenance capabilities, **Agilent Lab Advisor SW** helps you to keep your Agilent analytical instruments in top condition. Agilent Lab Advisor is independent of the chromatography software you are using.



Leaks

| Characteristics | |
|--------------------------|---|
| Characteristics | <ul style="list-style-type: none"> Lower pressure Potentially impacting retention times and peak shape |
| Parts affected | <ul style="list-style-type: none"> Potentially all parts in the flow path High potential at frequently operated fitting connections (for example, column inlet) and parts with high mechanical stress (rotor seal, needle, and needle seat) |
| Identification | <ul style="list-style-type: none"> Drops of solvent or residues of salt System diagnostic tests LA |
| Possible root cause | <ul style="list-style-type: none"> Loose or bad fitting connections Cracked capillaries Worn needle and needle seat |
| Instant action/first aid | <ul style="list-style-type: none"> Replace affected parts Renew or redo fitting connections |
| Preventive measures | <ul style="list-style-type: none"> Use proper fitting connections Replace fittings and wear parts in time |




Changes in System Pressure

Pressure fluctuations

| | Potential Cause | Recommended Action |
|---|---------------------------|--|
| ● | Air in the system | <ul style="list-style-type: none">• Prime and flush the instrument• Check for sufficient solvent supply• Check for correct plumbing (SSV/MCGV)• Check for correct degassing |
| ● | Malfunctions at pump head | <ul style="list-style-type: none">• Perform pump head diagnostic tests LA• Replace defective parts• Implement a proper maintenance schedule |
| ● | Cavitation effects | <ul style="list-style-type: none">• Check for flow restrictions (solvent bottle to pump head inlet)• Clean or replace parts• Verify that the solvent supply is positioned above the pump inlet |

In addition


Pressure fluctuations will typically also impact the UV signal due to refractive index effects.




Peak Shape Changes

Problems with Peak Shape


- Tailing
- Fronting
- Broadening
- Splitting/doubling




Tailing



Fronting



Broadening




Splitting/doubling

Changes in Peak Shape

Peak tailing

| If Applicable to Some Peaks | | Recommended Action |
|-----------------------------|---|--|
| | Secondary interactions | <ul style="list-style-type: none"> • Change pH • Change stationary phase |
| | Small peak eluting on the tail of a larger peak | <ul style="list-style-type: none"> • Change selectivity (column, mobile phase) • Switch to methods with higher resolution (UHPLC, 2D-LC) |


| If Applicable to All Peaks | | Recommended Action |
|----------------------------|---|---|
| | Silica-based – column degradation | <ul style="list-style-type: none"> • Use specialty, polymeric, or sterically-protected column |
| | Silica-based – basic interactions with stationary phase | <ul style="list-style-type: none"> • Use a stronger mobile phase or add an appropriate base (for example, TEA) |
| | Poor tubing connections; high dispersion volume | <ul style="list-style-type: none"> • Minimize the number of connections • Check the connections/fitting condition and proper seat of fittings • Use fittings with spring-load function |



Peak Tailing

Mobile phase-related factors – pH

The effect of pH on peak shape at or near the sample pKa



Column: ZORBAX SB-C8, 4.6 x 150 mm, 5 μ m
p/n: 883975-906



Mobile phase: 40% 5 mM KH_2PO_4 , 60% ACN

Flow rate: 1.0 mL/min

Temperature: Ambient

Peak Tailing

Mobile phase-related factors – mobile phase additives



Columns: Eclipse XDB-C8, 4.6 x 150 mm, 5 μ m, p/n: 993967-906 Mobile phase: 85% 25 mM Na_2HPO_4 : 15% ACN pH: 7 Flow rate: 1.0 mL/min
Temperature: 35 $^{\circ}\text{C}$ Sample: Amphetamines 1. Phenylpropanolamine 2. Ephedrine 3. Amphetamine 4. Methamphetamine 5. Phenteramine


Incorrect concentration of TEA in mobile phase negatively effects peak shape of basic compounds

Peak Tailing


Column contamination

Column: StableBond SB-C8, 4.6 x 250 mm, 5 µm, p/n: 880975-902 Mobile phase: 20% H₂O : 80% MeOH
Flow rate: 1.0 mL/min Temperature: RT Detection: UV 254 nm


QC test forward direction



QC test reverse direction



QC test after cleaning
100% IPA, 35 °C



Sample:
1. Uracil
2. Phenol
3. 4-Chloronitrobenzene
4. Toluene



Peak Tailing

Importance of having correct connections

Connection problems can lead to:

Poor chromatography

- Broad or tailing peaks
- Loss of resolution



InfinityLab Quick Connect and Quick Turn Fittings


- Spring-loaded design
- Easy-to-use
- Works for all column types
- Reusable
- Consistent ZDV connection

Quick Connect fitting


- Finger-tight up to 1300 bar
- Hand tighten the nut, then depress the lever

Quick Turn fitting

- Finger-tight up to 400 bar
- Up to 1300 bar with a wrench
- Compact design




Capillary for Quick Connect fitting



InfinityLab Quick Turn fitting

Capillary for Quick Turn fitting






Brochure: [5991-5164EN](#)

Changes in Peak Shape

Fronting


| Potential Cause | Recommended Action |
|----------------------|---|
| Channeling in column | <ul style="list-style-type: none"> • Replace the column • Use guard columns |
| Column overload | <ul style="list-style-type: none"> • Decrease sample amount • Use a higher capacity column (increase length, diameter, or change to high-capacity material) |



Changes in Peak Shape


Peak broadening

| Potential Cause | Recommended Action |
|---|---|
| Injection volume is too large | <ul style="list-style-type: none"> Decrease the injection volume |
| Long retention times | <ul style="list-style-type: none"> Use gradient elution or a stronger mobile phase |
| System settings | <ul style="list-style-type: none"> Check the data collection rate Adjust the detector setting or time constant to the fastest possible value without compromising signal-to-noise |
| Viscosity of the mobile phase is too high | <ul style="list-style-type: none"> Increase the column temperature |
| Detector cell volume is too large | <ul style="list-style-type: none"> Use the smallest possible cell volume |
| Improper fittings/connections | <ul style="list-style-type: none"> Ensure that your fitting connections are correct |
| Extra tubing volume on the system | <ul style="list-style-type: none"> Ensure that the tubing is narrow and as short as possible to avoid extra volume |
| Sample diluent is too strong | <ul style="list-style-type: none"> Reduce diluent strength |




Peak Broadening

Influence of data rate



Peak Broadening

Extracolumn dispersion (Volume)



For minimizing dispersion:

- Use short, small internal diameter tubing between the injector and the column and between the column and the detector.
- Make certain all tubing connections are made with matched fittings.
- Use a low-volume detector cell.
- Inject small sample volumes.

| Length | 10mm | 50mm | 100mm | 150mm |
|----------------|----------|--------|---------|---------|
| Tubing ID | Volume | Volume | Volume | Volume |
| 0.17mm (green) | 0.227 uL | 1.1uL | 2.27 uL | 3.3 uL |
| 0.12mm (red) | 0.113 uL | 0.55uL | 1.13 uL | 1.65 uL |

Retention

What Is the Specific Issue?

- Retention times of all peaks shift
- Retention time of only one peak shifts


Increase the detail to be more specific:

- Retention time of all peaks shift
- Retention time of all peaks shift earlier
- Retention time of all peaks shift to earlier times and the extent of the shift appears to be the same

Changes in Separation


Retention time changing

| Potential Cause | Recommended Action |
|---|--|
| Flow rate changing | Check "Pressure fluctuation", pump flow rate |
| Inconsistent online mobile phase mixing | Ensure gradient system is delivering constant composition check vs. manual preparation of mobile phase |
| Column temperature varying | Thermostat column and ensure constant lab temperature |
| Equilibration time insufficient with the gradient run or a change in isocratic mobile phase | Flush with at least 10 column volumes after solvent change or gradient conclusion |
| Selective evaporation of mobile phase component | Keep solvent reservoirs covered Prepare fresh mobile phase |
| Buffer capacity insufficient | Use >20 mM concentration of buffer |
| Contamination buildup | Occasionally flush the column with a strong solvent to remove contaminants |
| First few injections – adsorption on active sites | Condition the column using an initial injection of a concentrated sample |
| Column overloaded with sample | Decrease the injection volume or concentration |
| Active sites on silica packing | Add a competing base to mobile phase |
| Mobile phase composition is changing | Follow the 'best practices' |




Change in Volatile Buffer Concentration

Shift in Retention Time and Peak Shape



Mobile Phase Preparation

Agilent 1100 with quaternary pump
Column: Zorbax Eclipse XDB-C8 RR 3.5 μ m, 4.6 x 50mm
p/n 935967-906



Dial-a-Mix = A: water B: MeOH, pump 50% B

Vol/Vol = 250mL water + 250mL MeOH, pump 100%

Mix-to-Mark = 250mL MeOH, fill to 500 mL with water, pump 100%

Premixed (w/w) = 200g MeOH + 200g water, pump 100%

Detection: UV 254nm


Flow rate: 1 ml/min

Changes in Separation

Ghost peaks, carryover

| | Potential Cause | Recommended Action |
|---|---|--|
| ● | Peaks from previous injections | <ul style="list-style-type: none"> Flush the column to remove contaminants Check with blank injection |
| ● | Specific interaction with metal surfaces | <ul style="list-style-type: none"> Passivate instrument Use InfinityLab deactivator additive Use bio-inert LC equipment |
| ● | Contamination or unknown interferences in samples | <ul style="list-style-type: none"> Proper sample cleanup |
| ● | Ion pair – disequilibrium | <ul style="list-style-type: none"> Prepare sample in actual mobile phase to minimize disturbance |


**BIO
INERT**




[P/n 5191-4506 | Deactivator additive 50 mL](#)

[P/n 5191-3940 | Deactivator additive 25 mL](#)

Ghost Peaks



The LC system was equilibrated using starting conditions for 30 minutes, then a gradient run was made. Impurities were trapped and eluted out with the gradient.



Ghost Peaks


Where do they come from

- Organic
- Additives
 - TFA
 - Salts
- H₂O
- Sample or from a previous run
- Other
 - Glassware
 - pH meter
 - Filters


Unknown Phantom Peaks

Solvent contamination

Injections on Agilent 1100



Solvent: acetonitrile



Column Cleaning

Do what's recommended for your column

Tips for cleaning columns

- Flush with stronger solvents than your mobile phase
- Make sure the detector is taken out of the flow path
- Do not add your organic solvent directly to the buffer, as this may cause the buffer salts to precipitate out and lead to more backpressure



For reversed phase

Use at least 10 column volumes of each solvent for analytical columns

1. Start with your mobile phase, without buffer salts (water/organic)
2. 100% organic (MeOH or ACN)
3. Check the pressure to see if it has returned to normal; if not, then
4. Discard the column or consider more drastic conditions: 75% acetonitrile/25% isopropanol
5. 100% isopropanol
6. 100% methylene chloride, solvent wash for very nonpolar compounds
7. Hexane

*Always see your specific column user guide for instructions

[LC Column User Guides | Agilent](#)

LC Columns Are Not Indestructible


- Columns are packed using hydraulic pressure and can be damaged by excess pressure
- Silica dissolves (slowly) at higher pH
- Acid hydrolysis of bonded phase can occur at low pH
- Column failure
 - Void
 - Contamination

Columns must be stored properly

- Check your user guide

It's important to:

- Know the technical specifications for your column
- Choose a mobile phase that is right for your column
- Keep a record/history of your column



| | |
|--------------|---|
| Columns: | 4.6 x 150 mm, 5 µm |
| Purge: | 50% ACN/50% 0.02 M K ₂ HPO ₄ , pH 11 |
| Flow rate: | 1.5 mL / min |
| Temperature: | 25 °C |
| Detection: | Silicate concentration by silicomolybdate color reaction |

[LC Column User Guides | Agilent](#)

Choice of Your Column

Low and high pH can cause column failure

The InfinityLab Poroshell 120 portfolio offers choices for low and high pH

| Best All Around | Best for Low pH Mobile Phases | Best for High pH Mobile Phases | Best for Alternative Selectivity | Best for More Polar Analytes | HILIC for polar analytes | Chiral |
|--|---------------------------------------|--|-------------------------------------|---------------------------------------|--|----------------------------|
| EC-C18 1.9 µm, 2.7 µm, 4 µm | SB-C18 1.9 µm, 2.7 µm, 4 µm | HPH-C18 1.9 µm, 2.7 µm, 4 µm | Bonus-RP* 2.7 µm | Aq-C18* 2.7 µm | HILIC 1.9µm, 2.7 µm, 4 µm | Chiral-V 2.7 µm |
| EC-C8 1.9 µm, 2.7 µm, 4 µm | SB-C8 2.7 µm | HPH-C8 2.7 µm, 4 µm | PFP* 1.9 µm, 2.7 µm, 4 µm | SB-Aq* 1.9 µm, 2.7 µm, 4 µm | HILIC-Z 1.9 µm, 2.7 µm, 4 µm | Chiral-T 2.7 µm |
| Phenyl-Hexyl* 1.9 µm, 2.7 µm, 4 µm | | ← CS-C18 → 2.7 µm | | EC-CN* 2.7 µm | HILIC-OH5 2.7 µm | Chiral-CD 2.7 µm |
| | | | | | | Chiral-CF 2.7 µm |
| | | | | | | |


Every new column should be tested on your instrument

Performance verification based on Agilent checkout

- Run Agilent checkout before use
 - Record the difference between your instrument and the performance report (use as a base value)
- Run again if the column seems to lose performance
 - Compare with the results from first run

Performance verification based on in-house checkout

- Run in house checkout before use
 - Record key specifications, such as tailing factor, plates, and backpressure
- Run again if the column seems to lose performance
 - Compare with the results from first run





Baselines and Detection

Changes in Detection

Noisy baseline



| | Potential Cause | Recommended Action |
|---|---|--|
| ● | Gas bubbles in the mobile phase | <ul style="list-style-type: none"> • Apply degassing • Check the degasser performance |
| ● | Low difference between the sample and the mobile phase absorbance | <ul style="list-style-type: none"> • Check absorbance values of the sample versus the mobile phase |
| ● | Contamination | <ul style="list-style-type: none"> • Use degassed HPLC-grade solvents • Flush the system • Clean up the sample |
| ● | Detector optics | <ul style="list-style-type: none"> • Perform an intensity test • Check the signal with the flow cell removed if possible • Replace the lamp |
| | Pressure instability | <ul style="list-style-type: none"> • Check “Pressure fluctuation” |



Changes in Detection


Drifting baseline

| | Potential Cause | Recommended Action |
|---|---|--|
| ● | Contamination in the mobile phase | <ul style="list-style-type: none"> • Make up new mobile phase • If running a gradient, you might need to adjust the modifier |
| ● | Low difference between the sample and the mobile phase absorbance | <ul style="list-style-type: none"> • Check absorbance values of the sample versus the mobile phase |
| ● | Contamination | <ul style="list-style-type: none"> • Use degassed HPLC-grade solvents • Flush the system • Clean up the sample |
| ● | Detector | <ul style="list-style-type: none"> • Check the temperature stability • Check for leaks • Replace the lamp |
| | Pressure instability | <ul style="list-style-type: none"> • Check "Pressure fluctuation" |




Reference Wavelength

How settings can affect the baseline



Gradient: 0.02 m KH₂PO₄/ACN, from 12% ACN to 45% ACN in 12 min

Influence of Data Collection Rate on Noise



LC Troubleshooting Poster Available

LC Troubleshooting Guide

Your guide to solving common problems and staying productive

Agilent
InfinityLab

Places to Start

Solvents

- Use brown borosilicate bottles to avoid algae growth
- Prepare solvent volume to be used up within 1 to 2 days
- Use only HPLC-grade solvents filtered through 0.2 µm filters

Preparing and powering up the pump

- Inspect solvent bottles and inlet filters for damage or coloring
- Always use seal wash when installed and purge the pump
- Use the appropriate system conditioning method

Daily tasks

- Replace aqueous and organic mobile phases every second day
- Check seal wash solvent
- Flush the system with the composition of your application

Weekly tasks


- Change seal wash solvent and bottle and inspect solvent filters
- Check system backpressure and change filters if necessary

Pump shutdown

- Flush all channels to remove salt deposits and particulate matter
- Flush the system with appropriate storage solvent and power down the system

Handling of acetonitrile

- If possible, use 5 to 10% of water in your mobile phase
- Be sure to avoid ACN evaporation
- Don't leave ACN on the system for more than 2 to 3 days
- Perform a periodic warm water wash (60 to 70 °C) if you face problems



Maintenance

Agilent Lab Advisor software helps you manage your Agilent LC instruments to achieve high-quality chromatographic results in the most efficient way by ensuring high instrument performance, productivity, and reliability. It is available free-of-charge.

- Diagnostic tests to evaluate performance
- Easier maintenance of all Agilent LC modules
- Comprehensive reports generated to ease communication with Agilent service

Retention Time Drift



| Possible Cause | Solution |
|---|---|
| Inconsistent online mobile phase mixing | Ensure gradient system delivers constant composition; compare with manual preparation of mobile phase |

| Possible Cause | Solution |
|---------------------------------|--|
| Variation in column temperature | Thermostat or insulate column; ensure constant lab temperature |

| Possible Cause | Solution |
|---|---|
| Insufficient equilibration time with gradient run or change in isocratic mobile phase | Make sure at least 10 column volumes pass through column after sample run |

| Possible Cause | Solution |
|---|--|
| Selective evaporation of mobile phase component | Less vigorous helium sparging; keep solvent reservoirs covered; prepare fresh mobile phase |

| Possible Cause | Solution |
|-----------------------|---|
| Contamination buildup | Occasionally flush column with strong solvent |

| Possible Cause | Solution |
|-------------------------------|--|
| Column overloaded with sample | Decrease injection volume or concentration |

Pressure Fluctuation



| Possible Cause | Solution |
|--------------------|---|
| Leak in the system | Identify the channel and clean or replace check valve; replace pump seals |

| Possible Cause | Solution |
|-------------------------|--------------------------------|
| Buildup of particulates | Filter sample and mobile phase |

| Possible Cause | Solution |
|----------------|---|
| Bubble in pump | Perform solvent degassing; sparge solvent with helium |

Pressure Increase



| Possible Cause | Solution |
|-----------------|---|
| System blockage | Check flowpath (needle seat, capillaries, filter and frits) |

| Possible Cause | Solution |
|---|--|
| Water/organic systems: buffer precipitation | Test buffer organic mixtures to ensure compatibility |

High Column Backpressure



| Possible Cause | Solution |
|-----------------|---|
| Column blockage | Better sample cleanup; use guard column |

| Possible Cause | Solution |
|---------------------------------|--|
| Mobile phase viscosity too high | Use lower viscosity solvents or higher temperature |

| Possible Cause | Solution |
|-------------------------|-----------------------------------|
| Particle size too small | Use larger d _p packing |

| Possible Cause | Solution |
|--------------------|----------------|
| Plugged inlet frit | Replace column |

Drifting Baseline




| Possible Cause | Solution |
|--|--|
| Positive/negative direction: contaminant buildup/elution | Flush column; clean up sample; use pure solvents |

| Possible Cause | Solution |
|--|-------------------------------------|
| Positive/negative: difference in refractive index of injection solvent | Use mobile phase for sample solvent |

| Possible Cause | Solution |
|---------------------|---|
| Temperature changes | Insulate and thermostat column and tubing |

Noisy Baseline



| Possible Cause | Solution |
|----------------|---|
| Contamination | Use degassed HPLC-grade solvents; flush system; clean up sample |

| Possible Cause | Solution |
|-------------------|--|
| Detector problems | Check number of hours of UV lamp; replace UV lamp or flow cell |

Ghost Peaks



| Possible Cause | Solution |
|-------------------------------|---|
| Peaks from previous injection | Flush column to remove contaminants; check with blank injection |

| Possible Cause | Solution |
|---|-----------------------|
| Contamination; unknown interferences in samples | Proper sample cleanup |

| Possible Cause | Solution |
|--------------------------|---|
| Ion pair: disequilibrium | Prepare sample in actual mobile phase to minimize disturbance |

| Possible Cause | Solution |
|---------------------------|-------------------------|
| Contaminated mobile phase | Check your mobile phase |

| Possible Cause | Solution |
|--------------------|-------------------------------|
| Bubbles in solvent | Check and degas your solvents |

Peak Tailing



| Possible Cause | Solution |
|----------------------|--|
| Unwiped dead volumes | Minimize number of connections; ensure injector seal is tight; ensure fittings are properly sealed |

| Possible Cause | Solution |
|--------------------|-------------------------------------|
| Column performance | Change mobile phase; replace column |

| Possible Cause | Solution |
|---------------------------------|--|
| Silica-based column degradation | Use specialty, polymeric, or sterically protected column |

| Possible Cause | Solution |
|--|---|
| Silica-based: basic interactions with stationary phase | Use stronger mobile phase or add appropriate base (e.g., TEA) |

Peak Broadening



| Possible Cause | Solution |
|----------------------------|--|
| Injection volume too large | Decrease injection volume or solvent strength of injection solvent; use gradient methods |

| Possible Cause | Solution |
|----------------------------------|--------------------|
| Low sampling rate of data system | Increase data rate |

| Possible Cause | Solution |
|--------------------------------|-----------------------------------|
| Detector cell volume too large | Use smallest possible cell volume |

| Possible Cause | Solution |
|----------------------------|---------------------------|
| Injection volume too large | Decrease injection volume |

Sensitivity Problems



| Possible Cause | Solution |
|--|---|
| Peaks are outside of sensitivity range of detector | Dilute/concentrate sample to bring into linear region |

| Possible Cause | Solution |
|--|---|
| Sample-related losses during preparation | Use internal standard during sample preparation; optimize sample preparation method |

Leaks



| Possible Cause | Solution |
|--|---------------------------------------|
| White powder at fitting/ loose fitting | Tighten fittings; replace capillaries |

| Possible Cause | Solution |
|----------------|--|
| System leak | Identify location checking leak sensors/sensors; check flow cell |

Discover more best practices for using an Agilent LC system:
<https://www.agilent.com/chem/lc-best-practices>



Training courses are available at:
<https://www.agilent.com/crosslab/university>



Get answers. Share insights. Join the Agilent Community at:
<https://community.agilent.com>



For Lab Advisor software, please visit:
<https://www.agilent.com/chem/lab-advisor>



Resources for Support

- Column user guides: [LC Column User Guides | Agilent](#)
- LC Troubleshooting poster: [LC Troubleshooting Guide 5994-0709EN](#)
- Resource page: <http://www.agilent.com/chem/agilentresources>
 - Quick reference guides
 - Catalogs, column user guides
 - Online selection tools, how-to videos
- InfinityLab LC Supplies catalog: [5991-8031EN](#)
- LC Handbook: [5990-7595EN](#)
- YouTube – [Agilent channel](#) (maintenance videos)
- Consumables Community: [Agilent Collection of Columns, Supplies, and Standard Consumables - Agilent Community](#)
- App finder: [Application Finder | Agilent](#)
- Agilent University: [Agilent University](#)
- Your local product specialists
- Agilent Peak Tales podcasts: peaktales.libsyn.com
- Webinars, upcoming and recorded: [LC and LC/MS Column Webinars | Agilent](#)





AGILENT INFINITI/LR SUPPLIES FOR THE AGILENT 1260 INFINITI & LC WITH MULTISAMPLER

Quick reference guide

Agilent supplies for Agilent instruments

Agilent Technologies is committed to supporting our laboratory customers, as we have produced this list of the most commonly ordered Agilent supplies for the 1260 Infinity I/LC with Multisampler.

This Infiniti/LR supply list is an optional portfolio of LC components, columns, and supplies designed to work together seamlessly to maximize efficiency and performance.



| Infiniti/LR System Parts | Part # | MSRP |
|---|------------|-------------|
| Agilent 1260 Infinity I | G1315A-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC | G1315B-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315C-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315D-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315E-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315F-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315G-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315H-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315I-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315J-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315K-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315L-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315M-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315N-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315O-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315P-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315Q-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315R-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315S-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315T-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315U-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315V-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315W-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315X-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315Y-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315Z-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AA-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AB-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AC-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AD-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AE-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AF-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AG-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AH-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AI-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AJ-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AK-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AL-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AM-00 | \$11,400.00 |
| Agilent 1260 Infinity I, LC with Multisampler | G1315AN-00 | \$11,400.00 |



Contact Agilent Chemistries and Supplies Technical Support



Available in the U.S. and Canada, 8-5 all time zones

1-800-227-9770 option 3, option 3:

Option 1 for GC and GC/MS columns and supplies

Option 2 for LC and LC/MS columns and supplies

Option 3 for sample preparation, filtration, and QuEChERS

Option 4 for spectroscopy supplies

Option 5 for chemical standards

Option 6 for Prozyme products



gc-column-support@agilent.com

lc-column-support@agilent.com

spp-support@agilent.com

spectro-supplies-support@agilent.com

chem-standards-support@agilent.com

pzi.info@agilent.com



Agilent

Trusted Answers