

## Cannabidiol Oil Analysis with the Pyroprobe

Application Note  
Pharmaceutical

Author:


Karen Sam

### Abstract

This application note demonstrates multi-step analysis of CBD oil with RSDs of cannabidiol

With recent legislation paving the way for entertainment usage of marijuana, and the promise of cannabidiol (CBD) for treating a variety of ailments, CBD oil, derived from the cannabis plant, has a growing interest in nutraceutical and pharmaceutical industries. Similar in structure to psychoactive  $\delta$ 9-Tetrahydrocannabinol (THC) found in marijuana, this non-intoxicating extract is being credited with helping treat many medical issues. Diluted with hemp oil prior to use, CBD oil is a complex natural product, containing many volatile and non-volatile constituents. Analytical testing using the Pyroprobe can clarify ingredients in natural materials such as CBD oil by separating ingredients based on their volatility, then pyrolyzing the non-volatile portion, like the oil itself.

About 500 micrograms of CBD oil purchased from a therapeutic hemp company, was placed into a Drop-In-Sample-Chamber (DISC) tube and run using a multi-step sequence on a CDS 6000 Series Pyroprobe interfaced to a GC/MS. The Pyroprobe was programmed to heat first to 200°C, and then to 400°C to analyze semi-volatiles in the oil; then finally to 700°C to qualify pyrolysis products. Before each of these heating steps, the Pyroprobe waits for the GC to become ready, then starts the GC as the heating begins.



Chromatograms of all 3 multi-step runs are shown in Figure 1. At 200°C, some of the oil starts to vaporize, exhibiting as an unresolved mixture at the end of the chromatogram. Along with this, active compounds in the oil, like the sesquiterpene  $\alpha$ -Bisabolol (with anti-irritant, anti-inflammatory, and anti-microbial properties),  $\alpha$ -Caryophyllene (often found in aromatic plants), and the cannabinoid, CBD (active ingredient) vaporize. Then at 400°C, remaining CBD continues to vaporize, along with some fatty acids, alcohols like olivetol and phytol (a diterpene alcohol), and additional cannabinoids, including THC, whose amount must be small enough to pass regulations. Finally, at 700°C, the remainder of the oil pyrolyzes, breaking down the remaining triglycerides of the oils into long chain alkenes, alkanes, alkynes, aldehydes, and alcohols. A more detailed list of peak search results are found in Tables 1-3.



The worth of quantitative results is contingent on the reproducibility of the approach. This relies on instrument precision, as well as issues like consistent sample preparation and homogeneity. A series of 20 firings of a Pyroprobe filament at 1100°C produces an average measured temperature of 1100.15°C with a relative standard deviation of only 0.04%, establishing that Pyroprobe instruments perform with the highest precision. As previously observed, CBD oil contains a wide array of compounds. So fortunately, reproducibility seen with the Pyroprobe proceeds into complex samples such as CBD oil. As CBD's boiling point is between 160-180°C, a preliminary setpoint of 200°C for 20 minutes was used to vaporize cannabidiol while minimizing interference from other constituents in

### Experimental Parameters

Samples were pyrolyzed in a DISC tube, using a CDS Model 6200 Pyroprobe.

### Multi-step Analyses

Direct-Py Mode

Pyroprobe:

DISC Chamber: 200, 400, 700°C 30s

Interface: 300°C

Transfer Line: 325°C

Valve Oven: 300°C

### GC/MS

Column: 5% phenyl (30m x 0.25mm)

Carrier: Helium, 50:1 split

Injector: 320°C

Oven: 40°C for 2 minutes

10°C/min to 300°C

hold 15 minutes

Ion Source: 230°C

Mass Range: 35-600amu

### RSD Analyses

Trapping Mode

Pyroprobe:

DISC Chamber: 200°C 20min

Trap Rest: 50°C

Trap Final: 300°C 10min

Interface: 300°C

Transfer Line: 300°C

Valve Oven: 300°C

### GC/MS

Column: 5% phenyl (30m x 0.25mm)

Carrier: Helium 1.25mL/min, 16:1 split


Injector: 300°C

Oven: 80°C for 10 minutes

10°C/min to 300°C


Ion Source: 230°C

Mass Range: 35-600amu




the oil for RSD trials. One half a microliter of the CBD in hemp oil diluted in half with hexane and extracted to a 200°C setpoint produced an RSD for identified cannabidiol at around 3% (Figure 2).

Performing multi-step thermal analysis on both pharmaceuticals and natural products which make up pharmaceuticals, can help an analyst find distinct components, providing valuable information for competitive analysis and product development. Additionally, reliable results for thermal extraction of active and illegal ingredients like CBD and THC respectively can be achieved.



| Peak # | Identification       | 6  | Olivetol                |
|--------|----------------------|----|-------------------------|
| 1      | Caryophyllene        | 7  | Phytol                  |
| 2      | Humulene             | 8  | n-Hexadecanoic acid     |
| 3      | α-Bisabolol          | 9  | Linoleic acid           |
| 4      | Cannabidiol          | 10 | δ1-Tetrahydrocannabinol |
| 5      | 1,3,8-p-Menthatriene |    |                         |

Figure 1. CBD in hemp oil, 200°C (top), 400°C, and 700°C (bottom).



| Rep   | Area     | Counts | RSD |
|-------|----------|--------|-----|
| Rep 1 | 24649718 |        |     |
| Rep 2 | 23491890 |        |     |
| Rep 3 | 24151015 |        |     |
| Rep 4 | 24630748 |        |     |
| Rep 5 | 25671887 | 3.26%  |     |

Figure 2. Five analyses of CBD oil, 200°C, m/z 231.2.

Table 1: 200°C Library Search Results

| Ret.Time<br>(min) | Top Hit   |
|-------------------|---|
| 9.74              | Linalool  |
| 11.23             | a-Terpineol   |
| 11.32             | Decanal   |
| 13.71             | a-ylangene  |
| 14.40             | Caryophyllene   |
| 14.50             | a-Bergamotene   |
| 14.64             | a-ylangene  |
| 14.68             | cis-β-Farnesene   |
| 14.86             | Humulene  |
| 14.92             | Alloaromadendrene   |
| 15.05             | β-copaene   |
| 15.20             | Patchoulene   |
| 15.28             | g-Selinene  |
| 15.66             | Patchoulene   |
| 15.80             | a-Bisabolene  |
| 15.86             | β-Guaiene   |
| 15.95             | Columbin  |
| 16.02             | Formic acid, 3,7,11-trimethyl-1,6,10-dodecatrien-3-yl ester       |
| 16.08             | N,N'-Bis(Carbobenzyloxy)-lysine methyl(ester)                     |
| 16.42             | Caryophyllene oxide   |
| 16.54             | 5-Azulenemethanol, 1,2,3,4,5,6,7,8-octahydro-a,a,3,8-tetramethyl- |
| 16.69             | Ethyl iso-allocholate   |
| 16.89             | β-Guaiene   |
| 17.00             | a-acorenol  |
| 17.18             | Pregan-20-one, 2-hydroxy-5,6-epoxy-15-methyl-                     |
| 17.33             | Columbin  |
| 17.51             | a-Bisabolol   |
| 17.82             | 5,8,11,14-Eicosatetraynoic acid                                   |
| 18.66             | Phenanthrene  |
| 20.27             | n-Hexadecanoic acid   |
| 21.39             | Ethyl iso-allocholate   |
| 22.00             | 7,8-Epoxylostan-11-ol, 3-acetoxy-                                 |
| 22.26             | Astaxanthin   |
| 22.88             | Arenobufagin  |
| 24.45             | Cannabidiol   |
| 24.51             | 7,8-Epoxylostan-11-ol, 3-acetoxy-                                 |
| 24.58             | Cinobufotalin   |
| 32.21             | Lycoxanthin   |

Table 2: 400°C Library Search Results

| Ret.Time<br>(min) | Top Hit   |
|-------------------|---|
| 5.02              | Hexane, 2,4-dimethyl-                                 |
| 5.25              | 2-Octene  |
| 6.48              | 2-Methyl-1-hepten-3-yne                               |
| 6.78              | Heptanal  |
| 7.54              | 2-Heptenal, (Z)-                                      |
| 7.72              | 2,6-Dimethyl-1,3,5,7-octatetraene, E,E-               |
| 8.05              | Cyclohexene, 5-methyl-3-(1-methylethenyl)-, trans(-)- |
| 8.64              | 2,6-Dimethyl-1,3,5,7-octatetraene, E,E-               |
| 8.72              | D-Limonene  |
| 9.59              | Cyclohexene, 1-methyl-4-(1-methylethylidene)-         |
| 9.65              | Benzene, 1-methyl-4-(1-methylethenyl)-                |
| 9.85              | 13-Heptadecyn-1-ol                                    |
| 10.00             | p-Mentha-1,5,8-triene                                 |
| 10.18             | Cyclopentene,1-hexyl-                                 |
| 11.06             | 6-Dodecene, (Z)-                                      |
| 14.40             | Caryophyllene   |
| 14.96             | 5-Hexadecyne  |
| 15.06             | 9,17-Octadecadienal, (Z)-                             |
| 15.28             | b-Selinene  |
| 15.36             | a-Selinene  |
| 15.42             | β-Bisabolene  |
| 15.80             | a-Bisabolene  |
| 15.86             | β-Guaiene   |
| 15.95             | 3,5,11-Eudesmatriene-                                 |
| 16.42             | Caryophyllene oxide                                   |
| 17.00             | Aromadendrene oxide-(2)                               |

Table 2: 400°C Library Search Results -cont'd.

| Ret.Time<br>(min) | Top Hit                                       |
|-------------------|---|
| 17.05             | Caryophylla-4(12),8(13)-dien-5a-ol            |
| 17.11             | a-acorenol                                    |
| 17.33             | Isoaromadendrene epoxide                      |
| 17.40             | β-Guaiene                                     |
| 17.51             | a-Bisabolol                                   |
| 17.90             | 1-Hexadecanol, 2-methyl-                      |
| 18.21             | 1,3-Benzenediol, 5-pentyl-                    |
| 18.46             | Z,Z-3,15-Octadecadien-1-ol acetate            |
| 18.61             | Pregan-20-one, 2-hydroxy-5,6-epoxy-15-methyl- |
| 19.00             | 1-Heptatriacotanol                            |
| 19.06             | Neophytadiene                                 |
| 19.11             | 2-Pentadecanone, 6,10,14-trimethyl-           |

Table 3: 700°C Library Search Results

| Ret.Time<br>(min) | Top Hit   |
|-------------------|---|
| 1.71              | Carbon dioxide  |
| 2.12              | 1,3-Pentadiene  |
| 2.24              | Bicyclo[2.1.0]pentane                                       |
| 2.44              | 1-Hexene  |
| 2.71              | (Z),(Z)-2,4-Hexadiene                                       |
| 2.85              | 1,4-Cyclohexadiene  |
| 2.97              | (Z),(Z)-2,4-Hexadiene                                       |
| 3.07              | Benzene   |
| 3.16              | 1,3-Cyclopentadiene, 1-methyl-                              |
| 3.24              | trans-1,4-Hexadiene   |
| 3.29              | Bicyclo[3.1.0]hexane  |
| 3.39              | 1-Heptene   |
| 3.94              | 1-Methylcyclohexa-2,4-diene                                 |
| 4.31              | 3,4-Heptadiene  |
| 4.88              | 1-Octene  |
| 5.03              | Hexane, 2,4-dimethyl-                                       |
| 5.26              | 4-Octene, (E)-  |
| 5.38              | Cyclopentene, 3-propyl-                                     |
| 5.70              | Tricyclo[3.2.1.0(1,5)]octane                                |
| 5.89              | 1,4-Heptadiene, 3-methyl-                                   |
| 5.94              | 2,4-Octadiene   |
| 6.04              | E,Z-4-Ethylidenecyclohexene                                 |
| 6.07              | 2,4-Octadiene   |
| 6.10              | Ethylbenzene  |
| 6.19              | E,Z-4-Ethylidenecyclohexene                                 |
| 6.56              | 2-Nonene  |
| 6.65              | 3-Nonene, (E)-  |
| 6.74              | Cyclooctene, (Z)-   |
| 7.07              | 1,3-Nonadiene, (E)-   |
| 8.07              | 1-Decene  |
| 9.38              | 3-Heptenoic acid  |
| 9.61              | 1-Undecene  |
| 13.89             | 1-Hexadecanol   |
| 20.37             | n-Hexadecanoic acid   |
| 21.95             | (Z)-18-Octadec-9-enolide                                    |
| 22.19             | 9,12-Octadecadienoic acid (Z,Z)-                            |
| 24.40             | 9,12-Octadecadienoic acid (Z,Z)-, 2,3-dihydroxypropyl ester |

FOR MORE INFORMATION CONCERNING THIS APPLICATION,  
WE RECOMMEND THE FOLLOWING READINGS:

CDS Application Note #184a Pyrolysis GC-MS of Pharmaceutical Packaging

CDS Application Note #105a Cracking Products of Oleic Acid and Olive Oil

Sam, K. American Lab Multistep Thermal Characterization of Liquid-Filled Capsules and Medication Packaging Using GC/MS. Am. Lab, April 2018.