

# In-Depth Characterization of Monoclonal Antibodies Using Intact Mass Analysis and Middle-down Approaches on Orbitrap Ascend BioPharma Tribrid Mass Spectrometer

Jingjing Huang, Christopher Mullen, Kristina Srzentić, Graeme McAlister, Vlad Zabrouskov, Rafael Melani  
Thermo Fisher Scientific, San Jose, California, USA

## Abstract

**Purpose:** To assess the performance of the Thermo Scientific™ Orbitrap™ Ascend Tribrid™ mass spectrometer for comprehensive monoclonal antibody (mAb) characterization using intact and middle-down mass spectrometry (MS) approaches.

**Methods:** The middle-down methods were developed using a combination of ion activation techniques including electron transfer dissociation (ETD), electron transfer higher energy collision dissociation (EThcD), ultraviolet photodissociation (UVPD) and proton-transfer charge reduction (PTCR).

**Results:** In-depth sequence coverage of mAb subunits was achieved using the advanced MS/MS techniques on Orbitrap Ascend BioPharma Tribrid MS.

## Introduction

In-depth characterization of mAbs and their post-translational modifications (PTMs) is critical to ensure the safety and efficacy of biotherapeutics. Intact mass analysis can be used to obtain accurate masses of mAbs.<sup>1,2</sup> The utilization of the middle-down MS technique has gained traction as an encouraging method for the characterization of biotherapeutics.<sup>3-5</sup> In this work, we performed intact mass analysis and middle-down MS to characterize Trastuzumab and its subunits on Orbitrap Ascend BioPharma Tribrid mass spectrometer equipped with Native MS option, using different MS/MS fragmentation approaches.

## Materials and methods

### Sample Preparation


Trastuzumab was aliquoted into 10 µg/µL stock solution for native LC-MS experiments. For intact mass analysis under denaturing conditions with reverse phase column, the antibody was diluted with 0.1% formic acid to a series of concentrations at 10 ng/µL, 50 ng/µL, 100 ng/µL, 500 ng/µL and 1 µg/µL. For reverse phase LC-MS/MS middle-down experiments, Trastuzumab subunits were prepared at 1 µg/µL final concentration using DTT for reduction combined with or without previous IdeS digestion according to manufacturer's protocol.

### Methods

Native MS experiments were performed using a size exclusion column and an isocratic gradient with 50 mM ammonium acetate. Intact mass analysis (denatured) and middle-down MS experiments were conducted using a reverse-phase column. For intact mAb analysis under the native and denaturing conditions, an Orbitrap resolving power setting of 30,000 at 200 *m/z* was used. An Orbitrap resolving power setting of 240,000 at *m/z* 200 was used to acquire isotopically resolved full MS spectra of the Fc/2, light chain (LC) and Fd' subunits and setting of 7,500 at 200 *m/z* was used to acquire full MS spectra of heavy chain (HC) subunit. For middle-down analysis of all subunits, a quadrupole isolation window of 100 *m/z* centered around *m/z* 900 was used. Replicate ETD, EThcD and UVPD data were acquired using five different reaction/activation times.

### Data Analysis


All data analyses were performed using BioPharma Finder™ software version 5.2.



**Figure 1. Major LC-MS components used for biopharmaceutical characterization.** Experiments were using the Thermo Scientific™ Vanquish™ Flex UHPLC system and Orbitrap Ascend Tribrid mass spectrometer, together with the MAbPac™ reverse phase (P/N 088648) and SEC-1 size exclusion (P/N 077592) columns for LC separation, and BioPharma Finder software for data processing.


## Results

### Overview




**Figure 2. Illustration of trastuzumab sample preparation for intact mass and middle-down analyses.** Top – intact mAb was diluted and analyzed under native and denaturing conditions; Middle – IdeS digestion and TCEP/DTT reduction to yield Fc/2, LC and Fd' subunits; Bottom – TCEP/DTT reduction to generate LC and HC subunits.

### Intact mass analysis




**Figure 3. Intact mass analysis of trastuzumab under the denaturing and native conditions.** A) Full MS spectra acquired from intact trastuzumab under the denaturing (top) and native (bottom) conditions. The insets show an expanded view of the most abundant charge state with baseline resolved glycoforms. B) Deconvoluted masses measured for intact trastuzumab under the denaturing (top) and native (bottom) conditions using the ReSpect and Sliding Window algorithm




**Figure 4. Intact mass analysis of trastuzumab dilution series under the denaturing condition.** A) Full MS spectra of denatured trastuzumab with 10 ng to 1 µg loaded on column; B) Zoomed in spectra of the 49+ charge state with baseline resolved glycoforms; C) Deconvoluted results of trastuzumab obtained using BioPharma Finder ReSpect sliding window algorithm.

### Subunit mass analysis




**Figure 5. Subunit analysis of Fc/2, LC, Fd' subunits from IdeS digestion of trastuzumab.** A-C) Full MS spectra of the Fc/2, LC and Fd' subunits acquired at a resolution setting of 240,000 and 2 µscans. The insets show the isotopically resolved charge state at *m/z* ~900 Da (z=28 for Fc/2 and Fd', and z=26 for LC); D) Total ion chromatogram of Fc/2, LC and Fd' subunits; E-F) Deconvoluted results of three trastuzumab subunits obtained using the Xtract and sliding windows algorithm with BioPharma Finder software. Excellent mass accuracies (<3 ppm) were obtained for all three subunits.



**Figure 6. Subunit analysis of the LC and HC subunits from trastuzumab.** A) Full MS spectrum of the LC subunit acquired at a resolution setting of 240,000 and 2 µscans. The inset shows the isotopically resolved charge state (z=25) around *m/z* 939; B) Full MS spectrum of HC subunit acquired at a resolution setting of 7,500 and 10 µscan. The inset shows the isotopically unresolved charge state (z=55) around *m/z* 922; C) Total ion chromatogram of the LC and HC subunit showing baseline separation of two subunits; D) Excellent mass accuracy of <3 ppm was obtained for the LC subunit from the deconvolution of isotopically resolved data using the Xtract algorithm within BioPharma Finder software; E) Excellent mass accuracies of <3 ppm were obtained for three major glycoforms of the HC subunit from the deconvolution of the isotopically unresolved data using the ReSpect and sliding windows algorithm within BioPharma Finder software.


### Middle-down MS analysis




**Figure 7 (A-B). Middle-down analyses of trastuzumab subunits using ETD, EThcD, and UVPD.** A) Bar chart showing sequence coverage of the Fc/2, LC and Fd' subunits from IdeS digestion of trastuzumab using ETD, EThcD, and UVPD, with 10 files combined. B) Bar chart showing sequence coverage of the LC and HC subunits from reduction of trastuzumab using ETD, EThcD, and UVPD, with 10 files combined.



**Figure 7 (C-H). Middle-down analyses of trastuzumab subunits using ETD, EThcD, and UVPD.** C-E) Fragmentation maps of the Fc/2, LC and Fd' subunits obtained from five experiments using ETD, EThcD, and UVPD, with two replicates for each condition. F-H) Fragmentation maps of the LC and HC subunits obtained from five experiments using ETD, EThcD, and UVPD, with two replicates for each condition. High sequence coverage was obtained using all fragmentation techniques. EThcD and UVPD provided an improved sequence coverage than ETD due to the formation of additional fragment types.



**Figure 8. Combined PTCLR ion activation with other fragmentation methods benefits residue coverage by decluttering the congested MS2 spectra.** A) bar chart representing the improved residue coverage for Fc/2 subunit from 72% to 80% for EThcD fragmentations, and 66% to 75% for UVPD fragmentations, from two replicates of the selected experiment condition as labeled in the figure. When combined all the replicate results from experiments with and without PTCLR, EThcD experiments achieved 87% coverage for Fc/2 subunit and UVPD experiments achieved 89% coverage. B) fragmentation map of Fc/2 subunit for EThcD fragmentation from 5 ms reagent reaction time with two replicates (left). C) fragmentation map of Fc/2 subunit from EThcD 4ms plus PTCLR 30ms method with two replicates (middle) and D) cumulative fragmentation map of Fc/2 subunit from (B) and (C) (right). E) fragmentation map of Fc/2 subunit for UVPD fragmentation from 20 ms activation time with two replicates (left). F) fragmentation map of Fc/2 subunit from UVPD 15ms plus PTCLR 30ms method with two replicates (middle) and G) cumulative fragmentation map of Fc/2 subunit from (E) and (F) (right).



**Figure 9. Combined results of EThcD, UVPD, and PTCLR middle-down analyses provided excellent sequence coverage of trastuzumab subunits.** A.) Bar chart showing a nearly complete sequence coverage of Fc/2, LC and Fd' subunits obtained from the combined results of EThcD, UVPD and PTCLR raw files (10 in total); B) Bar chart showing a nearly complete sequence coverage of the LC subunit (97%) and high sequence coverage of the HC subunit (73%) when combining the raw files of EThcD, UVPD and PTCLR (10 in total); C) Fragmentation maps of the Fc/2, LC and Fd' subunits from the combined results of 10 raw files; D) Fragmentation maps of the LC and HC subunits from the combined results of 10 raw files.

## Conclusions

- The high mass accuracy and high resolution offered by the Orbitrap technology provides precise mass measurement for both intact mAbs (native and denaturing) and their subunits
- Compared to ETD, EThcD and UVPD provide an improved sequence coverage of mAb subunits due to the formation of additional fragment types. The combination of PTCLR with EThcD and UVPD led to further increase in sequence coverage of mAb.
- The combination of EThcD, UVPD and PTCLR offers a nearly complete sequence coverage (>90%) for subunits in the size of 20-25 kDa (Fc/2, LC and Fd' subunits) and high sequence coverage (>70%) for the HC subunit (~50 kDa)

## References

- Srzentic K, Damoc E, Scheffler K, Buchanan T, Radi K, Sutton J. Complete characterization of monoclonal antibodies under native and denaturing conditions. Thermo Fisher Scientific Application Note 73885 (2021).
- Scheffler K, Damoc E. Full characterization of heterogeneous antibody samples under denaturing and native conditions on the Q Exactive BioPharma mass spectrometer. Thermo Fisher Scientific Application Note 72348 (2017).
- Jake T, Kline, Christopher Mullen, Kenneth R. Durbin, Ryan N. Oates, Romain Huguet, John E. P. Syka, and Luca Fornelli. Sequential Ion-Ion Reactions for Enhanced Gas-Phase Sequencing of Large Intact Proteins in a Tribrid Orbitrap Mass Spectrometer. Journal of the American Society for Mass Spectrometry 2021 32 (9), 2334-2345.
- Romain Huguet, Christopher Mullen, Kristina Srzentić, Joseph B. Greer, Ryan T. Fellers, Vlad Zabrouskov, John E. P. Syka, Neil L. Kelleher, and Luca Fornelli. Proton Transfer Charge Reduction Enables High-Throughput Top-Down Analysis of Large Proteoforms. Analytical Chemistry 2019 91 (24), 15732-15739.
- Scheffler, K; Damoc E; Bailey, A; Josephs J. Enabling Mass Spectrometric Analysis of Intact Proteins in Native Conditions of A Hybrid Quadrupole-Orbitrap Mass Spectrometer, ASMS Poster 2016.

## Trademarks/licensing

© 2024 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. PO190-2024-EN