

Poster Reprint

ASMS 2026
WP 221

Rapid LC-MS/MS compound method development and sample analysis for DMPK workflows

Bianca Silva¹, Mike Knierman¹, Steve Ainley², Jeremy Lawton², Wayne Lootsma²

¹Agilent Technologies

²Sound Analytics LLC

Introduction

Bioanalysis – DMPK – studies

Drug metabolism and pharmacokinetics (DMPK) studies are essential to understanding the absorption, distribution, metabolism, and excretion (ADME) of drug candidates and play a critical role in drug discovery and development. The characterization of DMPK properties for modern therapeutics increasingly requires advanced analytical techniques due to long half-lives and low systemic concentrations. DMPK labs need a workflow that includes robust and easy to maintain instrumentation, easy to use and efficient software to enable the analysis of hundreds of compound studies a week. Rapid review of the results is necessary to deliver high quality results in a timely fashion. Training requirements of new lab personnel needs to be minimized on the operation of this workflow to keep lab operations running well.

Bioanalysis - Workflow

Bioanalysis are based on synthetic molecule testing *in vivo* and *in vitro* assays. Depending on the assay, more than 1000 of injections including ~ 700 compounds screening are performed on weekly basis. Liquid chromatography–tandem mass spectrometry (LC-MS/MS) has become a cornerstone of DMPK research, offering high sensitivity, selectivity, and throughput for the quantitative analysis of diverse compound classes. A simplified workflow is showed below.




Figure 1. General workflow for Molecules in Bioanalysis studies.

The Agilent solution for DMPK/Bioanalysis analysis

Agilent 6495D triple quadrupole LC/MS systems, combined with the **Sound Analytics LeadScape** solution, are designed to streamline method development and support robust, high-quality bioanalytical assays across multiple modality types. These capabilities enable efficient evaluation of DMPK properties and facilitate the optimization of drug candidates during early discovery. **LeadScape for Agilent** software delivers the speed of developing compound specified methods, rapid sample analysis and review, and enterprise method management in an easy – user friendly interface. Additionally, the 6495D LC/TQ does not require any special settings to switch between peptide and small molecules methods.

Experimental

Methods

An Agilent Infinity II Bio UHPLC coupled to an Agilent 6495D triple quadrupole LC/MS (LC/TQ) running MassHunter and LeadScape software was used for the analysis of five compounds and one internal standard. Compound specific methods were rapidly developed using LeadScape. Samples were prepared using eight calibration levels in extracted commercial biological matrix samples along with several spiked with a blinded concentration.



Advantages:

- ✓ Rapid MRM optimization
- ✓ Rapid methods and worklists creation
- ✓ Easy data review
- ✓ Enterprise-wide compound database

MRM Optimization

- Fast LC method template
- Compound name and formula/MW
- 2 min per compound
- Storage in local or global Database
- Only performed once

Text file with compound info

	A	B	C	D	E
1	Compound ID	Vial Position	Plate Code	Plate Location	Formula
2	Sulfadiazine	B1	*96Agilent*	P1	C10H10N4O2S
3	Sulfamerazine	B2	*96Agilent*	P1	C11H12N4O2S
4	Sulfamethazole	B3	*96Agilent*	P1	C10H11N3O3S

Fast LC method

Time [min]	A [%]	B [%]	Flow [mL/min]
0.00	30.00	70.00	0.800
0.03	30.00	70.00	0.800
0.04	30.00	70.00	0.100
0.32	30.00	70.00	0.100
0.33	30.00	70.00	0.800

Review MRM transitions



Automated Method Creation

- Uses MRMs from database to merge with LC-MS/MS method template
- Only need compound ID and positions
- Creates 100s of methods in < 1min from database.
- Creates automated worklists.

Setup
Tuning

Setup
Analyze

Queue

Review
MS Tune

Review
ChromaTune

Review
Analyze

Edit
Databases

Fast Chromatography

Time [min]	A [%]	B [%]	Flow [mL/min]
0.00	90.00	10.00	0.400
0.20	90.00	10.00	---
1.20	0.00	100.00	---
1.30	0.00	100.00	---
1.35	90.00	10.00	---

Text file with analysis info

Sample Name	SampleID	Acquisition Method	Plate Code	Plate Location	Val Position	Set Name	Data File Name	Quant Type	Compound ID	Concentration
ASMS_Sample01_Blank1	SampleID_01	ASMS_0514	*96Agilent*	P1	A1	Set1	Data_0514	Blank	Sulfadiazine	0
ASMS_Sample02_Std1	SampleID_02	ASMS_0514	*96Agilent*	P1	C1	Set1	Data_0514	Standard	Sulfadiazine	6.75
ASMS_Sample03_Std2	SampleID_03	ASMS_0514	*96Agilent*	P1	C2	Set1	Data_0514	Standard	Sulfadiazine	12.5
ASMS_Sample04_Std3	SampleID_04	ASMS_0514	*96Agilent*	P1	C3	Set1	Data_0514	Standard	Sulfadiazine	25
ASMS_Sample05_Std4	SampleID_05	ASMS_0514	*96Agilent*	P1	C4	Set1	Data_0514	Standard	Sulfadiazine	50
ASMS_Sample06_Std5	SampleID_06	ASMS_0514	*96Agilent*	P1	C5	Set1	Data_0514	Standard	Sulfadiazine	100
ASMS_Sample07_Unknown	SampleID_21	ASMS_0514	*96Agilent*	P1	A2	Set1	Data_0514	Unknown	Sulfadiazine	0
ASMS_Sample08_Blank2	SampleID_22	ASMS_0514	*96Agilent*	P1	A1	Set1	Data_0514	Blank	Sulfadiazine	0

Methods created

Method Information																																											
Method Name	ASMS_0514																																										
Start Method	D:\Projects\SoundAnalytics\Methods\ASMS_New_Starter_ASMS.m																																										
Compounds	<table><thead><tr><th>CompoundID</th><th>Compound List</th><th>IS Index</th></tr></thead><tbody><tr><td>Sulfadiazine</td><td></td><td>4</td></tr><tr><td>Sulfadiazine</td><td></td><td>4</td></tr><tr><td>Sulfamerazine</td><td></td><td>4</td></tr><tr><td>Sulfamerazine</td><td></td><td>4</td></tr><tr><td>Sulfamethazole</td><td></td><td>4</td></tr></tbody></table>	CompoundID	Compound List	IS Index	Sulfadiazine		4	Sulfadiazine		4	Sulfamerazine		4	Sulfamerazine		4	Sulfamethazole		4																								
CompoundID	Compound List	IS Index																																									
Sulfadiazine		4																																									
Sulfadiazine		4																																									
Sulfamerazine		4																																									
Sulfamerazine		4																																									
Sulfamethazole		4																																									
Conditions	<table><thead><tr><th>Polarity</th><th>Q1</th><th>Q3</th><th>EP</th><th>EP</th><th>CE</th><th>CEP</th></tr></thead><tbody><tr><td>Positive</td><td>251.1</td><td>156.0</td><td>140</td><td>10</td><td>15</td><td>11</td></tr><tr><td>Positive</td><td>251.1</td><td>92.0</td><td>140</td><td>10</td><td>25</td><td>11</td></tr><tr><td>Positive</td><td>251.1</td><td>65.0</td><td>140</td><td>10</td><td>60</td><td>11</td></tr><tr><td>Positive</td><td>251.1</td><td>92.1</td><td>140</td><td>10</td><td>35</td><td>11</td></tr><tr><td>Positive</td><td>254.1</td><td>156.0</td><td>140</td><td>10</td><td>15</td><td>11</td></tr></tbody></table>	Polarity	Q1	Q3	EP	EP	CE	CEP	Positive	251.1	156.0	140	10	15	11	Positive	251.1	92.0	140	10	25	11	Positive	251.1	65.0	140	10	60	11	Positive	251.1	92.1	140	10	35	11	Positive	254.1	156.0	140	10	15	11
Polarity	Q1	Q3	EP	EP	CE	CEP																																					
Positive	251.1	156.0	140	10	15	11																																					
Positive	251.1	92.0	140	10	25	11																																					
Positive	251.1	65.0	140	10	60	11																																					
Positive	251.1	92.1	140	10	35	11																																					
Positive	254.1	156.0	140	10	15	11																																					

Automated worklists created

Type	Batch Name	Settings	Completion	Live Plot
ASMS	ASMS_Sample01_Blank1	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample02_Std1	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample03_Std2	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample04_Std3	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample05_Std4	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample06_Std5	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample07_Unknown	ASMS_0514	ASMS_0514	ASMS_0514
ASMS	ASMS_Sample08_Blank2	ASMS_0514	ASMS_0514	ASMS_0514

Data Review

- Thousands of runs reviewed for less than 2 hours
- Screening review or quantitation analysis available.

Review results




Figure 2. LeadScape for Agilent workflow.

Results and Discussion

From Methods creation to Review data – Few steps: multiple runs

The combination of automated method generation and the high sensitivity of the 6495D LC/TQ instrument resulted in accurate and precise quantitative data at low concentration levels commonly encountered in DMPK studies. Sample acquisition and processing were completed quickly, demonstrating the efficiency of the software-driven workflow. Overall, these results highlight the advantages of combining advanced mass spectrometry instrumentation with intelligent software solutions to support fast, reliable, and high-quality data generation for DMPK applications.



Figure 3. Easy to review individual chromatograms vs entire batch and calibration curves.

Conclusions

- **User-friendly workflow:** A fast, reliable, and efficient bioanalysis workflow generated high-quality data using LeadScape for Agilent.
- **Excellent performance:** The 6495D LC/TQ demonstrated the sensitivity and the speed necessary for DMPK laboratories.
- **High Throughput:** The integration of LeadScape software with the Agilent LC/MS 6495D LC/TQ system enabled rapid optimization and acquisition of MRM transitions, significantly reducing method development time. Compound specific acquisition methods were generated quickly and applied to the study set efficiently, allowing multiple compounds to be analyzed within a single system setup. This approach eliminated time-consuming manual steps typically associated with method setup and sample organization.

<https://www.agilent.com/en/promotions/asms>

This information is subject to change without notice.

DE-014949

© Agilent Technologies, Inc. 2026

Published in USA, June 30, 2026