

Poster Reprint

ASMS 2024
Poster number TP 646

Flash Characterization of mAbs Using a Combination of Reagents via Automated On-Line Microdroplet Reaction

Jim Lau, Hui Zhao, [Mike Knierman](#)

Agilent Technologies, Inc. Santa Clara, CA USA

Introduction

Recently, an unmodified Agilent Jet Stream ESI source (AJS) was used to demonstrate “Flash Characterization” of mAbs by microdroplet facilitated reductions or enzyme digestions with good reproducibility and high yield ¹. These results motivated us to evaluate treatments that combined enzymatic digestion, reduction, and/or deglycosylation of mAb in a single microdroplet experiment. New optimization parameters for the Agilent Injection Program on Agilent 1290 Infinity II LC were presented. Microdroplet experiments, including GlycINATOR[®] deglycosylation, the combination of IdeS cleavage and DTT reduction and IdeS cleavage and GlycINATOR[®] deglycosylation were examined for reproducibility, intact mAb remaining after treatments and simplified subunit characterization. GlycINATOR is an endoglycosidase that cleaves N-linked glycans and leaves the core GlcNAc residue with or without a fucose on the asparagine. Additional antibodies of Herceptin and sigmamab are also studied with the combination treatment in the microdroplet reaction.

Experimental

The mAbs, enzymes, and DTT were automatically combined and injected in flow injection (FIA) mode using an Agilent Injection Program on Agilent 1290 Infinity II LC. The UHPLC system was coupled with an Agilent 6545XT AdvanceBio LC/Q-TOF system using the AJS in ESI positive mode for ultrafast microdroplet reaction. Spectra of the intact mAbs, digested, reduced or deglycosylated fragments were collected using the Q-TOF and analyzed with Agilent MassHunter Quantitative Analysis (for Q-TOF) and BioConfirm v12 software.



Agilent 1290
Infinity II LC
with 6545XT
AdvanceBio
LC/Q-TOF
System

Experimental

Antibody and Reagent Preparation

- All mAbs including NIST IgG1 were diluted to 0.5 mg/mL in 5 mM Ammonium Bicarbonate (ABC)
- IdeS (FabRICATOR[®]) and Glycinator[®] were diluted to 1 unit/ μ L in 5mM ABC
- DTT was prepared at 10 mg/mL in 5 mM ABC

Chromatographic Conditions

UHPLC: Agilent 1290 Infinity II

Flow injection analysis (FIA)

Column oven temperature: Ambient

Injection volume: 0.1 to 0.5 μ L

Autosampler: 5 \pm 2 $^{\circ}$ C

Mobile Phase A: 5 mM Ammonium Bicarbonate (ABC)

Gradient

| Time, min | Flow Rate, mL/min | %A |
|-----------|-------------------|-----|
| 0 | 0.3 | 100 |
| 0.1 | 0.3 | 100 |
| 0.2 | 0.025 | 100 |
| 1.9 | 0.025 | 100 |
| 2.0 | 0.3 | 100 |
| 3.5 | 0.3 | 100 |

Injector Program

| Function | Parameter |
|----------|--------------------------------------|
| Draw | Draw 1.0 or 2.0 μ L from reagent |
| Draw | Draw 0.1 to 0.5 μ L from mAb |
| Mix | Mix air and repeat five times |
| Wait | Wait 0.5 min |
| Draw | Draw 1.0 or 2.0 μ L from reagent |
| Mix | Mix air and repeat five times |
| Remote | Set remote line "Start" for 125 ms |
| Wait | Wait 0.5 min |
| Inject | Inject |


Optimized MS Conditions-Agilent 6545XT LC/Q-TOF

| Parameters | |
|------------------------|------------------|
| Drying gas temperature | 350 $^{\circ}$ C |
| Drying gas flow | 13 L/min |
| Nebulizer gas | 60 psi |
| Sheath gas temperature | 400 $^{\circ}$ C |
| Sheath gas flow | 12 L/min |
| Capillary voltage | 4500 to 5000 V |
| Nozzle voltage | 2000 V |
| Ion mode | AJS ESI Positive |
| Fragmentor | 380V |
| Skimmer | 45V |
| MS range | m/z 1350-10000 |
| Acquisition rate/Time | 2 spectra/s |

Microdroplet GlycINATOR® Deglycosylation²

- GlycINATOR® (EndoS2) is an endoglycosidase from *Streptococcus pyogenes* that specifically hydrolyzes glycans at the Fc N-glycosylation sites of IgG.

Microdroplet GlycINATOR® Deglycosylation on Herceptin




Reference


¹Lau, J, et al. Flash Characterization of Antibodies via Microdroplet Reactions in an Unmodified Jet Stream Source. Agilent Application Note 5994-6752EN, 2023

²Genovis, <https://www.genovis.com/smartenzymes/antibody-deglycosylation/glycinator/>


Microdroplet GlycINATOR® Deglycosylation on NIST IgG1




Microdroplet GlycINATOR® Deglycosylation on Sigmamab




Microdroplet IdeS GlycINATOR® Reaction on NIST IgG1




Microdroplet IdeS GlycINATOR® Reaction on Herceptin



Microdroplet FabRICATOR®(IdeS) / DTT Reaction on NIST IgG



Microdroplet FabRICATOR®(IdeS) / DTT Reaction on Herceptin



Conclusions

- Flash Characterization of mAbs with an automatic injector, combination treatments, and microdroplet reaction in the AJS source demonstrates simplification, reproducibility, and cost reduction of protein characterization.
- Combination treatments are viable with microdroplet reactions to accelerate the reduction, cleavage, and deglycosylation of native IgG1 significantly.