

Investigating cellular metabolism on Orbitrap Excedion Pro mass spectrometer by the simultaneous use of differently labeled elemental isotope substrates

Rahul Deshpande¹, Ayush Midha², Claire Dauly³, Bashar Amer¹, Isha Jain² and Susan S. Bird¹

¹Thermo Fisher Scientific, San Jose, CA, USA; ²Gladstone Institute, UCSF, CA, USA; ³Thermo Fisher Scientific, Bremen, DE

Abstract

Purpose:

In this study, a Thermo Scientific™ Orbitrap™ Excedion™ Pro hybrid mass spectrometer is accessed for use in stable isotope workflow.

Methods:

Pano2 cells were grown under different oxygen tensions (21% and 1.5%). Cells were unlabeled or labeled with ¹³C6 Glucose, ¹⁵N2 Glutamine or both ¹³C6 Glucose and ¹⁵N2 Glutamine. Metabolite separation was carried out on a Mixed Mode HILIC (HILIC-AEX) Column coupled to a Thermo Scientific™ Vanquish™ Horizon UHPLC system. Data acquisition was conducted on an Orbitrap Excedion Pro MS, designed for extended dynamic range, enhanced sensitivity, and minimal in-source fragmentation. Thermo Scientific™ Compound Discoverer™ 3.4 software was used for the data processing.

Result

The importance of resolution and mass accuracy is showcased for stable isotope workflows with Orbitrap Excedion Pro MS. Differences are shown in the cells grown at different oxygen conditions.

Introduction

Stable isotope labeling (SIL) in metabolomics is a notable tool to elucidate metabolic pathway associations (tracer analysis) and rate of change of metabolites (flux analysis). Analytical methods employing SIL typically focus on defined pathways of interest employing a targeted approach, yet the interconnected network of metabolic pathways may involve unappreciated pathways. Further, SIL is used to increase confidence in unknown identification and to determine true sample related mass spectral features in untargeted analysis. As the role of SIL in untargeted metabolomics and targeted analysis expands, we explore the importance of resolution and mass accuracy for these analysis.

Tissue Hypoxia




Figure 1. Tissue hypoxia refers to a condition where tissues in the body do not receive enough oxygen to function properly and is involved in various pathologies.

Materials and Method

Panc02 cells

Panc02 cells are derived from mouse pancreatic cancer cells and can be used to study the metabolic adaptability of cancer cells to hypoxia. The cell line was grown under normoxia (21% O₂) and hypoxia (1.5% O₂) with the following media:

1. Unlabeled Media
2. U-¹³C Glucose
3. ¹⁵N₂ Glutamine
4. U-¹³C Glucose and ¹⁵N₂ Glutamine

Metabolite Extraction and Data Acquisition

Metabolites were extracted from harvested cells using methanol followed by sonication. The methanol extract were dried and resuspended in 50% Acetonitrile,. Data acquisition was carried out using different mass resolution settings on the Orbitrap Excedion Pro MS.

Data analysis

Data processing was performed using Compound Discoverer 3.4 software.

Orbitrap Excedion Pro Mass Spectrometer

The new quadrupole-Orbitrap mass spectrometer.




Figure 2. Orbitrap Excedion Pro MS and its schematics.

Mass Resolution




Figure 3. Detection of isotopologue pattern and fine isotope structure is obtained from use of high resolution. Here Adenosine triphosphate (ATP) which is an important energy metabolite is shown. It has a molecular formula of C₁₀H₁₆N₅O₁₃P₃. Resolution of 240K or more is required for the determination of fine isotopic pattern which can be achieved on Orbitrap Excedion Pro MS.

Mass Accuracy




Figure 4. High Mass Accuracy is a prerequisite for accurate determination of isotopologues. Here ATP molecule labeled differentially with ¹³C and ¹⁵N is shown The high mass accuracy (<1ppm) of the instrument can predict the molecular formula of the isotopologues with high precision. 1. Samples labeled with ¹³C glucose. 2. Samples labeled with ¹⁵N glutamine. 3. Samples labeled with both ¹³C glucose and ¹⁵N glutamine. A i: Unlabeled sample. B: Normoxia (21% O₂) C: Hypoxia (1.5% O₂)

Compound Discoverer 3.4 Software




Figure 5. Workflow for SIL in Compound Discoverer 3.4 software




Figure 6. Compound Discoverer 3.4 software automatically detects isotopologues based on compound formulae (unlabeled sample) file




Figure 7. Compound Discoverer 3.4 software plots the ¹³C relative abundances of isotopologues of detected compounds. Here we see a difference in the labeling of proline in normoxia vs hypoxia., with the proline in hypoxia incorporating less ¹³C compared to normal conditions.

Conclusion

We have successfully showed the utility of Orbitrap Excedion Pro MS with its high resolution and accurate mass for stable isotope applications using Compound Discoverer 3.4 software. We applied this to study hypoxia in pancreatic cancer cell lines.

Trademarks/licensing

© 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage the use of these products in any manner that might infringe the intellectual property rights of others.