

High-sensitive multi-attribute analysis of ADCs under native conditions by using an online multiple heart-cut 2D-LC-HRAM mass spectrometry system

Frank Steiner¹, Maria Gruebner¹, Xuepu Li², Xiaoxi Zhang², Min Du³

¹Thermo Fisher Scientific, Germering, Germany, ²Thermo Fisher Scientific, Shanghai, China, ³Thermo Fisher Scientific, Lexington, MA, US

Abstract

Purpose: To develop a multi-attribute analysis for cysteine-linked antibody-drug conjugates (ADCs) using a two-dimensional liquid chromatography (2D-LC) system coupled with high-resolution accurate mass (HRAM) mass spectrometry (MS).

Methods: A Thermo Scientific™ Vanquish™ 2D UHPLC system, coupled with a Thermo Scientific™ Orbitrap Exploris™ 240 mass spectrometer with Biopharma option, were utilized to configure the multiple heart-cut 2D-LC-HRAM MS system. HIC-SEC and SCX-SEC methods were developed and optimized for the charge variants, drug-to-antibody ratio (DAR) and drug load distribution (DLD) characterization of ADCs under native conditions.

Results: The configured system successfully characterized multi-attributes, such as charge heterogeneity, DAR, and DLD of polatuzumab vedotin, a cysteine-linked ADC, under native conditions with high sensitivity.

Introduction

ADCs are a rapidly expanding class of biotherapeutics for cancer treatment. Critical quality attributes such as charge heterogeneity, DAR, and drug distribution are essential for ensuring the safety, stability, and efficacy of ADCs. Strong cation exchange (SCX) and hydrophobic interaction chromatography (HIC) are conventional methods employed for charge variant analysis and DAR determination of ADCs¹. However, these methods frequently utilize MS-incompatible salts, limiting their use for MS characterization. Online 2D-LC addresses this by enabling desalting in the second dimension (2D) and offering multi-attribute analysis in one method. Here, multiple heart-cut 2D-LC-HRAM MS methods were developed to achieve high sensitivity for the multi-attribute analysis of polatuzumab vedotin under native conditions.

Materials and methods

Sample Preparation

Commercial product polatuzumab vedotin, dissolved in deionized water with a concentration of 5.0 mg/mL, was provided by the customer and used directly in the experiment.

Test Methods

Figure 1 shows the configuration of the 2D-LC-HRAM MS system, Table 1 shows the UHPLC and MS conditions.

Figure 1. Overview of the experimental study design for polatuzumab vedotin multi-attribute analysis using the online multiple heart-cut 2D UHPLC-HRAM MS system.




Table 1: UHPLC and MS conditions.

1D SCX and HIC chromatographic conditions:	
Columns	SCX method: Thermo Scientific™ ProPac™ 3R SCX column, 100 × 4 mm, 3 μm (P/N 43103-104068) HIC method: Thermo Scientific™ MabPac™ HIC-Butyl, 100 × 4.6 mm, 5 μm (P/N 088558)
Mobile Phase	SCX Eluent A: Thermo Scientific™ CX-1 pH gradient buffer A (10-fold diluted) Eluent B: Thermo Scientific™ CX-1 pH gradient buffer B (10-fold diluted) HIC Eluent A: 1.5 M ammonium sulfate in 50 mM sodium phosphate buffer, pH 7.0 Eluent B: 50 mM sodium phosphate buffer/PA (80/20, v/v)
Column Temperature	40°C for SCX method, 30°C for HIC method
Autosampler Temperature	4°C
Flow Rate	0.3 mL/min for SCX, 0.5 mL/min for HIC
Autosampler Wash	10% MeOH in water
Injection Volume	10 μL for method optimization; 20 μL for HRAM MS characterization
Detector	280 nm, 0.2 Hz

2D SEC chromatographic and MS conditions:

Column	Thermo Scientific™ MabPac™ SEC-1, 300 × 4 mm, 5 μm (P/N 074696)	Parameters	Value/ settings	Parameters	Value/ settings
Mobile Phase	100 mM ammonium acetate	Spray Voltage: Positive Ion (V)	3800	Microscans	10
Flow Rate	0.12 mL/min	Sheath Gas (Arb)	35	Normalized AGC Target	300%
Column Temperature	30°C	Aux Gas (Arb)	10	RF Lens (%)	200
Autosampler Temperature	4°C	Ion Transfer Tube Temp (°C)	250	Maximum Injection Time Mode	custom
Autosampler Wash	10% MeOH in water	Vaporizer Temp (°C)	175	Data Type	Profile
Injection Volume	1D fraction volume	Orbitrap Resolution	60,000	Source Fragmentation	120 eV
Detector	280 nm, 0.2 Hz				

Gradient and column switching valve position:

Fraction collection from SCX and HIC						
1D Gradient	SCX method:			HIC method:		
	Time(min)	A%	B%	Time(min)	A%	B%
	0	100	0	0	95	5
	2	100	0	2	95	5
	2.1	95	5	25	0	100
	32	50	50	30	0	100
	33	0	100	31	95	5
	35	0	100	40	95	5
	36	100	0			
	41	100	0			
2D Gradient	Isocratic elution, 41 minutes					
Column switching valve position	CutValve: 1-2; DivertValve: 1-6; ArrayValves: Triggered by the valve switching time set in the '1D'					
SEC separation and HRAM MS data acquisition						
1D Gradient	Isocratic elution with 100% A, 30 minutes					
1D flow rate	0.10 mL/min					
2D Gradient	Isocratic elution, 30 minutes					
2D flow rate	0.12 mL/min					
Column switching valve position	CutValve: 1-6; DivertValve: 0-10 min (1-6); 10-25.5 min (1-2); 25.5-30 min (1-6) ArrayValves: Corresponding to the fractionation loop that was eluted into the 2D					

Data Analysis

Thermo Scientific™ Chromeleon™ Software Chromatography Data System (CDS) 5.3.2 was used for instrument control, data acquisition, and UHPLC data analysis. Thermo Scientific™ Biopharma Finder™ Software version 5.3 was used for HRAM MS data deconvolution and analysis.

Results

Multiple heart-cut 2D-LC-MS system setup

Five 250 μL fractionation loops were used in the 2D-LC system to offer flexible 1D fraction volumes, as shown in Figure 2. To minimize the fraction dispersion in the fractionation loops, the fluidic connections were set up to enable reverse flow directions in the loops for fraction transfer to the 2D (backflush, shown in Figure 2). Three 1D SCX peaks (main, acidic, and basic peaks) were utilized to demonstrate the benefit of the backflush fraction transfer; the results are listed in Figure 3.

Figure 2. Fluidic scheme of multiple heart-cut 2D-LC-HRAM MS setup.






Figure 3. The 1D SCX and 2D SEC UV chromatograms for three fractions, with a fraction volume of 105 μL from 1D, showed significant improvements in peak shape and sensitivity when fractions were transferred using the backflush system compared to forward flush.



SCX-SEC-HRAM MS analysis for polatuzumab vedotin


The SEC was used in the 2D method to keep the ADCs under the native state. In the 1D SCX UV chromatogram, 18 peaks were separated and transferred to the 2D system for intact mass analysis, as shown in Figure 4.

Figure 4. The 1D SCX-UV chromatogram and the relative abundance of each peak. The relative abundance of these charge variants was calculated using their 1D peak areas. The fraction volume was determined based on the peak width, ranging from 48 μL (0.16 min width) to 150 μL (0.5 min width).



After acquiring the MS data using Chromeleon CDS, the raw data were imported into the BioPharma Finder software for deconvolution and identification. Figure 5 shows the deconvolution results of peaks 11 and 12. The average mass of the main components detected in peaks 11 and 12 showed a mass shift of 128 Da, indicative of a lysine (Lys) truncation. This Lys truncation modification was also identified for the subsequent peaks 13, 15, 16, 17, and 18. And a D_n + 1 Lys and D_{n+2} + 0 Lys pattern was found in these peaks, indicating a drug payload and Lys truncation induced charge heterogeneity.


Figure 5. The full mass spectra (top) and deconvoluted spectra (bottom) of peak 11 and peak 12. The expanded view on the right of each spectra displays the details of the blue shaded area. The major payloads/glycoforms/PTMs were labelled in the deconvoluted spectra.



Besides the payload conjugation and Lys truncation, other modifications such as deamidation and oxidation can also affect the charge heterogeneity. The summary of the identified principal species for each 1D SCX peak is listed in Figure 5. The DAR₀, DAR₂, DAR₄, and DAR₈ species were eluted sequentially, which


suggests that the payload conjugation alters the surface charge of the mAbs, thereby increasing their retention on the SCX column. Additionally, deamidation and oxidation resulted in the formation of acidic variants, whereas Lys truncation leads to the formation of basic variants.

Figure 6: The summary of the main species for all charge variant peaks.



Owing to the high sensitivity of the Orbitrap Exploris 240 MS, it was feasible to identify charge variants even at very low abundance levels. As illustrated in Figure 7, peak 18, which has an abundance of merely 0.43% as detected by the 1D UV, still exhibits a high signal-to-noise (S/N) ratio in the raw mass spectra.


Figure 7. The 1D SCX-UV chromatogram, full mass spectrum, and deconvoluted spectrum for peak 18. The identification results were labelled on the deconvoluted spectrum.



HIC-SEC-HRAM MS analysis for polatuzumab vedotin


In the HIC-SEC-LC HRAM MS study, a total of 12 peaks were separated in the 1D HIC chromatogram. The five main peaks separated on the HIC chromatogram were putatively assigned with DAR values from 0 to 8 based on retention time.

Figure 8: The 1D HIC UV chromatogram showed a total of 12 separated peaks. The highlighted five main peaks were collected in the first cycle and transferred to SEC-HRAM MS for native MS characterization. The average DAR of polatuzumab vedotin was calculated to be 3.48 using the peak area in HIC chromatography, which closely matches its theoretical value (3.4 to 3.5)².




The deconvoluted mass analysis for these five main peaks revealed major species with DAR₀ (147891.86 Da), DAR₂ (150527.46 Da), DAR₄ (153162.74 Da), DAR₆ (155798.10 Da), and DAR₈ (158433.65 Da), respectively. Different glycoforms of the polatuzumab vedotin were also identified in this study. The identification results in Figure 9 reveal the glycosylation pattern for peak 4.

Figure 9. The full mass spectra (left) and deconvoluted spectra (right) of peak 4. The expanded view below displays the details of the blue shaded area. The predominant glycoform for DAR₂ was A2G0F+A2G0F, with a deconvoluted mass of 150527.46 Da, which is within 1 ppm deviation of the theoretical mass, demonstrating the high mass accuracy of the Orbitrap Exploris 240 HRAM MS.




The high sensitivity of the Orbitrap Exploris 240 HRAM MS effectively facilitates the identification of all the minor peaks in the 1D HIC, as shown in Figure 10.

Figure 10. The 1D HIC chromatogram (top), raw MS spectra (middle left), and deconvoluted mass spectra (middle right) of peak 5, which has a relative area of only 0.91% in the 1D HIC chromatogram, were successfully identified using the high-quality mass data. The results indicate that this peak is a DAR₂ isomer.




In addition to the expected positional isomers, ADCs with odd-number payloads such as DAR₁, DAR₃, and DAR₅ were also identified. The main species of each HIC peak in 1D were summarized in Figure 11.

Figure 11. The summary of the main species for all the 1D HIC peaks, the relative area was calculated based on the 1D UV peak area.



Notably, ADCs with linkers only, such as DAR₁+1 linker (peak 3) and DAR₃+1 linker (peak 7), were identified in our study. This indicates that the ADCs either lost one payload during production or storage, or that only one payload was linked to the free thiol after a disulfide bond was reduced. These are not the desired products in the ADC drug production process and can therefore be considered as impurities.

Figure 12: The identification results of peak 7, which was identified as DAR₃+1 linker.



Conclusions

In this study, an online multiple heart-cut 2D-LC-HRAM MS system was configured and optimized for the multi-attribute analysis of polatuzumab vedotin. This approach provides detailed and comprehensive insights into the charge heterogeneity analysis, DAR, and DLD analysis of polatuzumab vedotin under native conditions.

- The online multiple heart-cut 2D-LC system enables the direct transfer of SCX and HIC fractions to mass spectrometry analysis by utilizing SEC as a desalting step, thereby overcoming the MS-incompatibility buffer in the LC separation.
- The backflush transfer of fractions was employed to enhance peak shape and sensitivity in the 2D SEC separation.
- The 2D SEC-MS analysis enables the characterization of ADCs under native conditions.
- The high sensitivity of the Orbitrap Exploris 240 MS allows for the identification of fractions with low abundance in the 1D chromatography.

References:

- Bobály, B., Fleury-Souverain, S., Beck, A., Veuthey, J. L., Guillaume, D., & Fekete, S. Current possibilities of liquid chromatography for the characterization of antibody-drug conjugates. *Journal of pharmaceutical and biomedical analysis*, 147(2018): 493-505.
- Deeks, Emma D. Polatuzumab vedotin: first global approval. *Drugs*, 79.13 (2019): 1467-1475.

Trademarks/licensing

© 2025 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.