

Poster Reprint

**ASMS 2026**  
**WP 504**

# Lipid Insight Unblocked: Combining Nontargeted LC/MS Chemometrics With Automated Lipid Annotation

Mark Sartain<sup>1</sup>, Emma E. Rennie<sup>1</sup>, Sarah M. Stow<sup>1</sup>, Julie Horner-Buxton<sup>1</sup>, and Jeremy Koelmel<sup>2,3</sup>

<sup>1</sup>Agilent Technologies, Inc., Santa Clara, CA, USA

<sup>2</sup>Yale University, New Haven, CT

<sup>3</sup>Innovative Omics, Sarasota, FL

## Introduction

As adoption of non-targeted lipidomics expands, software tools must provide high confidence, high coverage annotations in a user-friendly framework.

MassHunter Explorer has a workflow-driven interface that enables MS and MS/MS feature extraction, statistical analysis, and MS/MS-based identification tools, while LipidMatch offers complementary deep lipid annotation and the ability for the user to review these annotations for confident lipid assignment.

We demonstrate a workflow that exchanges information between these software packages conveniently and show in-depth coverage with >800 lipid annotations from a porcine brain lipid sample.

These results highlight the importance of chemometrics, streamlined curation, and integrated visualization tools to support robust lipidomics interpretation.

## Experimental

### Sample Preparation

A simple experiment was conducted as a proof of principle for this workflow. A commercially available porcine brain lipid extract was spiked with Avanti Splash II Lipidomix standard. Three replicates were heated at 99 °C for 2 hours, and three replicates were left untreated (control) (Fig 1). The samples were dried under N<sub>2</sub> and resuspended in 1:1:0.1 butanol:methanol:10mM ammonium acetate.




Figure 1. Sample study design

## Experimental

### Method

A previously optimized 20-minute reverse-phase chromatography method<sup>1</sup> was applied with two modifications: an Agilent Altura ZORBAX Eclipse Plus C18 column was used and 10mM ammonium formate was replaced with 10mM ammonium acetate. Lipid extracts were analyzed in + and - ion modes with MS1-only and Auto MS/MS (Iterative) methods (Fig 2).



Figure 2. Agilent's standardized Omics LC configuration used a [1290 Infinity III Bio LC](#) coupled to a [Revident LC/Q-TOF](#) that provides robust and reliable lipidomics analysis.

### Software Workflow

Data analysis was performed with a combination of MassHunter Explorer 2.0, LipidMatch Suite 5.73 (Innovative Omics), and Mass Profiler Professional (MPP) 15.1 software. File format exchange is depicted in Fig 3.





Figure 3. Lipid annotation workflow with LipidMatch

## Results and Discussion



### MassHunter Explorer-Compound (Feature) Extraction and Chemometrics

MS-only and iterative MS/MS data files were acquired from control and heat-treated porcine brain lipid extracts and were processed in Explorer, where feature extraction, MS/MS extraction, normalization, filtering, and statistical analysis were performed (Fig 4).

All the extracted compound information and associated MS/MS in a project can be exported as a single .pfa file for easy import into other applications.

Figure 4. Explorer UI displaying results from the Statistics procedure. A volcano plot was used to find compounds that change in abundance with heat. The "Export to PFA" option is highlighted.

### LipidMatch Software for Automated Lipid Annotation

Explorer .pfa files from + and - ion mode projects were directly imported into the LipidMatch Flow interface, a software enabling deep annotation based on >500,000 lipids with fragmentation in libraries and a <5% false-positive annotation rate.<sup>2</sup> 827 and 256 high-scoring lipids (score=1+3) were annotated from + and - ion mode data, respectively.

Resulting lipid annotations were then reviewed with LipidMatch Visualizer (Fig 5). The lipid class phosphatidylserine (PS), known to be reduced with heat treatment, was inspected further with a Kendrick mass defect plot.

Reviewing homologous series patterns revealed further lipids as confirmed by carbon:double bond content even without MS/MS. In this manner, annotation evidence could be reviewed, filtered, and curated.




Figure 5. LipidMatch Visualizer zoomed view of the Kendrick Mass Defect Plot. The neutral loss (NL) corresponding to the PS headgroup was used to filter features indicating PS lipids.

## Results and Discussion

### Compound Review in Explorer and LipidMatch

LipidMatch annotations were mapped to the MassHunter Explorer project to enable easier review of compound behavior in statistical tests and visualization of abundance changes across sample groups (Fig 6).

A



B




Figure 6. Example of compound review and inspection of annotation results. (A) In Explorer, the selected compound group 301 annotated as lipid Coenzyme Q10 was found to be markedly reduced with heat treatment. To enable this annotation mapping a curated list was exported from the LM Visualizer tool and a compound database (.cdb) was created. This database was then leveraged in Explorer to support compound matching based on accurate mass with a tight retention time tolerance, resulting in 651 compounds in the positive-ion project with lipid annotations. (B) Inspection of LipidMatch results showed that Coenzyme Q10 was identified via spectral library matching, in contrast to most other lipids that were annotated via in silico library matching. LipidMatch supports the identification of lipids with both in silico approaches and classical spectral library matching..

<https://www.agilent.com/en/promotions/asms>

This information is subject to change without notice.

DE-014862

### Advanced Visualizations with MPP

MPP software can import Explorer results (.pfa) to enable advanced chemometrics and lipid visualizations (Fig 7.)




Figure 7. MPP lipid matrix plot of triacylglycerol (TG) compound abundances across control and heat-treated samples. Most TG lipids were reduced with heat. TGs with polyunsaturated fatty acids (PUFAs) were even more prone to degradation. The box plots highlight two TG lipids with different degrees of saturation.

## Conclusions

### Lipidomics workflow expands annotation through chemometrics, user curation, and visualization tools.

- Key enablement is .pfa file import to LipidMatch
- LipidMatch can be used to rapidly annotate lipids in an automatic fashion and confirm annotations and determine unknowns using an interactive visualizer
- Offers an alternative annotation approach, either independently or in conjunction with Agilent Lipid Annotator software.

## References

- <sup>1</sup> Hyunh, K et al. LC/MS dMRM Method Refinement Expands Targeted Lipidomics Studies from Plasma to Cells and Tissues. Agilent Application Note 5994-8365EN, 2026.
- <sup>2</sup> Koelmel, JP et al. LipidMatch: an automated workflow for rule-based lipid identification using untargeted high-resolution tandem mass spectrometry data. BMC Bioinformatics. 2017 Jul 10;18(1):331. doi: 10.1186/s12859-017-1744-3.