

Poster Reprint

ASMS 2026
WP 435

Unraveling Differential Lipids in Aging Tissues Using a Novel Data Analysis Workflow

Fernando “Ralph” Tobias¹, Mark Sartain¹, Karen Yannell¹, Brenna Keegan¹, Almudena Veiga-Lopez²

¹Agilent Technologies, Inc.

²University of Illinois Chicago

Introduction

Aging alters lipid metabolism across multiple tissues, driving metabolic inflexibility. While lipid remodeling in adipose and liver tissues is recognized, high-resolution approaches are required to characterize these dynamic shifts.

Here, nontargeted LC/Q-TOF lipidomics was performed on a multi-tissue sample set (WAT, BAT, liver) across three age cohorts (72 samples) using the 1290 Infinity III Bio LC coupled to the Revident LC/Q-TOF system.

To streamline discovery, we employed a novel workflow coupling MassHunter Explorer 2.0 (MH Explorer) for data extraction and chemometrics and MS-DIAL for complementary lipid annotation. This approach enabled the identification of 743 lipid species, using visualization tools in MH Explorer to quickly identify age-dependent changes.

These findings provide a comprehensive view of metabolic adaptation and underscore how integrated software solutions transform complex lipidomic datasets into actionable biological insights.



Figure 1. Agilent's omics standardized LC configuration used a [1290 Infinity III Bio LC](#) (left) which is iron free to give best peak shape of metal-sensitive analytes. Coupled to it is a [Revident LC/Q-TOF](#) (right) which provides robust and reliable lipidomics analysis.

Experimental

Sample Preparation

Age-matched tissue samples (Figure 2) were collected from each animal (10 – 15 mg) and was subjected to lipid extraction using a modified methanol/MTBE/water protocol. EquiSPLASH internal standard mix in methanol was added prior to tissue pulverization using a BeadBug homogenizer. Samples were resuspended in 9/1 methanol/chloroform for LC-MS/MS analysis.

Experimental




Figure 2. Study design for the aging lipidomics study.

Instrumentation

Lipid separation was performed using a 16-minute reversed-phase method on a ZORBAX Eclipse Plus C18 column with previously established chromatography¹. Samples were analyzed on the Agilent Revident LC/Q-TOF, tuned in m/z 1700 stable (+/-) ion modes.

A pooled QC sample was used for both system suitability and for iterative MS/MS acquisition ($n = 6$) per polarity for in-depth lipid identification coverage. MS1 acquisition was used for individual samples. The high resolution, isotopic fidelity and dynamic range of Revident enabled confident identification of 743 lipid species.

Data Analysis

Iterative MS/MS data files were searched in Lipid Annotator and in MS-DIAL 5.6². Prior to exporting as an SDF format, data was manually curated within the MS-DIAL software.

The Instrument type in the Property Settings of MS-DIAL was set to "LC-ESI-QTOF" for the MS/MS data to be annotated using this instrument type. Once the SDF file has been created, it was subsequently imported to ChemVista Library Manager. ChemVista can further refine the database or be exported as a CDB file.

Positive and Negative ion data (Iterative MS/MS and MS1 data) were imported to MassHunter Explorer 2.0 as two separate projects, in which peak picking, statistical analysis and identification using databases created from MS-DIAL data and/or Lipid Annotator.

Additional databases can be used for non-targeted lipidomics such as the Plasma Lipid PCDL (5994-7627EN).

Seamless Incorporation of MS-DIAL, Lipid Annotator, ChemVista and MassHunter Explorer




Figure 3. Lipidomics acquisition and data analysis workflow. MS1 data is collected for each individual sample and Iterative MS/MS is collected on pooled sample(s) to build a custom database to the study set. Iterative MS/MS data can be analyzed by Lipid Annotator or MS-DIAL and can be directly imported to Agilent ChemVista Library Manager. ChemVista can then export it as a .CDB file to be used with MassHunter Explorer 2.0.

Integrating MS-DIAL MS/MS Results into MassHunter Explorer 2.0

- MS-DIAL Setup and Export:** Set the instrument type to "LC-ESI-QTOF" in Project Properties and export curated MS/MS spectra data as an SDF file (Figure 4).
- Import to ChemVista:** Load the SDF file into ChemVista Library Manager. This also automatically transfers curated lipid data which includes retention times and other chemometric data (precursor mass, formula, InChiKey, SMILES).
- Export or Refine:** The database is ready for immediate export to MH Explorer, with the option for additional curation within ChemVista if required.
- Integrate with Explorer:** In the Identify Tab of MH Explorer, add the MS-DIAL database to finalize the workflow. Additional databases from Lipid Annotator and the Plasma Lipid PCDL can also be added.




Figure 4. Exporting from MS-DIAL. After curating the results in the software, export the alignment results by checking the highlighted boxes in the Alignment Result Window. This will generate an SDF file that contains MS/MS data of the filtered results, complete with retention time, precursor mass, ion type, chemical formula and structure.

Multivariate Lipidomics Analysis with MassHunter Explorer 2.0




Figure 5. Visualize data as you filter with publication-quality figures. A: 3D Principal Component Analysis (PCA) plot of multi-tissue sample set of liver, BAT, and WAT. B: Hierarchical clustering heatmap of 645 lipids identified by MS/MS and retention time in positive ion mode.

- Following sample normalization by tissue weight, a PCA plot verified tissue-specific clustering (Figure 5A). Furthermore, after filtering for features identified by MS/MS and RT, a heatmap was produced which validates the tissue specific separation and also lipid class-specific alterations (Figure 5B).
- Volcano plots (Figure 6) enable rapid identification of statistically significant features across experimental cohorts.
- Integration with MS-DIAL results supports high-confidence identification of additional lipid species. As shown, a TG-estolide was found elevated in the OLD liver cohort, supported by a robust MS/MS spectral match (Figure 7).

Conclusions

- Complementary and Flexible:** MassHunter Explorer 2.0 easily integrates with current MS-DIAL analysis pipelines to enable an end-to-end workflow for non-targeted lipidomics.
- Comprehensive:** ChemVista manages custom MS/MS libraries from MS-DIAL, Lipid Annotator and other databases sources to use for Identification in MH Explorer 2.0.
- Confident:** Statistically significant results can be obtained quickly and efficiently using powerful feature finding algorithms and statistical tools. Publication-quality visualizations can be created in the same workflow.




Figure 6. Volcano plot of lipid alteration liver in old and young mice and highlighting a TG-estolide, 2-fold difference and p-value cutoff of 0.05.

Compound Group	Mass(avg)	m/z	RT(avg)	Score	Hits	Formula	Name
342	952.8453	970.8793	11.747	87.75	1	C62 H11...	TG 59:4 TG 18:1_23:1_18:2
345	1139.0057	1157.0396	12.000	93.76	1	C73 H13...	TG 70:4 O2 TG 18:1_18:1_18:1;O(FA 16:0)
346	1125.9870	1126.9948					TG 68:5 O2 TG 16:0_18:1_18:2;O(FA 16:1)
347	704.5833	705.5904					SM 3:0;O2
348	882.7669	900.8015					TG 54:4 TG 18:1_18:1_18:2
351	818.7727	836.8067					TG 0-50:1 TG 0-16:0_16:0_18:1
355	858.7673	876.8014					TG 16:0_18:1_18:1
358	821.7479	822.7559					TG 16:0_16:0_16:1
363	796.6570	814.6924					TG 48:5 TG 12:0_18:2_18:3
369	1024.9385	1042.9723					TG 64:3 TG 28:0_18:1_18:2
373	1081.0016	1099.0356					TG 68:3 TG 18:1_18:1_32:1
374	778.7060	1575.4446					TG 46:0 TG 14:0_16:0_16:0
378	765.5310	766.5380					PE 18:1_20:4
379	984.9079	1002.9418					TG 61:2 TG 25:0_18:1_18:1
382	666.5795	684.6139					TG 38:0 TG 10:0_12:0_16:0
383	805.5623	806.5693					PC 18:2_20:4
385	723.5202	724.5276					PE P-36:4 PE P-16:0_20:4
386	1052.9700	1071.0036					TG 66:3 TG 30:0_18:1_18:2
391	743.5829	744.5901					PC O-34:2

Figure 7. Highlighting TG 18:1_18:1_18:1;O(FA 16:0) as a high confident ID in MH Explorer and the corresponding annotation in MS-DIAL.

References

- Huynh et al., A Comprehensive, Curated, High-Throughput Method for the Detailed Analysis of the Plasma Lipidome. Agilent Application Note 5994-3447EN, 2021.
- Takeda et al. Nature Communications, 2024