

Poster Reprint

ASMS 2026
MP 536

Comprehensive characterization of oligonucleotides and related impurities using advanced LC/Q-TOF analytical workflows

Saikat Banerjee, Preeti Bharatiya, Ashish Pargaonkar, Vivek Dhyani;
Agilent technologies, India

Introduction

Therapeutic oligonucleotides are synthetic, short nucleic acid sequences designed to modulate gene expression by targeting RNA or DNA with high sequence specificity. However, their inherent structural heterogeneity makes impurity profiling essential, as during GMP manufacturing, impurities like phosphorothioate to phosphodiester, 5' and 3' truncations or extensions, depurination, and abasic species can directly impact both product safety and efficacy. Antisense oligonucleotides (ASO) like Nusinersen (7127.171 g/mol) are indicated for the treatment of spinal muscular atrophy (SMA). Their closely related variants coelute with the main peak, complicating their characterization in conventional one-dimensional LC/MS workflows. To address this challenge, a 2D LC/MS strategy was implemented using the Agilent 1290 Infinity III Bio LC system with multiple heart cutting, coupled to the 6545XT AdvanceBio LC/Q-TOF and MassHunter BioConfirm software. The first-dimension separation uses an ion pairing buffer-based method to achieve higher chromatographic peak resolution, while the second dimension provides MS compatible desalting, enabling enhanced mass spectrometric sensitivity and improved detection of low-level impurities.

Instrumentation

1290 Infinity III Bio 2D-LC & 6545XT Q-TOF includes:

- 1290 Infinity III Bio Flexible Pump (G7131A)
- 1290 Infinity III Bio High-Speed Pump (G7132A)
- 1290 Infinity III Bio Multisampler (G7137A)
- 1290 Infinity III Multicolumn Thermostat (G7116B)
- 1290 Infinity III Variable Wavelength Detector (G7114B)
- 1290 Infinity III Diode Array Detector (G7117B)
- 2D-LC System with MHC set up



Figure 2. Agilent 1290 Infinity III Bio 2D-LC coupled to an Agilent 6545XT AdvanceBio LC/Q-TOF.

Experimental


Table 1: 1D and 2D Method Parameters

Parameters	Details									
First-dimension										
Mobile Phase	Buffer: 0.2% Dibutyl amine and 0.2% HFIPA into 1L water; Mobile phase A: Buffer: methanol (80:20; v/v) Mobile phase B: Methanol: buffer (80:20; v/v)									
Flowrate	0.4 ml/min									
Gradient	Time (min)	0	2	5	10	20	27	31	32	35
	B (%)	05	05	35	70	70	95	95	05	05
Injection Vol	10 µL									
Column	Altura Oligo HPH-C18, 4.6 × 150 mm, 2.7 µm									
Column Temp	30 °C									
Detector	UV 260 nm									
Second- dimension										
Mobile Phase	Buffer: 15 mM Ammonium acetate in water Mobile phase A: Acetonitrile: buffer (70:30; v/v) Mobile phase B: Buffer: Acetonitrile (70:30; v/v)									
Flowrate	0.4 ml/min									
Gradient	Time (min)	0	0.5	1.2	2.4	3.0	3.5	4.0		
	B (%)	15	15	90	90	15	15	15		
Column	Altura Poroshell HILIC-Z, 2.1*100mm, 2.7 µm									
Column Temp	30 °C									
Flowrate	0.3 ml/min									
2D Cycle Time	46mins									
ASM Setting	Factor: 3; Flush factor: 5									

Table 2: MS Parameters

System	6545XT AdvanceBio LC/Q-TOF
Gas Temp	350 °C
Drying Gas	13 L/min
Nebulizer	35 psi
Sheath Gas Temp	250°C
Sheath Gas Flow	12 L/min
VCap	4500 V
Nozzle Voltage	2000 V
Fragmentor	180 V
Mass Range	<i>m/z</i> 400 to 3,200
Auto MS/MS Range	<i>m/z</i> 50 to 3,200
Polarity	Negative
Sort Precursors	By abundance only

Optimized 2D-LC/MS Workflow



Results and Discussion

Nusinersen Characterization Using BioConfirm

Step 1: Acquisition Setup 6545XT LC/Q-TOF

- MS & MS/MS data collected
- Wide mass range & optimized Collision energy used

Step 2: Define Sequence and Modifications

- Input nusinersen full length product (FLP) sequence into BioConfirm with backbone as phosphorothioate (PS)
- Define 2' O methoxyethyl (2' MOE) ribose modifications to enable accurate identification of targeted main peak and impurities

Step 3: BioConfirm Target Plus Impurities Workflow

- Identify the FLP at MS level alongside isoforms, truncated variants

Step 4: BioConfirm Sequence Confirmation Workflow

- Tolerance set to MS: ± 10 ppm & MS/MS: ± 20 ppm
- Identify and confirm presence of FLP, truncated, variants, isoforms, degradants etc.

NusiN
Total monoisotopic mass: 7122.2762
Total average mass: 7127.2001
Sequence molecular formula: C234H340N61O128P17S17

Linker assumption
☒ If two building blocks are next to each other
S

Open definition

1	5'	/MOErU/	*/MOErC/	*/MOErA/	*/MOErC/	*/MOErU/	*/MOErU/	12
13		/MOErU/	*/MOErC/	*/MOErA/	*/MOErU/	*/MOErA/	*/MOErA/	24
25		/MOErU/	*/MOErG/	*/MOErC/	*/MOErU/	*/MOErG/	*/MOErG/	35

Definitions

Building blocks

Name	Code	Formula
2-Methoxy...	/MOErA/	C13H19N5O5
2-Methoxy...	/MOErC/	C13H21N3O6
2-Methoxy...	/MOErG/	C13H19N5O6
2-Methoxy...	/MOErT/	C13H20N2O7
2-Methoxy...	/MOErU/	C13H20N2O7
Deoxyaden...	A	C10H13N5O3
Deoxycytdi...	C	C9H13N3O4
Deoxyaden...	dA	C10H13N5O3
Deoxycytdi...	dC	C9H13N3O4
Deoxyguan...	dG	C10H13N5O4
Deoxythym...	dT	C10H14N2O5
2-Fluoroad...	fA	C10H12FN5...
2-Fluorocyti...	fC	C9H12FN3O4
2-Fluorogu...	fG	C10H12FN5...
2-Fluorouri...	fU	C9H11FN2O5

Linkers

Name	Code	Formula
Phosphorot...	*	HPOS
5' Phosphate	/5Phos/	H2PO3
Phosphodie...	p	HPO2
Phosphorot...	s	HPOS

Chemical data dictionary editor

Figure 3. Nusinersen backbone nucleotide sequence and chemical modifications 2' O methoxyethyl (2' MOE) ribose in sequence using Chemical data dictionary set up in BioConfirm.


2D-LC Separation and HiRes Sampling Strategy

- Orthogonal 2D-LC using Altura oligo HPH (1D) and Altura HILIC-Z (2D) enabled high-resolution separation and MS-compatible desalting.
- HiRes sampling with multi-inject allowed combined transfer of heart-cuts, improving throughput and efficiency and identification of trace level peaks.

Two targeted experiments:

- Experiment 1:
 - Focused on main peak at 16.0 min with the coeluting impurities peaks.
 - Targeted HiRes cuts captured main peak impurities
 - Enabled detailed profiling of FLP and coeluting impurities
- Experiment 2:
 - Focused on all the other peaks before 16.0 min
 - Multi-inject features transferred 10 heart-cuts per run, ensuring high-resolution separation and full MS data acquisition across all peaks.

BIO



- Multi-inject
Use multiple 40 μ l loops for emulating a large loop without hardware modifications.
- Safe time by analyzing a High-Resolution series at once

Figure 5. Multi-Inject with Bio MHC Valve & Biocompatible Loops




Figure 6. Overlay of 1D chromatograms with HiRes cut markers: Experiment 1 (Coeluting peak at 16 min) & Experiment 2 (Peaks prior to 16 min).

