



Background Reduction Protocol and Instrument Tuning

For detailed instructions please use the Simplified Background Reduction Protocol for Agilent Triple Quadrupole LC/MS Technical Overview.



Materials Needed


LC/MS grade water	Squirt bottles (2)	4000 grit paper
LC/MS grade methanol	Lint-free lab tissue	Sonicator
LC/MS grade isopropanol	Lint-free cloth	

Clean the source

- Disconnect the exhaust polyethylene tube
- Open the source and rinse with H₂O:Methanol and Isopropanol; Wipe with lint free tissue
- Remove and clean the spray shield and capillary cap with isopropanol; Use 4000 grit paper if needed for discoloration; Sonicate in isopropanol and rinse with methanol
- Inspect and clean the nebulizer if needed
- Reinstall the nebulizer, cap, and spray shield and reconnect the exhaust tube.

Use MassHunter Data Acquisition for tuning

- Checktune – Daily/Weekly
- Autotune – Monthly/Quarterly/As Needed



Agilent 6400 Series
LC/MS/MS Concepts
Guide (pages 67–70)



Data
Acquisition

Build a Background checkout method using MassHunter Data Acquisition

Configuration - LC		Configuration - MS	
Injection Volume	No injection	Ionization Mode	Positive
Analytical Column	Zero dead volume union	Scan Type	MS2 Scan <i>m/z</i> 40 – 1000
Column Temperature	Not controlled	Drying Gas Temperature	250 – 350 °C
Mobile Phase A	Water	Drying Gas Flow	11 L/min
Mobile Phase B	MeOH or AcN	Nebulizer Pressure	30 – 50 psi
Flow Rate	0.35-0.5 mL/min	Sheath Gas Temperature	300 – 400 °C
Isocratic	50:50	Sheath Gas Flow	11 L/min
Stop Time	1 minute	Nozzle Voltage	0 V
		Capillary Voltage	4000 V
		Delta EMV	0 V


Review Results

After acquisition of the MS2 Scan, open the datafile using MassHunter Qualitative Analysis software. Review the total ion chromatogram (TIC) and observe the Y-axis for total number of counts compared to the table below.

If you want to improve your instrument background counts, please use the LC/MS cleaning protocol from the Technical Overview linked above, before beginning method development.

TIC background reference ranges for Agilent LC/TQ instruments.

LC/TQ	Extra Clean	Clean	Working Range	Dirty
Ultivo	< 1M	1M – 2M	2M – 6M	> 6M
6470	< 1M	1M – 2M	2M – 6M	> 6M
6495	< 20M	20 – 30M	30 – 100M	> 100M



For up-to-date information about Agilent's clinical research solutions and technologies, visit our webpage.

www.agilent.com
DE22541104
This information is subject to change without notice.

© Agilent Technologies, Inc. 2022
Printed in the USA, May 13, 2022
5994-4797EN

Compound Optimization

Step 1

- Prepare your samples
- Set up LC/MS method

Materials Needed

Compound reference materials	LC/MS grade organic	Agilent LC Column
Internal Standard reference materials	LC/MS grade water	Zero dead volume union

Prepare standards


- Maximum number of compounds per mix – 10
- Compound concentrations between 1000 – 5000 ng/mL
- Ensure isobaric compounds are separated into different vials
- Prepare compound mixes in 50:50 Water:Organic

Build a short isocratic MS2 Scan "Optimizer Method"

- 50:50 Water:Organic
- Flow rate: 0.25 mL/min
- 1 – 10 μ L injection

Step 2

- Create a project in Optimizer
- Follow the Tabs



Optimizer Setup Precursor Ion Selection Product Ion Selection Compound Setup

Optimizer Setup

- Sample Introduction: Injection
- Fragmentor: Coarse Range 80-200 V; Fine Step 5 V
- Collision Energy: Range 5-40 V (Repeat any compound with 40 V CE at a higher range)
- Cell Accelerator Voltage(CAV): 4 V
- Direct acquisition method to the one built in Step 1

Precursor Ion Selection

- Select positive and/or negative ions (+H, -H, etc.)
- Can exclude masses or set a minimum abundance (Optional)

Product Ion Selection


- Set max number product ions, Recommended number: 6
- Set low mass cutoff: *m/z* 40
- Can exclude masses or losses (Optional)

Compound Setup

- Name your compound and group if desired (If re-running a compound input a unique name)
- Enter the formula to calculate the nominal mass
- Specify sample position

Step 3

- Start Optimizer – An automatic 4 Step Process
- MS2 SIM > Product Ion > MRM > Narrow Product Ion




Optimizer Report is generated; Results can be imported to MassHunter Data Acquisition

Source Optimization

Step 1

- Open Source Optimizer Program
- Set up Project and Instrument Parameters



Project Parameters

- Optimize Method: Load the acquisition method to be optimized
- Project Folder: Choose where to save the results
- Project Name: Name your project

Worklist Parameters

- Name your samples and specify sample position

Instrument Parameters

- Input all the information from the table below. Ranges may vary based on instrument model.

Project	Type	Pre-Wait (min)	Replicate	Step-Wait (min)	Start Value	End Value	Step Size
1	Capillary	0	1	0	1000	6000	500
1	Nozzle Voltage	0	1	0	0	2000	500
1	Drying Gas Temp	10	2	10	100	290	25
1	Drying Gas Flow	5	1	5	10	20	1
2	Sheath Gas Temp	10	2	10	100	400	25
2	Sheath Gas Flow	5	1	5	5	12	1
3	Nebulizer	0	1	0	15	60	5
4	High Pressure RF	0	1	0	70	210	20
4	Low Pressure RF	0	1	0	40	160	20

Step 2

- Group Parameters as shown above into 4 Projects
- Review results and update the method as you go

Start the First Project (Capillary/Nozzle Voltage/Drying Gas Temp & Flow)

- Check the boxes of the parameters to optimize, grouped by "Project" in the table above
- Calculate the "injections" multiplied by the injection volume to ensure enough sample volume to complete the project (e.g. 20 μ L x 144 = 2.88 mL)
- Click "Create Methods" and "Submit".

Create Methods
Submit
Close

Project created: D:\MassHunter\Data\Jennifer Hitchcock\E 122 methods 144 injections Estimated time: 1274 minutes

- The Study Manager application is now open.
- A script may be added to the end of the run if needed, for example to shutdown if running over night. Once everything is set, click, "Start"

Review Study Results in Quantitative Analysis

- Open Batch – The instrument parameters are listed as .s files – Select and open the study folder and subsequent pre-created batch file based on your project name
- Locate the Final Conc. Column in the Batch table – Find the highest value in that column and enter corresponding source condition into the acquisition method you are optimizing. Update and save the method before moving onto the next project.
- With multiple compounds, settings compromise may be needed. It is recommended to skew the settings toward less responsive compounds.

Next Project

- Redirect Optimizer to the updated acquisition method so it uses the updated parameters
- Check the boxes for the next project group and repeat the steps above until you have completed all 4 projects and an optimized method is achieved.
- iFunnel Optimizer will be performed for 6490 and 6495 users only.