


Advancing Low flow LC/MS for single cell proteomics with variable flow and 50 cm microfabricated pillar array columns










HPLC 2025 Bruges

Jeff Op de Beeck

 The world leader in serving science




Why single cell proteomics



Application	Bulk result	Single cell data
Identification of cell subpopulations based on protein expression or metabolic profiles (tumors, tissues, immune cells, cell cultures) 		
Detection and analysis of rare cells (i.e., CTCs from liquid biopsies) 		
Analysis of limited sample material (exosomes, needle aspirates, biopsies) 		

Single-cell analysis challenges all the steps in the proteomics workflow

Single cell proteomics workflow



Injection of single cell samples with Vanquish Neo UHPLC



Back Flush trap and elute workflow


- Single cell samples vary in volume (between and within experiments) – sample volume range 1-10 μL
- Sample uptake + loading can be reduced to 2.4min
- Additional air drawn into sample loop will not affect ionization / analytical column performance

Vial bottom detection

- Samples can be fully aspirated by setting injection volume with some margin (1-2 μL added)
- Compatible with conventional sample vials, 96 and 384 well plates

Sensitivity versus throughput

Impact on instrument productivity and sensitivity (in LC-ESI-MS)




Low flow LC = reduced instrument productivity & robustness?

- Sample loading (μL volumes) has a much bigger impact at low flow rates
- Small dead volumes can have a huge impact on performance

Why does miniaturization improve sensitivity?

- Increase of analyte concentration with decrease of column ID (lower sample dilution)
- Improved ionization efficiency at lower flow rates (ESI)


Ultra Low Flow sensitivity with increased instrument productivity



Gradient profile					
Duration	Time	%B	Flow rate	Volume	Total volume
min	min	/	µL/min	µL	µL
0	0	2	0.75	0	0
0.1	0.1	4	0.75	0.075	0.075
2.8	2.9	15.5	0.75	2.1	2.175
0.1	3	15.6	0.1	0.01	2.185
13	16	22.5	0.1	1.3	3.485
15	31	40	0.1	1.5	4.985
0.5	31.5	99	0.1	0.05	5.035
13.5	45	99	0.1	1.35	6.385

83% instrument productivity

(including sample uptake, sample loading and column equilibration)



- “High” Flow rate (750 nL/min) is used in initial part of gradient to push analytes forward
- At the point where analytes start eluting, the flow rate is reduced to “low” flow rate (100 nL/min) to boost sensitivity
- Caution as **2.9 min gradient at 750 nL/min will give 21.75 min elution at 100 nL/min**


Single cell by dilution – Bulk HeLa



Vanquish Neo UHPLC

- Precise injection of sample volumes down to 10 nL
- Full dilution series using concentrated sample stock
- Ideal sample conditions (minimal handling & adsorption)

- “true” 30SPD Trap & Elute – 100 nL/min elution
- Bulk HeLa dilution (from 5 ng/μL sample)
- Orbitrap Exploris 240 MS without FAIMS interface
- 120k MS2 – 246ms IT – 80Th windows – 400-800 m/z
- Processing in Spectronaut 19 – triplicates together (MBR)



Consistent LC-MS method and workflow optimization...

How about real single cell samples?


Effect of HeLa cell size on proteome depth

Cellenion cellenONE

- Single HeLa cells isolated within narrow size bins (1 - 1.5 – 3 μm)
- One-pot single cell digestion into LoBind Eppendorf 384 well plate – proper silicone seal needed!



- “true” 30SPD Trap & Elute – 100 nL/min elution
- Thermo Scientific™ Orbitrap Exploris™ 240 without FAIMS
- 120k MS2 – 246ms IT – 60Th windows – 400-800 m/z
- Processing in Spectronaut 19 – cells within size bins (n=10) together (MBR)



- Good correlation between cell size and proteome coverage observed
- Selection of single cells in narrow size bins increases data consistency


Ultimate sensitivity with Orbitrap Astral analyzer



- 72SPD – 200 nL/min elution
- Single HeLa cells
- Thermo Scientific™

- Orbitrap Astral™ mass spectrometer with FAIMS interface
- Processing in Spectronaut 19


Protein Groups Single Cells



2-fold Increase in throughput
2-fold Increase in depth

Summary single cell proteomics workflow optimization

- **Careful Selection of workflow components:**
Ensure minimal sample loss and maximum sensitivity in all steps.
- **Vanquish Neo UHPLC:**
Ideal for method optimization and true single cell experiments.
 - Precise injection range: 10 nL to 100 μ L.
 - Trap and elute with full aspiration of precious samples.
- **Ultra Low Flow NanoLC:**
Implement without impacting instrument productivity using variable flow rate methods and micro Pillar Array Column technology.
- **Comprehensive Coverage:**
Achieve study of biological processes (2k protein groups/cell) on Orbitrap Exploris MS platforms at 30 samples/day.
- **Ultimate Sensitivity:**
Orbitrap Astral MS provides deeper insights with increased throughput.



Acknowledgements



- Simon Devos
- An Staes



- Kris Gevaert
- Marcel Bühler
- Emin Araftpoor

ThermoFisher SCIENTIFIC

- Natalie Van Landuyt
- Fernanda Salvato
- Bernard Delanghe
- Tabiwang Arrey
- Julia Kragenbring
- Runsheng Zheng
- Eugen Damoc
- Alec Valenta
- Paul Jacobs

And of course the Thermo Fisher
Scientific team in Ghent



- Christof Völkle

Thank you

