

Poster Reprint

ASMS 2022
Poster number ThP156

It's time for LC/MS to Replace HPLC for Routine Cannabinoid Testing

Sue D'Antonio¹; Brian Whiteley²; Nikolas Lau³; Anthony Macherone⁴

¹Agilent Technologies, Cedar Creek, TX; ²Agilent Technologies, Little Falls, DE; ³Agilent Technologies, Chicago, IL; ⁴Agilent Technologies, Santa Clara, CA

Introduction

The most common analytical technique for the quantitative determination of cannabinoids in cannabis plants and other cannabinoid products, is HPLC with UV detection. However, given the complexity of the matrix in these products, interferences and co-eluting chemicals can often bias the results when using nonselective UV detectors. We compared HPLC-UV with electrospray LC/MS to illustrate the power and need of mass spectrometry for the selective and quantitative determination of 16 cannabinoids including Δ^9 -THC, Δ^8 -THC, THCA, CBD, and CBDA in food-grade hemp seed oil.

Experimental

Methods

Three HPLC-UV and LC/MS methods were developed that demonstrated baseline resolution for the commercially available cannabinoids. Simultaneous SIM/SCAN data were acquired with the LC/MS method, and full spectral (3D) data were acquired with the HPLC-UV method. A standard curve was created in hemp seed oil over the concentration range of 0.05 $\mu\text{g}/\text{mL}$ through 50.00 $\mu\text{g}/\text{mL}$ for both methods. A Hemp standard consisting of 500 $\mu\text{g}/\text{mL}$ of CBD and 1.5 $\mu\text{g}/\text{mL}$ of Δ^9 -THC was diluted 10:1 and analyzed. Limits of detection (LOD), limits of quantitation (LOQ), linearity, selectivity, accuracy, and precision were determined for both methodologies. Ten samples of commercially available hemp products were purchased and analyzed six times each to determine selectivity and precision. The samples were diluted at 1:10 and 1:100 ratios with a 60/40 Ethanol and Water mixture.



Experimental

Method 1

Column: Agilent Poroshell EC-C18, 3.0 x 150mm, 2.7 μm
Flow Rate: 0.5 mL/min

Mobile Phase: A) 0.1% Formic Acid in H₂O

B) 100% ACN

C) 100% MeOH

D) 10 mM Ammonium Formate in H₂O

Gradient: Time	%A	%B	%C	%D
0.00	29	70	0	1
3.20	29	70	0	1
7.20	12	0	87	1
10.00	0	0	95	5

Column Temperature: 30 °C

Post Time: 5 minutes

UV: 228 nm

Method 2

Column: Agilent Poroshell EC-C18, 3.0 x 150mm, 2.7 μm

Guard Column: Poroshell 120 EC-C18, 3.0 x 5mm, 2.7 μm

Flow Rate: 0.8 mL/min

Mobile Phase: A) 0.1% 5mM Ammonium Formate in Water

B) 0.1% Formic Acid in Acetonitrile

Gradient: Time	%A	%B
0.00	25	75
10.00	10	90
11.00	10	90

Column Temperature: 30 °C

Stop Time: 11 minutes

Post Time: 3.5 minutes

Injection Volume: 5 μL

Autosampler Temperature: Ambient

Peak Width >0.0063 min, 40 Hz

UV: 230 nm

Method 3

Column: Agilent ZORBAX Bonus RP 3.0 x 150 mm, 1.8 μm

Flow Rate: 0.5 mL/min

Mobile Phase: A) Water

B) Methanol

C) 0.1 % Formic Acid + 2.2 mL of 5mM Ammonium Formate in Water

Gradient: Time	%A	%B	%C
0.00	23	72	5
12.50	0	95	5

Column Temperature: 50 °C

Stop Time: 12.5 minutes

Post Time: 6.5 minutes

Injection Volume: 0.25 μL


Autosampler Temperature: 25 °C

UV: 230 nm


Overall run time: 20.0 minutes (including re-equilibration)

UV Chromatograms


Method 1



Method 2




Method 3




UV Chromatograms compared to LC/MS TIC for Hemp Standard

Data with Method 2



Hemp standard CBD concentration is 500 µg/mL which is typical for many of the Pet CBD oils. The THC concentration is 1.5 µg/mL which is less than 0.3% (wt/wt) of the CBD oil note the 10¹ of the UV axis compared to 10⁶ of the Y axis of the SIM at 315.2



Hemp oil results with both UV and LC MS data

The following is the results of Delta 9 - THC concentrations in 7 Samples of commercial Pet Hemp oil run with all three methods.

Sample	Δ9-THC UV 1	Δ9-THC MS 1	Δ9-THC UV 2	Δ9-THC MS 2	Δ9-THC UV 3	Δ9THC MS 3
1	620 µg/ml	712 µg/ml	650 µg/ml	740 µg/ml	600 µg/ml	690 µg/ml
2	Not detected	5.5 µg/ml	Not detected	5.5 µg/ml	Not detected	5.5 µg/ml
3	Not detected	13.5 µg/ml	Not detected	15.4 µg/ml	Not detected	12 µg/ml
4	Not detected	4.1 µg/ml	Not detected	3.6 µg/ml	Not detected	2.8 µg/ml
5	Not detected	14.3 µg/ml	Not detected	15.1 µg/ml	Not detected	13.9 µg/ml
6	158 µg/ml	165 µg/ml	148 µg/ml	167 µg/ml	141 µg/ml	163 µg/ml
7	45 µg/ml	58 µg/ml	50 µg/ml	70 µg/ml	Not detected	66 µg/ml
Hemp oil	Not detected	Not detected	Not detected	Not detected	Not detected	Not detected
Hemp Std	1.4 µg/ml	1.5 µg/ml	1.6 µg/ml	1.5 µg/ml	1.5 µg/ml	1.5 µg/ml

1:10 dilutions calculated for this table. Average of 6 replicates

