

Poster Reprint

**ASMS 2022**  
**Poster number MP334**

# Simultaneous Identification, Confirmation and Quantitation of Impurity in Antibiotics Formulations Using Mixed Scan Mode of Ultivo LC/TQ

Vikrant Goel<sup>1</sup>, Saikat Banerjee<sup>1</sup>, and Julie Horner Buxton<sup>2</sup>

<sup>1</sup>Agilent Technologies India Pvt Ltd.

<sup>2</sup>Agilent Technologies Inc., Santa Clara, CA

## Introduction

It is a vital requirement of pharma- and agro-chemical as well as other related bulk manufacturers to identify and to quantify impurities<sup>(1)</sup>. The challenging task of assessing non-chromophores in various formulations can be achieved using a combination of Mixed Scan Mode on the Agilent Ultivo LC/TQ (Fig 1) in addition to UV based detection in a single injection.



Figure 1. Agilent Ultivo LC/TQ Mass Spectrometer.

Sulbactam is a  $\beta$ -lactamase inhibitor which is given in combination with  $\beta$ -lactam antibiotics to inhibit  $\beta$ -lactamase, an enzyme produced by bacteria that destroys the antibiotics<sup>(2)</sup>. Here, an LC/TQ-based method is used to achieve  $\leq 50$  times the allowed limit (0.01%) for Impurity A (CAS 23315-18-6) in Sulbactam formulations. Data obtained using both UV detection and Mixed Scan Mode on an Agilent Ultivo LC/TQ provide identification and confirmation in a single injection, allowing this method to be used for purity and impurity workflows in both R&D and QA/QC labs.

## Experimental

### Sample Preparation:

All formulations and standards were prepared using water as a diluent. The dilution protocol for Impurity A standards is seen in Table 1. A spike study was performed at a concentration of 450 ppb in a 5mg/ml solution of formulation for recovery experiments.

| S. No | Stock Conc | Stock vol | Diluent | Final Conc |
|-------|------------|-----------|---------|------------|
| 1     | 3 mg/ml    | 10 ul     | 990 ul  | 30ppm      |
| 2     | 30 ppm     | 100 ul    | 900 ul  | 3 ppm      |
| 3     | 3 ppm      | 300 ul    | 700 ul  | 900ppb     |
| 4     | 900ppb     | 500 ul    | 500ul   | 450 ppb    |
| 5     | 450 ppb    | 500 ul    | 500 ul  | 225 ppb    |
| 6     | 225 ppb    | 400 ul    | 600 ul  | 90 ppb     |
| 7     | 90 ppb     | 500 ul    | 500 ul  | 45 ppb     |
| 8     | 45 ppb     | 500 ul    | 500 ul  | 22.5 ppb   |
| 9     | 22.5 ppb   | 400 ul    | 600 ul  | 9 ppb      |

Table 1. Dilution protocol for Sulbactam Impurity A

### Instrumentation:

Chromatographic separations were performed on an Agilent 1260 Prime LC system equipped with an Agilent C18 column using a water/methanol gradient. Mass detection was done on an Ultivo LC/TQ running MassHunter software for LC-UV-MS acquisition and analysis. The mass spectrometer was operated in Mixed Scan mode<sup>(3)</sup>; scan parameters are in Figure 2.

| Acquisition Parameters |          |                       |                 |         |            |               |                       |                     |         |                      |          |                          |        |           |
|------------------------|----------|-----------------------|-----------------|---------|------------|---------------|-----------------------|---------------------|---------|----------------------|----------|--------------------------|--------|-----------|
| Scan type              | Polarity | Compound/Segment name | Precursor (m/z) | MS1 res | Mass (m/z) | Product (m/z) | MS2 range start (m/z) | MS2 range end (m/z) | MS2 res | Scan/Dwell time (ms) | Frag (V) | ISTD?                    | CE (V) | Threshold |
| Product Ion            | Positive |                       | 182.1           | Unit    |            |               | 20                    | 200                 |         | 100                  | 75       | <input type="checkbox"/> | 16     | 0         |
| Product Ion            | Negative |                       | 180.1           | Unit    |            |               | 20                    | 200                 |         | 100                  | 75       | <input type="checkbox"/> | 24     | 0         |
| Scan                   | Positive |                       |                 |         |            |               | 140                   | 200                 |         | 50                   | 75       | <input type="checkbox"/> |        | 0         |
| Scan                   | Negative |                       |                 |         |            |               | 140                   | 200                 |         | 50                   | 75       | <input type="checkbox"/> |        | 0         |
| MRM                    | Negative | MRM_N                 | 180.1           | Unit    |            | 116.1         |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> | 12     |           |
| MRM                    | Negative | MRM_N                 | 180.1           | Unit    |            | 64.9          |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> | 16     |           |
| MRM                    | Positive | MRM_P                 | 182.1           | Unit    |            | 91            |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> | 8      |           |
| MRM                    | Positive | MRM_P                 | 182.1           | Unit    |            | 81.8          |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> | 24     |           |
| SIM                    | Negative | SIM_N                 |                 |         | 180.1      |               |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> |        |           |
| SIM                    | Positive | SIM_P                 |                 |         | 182.1      |               |                       |                     | Unit    | 20                   | 75       | <input type="checkbox"/> |        |           |

Figure 2. Screenshot of Acquisition Method Showcasing Full Scan, Product Ion, SIM and MRM Scan in +/-ve mode.

### LC-UV-MS Data:

LC-UV analysis of a 450-ppb standard of Impurity A furnished no identified peaks, whereas multiple high intensity peaks were seen for the 5mg/ml Sulbactam formulations indicating that Impurity A is not a detectable chromophore (Fig 3a).

LC-MS analysis of the same Impurity A standard using an Agilent Jet Stream ionization source shows signal in both positive and negative modes: negative ion SIM at m/z 180.1 (Figure 3b) gives a high intensity peak (RT 4.94 min) in the standard and a lower intensity peak (RT 5.36 min) in the formulations while positive ion SIM at m/z 182.1 (Figure 3c) shows an intense peak (RT 4.95 min) in the standard and in formulations (RT 5.36 min).




Figure 3. Data obtained in LC-UV (top), LC-MS SIM -ve ion (middle) and LC-MS SIM +ve ion (bottom) modes.

### LC-MSMS Data:

Ultivo LC/TQ Mixed Scan Mode enables simultaneous acquisition of data in Full Scan, SIM, Product Ion and MRM mode in a single run. Mixed Scan Mode data obtained for a 450-ppb standard (0.01%) of Impurity A are shown in Figure 4.




Figure 4. Mixed Scan Mode (Full Scan, Product Ion, SIM and MRM Scan in +/-ve modes).

### Calibration Data in SIM and MRM modes:

A 7-point calibration curve for 9ppb to 900ppb Impurity A was constructed for positive and negative ionization mode SIM data. An accuracy of  $\pm 20\%$  and  $r^2 \geq 0.995$  was obtained for both ionization modes (Figure 5).




Figure 5. Calibration plot in SIM +/- ve modes.

## Results and Discussion

A 7-point calibration curve for 9ppb to 900ppb Impurity A was also constructed for positive and negative ionization mode MRM data. An accuracy of  $\pm 20\%$  and  $r^2$  of  $\geq 0.995$  was obtained for both ionization modes (Figure 6).




Figure 6. Calibration plot in MRM +/- ve modes.

### Quantitation and Recovery Experiment:

Sulbactam-related compound A (Impurity A) is present at 150-250 ppb in 5mg/ml formulations; this range is less than the allowed limit of 0.01%.

| Batch No | Calc Conc in SIM -ve | Recovery %                       | Calc Conc in SIM +ve | Recovery %                      |
|----------|----------------------|----------------------------------|----------------------|---------------------------------|
| 019      | 223.5                | 100 (659.7-223.5)/450 = 96.93 %  | 241.2                | 100 (699.3-241.2)/450 = 101.8 % |
| 119      | 146.7                | 100 (617.8-146.7)/450 = 104.69 % | 166.0                | 100 (674.5-166.0)/450 = 113 %   |
| 219      | 233.4                | 100 (648-233.4)/450 = 92.13 %    | 267.3                | 100 (716.5-267.3)/450 = 99.82 % |

Table 2. Recovery Calculation for Sulbactam Impurity A in formulations at 0.01% in SIM mode

A recovery experiment, conducted by spiking 450 ng of Impurity A in 5mg/ml formulations, showed recoveries between -10% to +15%. A comparative analysis can be seen for data acquired in SIM mode (Table 2) versus data acquired in MRM mode (Table 3).

| Batch No | Calc Conc in MRM -ve | Recovery %                       | Calc Conc in MRM +ve | Recovery %                       |
|----------|----------------------|----------------------------------|----------------------|----------------------------------|
| 019      | 227.8                | 100 (654.4-227.8)/450 = 94.8 %   | 240.4                | 100 (678.2-240.4)/450 = 97.29 %  |
| 119      | 145.5                | 100 (614.6-145.5)/450 = 104.24 % | 156.2                | 100 (650.7-156.2)/450 = 109.89 % |
| 219      | 229.5                | 100 (648.4-229.5)/450 = 93.09 %  | 254.2                | 100 (697.4-254.2)/450 = 98.49 %  |

Table 3. Recovery Calculation for Sulbactam Impurity A in formulations at 0.01% in MRM mode

## Conclusions

- Simultaneous identification, confirmation and quantitation are established in single injection.
- UV, SIM and MRM scan modes in a single method are useful in impurity profiling workflows.
- Quantitative analysis is enabled by fast polarity switching on an Ultivo LC/TQ without loss of sensitivity.
- Ultivo LC/TQ is the suitable solution for purity and impurity analysis for agricultural and pharma industry.

## References

1. <https://www.agilent.com/cs/library/applications/5991-2796EN%20QC%20Applications%20Compendium.pdf>
2. <https://pubchem.ncbi.nlm.nih.gov/compound/Sulbactam>
3. <https://www.agilent.com/cs/library/brochures/5991-8146en.pdf>