

Determination of Ochratoxin A and Sterigmatocystin via DAD and FLD

Method
HPLC

RP Mode
VFD0142J

Column: Eurospher II 100-3 C18, 150 x 3 mm ID
with integrated precolumn

Order No. 15XE181E2G

Phase: Eurospher II 100-3 C18

Conditions: Eluent: A: Water
B: Acetonitrile gradient grade
Gradient: 0 – 20 min 10 – 90 % B
20 – 25 min 90 % B
25 – 30 min 90 – 10 % B
30 – 40 min 10 % B
Flow rate: 1.0 ml/min
Temperature: 40 °C
Volume: 5 µl (10 µl sample loop)

Detection: Fluorescence RF-20AXs, Ex 329 nm Em 460 nm, 5 Hz, 0.1 s, Sensitivity: High, Gain x16
UV DAD6.1L with High Sensitivity KNAUER LightGuide UV Flow Cell 329 nm, 20 Hz, 0.05 s

Substances: 0.01 mg/ml Ochratoxin A and 0.001 mg/ml Sterigmatocystin

Keywords: OTA, Ochratoxin A, Sterigmatocystin, Mycotoxins, Food analysis, Food safety

Chromatogram:

1. Ochratoxin A
(FLD - blue)
2. Sterigmatocystin
(DAD – green)

