

Improved Detection of Pesticide Residues in Botanicals by LCMS

ASMS 2018

Jeffrey H. Dahl and Tairo Ogura
Shimadzu Scientific Instruments, Columbia Maryland

Improved Detection of Pesticide Residues in Botanicals by LCMS

Introduction

Pesticides and other chemicals are used in the production of many agricultural products, including botanicals for use as dietary supplements. Supplements are widely used but their raw materials, often sourced from remote locations, are subjected to fewer regulatory controls than staple foods. To ensure quality the US FDA requires identity and

quality testing, but most botanicals do not have specific regulations. To analyze complex botanicals for residual chemicals such as pesticides, LC-MS-MS is needed for high sensitivity, high confidence results. We developed an LCMS method with improved detection sensitivity for chemical residues in botanicals.

Photo credits: Echinacea, Giancarlo Dessi; Cayenne, H. Zell; Valerian, Lairich Rig; Ginseng, National Institute of Korean Language; Tumeric, Simon A. Eugster; Passionflower, Bob Peterson; St Johns Wort, Glyn Baker. All photos obtained through wikimedia commons under creative commons attribution-share alike 2.0 or higher



Echinacea E. purpurea



Passionflower (tea) Passiflora sp.



Cayenne Capsicum annuum



Valerian Valeriana officinalis



Tumeric Curcuma Longa



Korean Ginseng Panax ginseng



St. Johns Wort Hypericum perforatum

Figure 1 A selection of popular dietary supplements tested. Various parts of the above-pictured plants may be used in actual dietary supplement formulations.

Improved Detection of Pesticide Residues in Botanicals by LCMS

Method

Representative samples of powdered botanicals were removed from their gelatin capsules, homogenized and extracted with acetonitrile accompanied by shaking and sonication. Samples were additionally cleaned up using dispersive solid phase extraction to remove unwanted matrix components. Analysis was carried out by LC-MS-MS

using a triple quadrupole mass spectrometer. The mass spectrometer interface parameters were carefully adjusted to improve the signal for the majority of the analytes. Spiking experiments were used to determine recovery and matrix-matched standards were used to prepare calibration curves.

Table 1 Instrument parameters used for analysis

LC Column	: Raptor ARC18 (2.1×150 mm, 2.7 μm)
Mobile Phase A	: 0.1% Formic Acid with 5 mM Am. Formate
Mobile Phase B	: Methanol
Flow Rate	: 0.5 mL/min
Probe Voltage	: +0.5 kV or -0.5kV
Interface Temp	: 100 °C
Nebulizing Gas	: 3 L/min
Drying Gas	: 10 L/min
DL Temp	: 100 °C
Heat Block Temp	: 100 °C




Figure 2 Sample preparation of botanicals for LC-MS analysis

Results and Discussion

Eight popular botanical supplements were selected for testing, including Cayenne, Valerian, Passionflower tea, Korean Ginseng, St. John's Wort, Tumeric and two varieties of Echinacea. The Tea and Echinacea variety 1 were labeled as organic, while the other supplements were not labeled as organic. For each sample, a single-point standard addition sample at 500 ng/g dried material was prepared in addition to check matrix-specific effects. Compared with a conventional method, we found significant improvement in instrument response for many

analytes by careful adjustment of interface temperature and spray voltage. For quantitation, matrix matched calibration curves were linear within the quantitation limits established for each compound, which was compound dependent. Detection limits and quantitation limits were required to have 3:1 and 10:1 signal to noise respectively, and quantitation limits were required to have less than 20% RSD in triplicate injections. Using our newly developed method, we are able to characterize the extent of residual pesticides t in popular botanicals.




Figure 3 Representative chromatogram of pesticides spiked into a sample of Korean Ginseng at the 500 ng/g of dried material level.




Figure 4 Acetonitrile extracts (before dSPE cleanup) of various botanicals. From left: Valerian, Cayenne, Echinacea-1, Tumeric, Echinacea-2, St. John's Wort, Korean Ginseng, and Passionflower Tea.

Improved Detection of Pesticide Residues in Botanicals by LCMS

Table 2 Chemical residues detected in botanical extracts. Values reported in ng/g dry material.

	Cayenne	Echinacea-1	Echinacea-2	Korean Ginseng	St. John's Wort	Passionflower tea	Tumeric	Valerian root
Azoxystrobin	15	ND	ND	7	ND	ND	ND	ND
Carbaryl	ND	ND	ND	ND	ND	35	ND	ND
Carbofuran	3.7	ND	ND	ND	ND	ND	ND	ND
Chlorpyrifos	ND	ND	ND	ND	ND	ND	9.4	ND
Cypermethrin	140	ND	ND	ND	ND	ND	ND	ND
Dimethomorph	ND	ND	ND	36	ND	ND	ND	ND
Imidacloprid	11	ND	ND	ND	ND	ND	ND	ND
Metalaxyl	ND	ND	ND	ND	ND	ND	8.3	ND
Methoprene	196	ND	ND	ND	ND	ND	ND	ND
Novaluron	24	ND	ND	ND	ND	ND	ND	ND
Propiconazole	ND	ND	ND	24	ND	ND	ND	ND
Pyraclastrobin	6.5	ND	ND	ND	ND	ND	ND	ND
Tebuconazole	59	ND	ND	4.8	ND	ND	ND	ND
Trifloxystrobin	10	ND	ND	ND	ND	ND	ND	ND

Table 3 List of compounds measured and limits of quantitation in ng/g.

Abamectin	30	Daminozide	15	Hexythiazox	15	Propiconazole	60
Acephate	5	Deltamethrin	30	Imazalil	10	Propoxur	2
Acequinocyl	60	Diazinon	<2	Imidacloprid	4	Pyraclastrobin	10
Acetamiprid	<2	Dichlorvos	15	Kresoxim-methyl	4	Pyrethrin I	100
Aldicarb	<2	Dimethoate	<2	Malathion	2	Pyridaben	2
Allethrin	50	Dimethomorph	5	Metalaxyl	2	Resmethrin	35
Azoxystrobin	4	Dinotefuran	2	Methiocarb	4	Spinetoram	2
Bifenazate	2	Dodemorph	4	Methomyl	<2	Spinosad	<2
Bifenthrin	4	Endosulfan-sulfate	4	Methoprene	50	Spirodiclofen	10
Boscalid	4	Ethoprophos	2	Mevinphos	4	Spiromesifen	20
Buprofezin	<2	Etofenprox	4	MGK-264	500	Spirotetramat	2
Carbaryl	10	Etoxazole	<2	Myclobutanil	10	Spiroxamine	2
Carbofuran	<2	Fenhexamid	20	Naled	2	Tebuconazole	2
Chlorantraniliprole	2	Fenoxycarb	2	Novaluron	15	Tebufenozide	5
Chlorpyrifos	10	Fenpyroximate	10	Oxamyl	2	Teflubenzuron	15
Clofentazine	4	Fensulfothion	5	Paclobutrazol	2	Tetrachlorvinphos	4
Clothianidin	4	Fenthion	10	Permethrin	10	Tetramethrin	4
Coumaphos	4	Fenvalerate	100	Phenothrin	10	Thiacloprid	<2
Cyantraniliprole	2	Fipronil	2	Phosmet	10	Thiamethoxam	<2
Cyfluthrin	500	Fonicamid	25	Piperonyl butoxide	5	Thiophanate-methyl	5
Cypermethrin	60	Fludioxonil	2	Pirimicarb	2	Trifloxystrobin	<2
Cyprodinil	10	Fluopyram	2	Prallethrin	10		

Improved Detection of Pesticide Residues in Botanicals by LCMS




Figure 5 Recovery for each analyte in a 500 ng/g spike of each sample. Tumeric and Cayenne had the greatest number of analytes with low recovery, due to signal suppression.




Figure 6 (Left) Cypermethrin (four isomers) detected in Cayenne sample.
(Right) Carbaryl detected in organic passionflower tea sample.

Improved Detection of Pesticide Residues in Botanicals by LCMS

Conclusion

We developed a high performance method for sensitive detection of pesticides in popular botanical supplements with a simple sample preparation and applied the method to measure pesticides in selected botanical products offered for retail sale. We found matrix effects to be minimal with the exception of Tumeric and Cayenne. For these matrices, additional sample cleanup may be useful to minimize signal suppression by the matrix.

Several pesticides were detected in some of the dietary supplements. Cayenne had the greatest number of detections and with the highest levels approaching 200 ng/g. Significantly, the Passionflower tea, which was labeled as organic, was found to contain 35 ng/g carbaryl. Our rapid, sensitive, and selective method is well-suited to high throughput detection of pesticide residues in popular botanical products.

First Edition: June, 2018



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".

Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".

Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.