



Fight Smuggling and Distribution of Narcotics

Illicit narcotics move through a complex chain of dealers to reach the end user. During this time, the pureproduct is “cut” to make it more profitable. Identification – of both bulk and cut narcotics – is essential in the fight against smuggling and distribution. Ideally, in the fast-paced world of street drugs, identification can occur on-site and in real time, with no need for a laboratory.

PEOPLE
YOU
CAN
TRUST

 **Metrohm**

Metrohm provides the solution for tracing drugs throughout the supply chain. Metrohm products support local, state, and federal agents in safe, fast, on-the-street identification of narcotics and cutting agents.

Bulk Trafficking Samples

Fentanyl and cocaine bricks, white powders, synthetic party drugs i.e., MDMA and Ecstasy



Street Level Drugs

Fentanyl, cocaine, methamphetamine






Trace Level Drugs

Fentanyl, complex mixtures, counterfeit drugs



CONCENTRATION



FIND



IDENTIFY

Measure through plastic wrap, baggies, and most containers

Test from a safe distance of 2 meters away

Measure liquid samples

Extensive narcotics libraries, including fentanyl derivatives

Mixture matching detects and IDs major components in a sample

Use ID Kit to detect trace components in complex mixtures or very dilute samples

INCARCERATE



FOR MORE INFORMATION, PLEASE VISIT: TACTICALRAMAN.COM