



Cannabis Questions? More Answers in One Injection

As laws and regulations change, cannabis and all its related compounds of interest are showing up in a wide variety of consumer products. Instruments that combine powerful Time-of-Flight (TOF) and GCxGC technologies can make your analyses more efficient. Save time and money by quantifying and identifying both the terpene and cannabinoids in just one injection, as shown at right. That same injection can even reveal high concentration pesticides!



How are we able to do this? After GCxGC chromatographically resolves all of these species, TOF allows you to spectrally match the species against the libraries—giving you comprehensive cannabis analysis all from one injection.

Ask us how GCxGC and TOF can work for your analysis!

3000 Lakeview Avenue | St. Joseph, MI 49085 | 800-292-6141 | Phone: 269-985-5496
info@leco.com | www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.

LECO Corporation

LECO

EMPOWERING RESULTS

© 2019 LECO Corporation