

Determination of Melamine

Method VFD0125J
HPLC HILIC mode

Column: Eurospher II 100-3 HILIC, 150 x 3 mm

Order No. 15CE120E2G

Phase: Eurospher II 100-3 HILIC

Conditions: Eluent: A: 10 mM NH₄-Formate pH 4
B: ACN
Gradient: Isocratic 90 % B
Flow rate: 0.8 ml/min
Temperature: 25 °C
Volume: 1 µL

Detection: UVD 2.1L, 210 nm, 20 Hz, 0.05 sec

Substances: Melamine

Keywords: HILIC, AZURA Compact, Eurospher II HILIC

Chromatogram:

1 Melamine

