

Fast Separation of DMSO and NHS

Method
HPLC

Matrix

Column: Eurokat H, 10 µm, 30 x 8 mm

Order No. 03GX340EKN

Phase: Eurokat H, 10 µm

Conditions: Eluent: A: 0.01 N Sulfuric acid
Flow rate: 1.5 ml/min
Temperature: 75 °C
Volume: 10 µl

Detection: RI

Substances: Dimethylsulfoxide, N-Hydroxysuccinimide

Keywords: Organic solvents, Eurokat H

Chromatogram:

- 1 N-Hydroxysuccinimide
- 2 Dimethylsulfoxide

