

Fast Separation of Carbonic acids III

Method VFD0068J
HPLC

Column: Eurokat H, 10 µm, 300 x 8 mm

Order No. 30GX340EKN

Phase: Eurokat H, 10 µm

Conditions: Eluent: A: 0.01 N Sulfuric acid

Flow rate: 1.0 ml/min

Temperature: 75 °C

Volume: 10 µl

Detection: RI

Substances: Oxalic acid, Citric acid, Tataric acid, Fumaric acid, Succinic acid, Adipic acid

Keywords: Carboxylic organic acids, Eurokat H

Chromatogram:

- 1 Oxalic acid
- 2 Citric acid
- 3 Tataric acid
- 4 Fumaric acid
- 5 Succinic acid
- 6 Adipic acid

