

Separation of Dimethylsulfoxide and N-Hydroxysuccinimide

Method VPH0020J
HPLC RP Mode

Column: Eurospher 100-5 C18, 150 x 4.6 mm ID

Order No. 15EE181ESJ

Phase: Eurospher 100-5 C18

Conditions: Eluent: Acetonitrile / Phosphate-Buffer pH 4.0 1:99 (v/v)
Gradient: isocratic
Flow rate: 1.0 ml/min
Temperature: 25 °C
Volume: 10 µl

Detection: UV at 207 nm

Substances: Dimethylsulfoxide (DMSO) , N-Hydroxysuccinimide (NHS)

Keywords: DMSO , NHS

Chromatogram:

11. NHS
12. DMSO

