


Tetrachlorodibenzofuran Congeners on Rtx®-Dioxin2



GC_EV00949

Column	Rtx®-Dioxin2, 60 m, 0.25 mm ID, 0.25 µm (cat.# 10758)
Sample	WMS-01 Reference Material, Wellington Laboratories
Injection	splitless
Inj. Temp.:	250 °C
Oven	
Oven Temp:	130 °C (hold 1.0 min.) to 200 °C at 40 °C/min. to 235 °C at 3.0 °C/min. to 300 °C at 5 °C/min. (hold 10 min.)
Carrier Gas	He, constant flow
Flow Rate:	1.5 mL/min.
Detector	Micromass Ultima high-resolution mass spectrometer
Ionization Mode:	El
Notes	Detector mode: SIR
Acknowledgement	Chromatograms courtesy of Terry Kolic, Karen MacPherson, Eric Reiner, Ontario Ministry of Environment, Toronto, Ontario, Canada